

# Usability of Amide and C28 Core Shell and Fully Porous Column for Separation of Hydrophilic Compounds

Tomoyasu Tsukamoto<sup>1)</sup>, Norikazu Nagae<sup>1)</sup>, Anders Grahn<sup>2)</sup>, Ralf Jutvik<sup>2)</sup> and Scott Silver<sup>3)</sup>



- 1) ChromaNik Technologies Inc., 6-3-1 Namiyoko, Minato-ku, Osaka, Japan
- 2) Biotech AB, Box 133, 439 23, Onsala, Sweden
- 3) Innovations United, 300 East 57th Street, Suite 11J New York, NY 10022 USA

## Abstract

Hydrophilic Interaction liquid Chromatography (HILIC) proposed by Alpert in 1990 has been applied for analysis of hydrophilic compounds. Amide, diol, polyol, bare silica, ion exchange and zwitter ion phases have been used as a hydrophilic stationary phase along with an organic solvent rich mobile phase. A polar group embedded C18 or a long alkyl chain phase such as C30 or C28 also have been used to separate hydrophilic compounds without change in retention using an aqueous mobile phase on a reversed-phase mode. The reason why these reversed-phases can showed no change in retention under an aqueous condition is that a low contact angle of water on the surface of the pore of these reversed-phase packing materials makes an aqueous mobile phase keep in the pore because pressure yielded by capillarity is less atmospheric pressure, so that retention doesn't change. Both HILIC stationary phases and reversed-phases have completely opposite characteristics each other. Therefore both HILIC and reversed-phase modes are useful for separation of hydrophilic compounds. It is important to understand separation behavior of each mode. In this study, an amide column and a C28 column were compared and evaluated to separate hydrophilic compounds. SunShell HILIC-Amide and Sunniest RP-AQUA (C28) and SunShell RP-AQUA (C28) were used to separate nucleobases, amino acids and hydrophilic vitamins. When nucleobases were separated on HILIC and reversed-phase modes using an amide column and a C28 column, each elution order of samples is said to be opposite. Only uracil, however, showed a specific elution. It was considered that the polarity of uracil under an organic solvent rich condition was different from that on water rich condition to be due to keto-enol tautomerization. LC/MS analysis of amino acids was achieved using C28 column and a mobile phase added 5 mM heptafluorobutyric acid under gradient elution conditions.

## Collapse or Depermeating

**C18 phases exhibit decreased and poorly reproducible retention under more than 98% aqueous conditions. This problem traditionally has been explained as being the result of ligand collapse. Nagae<sup>1-3</sup> ascertained, however, that the mobile phase was being expelled from the pores of the packing material.**

**When the surface of packing materials isn't wet by water, water used as a mobile phase expels from the pore of the packing material by capillarity. This is a reason why reproducibility in retention is low under 100% aqueous conditions. Reversely pressure around the pore of the packing material to overcome a force worked by capillarity.**

- 1) N. Nagae, T. Enami and S. Doshi, LC/GC North America October 2002.
- 2) T. Enami and N. Nagae, American Laboratory October 2004.
- 3) T. Enami and N. Nagae, BUNSEKI KAGAKU, 53 (2004) 1309.

What does "Dewetting" mean?

A surface state changes from wetting to un-wetting? The surface of C18 is always un-wetting even if water exists in the pore, so that expression of "dewetting" is not scientific. Depermeating is a scientific expression!




Figure 6. Schematic diagram of C18 particle

Formula of Capillarity:  
 $h = \frac{2\gamma \cos\theta}{(\rho g)}$   
 $\gamma$ : Surface tension  
 $\rho$ : Density of liquid





Figure 7. Schematic of capillarity


## Amide vs C28



Comparing between amide and C28 column ...  
 • Each elution order of samples was almost opposite.  
 • Only elution order of uracil did not change

Figure 1. Comparison between reversed phase and HILIC

## Amide vs other HILIC phase



Comparing between SunShell HILIC-Amide and other HILIC columns ...  
 • SunShell HILIC-Amide showed stronger retention than others.  
 • Amide showed different selectivity comparing with bare silica HILIC column.

Figure 2. Comparison of three kinds of HILIC column

## Applications of HILIC




Figure 3. Separation of melamine and cyanuric acid




Figure 4. Separation of glycosides




Figure 5. Separation of synthetic sweeteners

## Conclusions


- Nucleobases were separated using an amide column and a C28 column, each elution order of samples is said to be opposite.
- Only uracil showed a specific elution. It was considered that the polarity of uracil under an organic solvent rich condition was different from that on water rich condition to be due to keto-enol tautomerization.
- LC/MS analysis of amino acids was achieved using C28 column and a mobile phase added 5 mM heptafluorobutyric acid under gradient elution conditions.
- Both amide and C28 column were useful for analysis of hydrophilic compounds

## Repellency and Hydrophobicity

### Repellency

#### Water-shedding property


- Repellency is expressed as a contact angle of water on a material.
- The larger a contact angle, the stronger repellency, if the contact angle is more than 90 degree.



### Hydrophobicity

#### Difficult to mixing with water

- Hydrophobicity is expressed as the ratio of concentrations of a compound between water and n-octanol using a mixture of both solvents.
- This value is well known as LogPow.



### Repellency $\propto$ Hydrophobicity

- Repellency and hydrophobicity are independent each other.
- Those two parameters are out of proportional.
- When hydrophobicity is high, it doesn't mean that repellency is always high.
- Capillarity depends on a contact degree.



- C28 has the smallest contact degree comparing with C18 and C8.
- In C28, the pressure for a mobile phase to go out from pore is smaller than atmospheric pressure.
- Aqueous mobile phase isn't expelled from the pore of C28 phase.

Table 4. Physical property of each compounds

	Trifluoromethane	Octadecane (C18)	Octane (C8)	Octacosane (C28)
Contact angle( $\theta$ )	120°	126°	140°	108°
Partition coefficient (LogP)	0.64	9.18	5.18	14.09
Solubility(mg/L)	4090	0.006	0.66	$8.84 \times 10^{-10}$

T. Enami and N. Nagae, BUNSEKI KAGAKU, 53 (2004) 1309.

## Stability of C28




Figure 8. Separation of nucleobases and retention time stability




Figure 9. Stability of C28 under acidic and basic conditions

## Separation of amino acids with C28




Figure 10. Separation of amino acids using UV and MS detections

## Applications of C28




Figure 11. Separation of organic acids