

iSee370

EYE MASSAGER

USER MANUAL

OBM Distribution / Breo-USA
13900 Sycamore Way
Chino, CA 91710
<http://www.breousa.com>
info@breousa.com

breo[®]

**Thank you for purchasing
BRE0 iSee370.**

Please read this manual carefully before using
the product, and keep it in a safe place for future
reference.

Contents

Introduction	3
Inside the Package	3
Function Deployment	4
Preparation	5
Operation	6
Precautions	7-8
Troubleshooting & Technical Parameters	9


Introduction

iSee370 adopts air pressure and hot compress massaging technology. It rubs, presses, hot compresses, and massages the ocular region in order to relax it. You can enjoy the soothing music while massaging to completely relax your brain.

Inside the Package


The package includes: Main Unit, Power Charger, Batteries, User Manual and Earphones.

Massager



Master Unit

Basic components



Power Charger




Earphones




User Manual

Function Deployment

Main Unit



Controller



- 1. Elastic band Button stand
- 2. LED light (orange)
- 3. Fabric covering
- 4. On / Off button
- 5. Mode Select Button
- 6. Quintuplet button
- 7. Pressure Reduction Valve


Preparation

Charging


1. Please use power adapter provided by our company only.
2. Open the USB cover on the Controller.
3. Insert the plug of the power charger to the DC port.
4. During charging, the light of the power charger is red. The light turns green when fully charged.

Wear

1. Please hold the fabric covering to the position directed to the eyes and put on the host machine.
2. Adjusted the elastic band button to a comfortable position



Front View



Side View

Operation

Starting Up Default

circulate air pressure + low temperature hot compress, timing 5 minutes, piano music is defaulted as the mode of music.

If you need to change the mode of music, time, vibration, hot compress and air pressure, please manually operate related function keys.

Manual mode

ON/OFF 

Switch — Press  Button — starting up (default mode)

Press and hold  Button — shutdown

Mode selection 

press  Button to switch among the following modes:

press  Button — Air pressure + Hot press + music / air pressure + music / hot press + music

Air pressure+ hot compress+ music / air pressure+ music / hot compress+ music

Air pressure

Open Air pressure — Press  Button — "air" character flashes to switch among the following modes:

Hot compress

Open hot compress — Press  Button, "heat" character flashes to switch between the following modes:

Low temperature/high temperature

Time

press  Button — Switching among 5'/10'/15'

Music

press  Button — music plays in order: Piano music / running water / warble / sea wave / cricket / water-drop

Press and hold  Button — turn off music

Precautions

1. Please read the manual carefully and operate accordingly.
2. Children under 7 years old cannot use the massager.
3. This massager should not be used by people who suffer from brain disease, eye disease, or another physical or mental disorder.
4. Take off contact lenses / eye glasses before use.
5. Put on / wear the unit first BEFORE turning the power on.
6. Do not use while walking or driving.
7. Keep your eyes closed for a few seconds after use, then open them slowly.
8. Keep the massager away from humidity, heat, and/or corrosive gas.
9. If the massager will not be used for a long time, remove the batteries to avoid battery leakage.
10. This massager is designed as a personal care product. It is not for diagnosis/treatment of any health disorder, nor meant for commercial use.
11. To clean, wipe with aseptic cotton and alcohol. Do not wash or immerse in water.
12. This product is not a substitute for healthcare professionals' advice and treatment. In case of doubt or if necessary, please consult your healthcare professional.

13. This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.
14. This massager features heat compression. Individuals who are incapable of or do not ordinarily feel heat/cold sensation should use this device with caution.
15. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:
 - Reorient or relocate the receiving antenna.
 - Increase the separation between the equipment and receiver.
 - Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
 - Consult the dealer or an experienced radio/TV technician for help

Troubleshooting & Technical Parameters

When operational problems occur, please refer to the chart below.
Do not disassemble the machine or try to repair it by yourself.

Problem	Possible Reason	Possible Solution
Power-ON failure	Cable is broken / in bad condition	Check the condition of the cable
	Batteries are low / out	Change batteries
	Batteries are not installed correctly	Check if positive/negative poles of batteries are positioned right Make sure the batteries are snug in place
	AC adapter is not connected properly	Check the connection of the AC adapter to the power socket .
Weak heat or air pressure	Batteries are low / out	Change batteries Make sure the batteries are installed correctly and snug in place

Technical Parameters

Power	DC: 6V
Ideal environment	Temperature: 5°C-40°C
	Relative Humidity: ≤80%
Dimension	Main Unit: 206mm ¼108mm ¼86mm
	Controller: 141mm × 55.5mm × 33mm
Net Weight	458g
Power	6W
Magnetic Strength	50~250 mT
Formentation	30~45 °C