ENGLISH

1 Loosen Straps

- a.** Release the two over-center buckles by pinching the two gray buttons on each of the buckles.
b. Loosen the two remaining straps.

2 Position the Brace

With the patient in a seated position and leg bent at 90°, place the brace on the leg aligning the center of the hinge with the top of the knee cap. **Note:** Make sure to push the brace posterior so that the hinges run along the mid-line of the leg.

3 Secure the Top Calf Strap

Wrap the top calf strap around the calf as high as possible and loop through the D-ring. Pull the tab until a snug and comfortable fit is achieved and secure.

4 Secure the Lower Thigh Strap

Wrap the lower thigh strap around the thigh and loop through the D-ring. Pull tab until a snug and comfortable fit is achieved and secure.

5 Adjust the Crescent Finger Pad

Remove the Crescent Finger Pad from underneath the Y-Pad. Relocate the pad to the desired position along the edge of the patient's knee cap.

6 Secure Lower Patello-Force Strap

- a.** With the patient's knee at 20° of flexion, wrap the lower Patello-Force Strap around the leg and secure the over-center buckle.
b. Adjust the strap using the fastener at the strap end until a snug and comfortable fit is achieved.

7 Secure Upper Patello-Force Strap

Repeat Step 6 for the upper Patello-Force Strap.

8 Adjust the Patello-Force Strap

Use the fasteners for additional adjustment.

9 Inflate the Patello-Force Strap

Insert the clear end of the squeeze bulb into the pad. Slowly inflate the pad by squeezing the bulb. The patient will know how much pressure is needed. **Note:** To deflate the pad use the black end of the squeeze bulb.

10 Final Adjustments

Have the patient stand and walk around. Make any final strap and/or air pressure adjustments needed.

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