


Z-12 Knee Brace
Z-12-D Knee Brace
Z-12 OA Knee Brace
Z-12 Adjustable OA Knee Brace*Application Instructions CP020179 Rev E 09/11***A.****B.****C.****ENGLISH****1. Unfasten Buckles**

Release all four straps on the back side of the brace by squeezing the gray levers on either side of the buckle and withdrawing the over-center buckle. Loosen the front strap(s).

2. Position the Brace

With the knee bent at 45°, lay the brace on the leg and align the hinge as shown above. Position the center of the hinge on an imaginary line that goes across the top of the knee cap. Push the hinges posteriorly until the hinges are in the center of the knee joint when viewed from the side. (A)

3. Anchor the Brace

Wrap the upper rear calf strap around the calf as high as possible over the calf muscle.

Insert the over-center tab into the base. Push the tab until you hear the hinge click on both sides. (B)

Adjust the strap by releasing the tab and pull the strap until a snug and comfortable fit is achieved. All straps have additional length adjustment at the opposite end of each strap. See step 6. (C)

4. Fasten Lower Calf Strap

Insert the over-center tab into the buckle base on the lower calf shell of the brace. Push the tab in until both sides click locked.

Adjust strap until a comfortable fit is achieved. (If the strap tab reaches its limit without achieving the desired fit, see step 6 for adjustment instructions.)

Center the strap pad.

5. Adjust Strap Length

All straps have additional length adjustment at the end opposite the tab. Release the strap and adjust the length of the strap using the hook and pile fastener. (D)

6. Fasten Rear Thigh Straps

Following the same instructions as in Step 5, fasten the lower and then upper rear thigh strap.

7. Adjust Front Tibia Strap

Tighten the front tibia strap. If the strap overlaps too far, see step 6 for adjustment instructions.

Center the strap pad. (E)

8. Adjust Front Thigh Strap (PCL/Combined Instability Versions)



Tighten the front thigh strap. Note: Overtightening the front thigh strap can cause the brace not to function properly.

Stand up and walk or jog around to see how your new AXIOM or AXIOM-D brace feels. Straps should not be so tight that blood circulation is restricted. Re-adjust as needed.

If this is your first time to wear a functional knee brace, it might take a few hours of wear to get used to how it feels. (F)

9. Adjust OA Condyle Pad

Press down on grey tab and slide into position for desired unloading. Adjust for comfort. (G)

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