



MedicShield

Hospital Grade Surface Spray Disinfectant/Sanitiser

Ready to use



MedicShield is a hard surface cleaner and sanitiser that is proven to be effective against a range of bacteria, fungi and viruses that can cause disease especially in those that are in a weakened state. MedicShield has application in Medical, Dental, Aged Care and Laboratory.

MedicShield is Suitable for the following:

In Intensive Care Units - Cleaning and disinfection of all hard surfaces.

In Operating Theatres - Cleaning and disinfecting tables disinfect equipment, cables and will not adversely affect surfaces.

For Blood spills - MedicShield is ideal for cleaning and disinfecting blood spills on all hard surfaces.

Telephones - MedicShield is ideal for disinfecting phones where large numbers of bacteria are found.

In Dental and GP Surgeries - MedicShield can be used in Dental and GP surgeries for cleaning and disinfecting surfaces including chairs, stools, tables benches and light fittings.

In Laboratories & Pathology Clinics - MedicShield is excellent for wiping down benches, disinfection of white rooms, cleaning and disinfection of fume cupboards, etc.

In the Nursery - Cleaning cribs, cots and change tables.

CSSD - Cleaning bench tops, stainless steel surfaces, sinks, trolleys, tap fittings.

In Toilets Showers and Sinks - Tile Surfaces, and all hard surfaces.

Directions

- 1) Remove all equipment, and any gross soiling by manual scraping and discard hygienically.
- 2) Shake well. Hold container upright 15cm to 20cm from Surface, Allow to contact the surface for at least five minutes.
- 3) Wipe out the area or rinse clean the area to remove the sanitising solution
- 4) Leave the area to dry thoroughly or dry with a clean towel.

Integra Industries

North Island: 42 Riverbank Rd, Otaki • South Island: 21 Glasgow St, Dunedin
Head Office: 149 King Edward St, South Dunedin. Phone: 0800 667 843

www.geller.co.nz

Proven effective on the following list of pathogens

Virucidal Efficacy

- Adenovirus type 4
- HIV-1 (AIDS Virus)
- Hepatitis B Virus (HBV)
- Hepatitis C Virus (HCV)
- Herpes Simplex 1 & 2
- Human Coronavirus
- Influenza A / Hong Kong
- Norovirus
- Respiratory Syncytial Virus (RSV)
- Rubella (German Measles)
- Vaccinia (POX Virus)

Fungicidal Efficacy

- Aspergillus niger (Mold & Mildew Claim)
- Candida albicans
- Trichophyton mentagrophytes (Athlete's Foot Fungus)

Bactericidal Efficacy

- Acinetobacter calcoaceticus
- Bordetella bronchiseptica
- Chlamydia psittaci
- Enterobacter aerogenes
- Enterobacter cloacae
- Enterococcus faecalis (Vancomycin Resistant)
- Escherichia coli (E.coli)
- Fusobacterium necrophorum
- Klebsiella pneumoniae
- Legionella pneumophila
- Listeria monocytogenes
- Pasteurella multocida
- Proteus mirabilis
- Proteus vulgaris
- Pseudomonas aeruginosa
- Staphylococcus aureus
- Staphylococcus aureus (MRSA)
- Staphylococcus aureus (VISA)
- Salmonella choleraesuis
- Salmonella enteritidis
- Salmonella typhi
- Salmonella typhimurium
- Serratia marcescens
- Shigella flexneri
- Shigella sonnei
- Staphylococcus epidermidis
- Streptococcus faecalis
- Streptococcus pyogenes (cause of Scarlet Fever)
- Xanthomonas axonopodis pv. (Citrus canker)

Available sizes with codes



750ml
MEDICHGSS750TS



5L
MEDICSHIELD05RTU