


# T-cell and B-cell receptor repertoire profiling for biomarker discovery Alex Chenchik, Mikhail Makhanov, Tianbing Liu, Dongfang Hu, Khadija Ghias, Paul Diehl, Lester Kobzik; Cellecta, Inc., Mountain View, CA

# Abstract


T-cell receptor (TCR) and B-cell receptor (BCR) repertoire profiling, also referred to as adaptive immune receptor repertoire (AIRR) profiling, holds great potential for the understanding of disease mechanisms and for the development of new treatments in infectious disease, autoimmunity, and immuno-oncology. This potential could be greatly improved by combining information about receptor clonotypes with immunophenotypes of T- and B-cells. We developed a new technology for combined profiling of all human TCR and BCR variable regions with phenotypic characterization of immune cells in the same workflow. The TCR and BCR immunophenotyping method proposed involves RT-PCR amplification and sequencing of the CDR3 regions of the TCR and BCR genes, as well as determining the gene expression levels of the most informative T- and B-cell phenotyping genes. Results show that this method allows for comprehensive profiling of all seven TCR and BCR chains from a single sample, in a highly reproducible manner, directly from microsamples including cancer tissue, whole blood, sorted cells and more. Bioinformatic analysis of the next-generation sequencing (NGS) data allows profiling of the TCR and BCR clonotypes and identification of major immune cell subtypes and their activation status. Preliminary data from rheumatoid arthritis blood samples and others will be presented.

# Introduction




- Genetic recombination in Tand B-cells generates diverse repertoires of TCR, BCR, and antibodies.
- Variable part (CDR3) of TCR and BCR recognizes foreign antigens presented by Major Histocompatibility Complex (MHC).
- Millions of different T- and B-cells with unique TCRs and BCRs define differences in our immune responses.
- Understanding the complex TCR-BCR repertoire can provide insights into disease mechanisms and suggest strategies for effective immunotherapies.

# Method



- Universal assay for targeted RNA-seq profiling of all TCR/BCR and key T/B biomarker
- genes Single-cell sensitivity, 10fold increase in sensitivity versus RNA-seq and SMART technology (Takara Bio)
- Could be run directly on cell lysate (single cells, sorted cells)
- Doesn't require rRNA, mitochondrial, or globin RNA depletion when starting with total RNA


### Full-length or CDR3 receptor coverage



- Comprehensive AIR repertoire coverage for all seven TCR and BCR chains in a single multiplex RT-PCR reaction from total RNA or DNA
- Detection of only functional AIR clonotypes without pseudogenes or nonrearranged genes
- Improved coverage and unbiased amplification of full-length or CDR3 regions with a highly validated primer set based on DriverMap™ technology
- Quantitative clonotype analysis with AIR RNA or DNA calibration standards and UMI (unique molecular identifiers)
- Integrated with MiXCR software package (MiLaboratories) for immune repertoire analysis


#### Integrated AIR profiling and immunophenotyping

Integrated AIR profiling and Immunophenotyping



- **Fig 1:** DriverMap<sup>™</sup> AIR and DriverMap™ IMP assays together allow in-depth characterization of TCR repertoire in CD8 naive, CD8 effector, CD4 naive, CD4 effector, and regular T-cell fractions.
- Integrated AIR profiling and Immunophenotyping directly in sorted cells without RNA purification is possible using the DriverMap™ technology.
- High-resolution immunophenotyping (matching) data from top TCR/BCR clonotypes based on the expression of 300 key cell typing and activation T/B markers.
- Candidates selected from a set of 3000 candidate genes described in >100 public databases, commercial assays, and publications.

#### View all 7 TCR/BCR chains in a single reaction




in 50ng of normal PBMC or whole blood RNA (10x10<sup>6</sup> reads per sample, triplicates).

Fig 2: Number of

clonotypes for 7 TCR/

BCR chains identified



- Fig 3 shows high reproducibility of TRB repertoire profiling for top 500-1,000 clonotypes in RNA samples with at least 5-10 TRB mRNA molecules
- Stochastic, unreproducible profiling of rare clonotypes (100,000s) present in whole/PBMC RNA samples at the single-molecule
- Fig 3: TRB clonotype repertoire analysis in 50 ng of whole blood in triplicate.

#### Reproducible clonotype repertoire analysis




Fig 4: Reproducible clonotype distribution of TCR and BCR chains across different input amounts ranging from 1.5 ng to 100 ng of starting total RNA of all seven chains in a single reaction. (5 x 10<sup>6</sup> reads/sample) The recommended optimal starting amount is 50-100 ng of total RNA.

Fig 5: Similar V gene usage for TRB genes in whole blood, whole blood micro-samples (30 ul dried blood), and PBMC samples. log<sup>10</sup> (V gene usage percentage) in triplicates.

# Sensitive clonotype detection





Fig 6: Comparison of TCR clonotypes detected by DriverMap™ AIR RNA vs SMART assay. Both assays were run with 50 ng of total RNA isolated from PBMCs. The DriverMap™ AIR RNA assay detects ~ 3X more TCR clonotypes than the SMART assay. (Barennes et al., 2020)

# Identify cancer-activated CDR3 clones





- BCR > up-to 1,000-fold activation
- TCR > up to 50-fold activation

Fig 7: Detection of Cancer-Activated CDR3 Clones in mRNA based DriverMap™ AIR assay normalized to gDNA based ImmunoSeq assay (Adaptive Biotechnologies)—in metastatic tumor samples. Green brackets indicate clonotypes highly expressed in RNA but not in DNA, indicating highly upregulated receptor clones in RNA vs. DNA. However, some DNA regions are also upregulated, indicating a higher copy number of cell (i.e., DNA) clonotypes in the tumor samples. (Sudmeier et al., 2021)

# Discussion

- Adaptive Immune Repertoire (AIR) Profiling assay: Quantitative, and comprehensive TCR/BCR repertoire analysis (all seven chains) in one or two multiplex RT-PCR reaction in bulk RNA or DNA samples (PBMC, whole blood, cancer tissue samples).
- **Direct AIR Profiling:** High sensitivity with minimum background detection of TCR/BCR clonotypes directly in micro samples (cancer tissue, whole blood), sorted cells, and single cells using DriverMap™ technology.
- T/B Immunophenotyping: Integrated analysis of top TCR/BCR clonotypes and expression profiling of cell typing, activation markers in sorted T and B cell subfractions, and single cells.

#### Cellecta offers AIR RNA or DNA assays as kits and custom services

#### **References:**

• Barennes, P., Quiniou, V., Shugay, M., Egorov, E., Davydov, A., & Chudakov, D. et al. (2020). Benchmarking of T cell receptor repertoire profiling methods reveals large systematic biases. Nature Biotechnology, 39(2), 236-245. doi: 10.1038/s41587-020-0656-3

• Sudmeier, L., Hoang, K., Nduom, E., Wieland, A., Neill, S., & Schniederjan, M. et al. (2021). The CD8<sup>+</sup> T cell landscape of human brain metastases. doi: 10.1101/2021.08.03.455000