



# Two Vector CRISPR/Cas9 System

CELLECTA


The CRISPR/Cas9 system can be used for knocking out gene expression by using a combination of sgRNA and the Cas9 nuclease. The Two Vector CRISPR/Cas9 system allows for faster knockout of the target gene in cells by first selecting cells expressing high levels of Cas9. The Two Vector CRISPR/Cas9 system is also great for creating custom sgRNA libraries for screening assays.

## Expression of Cas9 and sgRNA from Separate Vectors to Optimize Knockouts


- Obtain higher titers for Cas9 for difficult-to-transduce cells
- Reduce noise in screening experiments by standardizing expression of Cas9
- Cas9 and sgRNA vectors contain different antibiotic resistance genes for easy selection
- Tet-inducible and constitutive sgRNA expression available

## Knock Out Targets with Lentiviral sgRNA Constructs


- You provide the RefSeq number or gene ID
- We design 3-5 sgRNA constructs that target the transcript



pRSG16-U6-sg-UbiC-TagRFP-2A-Puro  
8.0 kb




pR-CMV-Cas9-2A-Hygro  
11.6 kb



When transducing many thousands of constructs into large populations of cells, as is required for loss of function screens with pooled sgRNA libraries, high titers are a necessity. The 2-Vector CRISPR system provides higher titers than the 1-vector system and allows selection of cells expressing a high level of Cas9 before introducing the sgRNA, which leads to more efficient knockouts with less variability for screens.

### High Expression of Cas9 Increases Knockout Rate



Cells expressing GFP were transduced with Cas9 at high MOI producing a population with approximately 3 Cas9 per cell on average, or a low MOI generating a population of approximately 1 per cell on average. After selecting the Cas9 transductants with hygromycin, each population of cells was then transduced with the same sgRNA to GFP and grown for 9 more days in media containing puromycin. Cells with a higher number of integrated copies of Cas9 have 5-fold fewer GFP-positive cells.

## Related CRISPR Products & Services

### Human CRISPR Genome-Wide sgRNA Library

- Built using solid support oligonucleotide synthesis method
- Provided as either plasmid or packaged pooled sgRNA library

### Custom CRISPR / sgRNA Libraries

- You decide which genes you want to knock out
- Cellecta designs 3-5 sgRNAs per gene (or to your specifications)

### Create Knockout Cell Lines with CRISPR

- Design and cloning of sgRNA constructs
- Viral packaging for sgRNA constructs and transduction into desired cell type
- Confirmation by PCR that both alleles are out-of-frame
- Functional assays such as proliferation, viability screens, or pathway activation assays available

