
**LOOKEE® TENS UNIT EMS
MESSAGE MUSCLE STIMULATOR**
Model : LK113

Instruction Manual

Product is subjected to the real object
Thanks for using our product.
Please keep this manual safe for future reference.
Please read this manual carefully before using it.

Precaution:


Please read the instructions carefully before use. If you have any questions, please contact the retailer or consult your doctor.

1. Product Profile

Working principle

It's through electricity pulse wave after acting on the human body, the electric energy is sucked by the human body and the energy form is changed, which causes a series of physical and chemical changes, produces local or systemic physiological reactions, and finally affects the pathological process and plays a role Corresponding Role.

List of accessories

- Please use a 5V/1A adapter (input: AC100-240V/50/60Hz /0.2A). Adapters that do not meet the requirements may reduce the safety of the product. Please make sure that the selected adapters meet the relevant safety medical level safety certification requirements (example: IEC60601-1).
- After the connection, it will be charged automatically.
- Note:**
- This product does not come with an adapter, please purchase separately.
- If you have any problems with the adapters, consult your dealer at the place of purchase.

2. Intended use


For TENS mode: To be used for temporary relief of pain associated with sore and aching muscles due to strain from exercise or normal household work activities.
For EMS mode: To be used to stimulate healthy muscles in order to improve and facilitate muscle performance.
For MASSAGE mode: To be used for recovery of fatigue, improves blood circulation neuralgia.

3. Contraindications

- Acute disease, especially in patients suspected or suffering from hypertension, thrombosis or malignant mass.
- Patients with various skin diseases.
- Patients with chronic pain disorders (e.g. metabolic disorders) with reduced pain.
- Patients undergoing medication.


4. Product structure description

4.1 Low frequency therapeutic apparatus

1 AB channel electrode connection port	2 CD channel electrode connection port
3 Type-C charging interface	4 Power/Set Time
5 AB channel intensity adjustment button	6 CD channel intensity adjustment button
7 Program	8 Mode

4.2 Display

4.3 Warnings and Notes

Warning

We recommend not to use this therapeutic apparatus in the following cases:

- 1) Do not use this device if the patient has any implantable electronic or metal equipment. (e.g. pacemaker).
- 2) Has epilepsy.
- 3) Severe lower extremity arterial circulation problems.
- 4) Abdominal or inguinal hernia site.
- 5) Do not use this device if you have a heart disease without consulting a doctor.
- 6) During pregnant.
- 7) Has cancer.
- 8) Be Fever.
- 9) Do not use therapeutic instruments in high humidity environments, such as bathing or showering in a bathroom bathtub.
- 10) Can not be used after drinking.
- 11) Low-frequency therapeutic instruments can not be used when using high-frequency surgical instruments.
- 12) Near the genitals.
- 13) Recommendation: do not use low-frequency therapeutic instruments near shortwave or microwave therapeutic equipment. (e.g .1 M), where the output of the low frequency therapeutic instrument used is unstable.
- 14) The high frequency surgical equipment and stimulator. When connected to a patient, burns may occur at the stimulator electrode and may damage the stimulator.
- 15) Patients undergoing medical treatment and feeling abnormal.
- 16) Heart patients with neurological abnormalities in the brain.
- 17) Patients whose body temperature exceeds 38°C degrees.
- 18) Infectious diseases Patients.
- 19) Patients who can not express their will.
- 20) Do not use this product while charging.

Note:

- 1) Make sure you are always aware of all warning items.
- 2) Stimulators should not be used in the throat and mouth. This may result in severe throat muscle spasms and muscle contraction, which is sufficient to block the airway and cause dyspnea.
- 3) Stimulator should not be used in front of the chest, the current introduced into the heart may cause arrhythmia To cause the heartbeat to stop.
- 4) This stimulator should not treat unattended patients, or suffer from emotional distress, dementia, and low IQ Personnel use.
- 5) Use guidance has been detailed; any improper use may lead to danger.
- 6) Do not use this device for any undiagnosed pain symptoms, consult a doctor for diagnosis before use.
- 7) Electrodes can not be placed in areas that may cause current to pass through the brain. (i.e. electrodes. One-sided can not be placed on the front and back of the head at the same time, or relative to the left and right sides of the head).
- 8) Patients suffer from heart disease, cancer or other health-threatening conditions, they should not consult a doctor in advance use this device

under the condition of birth.

- 9) If the intensity of stimulation makes you feel uncomfortable, reduce the intensity of stimulation to a comfortable level if the problem remains. If so, please contact your doctor.
- 10) Do not use this equipment in rooms where aerosol (spray) and pure oxygen are used.
- 11) Please do not use this equipment near any highly flammable substances, gases or explosives.
- 12) Please check that the cable and related connections are in good condition before each use.
- 13) Turn off the device before using or removing the electrode pad.
- 14) Electronic stimulators can only use electrode wires and electrodes recommended by manufacturers.

Consult your doctor before using the device, especially if:

- Acute disease, especially in patients suspected or suffering from hypertension, thrombosis or malignant mass.
- Patients with various skin diseases.
- Patients with chronic pain disorders (e.g. metabolic disorders) with reduced pain.
- Patients undergoing medication.

Adverse reactions

- Skin allergies caused by electrode sheathpad gels and electrode irritation burns are potential adverse reactions. If the skin appears sensitive response, should be discontinued and consult your doctor.
- If the intensity of stimulation makes you feel uncomfortable, reduce the intensity of stimulation to a comfortable level, and if the problem persists, please contact your doctor.


5. Product description

5.1 Instructions

Please make sure the battery is sufficient before use, if the battery is insufficient, please charge it in time.

Do not use this product while charging.

- 1) The electrode pad is inserted into the electrode line.
- 2) Insert the electrode plug into the treatment instrument.
- 3) Attach electrodes pad to the site to be treated.
- 4) Press the Power button and select the program (Refer to 7).

Attention

About electrode

- Don't use it by overlapping.
- Do not bend the electrode pad.
- Please try not to touch the adhesive surface with your fingers.
- After bathing or sweating after exercise, please wipe the parts to be treated before use.
- Don't pull the electrode wire. This will damage the wires and electrode pad.
- Do not apply to damaged skin.
- If the electrode pad is not sticky, please replace it in time.
- The electrode pad is only for one patient.
- If allergy is found, stop using immediately and consult your doctor.
- Please read the instruction of self-adhesive electrode before use.
- When the electrode pad is stuck to the clothes, don't tear it off heavily. As long as the inside of the side with the electrode pad is covered with water, it can be removed easily.

5.4 TENS/MASS/EMS

TENS

Percutaneous nerve electrical stimulation, meaning electrical stimulation of the nerve through the skin. It's a non-drug A method of treating pain caused by certain causes. TENS electrical stimulation from the skin by small electrical pulses skin is sent to the nerve to modify your perception of pain. Electrical stimulation doesn't cure any physical problems, it's just help control pain. Electrical stimulation doesn't work for everyone, but for most patients it's reducing or eliminating in addition to pain is effective, so that patients can return to normal activity.

MASS


Suitable for relaxing body, relieving fatigue, reducing stress and improving sleep.

EMS

Muscle electrical stimulation It has been widely used in sports and rehabilitation for many years and is widely used approved methods. Among the EMS uses in sports and fitness is to complement traditional muscle training, increase the effectiveness of muscle groups and adapt them to the body ratio to achieve the desired aesthetic effect. EMS applications two aspects of the effect. On the one hand to achieve the effect of muscle tissue enhancement (activation application); on the other hand to achieve the effect of relaxation and rest (relaxation application).

5.5 Paste the electrode pad to the area to be treated

- 1) Shoulder: The shoulder may be hard because of moving shoulder, or because the left shoulder is under heavy pressure, there will also be left shoulder pain, such as sometimes with the left shoulder to lift the heavier things. If take the spine as the center and paste it symmetrically on the hard part of the shoulder.

- 2) There are many possible causes of low back pain. Low back pain, which

is suitable for low frequency therapy, mainly includes lumbar acid or muscle pain, nerve pain and so on.

Take the spine as the center, paste on one side of the pain or left and right.

5.6 Adhesive joint: excessive movement leads to increased load and joint pain.

Please paste two places in the pain joint, please paste two places in the pain joint.

5.7 Adhesive upper arm: arm pain may be due to congestion or fatigue in the muscle group, blood circulation is not smooth.

Please paste the upper and lower parts of the arm sore, please paste the upper and lower parts of the arm sore.

5.8 Stickiness: cold feet may be caused by the destruction of the balance of the nervous system and poor blood circulation in the hands and feet. The swelling or tiredness of the feet may be due to the deterioration of blood and fluid reflux.

Please paste one with the heart of the foot and the other with the calf or bare foot.

Please paste one with the heart of the foot and the other with the calf or bare foot.

5.9 Sticker leg: leg pain or fatigue may be due to congestion or fatigue in the muscle group.

Please test and paste on your legs, please test and paste on your legs.

5.10 Adhesive hip: hip pain may be caused by a long sitting.

Please test and paste around the buttocks.

6. Treatment procedures

6.1 choice of treatment modality.

TENS	MASS	EMS
It is suitable for relieving all kinds of acute and chronic pain, intractable pain or pain caused by daily work related activities.	Promote blood circulation and eliminate fatigue	Suitable for strengthening muscle, preventing muscle atrophy and promoting muscle regeneration.

6.2 Choice of treatment procedure, select from 12 modes.

cervical spine	shoulder	waist
upper arm	joint	arm
hip	thigh	calf
knee	heel	sole of foot

6.3 strength settings

40 intensity levels


Choose the appropriate intensity between 0~40 according to your personal preference.

7. Operation mode

7.1 Shutdown mode

In shutdown mode, all outputs are off.


In shutdown mode, press button [] to enter standby mode, after 1 second full display, then enter normal display:


7.2 Standby mode

Standby mode panel displays standby information, when press any keys, it will generate a sound prompt.


- 1) press the mode key [] to select treatment mode TENS-MASS-EMS:


- 2) Press the program key [] to select the treatment program:

- 3) On standby mode, press +A [] or D+ []], the panel shows the current treatment channel symbol and current treatment intensity, the treatment intensity value is greater than 0, and then enter the treatment mode:


- 4) set treatment time
Press [] button to set the treatment time (10-60) min, the default time is 30min, press [] button to add 10min


7.3 Treatment mode (non-locking condition).


- 1) In the treatment mode, long press the [] button to enter the shutdown mode:


- 2) In the treatment mode, press +A [] or D+ [] strength value plus 1, maximum strength level 40, Press [] or [] to directly return to the standby state.


- 3) In the treatment mode, when press B- [] or -C [], the intensity value will decrease by 1 each time, and return to standby mode when 0 intensity.


- 4) In the treatment mode, if the electrode patch is not adhered correctly, it will return to the standby.


- 5) In the treatment mode, press the [] key to select the treatment mode, and the treatment intensity returns to 0 at the same time.

- 6) In the treatment mode, press the [] key to select the treatment program, and the treatment intensity returns to 0 at the same time.


- 7) In treatment mode when no button press will automatically lock, treatment time start countdown only press B- [] or -C [] key to unlock:


7.4 Charging mode

In any mode, after plugging in the charging cable, it will exit the current mode and enter the charging mode. When charging mode panel only displays battery symbol, when it's under-charged battery symbol flashes, when it's fully-charged battery symbol is always lit.


7.5 Low-power display

When the battery is in low power state, a low power alarm will appear, and the battery symbol flashes.

7.6 Shutdown mode

In non-charge and standby mode, when no key is pressed within 30 seconds, the shutdown mode will automatically enter:


8. After the treatment is completed

8.1 Arrangement

The electrode pad is exposed and easy to become dirty and make the viscosity decrease, so please stick the electrode pad to the original place after treatment.

- 1) Turn off the power and unplug the electrode from the body. pulling the plug from the body, do not pull the wire, but hold the plug out.
- 2) Remove electrodes from skin.
 - Grab a corner and remove the electrode pad. Please do not contact the adhesive surface in a large area.
 - Do not pull the wire. Otherwise, it may lead to disconnection. (with pictures of pulling electrodes from the skin) Stick the electrode pad gently to the original electrode pad sticker.
 - Do not press the electrode hard, will cause paste surface damage.
 - Do not pull the electrode sheathpad or wire connection.
 - Please operate patiently to prevent excessive bending of the electrode pad.
 - Please keep the electrode pad sticker clean.
- 3) Store the body, electrode sheathpad and wire into the original packing box to prevent contamination.


Attention

keeping the body. Do not roll the wire on the body, wire may be excessively bent near the plug, resulting in breakage.

8.2 Maintenance and custody

- Maintenance methods
Please clean the body regularly.
- After power off, be sure to unplug the electrode sheathpad from the body and maintain it.
 - If there is dirt in the body, please wipe it with dry soft cloth.
 - If the is particularly dirty, use soft cloth with some water or neutral lotion, fully wring dry after wiping the body, and then wipe with dry soft cloth.
- ※ Do not let water and other liquids flow into the body.

Maintenance of electrodes

When the paste surface of the electrode sheetpad is dirty and difficult to paste, the soft cloth with water can be used to wet the adhesive surface and wipe off the dirt, which can temporarily restore the viscosity. But if the moisture is too much, it will lead to a decrease in viscosity.

※ Electrodes not washable,Otherwise, it may cause electric shock or accident.)

- Replacement cycle
The electrode sheetpad is consumable. When paste surface damage, please stop using and buy other products.
• Do not put electrodes in high temperature, humidity, direct sunlight.
• Maintenance and custody precautions
• Don't wipe dirt with thinner, gasoline, etc
• During maintenance, please handle the electrodes carefully
- Not bendable
- Do not force or drag.
- Do not touch water
• Do not keep it in the following places
- Places with water.
- High temperature, humidity, direct sunlight, dust and salt in the air.
- Places of inclination, vibration, shock.
- Places where chemicals are kept and where corrosive gases are produced.
• Please place where the baby can not touch.
If you do not comply with the above notice and the correct method of use, we will not be responsible for the quality of the product.
※ Do not carry out any repairs privately.

8.3 Application program

Table with 7 columns: Pro. Name, Time min., Frequency (Hz), Pulse width (uS), Wave (EMS), Ramp up and ramp down time (s), Keep Time(s), Release Time(s). It lists various programs like EMS, TENS, and MASSAGE with their respective parameters.

Table with 5 columns: Program Name, Time min., Frequency(Hz), Pulse width(uS), Wave(TENS). It lists programs like EMS, TENS, MASSAGE, and various intensity levels.

Table with 5 columns: Icon, Time range, Frequency range, Pulse width range, Modulation type. It lists different modulation modes like TENS, EMS, and MASSAGE.

Table with 7 columns: Pro. Name, Time min., Frequency (Hz), Pulse width (uS), Ramp up and ramp down time (s), Keep Time(s), Release Time(s). It lists various programs like EMS, TENS, and MASSAGE with their respective parameters.

9. Anomalies (failures and solutions)

Table with 3 columns: Issues, Possible causes, Solution. It lists common problems like 'Unable to boot', 'Automatic intensity zero', and 'Feeling unwell' with their causes and solutions.

10. Product specifications

Table with 2 columns: Product name, Model, Charge time, Output current, Output voltage (peak). It lists the product name LOOKEE® TENS UNIT EMS MASSAGE MUSCLE STIMULATOR and its specifications.

Table with 2 columns: Parameter, Value. It lists technical specifications like Power supply, Life time, Weight, Size, and Operating environment.

11. Final disposal

Please do not dispose of the product in the household waste at the end of it useful life. To protect the environment, dispose of empty batteries at appropriate collection sites according to national or local regulations.

12. Standardized symbolic description

Table with 2 columns: Symbol, Description. It lists various symbols used in the manual, such as warning, attention, electric shock, disposal, and CE marking.

13. Electromagnetic compatibility information

FCC STATEMENT

FCC CAUTION
Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.
This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:
(1) this device may not cause harmful interference, and
(2) this device must accept any interference received, including interference that may cause undesired operation

Note:
This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning

the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
• Increase the separation between the equipment and receiver.
• Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
• Consult the dealer or an experienced radio/TV technician for help.

! Note:

- Important information regarding Electro Magnetic Compatibility (EMC).
• With the increased number of electronic devices such as PC's and mobile (cellular) telephones, medical devices in use may be susceptible to electromagnetic interference from other devices.
• Do not use mobile (cellular) telephones and other devices, which generate strong electrical or electromagnetic fields, near the medical device.

! Warning:

- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation.
• Portable RF communications equipment (including peripherals such as an-tenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the LK113, including cables specified by the manufacturer.

Table with 3 columns: Emissions test, Compliance, Electromagnetic environment-guidance. It lists tests like RF emissions, Harmonic emissions, and Voltage fluctuations with their compliance levels and guidance.

Guidance and manufacturer's declaration - electromagnetic immunity - for all EQUIPMENT and SYSTEMS

Table with 4 columns: Immunity test, IEC 60601 test level, Compliance level, Electromagnetic environment-guidance. It lists immunity tests like Electrostatic discharge (ESD), Electrostatic transient / burst, surge, Voltage dips, and Proximity magnetic fields.

Table with 4 columns: Power frequency magnetic field IEC 61000-4-8, 30A/m, 30A/m, Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Note: U1 is the a. c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration - electromagnetic immunity - for EQUIPMENT and SYSTEM

The LK113 is intended for use in the electromagnetic environment specified below. The customer or the user of the LK113 should assure that it is used in such an environment.

Table with 4 columns: Immunity test, IEC 60601 test level, Compliance level, Electromagnetic environment-guidance. It lists tests like Conducted RF and Radiated RF with their compliance levels and guidance.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies. Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies. Note 2: These guidelines may not apply in all situations. Electromagnetic is affected by absorption and reflection from structures, objects and people.

Table with 4 columns: Immunity test, IEC 60601 test level, Compliance level, Electromagnetic environment-guidance. It lists tests like Proximity magnetic fields.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies. Note 2: These guidelines may not apply in all situations. Electromagnetic is affected by absorption and reflection from structures, objects and people.

a The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz, 3,5 MHz to 4,0 MHz, 5,3 MHz to 5,4 MHz, 7 MHz to 7,3 MHz, 10,1 MHz to 10,15 MHz, 14 MHz to 14,2 MHz, 18,07 MHz to 18,17 MHz, 21,0 MHz to 21,4 MHz, 24,89 MHz to 24,99 MHz, 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.

b Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the LK113 is used exceeds the applicable RF compliance level above, the LK113 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the LK113.

c Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM - for EQUIPMENT and SYSTEMS

Recommended separation distances between portable and mobile RF communications equipment and the LK113.

The LK113 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the LK113 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the LK113 as recommended below, according to the maximum output power of the communications equipment.

Table with 6 columns: Rated maximum output of transmitter (W), Separation distance according to frequency of transmitter (150 kHz to 80 MHz, 150 kHz to 80 MHz in ISM and amateur radio bands, 80 MHz to 800 MHz, 800 MHz to 2.7 GHz). It lists separation distances for different transmitter power levels.

For transmitters rated at a maximum output power not listed above the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

14. GUARANTEE

- 1. The warranty period for device is one year from date of delivery. In case of a warranty claim, the date of delivery has to be proven by means of the sales receipt or invoice.
2. Repairs under warranty do not extend the warranty period.
3. The following cases are excluded under the warranty
• All damage which has arisen due to improper treatment, e.g. nonobservance of the user instruction.
• All damage which is due to repairs or tampering by the customer or unauthorized third parties.
• Damage which has arisen during transport from the manufacturer to the consumer or during transport to the service centre.
• Accessories which are subject to normal wear and tear.

4. Liability for direct or indirect consequential losses caused by the unit is excluded even if the damage to the unit is accepted as a warranty claim.

Customer Support:

LOOKEE® Tech

24/7 Customer Support: Support@LookeeTech.com

Live Chat with us at www.LookeeTech.com

BEST-SUPPLIER INTERNATIONAL INC. (LOOKEE® Tech)
999 W Broadway, #720, Vancouver, BC V5Z 1K5, Canada
US Office: 99 Wall Street #1698, New York, NY 10005

Name: Shanghai International Holding Corp. GmbH (Europe)
Add: Eiffestrasse 80, 20537 Hamburg, Germany

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11. Final disposal

⚠ Please do not dispose of the product in the household waste at the end of it useful life. To protect the environment, dispose of empty batteries at appropriate collection sites according to national or local regulations.

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	Note: please refer to the instruction before using
	Attention: see Instructions for use!
	Type BF applied parts
	Electric shock protection type :Class II equipment
	Disposal in accordance with Directive 2002/96/EC (WEEE)
	Complies with the European Medical Device Directive (93/42/EEC), Notified Body is SGS Belgium NV
	Manufacturer: BEST-SUPPLIER INTERNATIONAL INC. (LOOKEE® Tech)
	Authorized representative in the European Community
	Serial number
	Temperature limitation
	Humidity limitation
	Atmospheric pressure limitation
	Level of protection for ingress of water or particulate matter into ME EQUIPMENT

13. Electromagnetic compatibility information

- ⚠ **Note:**
- Important information regarding Electro Magnetic Compatibility (EMC).
- With the increased number of electronic devices such as PC' s and mobile (cellular) telephones, medical devices in use may be susceptible to electromagnetic interference from other devices. Electromagnetic interference may result in incorrect operation of the medical device and create a potentially unsafe situation. Medical devices should also not interfere with other devices. In order to regulate the requirements for EMC

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Note: U_r is the a. c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration – electromagnetic immunity – for EQUIPMENT and SYSTEM

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Nom : Shanghai International Holding Corp. GmbH (Europe) Add : EIFFESSRASSE 80,