
stryker[®]

4K Coupler, C-Mount, with Advanced Imaging Modality

REF 1688-020-122

CE **R_x** ONLY

English	EN-1
Français	FR-13
Deutsch	DE-25
Italiano	IT-37
Português	PT-49
Español	ES-61
Nederlands	NL-73
Dansk	DA-85
Suomi	FI-97
Norsk	NO-109
Svenska	SV-121
Polski	PL-133
Ελληνικά	EL-145
Türkçe	TR-157
Русский	RU-169
简体中文	CH-183
日本語	JA-193
한국어	KO-205

English

4K Coupler, C-Mount, with Advanced Imaging Modality


Product Description and Intended Use

The Stryker 4K Coupler, C-Mount, with Advanced Imaging Modality (or “4K Coupler”) is used to attach endoscopes to Stryker surgical camera heads that do not have an integrated coupler. The coupler facilitates the easy exchange of different-sized endoscopes during a surgical procedure and allows users to focus the image captured through the endoscope.

When used in a system with a Stryker light source with Advanced Imaging Modality (e.g., 0220-220-300, 0220-230-300), the coupler can visualize the near-infrared light produced by the light source.¹ Refer to the light source user manual for complete system requirements.

¹ To verify if the 4K Coupler can visualize near-infrared light produced by other Stryker light sources, refer to the light source user manual.

Note: the 4K Coupler is not intended for use with the InfraVision™ Illuminator (0220-180-521).

1. Endobody clamp
2. Endobody (endoscope end)
3. Rear adapter (camera end)
4. Focusing ring

Indications

The 4K Coupler is indicated for use in general laparoscopy, general arthroscopy, nasopharyngoscopy, ear endoscopy, sinuscopy, and plastic surgery wherever a laparoscope/endoscope/arthroscope/sinuscope and the Stryker 1688 Video Camera are indicated for use.


Warnings and Cautions

1. Read these instructions thoroughly before using the device.
2. Before using this device, read the 1688 Video Camera user manual (P38139 English, P38888 multilingual) for warnings, indications, intended use, and other information about using the camera system.
3. Federal law (USA) restricts this device to use by, or on order of, a physician.
4. Inspect the coupler for any damage that might have occurred during shipment or prior use.

5. The coupler contains optomechanical parts that are precisely aligned. Dropping the coupler can cause mechanical damage and/or misalignment.
6. Test the coupler for functionality before using it in a surgical procedure.
7. This coupler is shipped non-sterile. Clean and sterilize the coupler prior to the first use and after every subsequent use. Follow the cleaning, disinfection, and sterilization instructions provided in this user manual.
8. Do not use the coupler and the camera system in the presence of flammable or explosive gases or in an oxygen-rich environment.
9. No modification of this device is authorized by Stryker.

Operating the Endobody Clamp

To perform setup, disassembly, and cleaning of the coupler, the endobody clamp (or “endoclamp”) must be depressed and/or released.

- To depress the endoclamp, push down as illustrated by the down arrow.
- To release the endobody clamp, release pressure on the clamp and allow the spring to push the clamp into its original position, as illustrated by the up arrow.


Setting Up the Coupler

When attaching or removing the coupler, grip only the rear adapter, as twisting other parts of the coupler with force can result in mechanical damage.

Do not overtighten the coupler, as this can damage the front window of the camera.


1. Gripping the rear adapter, screw the coupler clockwise onto the camera head until it forms a tight seal.

2. Depress the endobody clamp **1** and insert an endoscope into the endobody **2**.

Before each use, check the outer surface of the endoscope to ensure there are no rough surfaces, sharp edges, or protrusions.

3. Release the endobody clamp to secure the endoscope.
4. Attach a light cable from the light source to the light post on the endoscope 3.
5. With the camera console turned on, rotate the coupler's focusing ring 4 to establish a sharp focus on the display monitor.

Note: to remove the coupler, grip the rear adapter and unscrew the coupler counterclockwise from the camera head.

Reprocessing

These reprocessing instructions are provided in accordance with ISO 17664, ISO 15883 AAMI TIR12, and AAMI TIR30. The instructions have been validated by Stryker as being capable of preparing the device for re-use. To achieve the desired result, the processor shall ensure that the following instructions are performed as written in their entirety and as appropriate in the processor's facility. This normally requires routine monitoring and validation of the facility's reprocessing procedures. Stryker recommends users observe these standards when reprocessing medical devices.

Overview

Reprocessing the device involves manual or automated cleaning with either an enzymatic or a non-enzymatic detergent, optional disinfection, and sterilization.

- **Step 1** (required): Cleaning with Enzymatic or Non-Enzymatic Detergent
- **Step 2** (optional): Disinfection
- **Step 3** (required): Sterilization

Warnings

- This device must be cleaned and sterilized prior to the first use and after every subsequent use.
- Use only the sterilization cycles outlined in this document. Using unspecified sterilization cycles can damage the device or result in incomplete sterilization.
- Separate the camera head, coupler, and endoscope prior to cleaning, disinfection, and sterilization. Failure to follow this instruction will render the devices non-sterile. (Refer to the camera head and endoscope product manuals for reprocessing instructions of those devices.)

- The sterilization parameters presented in this document apply only when the device is sterilized outside of a sterilization tray. When using a sterilization tray, consult the instructions provided with the tray for proper sterilization parameters. Stryker recommends sterilizing the device inside of a sterilization tray.
- To avoid health risks from aerosol contamination, brush the device only when it is submerged in liquid.
- Wear appropriate protective equipment: gloves, eye protection, etc.
- Devices repaired by or purchased from third-party service organizations could expose patients to significant risk. These devices are no longer validated by Stryker for cleanliness, disinfection, and sterilization, or for safety and efficacy.
- The user shall defer to the facility's procedures regarding occupational exposure to bloodborne pathogens.

Cautions

- Do not use brushes or pads with metal or abrasive tips during manual cleaning, as permanent scoring or damage could result.
- The device cannot withstand automated disinfection.
- The 4K Coupler is not autoclavable. Steam sterilizing couplers that are not marked **AUTOCLAVE** will result in product damage.
- To minimize galvanic corrosion, avoid soaking dissimilar metals in close proximity.

Limitations on Reprocessing

- Do not cross-sterilize the device. Using multiple sterilization methods can significantly reduce the performance of the device.
- Repeated automated cleaning can degrade the product's cosmetic appearance.
- Damage incurred by improper processing will not be covered by the warranty.

Materials and Equipment

All materials and equipment required to reprocess the coupler shall be supplied by the user unless otherwise noted.

Item	Description
All phases	
Gloves, eye protection, etc.	Wear protective equipment as required by the medical facility and procedure
Cleaning	
Water basin	Large enough to accommodate the device
Lukewarm water	To prepare cleaning solutions
Detergent ¹	Used in cleaning solution to remove surgical debris
Soft-bristle brush ²	To clean exterior or hard-to-reach areas of the device
Reverse osmosis/deionized water ³	To rinse the device

Clean lint-free cloth or filtered pressurized air (≤ 40 psi)	To assist with drying
Automated washer	For using the automated cleaning procedure
Disinfection	
Water basin	Large enough to accommodate the device
Disinfecting solution ⁴	$\geq 2.4\%$ glutaraldehyde
Lukewarm water	To prepare disinfecting solution
Reverse osmosis/deionized water ³	To rinse the device
Clean lint-free cloth or filtered pressurized air (≤ 40 psi)	To assist with drying
Sterilization	
Sterilization system	<ul style="list-style-type: none"> • STERRAD[®] 100S, NX[®], 100NX[®], NX ALLClear[®], or 100NX ALLClear • STERIS/AMSCO[®] V-PRO[®] 1, V-PRO 1 Plus, V-PRO maX, V-PRO maX 2, or V-PRO 60
Sterilization wrap ^{5,6}	To maintain sterile barrier
Sterilization tray ^{6,7}	Optional. Must be compatible with sterilization method.

¹ The following detergents were validated for cleaning efficacy according to the manufacturer's instructions. Choose one of the detergents listed below or a substantially equivalent detergent:

Detergent	Type	Minimum Concentration	Minimum Soak Time
ENZOL [®] Enzymatic Detergent	Enzymatic	1 oz/gal.	1 minute
Prolystica [®] 2X Neutral Detergent	Non-enzymatic	1/8 oz/gal. (1 ml/L)	5 minutes
Neodisher [®] MediClean Forte	Enzymatic	5/8 oz/gal. (5 ml/L)	5 minutes

² Cleaning was validated with an M16 soft-bristle brush.

³ Rinsing the device during cleaning and disinfection was validated using reverse osmosis/deionized (RO/DI) water at $\leq 30^{\circ}\text{C}$.

⁴ Disinfection was validated using CIDEX[®] Activated at 25°C with a soaking time of 45 minutes.

⁵ Sterilization was validated using Kimberly-Clark[®] KC600 KIMGUARD sterilization wrap.

⁶ For United States users: when sterilizing the device, use only sterilization wraps and sterilization trays that have been cleared by the FDA to use with the selected sterilization cycle.

⁷ Stryker sterilization trays 0233-032-301, 0233-032-302, and 0233-410-002 are validated as compatible with the 4K Coupler.

Instructions for Reprocessing

Point of Use

- Disassemble the coupler from the endoscope and camera head. To disconnect the endoscope, depress the endobody clamp on the coupler and remove the endoscope from the coupler. To disconnect the coupler, grip the rear adapter of the coupler and unscrew it counterclockwise from the camera head.
- Wipe excess soil from the device.

Containment and Transportation

- Reprocess the device as soon as possible following use.
- Transport the device in a tray to avoid damage. Follow the facility's internal procedures for the transportation of contaminated surgical instruments and devices.

Cleaning

- 1 Follow the Preparation for Cleaning instructions below.
- 2 Then clean the device using either the Manual Cleaning or Automated Cleaning instructions below.

Note: for necessary materials and equipment, see the Materials and Equipment table.

Preparation for Cleaning

1. Fill a wash basin with lukewarm water.
2. Measure and dispense the desired amount of detergent into the water.
Note: see the Materials and Equipment table for validated detergents with their minimum concentration.
3. Gently mix the detergent into the water by hand.
4. Submerge the device into the prepared wash basin.
5. With the device immersed in the solution, thoroughly brush the exterior with a soft-bristle brush, focusing on any mated or rough surfaces.
6. Actuate and brush any movable parts in all extreme positions.
7. Rinse the device with water until all detergent residue is removed.
8. Once all detergent residue is removed, continue to rinse for 30 seconds.
9. Drain excess water from the device and dry it using a clean cloth or pressurized air.
10. Visually inspect the device for cleanliness, paying close attention to hard-to-reach areas. If visible soil remains, repeat steps 1– 10.

Manual Cleaning

1. Prepare a fresh solution of detergent with lukewarm water.
Note: see the Materials and Equipment table for validated detergents with their minimum concentration and soak time.
2. Wipe the entire surface of the device using a soft clean cloth dipped in the detergent solution.
3. Immerse the device in the detergent solution, ensuring the solution contacts all inner and outer surfaces.
4. Soak the device in the solution according to the manufacturer's recommendations.
5. With the device immersed in the solution, thoroughly brush the exterior with a soft-bristle brush, focusing on any mated or rough surfaces.
6. Actuate and brush any movable parts in all extreme positions.
7. Rinse the device with water until all detergent residue is removed.
8. Once all detergent residue is removed, continue to rinse for 30 seconds.
9. Drain excess water from the device and dry it using a clean cloth or pressurized air.
10. Visually inspect the device for cleanliness, paying close attention to hard-to-reach areas. If visible soil remains, repeat steps 1– 10.

Automated Cleaning

The device cannot withstand an automated disinfection method. When programming the washer, do not include a thermal rinse cycle or the device will be damaged.

1. Place the device in the automated washer on an incline to facilitate drainage.
2. Program the washer using the following parameters:

Phase	Recirculation Time	Temperature	Detergent Type
Pre-Wash	2 minutes	Cold Water	N/A
Wash 1	5 minutes	Set Point 60 °C (140 °F)	See Materials and Equipment table
Rinse 1	2 minutes	Hot Water	N/A
Dry Phase	2 minutes	115 °C (239 °F)	N/A

3. Drain excess water from the device and dry it using a clean cloth or pressurized air.
4. Visually inspect each device for cleanliness, paying close attention to hard-to-reach areas. If visible soil remains, repeat steps 1– 4.

High-Level Disinfection (Optional)

The device must be sterilized after disinfection. Failure to sterilize the device before reuse presents an acute infection control risk to the patient and operator.

Note: for necessary materials and equipment, see the Materials and Equipment table.

The device can be disinfected using a disinfecting solution that has the following active ingredient: $\geq 2.4\%$ glutaraldehyde.

1. Clean and prepare the device as recommended in this user manual.
2. Prepare the disinfecting solution and verify the minimum effective concentration according to the manufacturer instructions.
3. Immerse the device in the solution, filling all mated surfaces and crevices. Ensure air bubbles are removed from the surface of the device.
4. Allow the device to remain in contact with the disinfecting solution according to the manufacturer's recommended soak time.
5. Thoroughly rinse and flush the device with running, reverse osmosis/deionized water to remove the disinfectant.
6. Dry the device with a sterile, lint-free cloth immediately after rinsing. Filtered pressurized air may be used to assist in drying.

Sterilization

After performing the cleaning instructions specified above, perform one of the following sterilization cycles.

Note: for necessary materials and equipment, see the Materials and Equipment table.

STERRAD

1. Clean and prepare the device as recommended in this user manual.
2. If using a sterilization tray (optional), follow any additional instructions provided with the tray. Use only trays that are compatible with STERRAD.
3. Double wrap the device (or tray) prior to sterilization.
4. Sterilize the device using one of the following STERRAD sterilization systems:
 - 100S (Short or Long cycle)
 - NX (Standard or Advanced cycle)
 - 100NX (Standard or DUO cycle)
 - NX ALLClear (Standard or Advanced cycle)
 - 100NX ALLClear (Standard or DUO cycle)
5. Allow the device to cool to room temperature before reassembling it to an endoscope or camera head. Otherwise, the lens can fog during use.

STERIS/AMSCO V-PRO

1. Clean and prepare the device as recommended in this user manual.
2. If using a sterilization tray (optional), follow any additional instructions provided with the tray. Use only trays that are approved for sterilization with V-PRO.
3. Double wrap the device (or tray) prior to sterilization.
4. Sterilize the device using one of the following V-PRO sterilization systems:
 - V-PRO 1 (Standard cycle)
 - V-PRO 1 Plus (Non-Lumen or Lumen cycle)
 - V-PRO 60 (Non-Lumen, Lumen, or Flexible cycle)
 - V-PRO maX (Non-Lumen, Lumen, or Flexible cycle)
 - V-PRO maX 2 (Non-Lumen, Lumen, Flexible, or Fast Non-Lumen¹ cycle)

¹ **Warning:** When using the Fast Non-Lumen cycle, follow the instructions provided with the V-PRO maX 2 system regarding use of rigid containers or sterilization wrap. The cycle is intended to be used with sterilization pouches. Also observe the Fast Non-Lumen cycle weight limit of up to 11 pounds (4.99 kg) of instruments.
5. Allow the device to cool to room temperature before reassembling it to an endoscope or camera head. Otherwise, the lens can fog during use.

Maintenance

Inspection

Inspect the coupler for cleanliness before each use. If fluid or tissue buildup is present, repeat the above cleaning and sterilization procedures.

Inspect the coupler before each use. If a problem listed below is observed or suspected, contact your Stryker representative or return the device to Stryker for service.

- Ensure the endoclamp can be depressed fully and it returns to its original position when released. Ensure the endoclamp can securely hold the endoscope.
- Ensure the focusing ring is able to rotate and provide focus.
- The image quality performance of the coupler must be deemed acceptable. The coupler shall be able to provide an image that is clear and focusable.
- Observe any problems such as cloudiness, particles in the system, scratches, cracks, sharp edges, etc.

Transport and Storage Conditions

Observe the following environmental conditions for transport and storage:

- Temperature: -18–60 °C
- Relative Humidity: 15–90%

Storage

Store the coupler in a dry, clean, and dust-free environment at room temperatures.

Service Life

The expected service life of the 4K Coupler is 2 years.

Expected service life is determined by the number of times the device can be expected to be reused and/or reprocessed before it may require repair. In addition, the device's service life is largely determined by wear, reprocessing methods, and any damage resulting from use. To extend the time between device servicing, always follow the care and handling instructions in this user manual.

Before each use, test the device functionality and inspect it for any sign of damage per the Inspection section. If the device does not properly function or appears to be damaged, return it to Stryker for evaluation and/or repair. Repair through Stryker, the equipment manufacturer, brings the device back to manufacturer specifications. Clean and (when applicable) sterilize all potentially contaminated devices before returning them to Stryker.

Disposal

The device must be disposed of according to local laws and hospital practices. The device does not contain any hazardous materials.

Symbols

This device and its labeling contain symbols that provide important information for the safe and proper use of the device. These symbols are defined below.

Federal law (USA) restricts this device to use by, or on order of, a physician

Caution (consult instructions for use)

Consult instructions for use

Date of manufacture

Legal manufacturer

Device is shipped non-sterile and must be sterilized before use

Product catalog number

Serial number

Quantity

The device meets European Union medical device requirements

Stryker European representative

Product is manufactured in Singapore

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Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: the Stryker logo. All other trademarks are trademarks of their respective owners or holders.