
stryker®

Precision Ideal Eyes™ HD Autoclavable Laparoscopes

User Guide

CE₀₁₉₇ Rx ONLY

Contents

English	1
Français.....	9
Deutsch	19
Italiano.....	29
Português.....	39
Español	49
Nederlands	59
Dansk	69
Suomi.....	79
Norsk.....	87
Svenska.....	97
Polski.....	105
Ελληνικά.....	115
Türkçe	125
Русский	135
简体中文.....	145
日本語.....	153
한국어.....	161

Product Description

Stryker Precision Ideal Eyes™ (Precision IE) HD autoclavable laparoscopes are tubular optical instruments used to provide a view of internal patient anatomy for examination, diagnosis, and therapy during laparoscopic procedures.

Precision Ideal Eyes™ HD autoclavable laparoscopes are offered in a variety of sizes, working lengths, and directions of view.

Precision Ideal Eyes™ HD Autoclavable Laparoscopes

- Visualize with high image definition
- Visualize using the near-infrared light produced by the Stryker InfraVision™ Illuminator (220-180-521)
- Autoclavable
- Precision Ideal Eyes™ laparoscopes are easily identifiable by their distinct purple eye piece and are classified with the following part numbers:

Precision Ideal Eyes™ HD Autoclavable Laparoscopes			
Angle of view	0°	30°	45°
5 mm	502-503-010	502-503-030	502-503-045
10 mm	502-103-010	502-103-030	502-103-045
Bariatric 5 mm	502-513-010	502-513-030	502-513-045
Bariatric 10 mm	502-113-010	502-113-030	502-113-045

All Stryker laparoscopes are intended for use by qualified persons during laparoscopic procedures to provide an internal view or image of the patient. In all cases, the surgeon is best advised to use a method in which his/her own practice and discretion dictate to be the best for the patient. These instructions are recommended for the proper function of the device. This material is not intended as a reference for laparoscopic examination or surgery, nor is it intended to replace training, including relevant preceptorship(s).

1. Distal End
2. Light Post
3. Proximal End

Indications

Standard Laparoscopes

The Stryker Laparoscope is intended to be used for gynecological and general procedures that clinicians deem appropriate for the patient. Laparoscopes have been available for many years and laparoscopy is considered a mature surgical technique, which has been established in general and gynecological procedures.

Bariatric Laparoscopes

Laparoscopes are intended to be used by surgeons in diagnostic and therapeutic procedures. Laparoscopic minimally invasive procedures are performed in the abdominal cavity by means of small skin punctures that allow the insertion of the laparoscope and laparoscopic instruments. This includes, but is not limited to such uses as gallbladder and appendix removal, hernia repair, gastric bypass, laparoscopic Nissen and examination of the abdominal cavity, appendix, gallbladder and liver. Bariatric laparoscopes allow surgeons to perform procedures on morbidly obese segments of their patient population.

Contraindications

There are no known contraindications.

Warnings

1. **Rx ONLY** Federal law (USA) restricts this device to sale by or on the order of a physician.
2.  This laparoscope is shipped non-sterile. To prevent infection, clean and sterilize the laparoscope prior to the first and every subsequent use. Prior to cleaning or sterilization, remove any protective sheathing that was used to protect the laparoscope during shipment.
3. Laparoscopic procedures should be performed only by persons having adequate training and familiarity with laparoscopic techniques. Consult medical literature relative to techniques, complications and hazards prior to performing any laparoscopic procedure.
4. Failure to use protective filters or suitable filtering spectacles during the activation of a surgical laser beam may cause eye damage to the user.
5. Activation of a 600 – 1100 nm wavelength surgical laser during a laparoscopic procedure can cause the image to become saturated, and therefore unclear. Verify compatibility prior to surgery.
6. When used with a light source, the laparoscope can reach temperatures exceeding 41°C (106°F) near the light cable inlet and the distal tip, creating a risk for patient burns and operating room fires. To reduce these risks:
 - Never leave the laparoscope in direct contact with patient tissue or combustible materials, such as drapes or gauze.
 - Always start the procedure with the light source at its lowest possible setting. As needed, gradually increase the light output to the minimum level where adequate illumination is achieved.
 - Before disconnecting the light cable or connectors, turn off

the light source and allow the laparoscope to cool.

7. Do not use the laparoscope as a probe or apply pressure with the laparoscope tip to surrounding tissues. The temperature and construction of the laparoscope can damage patient tissue.
8. Before each use, the outer portion of the laparoscope intended to be inserted into the patient should be checked to ensure there are no rough edges or protrusions that may accidentally cut or damage patient tissue.
9. Contact with a rotating or cutting instrument may cause the surface of the laparoscope to become sharp or jagged or may result in portions of the laparoscope being broken into the patient.
10. Failure to remove the laparoscope and laparoscopic accessories during a cardiac defibrillator discharge may lead to punctured or injured patient tissue.
11. To ensure a clear image, clean the window at the laparoscope's distal tip with sterile alcohol wipes or cotton tip applicators and isopropyl alcohol.

Cautions


1. Keep the laparoscope clean and dry when not in use. The stainless steel exterior of the laparoscope is rust-resistant, not rust-proof.
2. Bending the laparoscope, or using the laparoscope as a lever or pry bar, may result in lens damage and may render the laparoscope unusable.
3. Contact with a surgical laser beam may damage the laparoscope surface and internal optics.
4. Third party repair is not recommended for Precision Ideal Eyes™ Laparoscopes. Repair by any party other than Stryker Endoscopy could compromise InfraVision™ IR Illuminator compatibility.
5. Some surgical lasers may not be visible when used with Precision Ideal Eyes™ laparoscopes. Prior to surgery, verify that the selected laser can be viewed through the laparoscope.
6. Only Infrared laparoscopes can be used in conjunction with the InfraVision™ IR Illuminator (REF P/N 220180521). Using non-infrared laparoscopes with the InfraVision™ IR Illuminator can lead to decreased visibility of infrared wavelengths. Precision Ideal Eyes™ laparoscopes are easily identifiable by their distinct purple eye piece.

Electrical Safety

1. Stryker laparoscopes should only be connected to Type F devices. Stryker laparoscopes are classified as Type CF when attached to a Type CF device, and Type BF when attached to a Type BF device.
2. When laparoscopes are used with energized endoscopically-used accessories, patient leakage currents may be additive. To minimize total patient leakage current, use only Type CF or Type BF endoscopic accessories. Type CF applied parts should always be used together with other Type CF applied parts.
3. Consult literature for the safe use of electrosurgical equipment. Inadvertent burns can occur when the patient return path is obstructed.


Connecting a Light Cable

Three adapters (Wolf, Storz and ACMI) are included to enable connection of the laparoscope to different styles of light cables. SafeLight™ adapters may be purchased separately.

Notes:

- For best performance, ensure the light cable is not worn or damaged.
- Remove the adapter from the laparoscope and light cable before cleaning and sterilization.

- 1 Connect the light post to the appropriate adapter.
- 2 Connect the adapter to the appropriate cable.

Using Polishing Paste

Polishing paste is included with the laparoscope. If any of the three optical surfaces (distal end, light post or proximal end) has deposits that cannot be removed in normal cleaning practices, they can be removed with the paste prior to cleaning and sterilization.

1. Apply a small amount of polishing paste to a clean cotton swab.
2. Gently rub the cotton swab on the optical surface, cleaning in a circular motion.
3. Use acetone or alcohol on a clean cotton swab to remove the paste residue.
4. Clean and sterilize the laparoscope before its next use.

Caution: Polishing optical surfaces should not be part of routine cleaning procedures. Repeated application of polishing paste to the laparoscope lens can cause damage. It should only be performed when the laparoscope image is cloudy.

Reprocessing

These reprocessing instructions are provided in accordance with ISO 17664, AAMI TIR 12, AAMI ST79, and AAMI ST81. While they have been validated by Stryker as being capable of preparing the device for re-use, it remains the responsibility of the processor to ensure that the reprocessing, as actually performed (using equipment, materials, and personnel in the reprocessing facility), achieves the desired result. This normally requires validation and routine monitoring of the process. Stryker recommends users observe these standards when reprocessing medical devices.

Warnings

- This device must be cleaned and sterilized prior to the first use and after every subsequent use.
- The sterilization parameters presented in this document apply only to laparoscopes sterilized outside of a sterilization tray. When using a sterilization tray, consult the instructions provided with the tray for proper sterilization parameters. Stryker recommends sterilizing device inside of a sterilization tray.
- Steam sterilization is the recommended method for Stryker laparoscopes.
- Use only the sterilization cycles outlined in this document. Using unspecified sterilization cycles may damage the device or result in incomplete sterilization.
- Failure to follow cleaning instructions may result in inadequate sterilization.
- Prior to cleaning, disinfection, or sterilization, remove any protective sheathing that was used to protect the laparoscope during shipment. These include rubber caps on the distal and proximal end of the laparoscope as well as a plastic cover on the light post.
- Prior to cleaning, disinfection, or sterilization, separate the laparoscope from the camera head, coupler, and adapter, which may come pre-attached to the laparoscope. If any of these components are cleaned, disinfected, or sterilized as a single unit, disconnecting the devices during use will compromise the sterility of the products. (Refer to the camera head and coupler product manuals for reprocessing instructions.)
- Unwrapped immediate-use (“flash”) steam sterilization is not recommended because water source impurities may cause hard films to build up on the outer optical windows of the laparoscope, reducing optical performance.
- Wear appropriate protective equipment: gloves, eye protection, etc.

Cautions

- Do not use brushes or pads with metal or abrasive tips during manual cleaning, as permanent scoring or damage could result.
- To minimize galvanic corrosion, avoid soaking dissimilar metals in close proximity.
- Allow the device to air cool following steam sterilization. Rapid cooling or “quenching” in a liquid will damage the device and void the warranty.

Limitations on Reprocessing

- Using multiple sterilization methods may significantly reduce the performance of the device and is not recommended.
- Proper reprocessing has a minimal effect on this device. End of life is normally determined by wear and damage due to use.
- Damage incurred by improper reprocessing will not be covered by the warranty.

Instructions for Reprocessing

Point of Use

- Wipe excess soil from the device using disposable paper towels or lint free cloth.
- If an automated reprocessing method will be used, rinse any channels in the device with 50 mL of sterile distilled water immediately after use.

Containment and Transportation

- Reprocess the device as soon as reasonably practical following use.
- Transport the device in a sterilization tray to avoid damage.

Preparation for Cleaning

1. Remove any light post adapters from the scope prior to cleaning.
2. Prepare an enzymatic detergent¹ according to the manufacturer's recommendations.
3. Wipe the entire device with the detergent, using a clean cloth.
4. Immerse the device in the detergent. Inject any inside regions of the device with 50 mL of the detergent solution to ensure all parts of the device are reached.
5. Soak the device in the detergent for 15 minutes.

Manual Cleaning

1. Brush

- Prepare a fresh solution of enzymatic detergent¹ according to the manufacturer's recommendations.
- Thoroughly brush the exterior of the device with a soft-bristled brush, focusing on any mated or rough surfaces.
- Using a syringe, inject any mated surface 5 times with 50 mL of the detergent.

2. Rinse

- Rinse the device with water² at ambient temperature until all detergent residue is removed.
- Flush any mated surfaces 5 times.
- Once all detergent residues is removed, continue to rinse for a minimum of 30 seconds.
- Drain excess water from the device and dry it using a clean cloth or pressurized air.
- Visually inspect the device for cleanliness, paying close attention to hard-to-reach areas. If visible soil remains, repeat steps 1 and 2.

3. Soak

- Prepare a non-enzymatic detergent³, according to the manufacturer's recommendations.
- Fully immerse the device and inject any mated surfaces with 50 mL of the detergent.
- Soak the device for 15 minutes.

4. Brush

- Thoroughly brush the exterior of the device using a soft-bristled brush.
- Inject the prepared detergent into any mated surfaces 5 times.

5. Rinse

- Thoroughly rinse the device with water² until all detergent residue is removed. Flush any crevices 5 times. Once all detergent residue is removed, continue to rinse for 30 seconds.
- Drain the excess water from the device.

¹ Enzol® at 1 oz/gal (US) at 35 – 40°C was validated for cleaning efficacy.

² Reverse Osmosis/Deionized (RO/DI) water was validated for cleaning efficacy.

³ Prolystica 2x Neutral Detergent at 1/8 oz/gal was validated for cleaning efficacy.

Automated Cleaning

1. Brush

- Thoroughly brush the device with an appropriately sized soft-bristled brush, focusing on hard-to-clean areas.

2. Rinse

- Rinse the device with water² at ambient temperature until all visible detergent residues have been removed.
- Flush all crevices and mated surfaces at least five times.
- Continue to rinse the device for a minimum of 30 seconds after all detergent residue has been removed.

3. Automated Wash

- Place the device in an automated washer on an incline to facilitate drainage.
- Program the washer with the following parameters, then activate the wash:

Pre-Wash	
Recirculation Time	2 Minutes
Water Temperature	Cold
Detergent Type	N/A

Enzyme Wash	
Recirculation Time	2 Minutes
Water Temperature	Hot
Detergent Type	Enzymatic Detergent ⁴

Wash 1	
Recirculation Time	2 Minutes
Water Temperature	66°C Set Point
Detergent Type	Non-Enzymatic Detergent ⁵

Rinse 1	
Recirculation Time	2 Minutes
Water Temperature	Hot (60°C)
Detergent Type	N/A

Dry Phase	
Recirculation Time	7 Minutes
Water Temperature	115°C
Detergent Type	N/A

⁴ Enzol® at 1 oz/gal (US) was validated for cleaning efficacy.

⁵ Prolystica 2x Neutral Detergent at 1/8 oz/gal was validated for cleaning efficacy.

Drying

- Dry the device using a cleaning cloth.
- Filtered pressurized air can be used to assist in drying.

Inspecting

1. Visually inspect the device, including all internal surfaces, for remaining soil.
2. If soil remains, repeat manual or automated cleaning procedure, focusing on those areas.

Disinfection (optional)

Warning: Not to be used in place of sterilization.

1. Remove any light post adapters from the scope prior to disinfection.
2. Disinfect the device in a disinfecting solution that has the following active ingredients:
 - ≥ 2.4% glutaraldehyde (with minimum immersion time of 45 minutes at 25°C)⁶
 - ≥ 0.55% ortho-phthalaldehyde (with a minimum soaking time of 12 minutes at 25°C)⁷
3. Prepare the disinfecting solution according to the manufacturer's instructions.
4. Per the manufacturer's recommendations, immerse the device in the disinfecting solution for the required time at the appropriate temperature.
5. Using a syringe, flush any openings a minimum of 10 times each with the same disinfecting solution.
6. Thoroughly rinse and flush all interior and exterior surfaces with running deionized water until the device is visibly clean and all disinfecting solution and/or residue is removed.
7. Dry all parts with a lint-free towel immediately after rinsing.

⁶ Cidex® Activated was validated for disinfection efficacy.

⁷ Cidex® OPA was validated for disinfection efficacy.

Thermal Disinfection (optional)

Warning: Not to be used in place of sterilization.

1. Disassemble any adapters from the scope prior to disinfection.
2. Place the scope into the washer disinfectant.
3. Program the washer using the following parameters:

Thermal Disinfection	
Recirculation Time	1 Minute
Water Temperature	90°C

Sterilization

After performing the cleaning instructions specified above, perform one of the following sterilization cycles:

Steam Sterilization for Precision Ideal Eyes™

Laparoscopes

1. Clean and prepare the device as recommended in the Cleaning and Disinfection sections.
2. Disassemble any adapters from the scope prior to sterilization.
3. If using a sterilization tray, place the scope and adapters in the tray and follow instructions provided with the tray. Use only trays that are approved for sterilization with steam.
4. Double wrap the device (or tray, if being used) prior to sterilization.
5. If sterilizing the device without a tray, use the following parameters for the appropriate length:

 **Warning:** Drying time depends on several variables, including altitude, humidity, type of wrap, preconditioning, size of chamber, mass of load, material of load, and placement in chamber. Users must verify that drying time set in their autoclave yields dry surgical equipment.

Standard Length Laparoscopes

Pre-Vacuum (U.S.)	
Wrapping	Double
Temperature	132°C (270°F)
Sterilization Time	4 minutes
Dry Time	30 minutes

Pre-Vacuum (E.U.)	
Wrapping	Double
Temperature	134°C (273°F)
Sterilization Time	3 minutes
Dry Time	30 minutes

Gravity	
Wrapping	Double
Temperature	132°C (270°F)
Sterilization Time	15 minutes
Dry Time	35 minutes

"Flash" Pre-Vacuum	
Wrapping	-
Temperature	132°C (270°F)
Sterilization Time	4 minutes
Dry Time	-

Bariatric Length Laparoscopes

Pre-Vacuum (U.S.)	
Wrapping	Double
Temperature	132°C (270°F)
Sterilization Time	4 minutes
Dry Time	45 minutes

Pre-Vacuum (E.U.)	
Wrapping	Double
Temperature	134°C (273°F)
Sterilization Time	3 minutes
Dry Time	45 minutes

"Flash" Pre-Vacuum	
Wrapping	-
Temperature	132°C (270°F)
Sterilization Time	4 minutes
Dry Time	-

6. Allow the laparoscope to completely cool before reassembling it to a coupler or camera head. Otherwise, the lens will fog during use.

STERRAD®

Warning: Not all sterilization trays are compatible with STERRAD® systems. Using an incompatible tray may result in incomplete device sterilization. Consult the instructions that came with your sterilization tray to determine which sterilization method is compatible with your tray and devices. If a compatible tray is not available, the devices can be double-wrapped for use with the STERRAD® system.

1. Clean and prepare the endoscope as recommended in the Cleaning and Disinfection sections.
2. Disassemble any adapters from the scope prior to sterilization.
3. Double wrap the device (or tray, if being used) prior to sterilization.
4. Sterilize the endoscope following the instructions of the manufacturer, using the STERRAD® 100S, 200, NX™ or 100NX™ Sterilization System. Select the standard cycle.
5. Allow the endoscope to completely cool before reassembling it to a coupler or camera head. Otherwise, the lens will fog during use.

STERIS®/Amsco® V-PRO

1. Clean and prepare the endoscope as recommended in the Cleaning and Disinfection sections.
2. Disassemble any adapters from the scope prior to sterilization.
3. Double wrap the device (or tray, if being used) prior to sterilization.
4. Sterilize the endoscope using the V-PRO maX Sterilizer (Non-Lumen or Standard Cycle), the V-PRO 1 Plus Sterilizer (Non-Lumen or Standard Cycle), or the V-PRO 1 Sterilizer (Standard Cycle).
5. Allow the endoscope to completely cool before reassembling it to a coupler or camera head. Otherwise, the lens will fog during use.

Maintenance, Inspection, and Testing

- Inspect the device on a continual basis. If a problem is observed or suspected, the device should be returned for repair.
- Inspect all components for cleanliness. If fluid or tissue buildup is present, repeat the above cleaning and disinfection procedures.
- Inspect the optical surfaces of the laparoscope for clarity. If there are deposits that need to be removed, follow the instructions in the "Using Polishing Paste" section.

Caution: Polishing optical surfaces should not be part of routine cleaning procedures. Repeated application of polishing paste to the laparoscope lens can cause damage. It should only be performed when the laparoscope image is cloudy.

Expected Service Lifetime

The recommended service interval for Stryker Precision Ideal Eyes™ Laparoscopes is one year.

Technical Specifications

Transport and Storage Conditions:

Temperature: -18 – 60°C (0 – 140°F)

Relative Humidity: 15 – 90%

Insertion portion width of 10 mm Precision Ideal Eyes™ laparoscopes:

10.00 ±.05 mm (standard)

10.00 +.02/-.05 mm (bariatric)

Insertion portion width of 5 mm* Precision Ideal Eyes™ laparoscopes:

5.5 +.02/-.05 mm (standard and bariatric)

*Precision Ideal Eyes™ Laparoscopes marked on device as 5.5 mm are part of and compatible with the 5 mm family of products.


Laparoscope Disposal

Laparoscopes must be disposed of according to local laws and hospital practices.

Warranty and Return Policy

Go to <http://stryker.com/EndoSalesTerms>.

Symbol Definitions

	Consult instructions for use
	Caution; see instructions for use
	Non-Sterile
	Manufacturer
	Date of Manufacture
	Catalog Number
	Serial Number
	Quantity
	Authorized Representative in the European Community
	Conformance to 93/42/EEC
	Made in Germany
	Temperature Range
	Humidity Range

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Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: **Precision Ideal Eyes™**, **InfraVision™**, **SafeLight™** and the **Stryker logo**. All other trademarks are trademarks of their respective owners or holders.