

General Instructions for Using the iNPLANT!TM

- 1 Open the sterile barrier to release the inner bag.
- 2 Holding the inner bag level, apply pressure towards the center to burst the saline packet located within the funnel. Keep applying pressure up and down the bag for at least 60 seconds to allow the wicking fabric to evenly distribute the saline.
- 3 Peel open the inner sterile packaging and transfer the funnel to the sterile field.
- 4 With one hand holding the funnel flat, remove the wicking fabric by pulling on the tab at the distal end. This will distribute the saline throughout the inner surface of the funnel, activating the lubricious coating.
- 5 Fully open the Glide Track on the back of the funnel. Pour the implant from the basin into the funnel. We advise adding additional antibiotic impregnated saline into the funnel. This will not escape as the end is sealed.
- 6 Close the Glide Track by forcefully engaging both tracks right over left. Confirm the closure in both directions with auditory and tactile feedback.
- 7 To ensure the funnel is fully lubricated, we suggest sliding the implant back and forth within the funnel.
- 8 Cut the tip of the funnel to the desired size according to the guidelines printed on the funnel. We recommend using serrated or curved Mayo scissors. Deliver the implant into the breast by squeezing the implant towards the funnel opening.



Sterilized by
Gamma



Lower Limit Temperature



Use By Date



Non-Phthalate



Do Not Reuse



Do Not Resterilize



Date
Manufactured



Consult Instructions
for Use



Manufacturer



Caution



Keep Dry



Catalog Number



Batch Code



Not made with Natural
Rubber Latex

DO NOT use the product if the package has been opened or damaged, the tamper-proof seal has been violated or the water pouch has been ruptured or deflated



U.S. Federal Law restricts this device to sale by or on the order of a licensed physician

Description

The iNPLANT!TM is a sterile, transparent, cone-shaped sleeve constructed of a flexible, Pthalate-free polymeric film and has a lubricious hydrophilic coating on the inside surface.

Intended Use

The iNPLANT!TM is intended to facilitate the delivery of silicone gel implants by providing a shell-tissue interface with less friction during insertion of the implant. It is also intended to facilitate implementation of a "no-touch technique" in implant insertion.

Product Warranty

Proximate Concepts, LLC stands behind the iNPLANT (TM) funnel. If you are not completely satisfied after the use of one iNPLANT Funnel, contact our Customer Service department within 30 days to receive a refund. Once the Return Form and unwanted items are received and processed, we will issue a refund to your original form of payment. Please note that Shipping Fees and interest fees are non-refundable.

If you experience difficulty in its use resulting in the need for a second funnel being used per patient, we will replace it at no charge.

Prior to 1st Time Use

We **strongly** advise you to watch the instructional video available on our website.

[Inplantfunnel.com/instructions](https://inplantfunnel.com/instructions)



24/7
Instructional Hotline

201-961-5084



Address

1566 Lemoine Avenue,
Fort Lee, NJ 07024



E-mail

orders@inplantfunnel.com



Customer Service

201-913-8833
Fax: 201-351-3593