

Random Colors

The Multi-Color Eyeset can randomly change colors for a unique experience.

Press and hold the color select button(s) for a few seconds until the lights stop flashing. Release the button. Repeat for each side. Disconnect the eyeset to reset it. The eyeset does not stay in random mode between uses.

Color Lock

To prevent accidental color changes, the eyesets may be locked.

Set the desired color first. Press and hold the color select button(s). After several seconds the lights will go out; **do not release** the button. Continue to hold the button until the lights come back on, then release. Repeat for each side. Repeat the process to unlock.

Note: the eyeset will remain locked even if disconnected or turned off.



Cleaning and Care

Wipe with a clean, soft cloth and a mild detergent to remove smudges from the eyeset. Use alcohol wipes to disinfect the eyeset between uses.

Disconnect the eyeset cable from the DAVID device and the eyeset before storage.

Mind Alive Inc.

6716 75 St NW
Edmonton, Alberta
Canada T6E 6T9

Toll Free: 800.661.MIND
Phone: 780.465.MIND
E-mail: info@mindalive.com
www.mindalive.com

Tru-Vu Omniscreen Multi-Color Eyeset



For DAVID Audio-Visual
Entrainment Devices



***Enhance your AVE experience
with the benefits of color therapy***

Color Therapy

The use of colors in meditation and spirituality has a long history. While there has been little research done on the use of color or its application with Audio-Visual Entrainment, there are many anecdotal accounts of the positive effects of color therapy.

Different colors are thought to influence our moods and emotions. Colors affect our perceptions of space. Often certain colors are associated with “calming” effects or a sense of well-being

The Multi-Color Eyeset uses red, green, and blue emitted light to create a range of colors that can be used for basic color therapy, or just to enhance the enjoyment of your AVE experience. Feel free to experiment with different colors to find your favorites.

Color Preferences
Everyone has colors they prefer. Try your favorite AVE session with your favorite color (or colors). Try mixing two colors for different results.

Color Suggestions

White

All of our studies are based on white light stimulation. Use with any session.

Blue

Calming and relaxing. Best for deep meditation or sleep. Use with Meditation or Sleep categories of sessions.

Cyan

Blue tinted light has been shown to reduce the risk of seizure (Takahashi & Tsukahara). Use with any session.

Green

Effective for pain reduction and meditation. Use with Alpha and Schumann sessions (7-10Hz).

Yellow

Improve cognition and focus. Use with Beta, ADD and Brain Brightener sessions (14-20Hz).

Red

Increase physical arousal and energize yourself. Caution: may cause anxiety. Use with Beta and Energizer sessions (14-20Hz).

Magenta

Best for contemplative meditation and building creativity. Use with Theta sessions (4-7Hz)

How to Use

Connect the Eyeset

The Multi-Color Eyeset connects to your DAVID device using the same 2.5mm patch cord as a white eyeset. Plug into the Omniscreen () jack.

The Multi-Color Eyeset remembers the last color used for each field.

Start a Session

Start a session as directed by your device's manual. This eyeset is

designed to be used with your eyes closed. The Multi-Color Eyeset may take a few more seconds to turn on once the session is started.

Select a Color

Gently press the buttons underneath the white screen to cycle through the colors.



Each visual field (left and right) must be set independently. Adjust the intensity normally, from your DAVID unit.

More Information

For more information on AVE or the DAVID devices, please visit our website at www.mindalive.com