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
DE

USAGE

The FUCHSI Kegel Exercise Weights helps you to develop a strong pelvic floor musculature easily and regularly:

- Start training with the pink Kegel, the lightest one (45g). After washing, carefully insert it about 3-4cm deep into the vagina, as you would introduce a tampon.
- Practice with the pink Kegel daily and hold it for 10 minutes when standing or walking. Increase the hold time slowly over the week. (Caution: sore muscles possible,so please pause.)
- As soon as you can easily hold pink Kegel for 30 minutes, increase the training weight with the red Kegel (75g).
- And this is how it goes on until you can also hold the red Kegel for 30 minutes while standing and walking. Then you can grab the last weight of our set and continue the exercises with the yellow Kegel,the 105g.

EN

SAFETY INSTRUCTIONS

- Thoroughly clean your FUCHSI Kegel Exercise Weights prior to initial use and after each use with warm water and and ensure clean conditions (such as washing hands).
- Check the Kegel weights for scratches or other damage to the surface and only use flawless Kegel weights. A Kegel weight scratched or otherwise damaged surface can lead to injuries.
- In case of vaginal infections, inflammation, injuries or painful discomfort in the genital area as well as during menstruation and in case of a problem pregnancy, the application of the Kegel Weights should be avoided.
- In addition, use should be avoided in the following cases: known hypersensitivity, virginity.
- Cervical capsule, contraceptive ring, diaphragm or intrauterine device may be worn during exercise. As a precaution, contact your doctor first.
- Use the FUCHSI Kegel weights after pregnancy only after the puerperal period (6-8 weeks after birth).

EN

CLEANING INSTRUCTIONS

- Clean the FUCHSI Kegel weights after each use with warm water and soap.
- Cleaning by boiling or steam and dry sterilization is not suitable as it leads to deformations of the material.
- Keep the FUCHSI Kegel weights set in clean conditions when not in use.

For any further questions please contact info@fuchsi.shop.

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GUTSCHEINCODE: FUCHSI20

PROMOTE TO A FRIEND

As a thank you for your purchase we offer you and your friends 20% discount on your next Fuchsi product!The coupon you can pass on to any number of your friends you want and it can be used on every Amazon account once.

DISCOUNT CODE: FUCHSI20

DE/EN