

Cat. #	Description	Product Type	Size	Applications	Species Reactivity	ICC Dilution
AP1651	Actin (N-terminal region)	Rabbit pAb	50 µl	WB, ICC, E, IHC, IP	Hu, Rt, Ms, Ck	1:50
AM2021	Actin (C-terminal region)	Mouse mAb	50 µl	WB, E, ICC, IHC	Hu, Rt, Ms, Ck	1:50
PF7501	Phalloidin:FITC	Reagent	100 µl	ICC, IHC	Hu, Rt, Ms, Ck, F	1:500
MS3031	Anti-Mouse Ig:DyLight® 594	Goat pAb	100 µl	ICC, IHC	Ms	1:200
RS3271	Anti-Rabbit Ig:DyLight® 594	Goat pAb	100 µl	ICC, IHC	Rb	1:200

Applications: WB = Western blot, E = ELISA, ICC = Immunocytochemistry, IP = Immunoprecipitation, IHC = Immunohistochemistry, FC = Flow Cytometry  
 Species: H = Human, R = Rat, M = Mouse, C = Chicken, F = Fish, Fr = Frog, Rb = Rabbit

### Kit Summary


The actin filament regulation kit can be used to compare labeling between actin filaments and total actin. The kit contains two antibodies for detecting total actin, and Phalloidin:FITC for labeling actin filaments. In addition, secondary reagents conjugated to DyLight® 594 are included for dual labeling experiments with Phalloidin:FITC.

### Buffers and Storage

Mouse monoclonal, rabbit polyclonal, and secondary reagents are supplied in phosphate-buffered saline, 50% glycerol, 1 mg/ml BSA, and 0.05% sodium azide. Store at –20°C. Stable for 1 year.

### Background

Actin is a major cytoskeletal protein involved in diverse cellular functions including cell motility, adhesion, and morphology. Six different actin isoforms have been identified in vertebrates. There are four  $\alpha$  isoforms: skeletal, cardiac, and two smooth muscle (enteric and aortic) actins, along with two cytoplasmic actins ( $\beta$  and  $\gamma$ ). Actin exists in two principal forms, globular monomeric (G) actin, and filamentous polymeric (F) actin. The assembly and disassembly of actin filaments, and also their organization into functional networks, is regulated by a variety of actin-binding proteins (ABPs). Phosphorylation may also be important for regulating actin assembly and interaction with ABPs. In Dictyostelium, phosphorylation of Tyr-53 occurs in response to cell stress and this phosphorylation may alter actin polymerization. In B cells, SHP-1 tyrosine dephosphorylation of actin leads to actin filament depolymerization following BCR stimulation.



FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

### Adherent Cell Fixation

1. Remove cell growth medium from culture plate containing cells, and rinse cells once with Hank's buffered saline solution (HBSS) or other rinse buffer acceptable for your cell type.
2. Fix cells with 4% Paraformaldehyde/0.2% NP-40 in HBSS for 30 minutes at room temperature.

**Note:** Some antibodies work better for immunocytochemistry using one of the following methods:

A. Methanol/Acetone fixation: Fix and permeabilize in 1:1 Methanol/Acetone at -20°C for 10 min.

B. Aldehyde/Acetone fixation: Fix cells with 4% Paraformaldehyde in HBSS for 30 minutes at room temperature, then permeabilize for 15 min. with 100% Acetone for at -20°C.

3. Remove fixation solution and rinse cells two times with phosphate buffered saline solution (PBS).
4. Block non-specific binding sites with 1% bovine serum albumin (BSA) in PBS for 30 minutes at room temperature.

**Note:** Normal animal serum (e.g. horse, goat) that matches the species of the secondary antibody can be substituted for BSA, if non-specific labeling occurs with certain secondary reagents.

### Primary Antibody Labeling

5. Make primary antibody dilutions in 1% BSA in PBS, using the recommended dilution described in the table above. For some cell types, the optimal antibody dilution may need to be empirically determined. Titrations of 1:50 to 1:500 can be useful to determine the optimal dilution for each primary antibody.
6. Remove the blocking solution from step #4, then add primary antibody dilutions and incubate for 1-2 hours at room temperature.
7. After primary antibody probing, rinse cells three times with PBS.

### Secondary Antibody Labeling

8. Make secondary dilutions in 1% BSA (or normal serum) in PBS.
9. Suggested dilutions for secondary antibodies used at ECM Biosciences:

RS3261	Goat anti-Rabbit Ig:DyLight® 488	(Green; Abs./Em. = 493/518)	1:200
MS3011	Goat anti-Mouse Ig:DyLight® 488	(Green; Abs./Em. = 493/518)	1:200
RS3271	Goat anti-Rabbit Ig:DyLight® 594	(Red; Abs./Em. = 593/618)	1:200
MS3031	Goat anti-Mouse Ig:DyLight® 594	(Red; Abs./Em. = 593/618)	1:200

10. Add secondary antibody to cells for 30 minutes at room temperature.  
**Note:** Fluorescent secondary antibody labeling should be performed in the dark.
11. For long-term storage (months at 4°C), remove PBS and add SlowFade Gold (Invitrogen) to the cells and seal the slides or plates.

DyLight® is a trademark of Thermo Scientific, Inc. and its subsidiaries.

FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.