

— PROTOCOL —

OneMARK 100

Catalog Number	Unit Size	Reactions
DM101-0100	600 ul	

Storage : Store at RT and 4°C up to 6 months. Store at -20°C up to 1 year.

Description

OneMARK 100 with the Novel Green was designed to show virtually uniform spacing over a wide fragment range. The ladder is supplied in a ready-to-use format containing the fluorescent DNA stain and tracking dyes. High quantum yield and excellent stability make the fluorescence dye the ideal fluorophore for DNA staining applications and a superior replacement for the widely used dyes, ethidium bromide or SYBR® Green I.

The OneMARK 100 with the Novel Green was optimized for direct loading onto unstained agarose gels.

The ladders provide highest level of convenience during the routine handling and avoid commonly used gel staining procedures with ethidium bromide or SYBR® Green I.

The OneMARK 100 includes fragments ranging from 100-3,000 base pairs. The 500 and 1,500 base pair bands have increased intensity to serve as reference points. The approximate mass of DNA in each band is provided (0.54 ug per loading) for approximating the mass of DNA in comparably intense samples of similar size.

Application

- No-post-staining procession
- Direct loading onto your agarose gel for analysis

Source

PCR products and double-stranded DNA digested with appropriate restriction enzymes are phenol-extracted and equilibrated to 10 mM Tris-HCl (pH 8.0) and 1 mM EDTA.

Note:

OneMARK 100 is light sensitive and should be stored and protected from light.

Range: 100-3,000 bp


Number of bands: 12

Concentration: 90 ug/ml

Package: 54 ug /600 ul

Recommended Load: 6 ul / well

Containing orange G, xylene cyanol FF as the tracking dyes.



1.5 % TBE agarose gel

The gel was observed with the blue-light transilluminator.