Symphony®

EN Instructions for use | **ES** Instrucciones de uso | **FR** Mode d'emploi

Read all instructions before using this product.

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silicone that touch breastmilk
made without BPA

silicone
sin BPA

silicone
de BPA

Table of Contents

1. Important Safeguards	3-6
2. Product Description	7-11
3. Getting Started	12
4. Cleaning Overview.....	13
5. Wash – Before first use and after each use*	14-15
6. Sanitize – Before first use and once per day*	16-17
7. Tubing & Breast Pump Care	18-19
7A. Tubing care.....	18-19
7B. Cleaning and disinfecting the breast pump	19
8. Powering Your Pump.....	20-21
8A. Card and Cord Protector	20
8B. Powering On	20-21
8C. Changing the Program Card.....	21
8D. Traveling Outside the U.S.	21
9. Assembly of Pump Kit.....	22-24
10. Breast Shield Sizing	25
11. Determine the Program Card	26-27
11A. Determining Program Card by LCD Screen	27
11B. Determining Program Card by Card Image	27
12. Pumping	28-30
12A. Before pumping.....	28
12B. Preparing for pumping.....	29
12C. Finish Pumping.....	30
13. Program Cards	31-35
13A. Symphony® Standard Program Card	31
13B. Symphony® Preemie+™ Program Card.....	32-33
13C. Symphony® PLUS™ Program Card	34-35
14. Storing, Preparing & Feeding Breast Milk.....	36
14A. Storing breast milk.....	36
14B. Preparing breast milk.....	36
14C. Feeding breast milk	36
15. Troubleshooting.....	37-38
16. Error Codes	39
17. Supplemental Information	40
18. Warranty	41
19. Meaning of Symbols.....	42-43
20. EMC / Technical Description	44-49
21. Technical Specifications	50
22. Disposal.....	51

1. Important Safeguards

IMPORTANT

- Plastic bottles and component parts become brittle when frozen and may break when dropped.
- Bottles and component parts may become damaged if mishandled, e.g. dropped, over-tightened, or knocked over.
- Take appropriate care in handling bottles and components.
- Do not use the breast milk if bottles or components become damaged.

Intended Use

The Symphony® Breast Pump is a powered breast pump to be used by lactating women to express and collect milk from their breast.

Product Description

The Symphony is a hospital grade (multi-user), electric breast pump that is capable of single and double pumping. There are three different programs cards available for use with the Symphony: Standard, Preemie+™, and PLUS™. The program cards are utilized to generate the pump suction patterns.

The Standard program card is designed to build and maintain milk production with 2-Phase Expression® technology. The 2-Phase Expression technology mimics the baby's natural sucking rhythm. Short, stimulation cycles in the stimulation phase are followed by longer ones in the expression phase.

The Symphony Preemie+ program card was developed to help mothers of preterm and term infants initiate their milk supply. Preemie+ program card contains two pump programs, Preemie+ and Standard. The Preemie+ program is designed to be used during the immediate postnatal period, to initiate milk production. The Standard program with 2-Phase Expression technology is utilized after initiation of milk production.

Symphony PLUS program card has the same technology as Preemie+, with INITIATE and MAINTAIN as the new interface descriptions being inclusive of the preterm and term population. Research has shown that Initiation Technology™ is beneficial for mothers of both preterm and term infants. The goal of the INITIATE program is to initiate milk production with Medela's Initiation Technology for both preterm and term infants. The MAINTAIN program is designed to build and maintain milk production with the 2-Phase Expression technology.

The warning symbol identifies all instructions that are important to safety. Failure to observe these instructions can lead to injury or damage to the breast pump. When used in conjunction with the following words, the warning symbols stand for:

WARNING Can lead to serious injury or death.

CAUTION Can lead to minor injury.

NOTE Can lead to material damage.

TIP Useful or important information that is not related to safety.

PLEASE SAVE THESE INSTRUCTIONS.

Questions? Visit www.medela.com or call us at 1-800-435-8316.

1. Important Safeguards (cont.)

When using electrical devices, basic safety precautions shall always be followed, including the following:

READ ALL INSTRUCTIONS BEFORE USING THIS PRODUCT.

⚠️ WARNING: To avoid fire, electrocution, or serious burns:

- Always unplug electrical product immediately after use.
- Do not use while bathing or showering.
- Do not place or store product where it can fall or be pulled into a tub or sink.
- Do not submerge pump unit into water or other liquid.
- Do not reach for a product if it has fallen into water or other liquid. Immediately unplug from electrical outlet and contact manufacturer.
- Do not spray or pour liquid directly onto the pump unit.
- Do not leave product unattended when plugged into an electrical outlet.
- Close supervision is necessary when this product is used near children.
- Use the product only for its intended use as described in this manual.
- Do not use attachments not recommended by the manufacturer, as this may cause hazards.
- Never operate this product if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water.
- The Symphony breast pump is not heat resistant. Keep away from radiators and open flames.
- Never use while sleeping or when drowsy.
- Never drop or insert any object into any opening or tubing.
- Do not use outdoors.
- Keep appliance dry and do not place in direct sunlight.
- Do not operate where aerosol (spray) products are being used or where oxygen is being administered.
- Equipment not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or nitrous oxide.
- Separation from the mains is only assured through the disconnection of the power cord from the wall socket outlet.
- Position the device to ensure easy disconnection of the power cord from the supply mains.
- Use Medela original authentic accessories only.
- There are no user-serviceable parts inside the Symphony. DO NOT ATTEMPT TO REPAIR THE Symphony breast pump YOURSELF. No modifications to the device are permitted.
- Always inspect power cord wires prior to use for damage or exposed wire. If damage is found, immediately discontinue use of power cord and call Medela Customer Service at 1-800-435-8316.
- Only use the power cord that comes with the Symphony. Make sure the voltage of the power cord is compatible with the power source.

WARNING: To avoid health risk and reduce the risk of injury:

- Breast pump kits are personal care items intended for use by a single user only and should not be shared between users.
- Inspect the breast pump kit components prior to each use.
 - Do not use the components if they are damaged.
 - Do not use the components if dirt, mold or other contamination is visible.
- Inspect all appropriate pump components before each use.
- Clean and sanitize all parts that come in contact with your breast and breast milk prior to first use.
- Do not modify any breast pump kit components.
- Do not drive while pumping.
- Do not use the Symphony breast pump while operating a motor vehicle.
- Never use while pregnant, as pumping can induce labor.
- Wash all parts that come in contact with your breast and breast milk after every use.
- Do NOT continue pumping for more than 2 consecutive pumping sessions if no results are achieved.
- Do not thaw frozen breast milk in a microwave or in a pan of boiling water. Do not microwave breast milk. Microwaving can cause severe burns to baby's mouth from hot spots that develop in the milk during microwaving. (Microwaving can also change the composition of breast milk.)
- If you are a mother who is infected with Hepatitis B, Hepatitis C or Human Immunodeficiency Virus (HIV), pumping breast milk will not reduce or remove the risk of transmitting the virus to your baby through your breast milk.
- If problems or pain occurs, consult your lactation consultant or physician.

1. Important Safeguards (cont.)


CAUTION: Can lead to minor injury:

- Pump only with the breast pump in an upright position.
- Do NOT shorten tubing between breast shield and pump.
- Make sure the voltage of the device is compatible with the power source.
- Plug the power cord into the breast pump first and then into the wall socket.
- Do NOT use abrasive cleaners/detergents when cleaning breast pump or breast pump parts.
- Never put breast pump in water or in a sterilizer, as you can cause permanent damage to the breast pump.
- Do not attempt to remove the breast shield from your breast while pumping. Turn the breast pump off and break the seal between your breast and breast shield with your finger, then remove breast shield from your breast.
- If pumping is uncomfortable or causing pain, turn the unit off, break the seal between the breast and the breast shield with your finger and remove the breast shield from your breast.
- Do not store wet or damp parts as mold may develop.
- Inspect the tubing for condensation or milk after each use. Wash or replace the tubing immediately after use, if milk is observed in the tubing.
- If tubing becomes moldy, discontinue use and replace tubing.
- Do not use the Medela Symphony breast pump kit with a non-compatible breast pump.
- Do not use the breast pump kit for any other purpose than expressing breast milk.
- Separate and wash all parts that come in contact with breast milk immediately after use to avoid dry up of breast milk residue and prevent growth of bacteria.
- Portable and mobile radio frequency communications equipment can affect the Symphony breast pump.
- The Symphony breast pump needs special precautions regarding EMC (Electromagnetic compatibility) and needs to be installed and put into service according to the EMC information provided in these instructions for use.
- Wash hands thoroughly with soap and water before touching breast pump, kit and breasts and avoid touching the inside of containers or lids.
- Pressing shields too firmly onto breasts may affect milk flow.
- Only use drinking-quality tap or bottled water for cleaning.
- Unplug breast pump before cleaning.
- Contact your healthcare professional or lactation professional if you can only express only minimal or no milk or if expression is painful.
- Do not try and express with vacuum that is too high and uncomfortable (painful). The pain, along with potential breast and nipple trauma, may decrease milk output.

2. Product Description

2. Product Description (cont.)

Authentic Medela Spare Parts are made with the quality materials and workmanship you expect and trust, only from Medela.

Non-Medela spare parts can vary significantly in terms of product design, materials, and workmanship; all of which may affect the performance of your Medela breast pump.

Quality & performance
you can trust

Optional Accessories

Symphony® Hard Case	6007080
Symphony Trolley	6007092
Symphony Vehicle Lighter Adaptor	67173
Symphony Cooler Carrier - Black	37002
Breast Milk Transport Bag™ - NICU	67235-06

Breast Pump Kits

Symphony Double Pumping Kit - Retail	67099-06
Symphony Double Pumping Kit - Retail - (Not Case Qty)	67099
Symphony & Lactina® Double Pumping Kit - Sterile	67316S
Symphony Double Pumping Kit - Sterile	67399S
Symphony Breast Milk Initiation Kit - Sterile	67350S
Symphony Harmony™ Double Duet Breast Pump Kit - Sterile	67115S
Symphony Harmony Double Duet Breast Milk Initiation Kit - Sterile	67355S
Symphony Lactina Breast Milk Initiation Kit - Sterile	67340S
Symphony Single Pumping Kit - Sterile	67026S
Symphony Harmony Duet Breast Pump Kit - Sterile	67114S
Symphony Lactina Double Pumping Kit - Retail	67116-06
Symphony Lactina Double Pumping Kit - Retail - (Not Case Qty)	67116
Symphony & Lactina Carrier Kit	67206-03
Symphony & Lactina Carrier Kit - (Not Case Qty)	67206
Symphony to Lactina Conversion Kit	67090
Lactina to Symphony Conversion Kit	67091
Harmony to Symphony Conversion Kit	67023

2. Product Description (cont.)

Replacement Parts

Symphony® Diaphragm Kit	8007199
Symphony Container Stand	8100552
Symphony Membrane Cap	8100769
Symphony Protective Membrane	8200094
Symphony Power Cord	9280043
Battery Cover	8100551
Slotted Screw	4120017
Rubber Foot	8120025
Fuse	9310035
Trolley Repair Kit	6007116
Strap for Symphony Carry Case	8417009
Symphony Standard Program Card - English	67208
Symphony Standard Program Card - French	67209
Symphony Standard Program Card - Spanish	67210
Symphony Preemie™ Program Card - English	67227
Symphony PLUS™ Program Card - English	67242
Symphony PLUS Program Card - Spanish	2009190
Symphony Program Card/Power Cord Protector	8100589
Symphony Program Card/Power Cord Protector Replacement Screw	4130003
PersonalFit™ Connector	87071
PersonalFit 21mm Breast Shield	87072
PersonalFit 24mm Breast Shield	87073
PersonalFit 27mm Breast Shield	87274
PersonalFit 30mm Breast Shield	87075
PersonalFit 36mm Breast Shield	87084
Replacement White Membranes	87088
Replacement Yellow Valve Heads & White Membranes	87089

Replacement Parts cont.

Symphony PLUS Sticker - English	2009092
Symphony Standard Quick Start Card (cable tie sold separate)	1907703
Symphony Preemie ⁺ Quick Start Card (cable tie sold separate)	1908078
Symphony PLUS Quick Start Card (cable tie sold separate)	1908524
Cable Tie for Quick Start Cards	6997004
Symphony Maintenance & Technical Guide	1907555

3. Getting Started

Note

It's important for you to do the following before using Symphony® for the first time:

1. For models that include a rechargeable battery, fully charge the battery 12 hours uninterrupted. For further battery care information see **Section 8**.
2. Separate all parts
3. Clean – see instructions:
 - Wash – **Section 5**
 - Sanitize – **Section 6**

Breast shield assemblies come assembled. Separate all parts that will come into contact with the breast and breast milk before cleaning.

Parts to clean:

- Breast shields
- Breast milk bottles
- Lids
- Connectors
- Yellow Valve Heads
- White Membranes

4. Cleaning Overview

	Breast pump kit & bottles 	Tubing 	Symphony® Breast Pump and case 
Before 1st use	✓ wash ✓ sanitize		
After each use	✓ wash		
Once per day	✓ sanitize		
As needed		✓ wash only if milk is in tubing or tubing is soiled	✓ wipe with clean, damp cloth

Detailed cleaning instructions on pages 14-19:

- Wash – page 14-15
- Sanitize – page 16-17
- Tubing & Breast Pump Care – pages 18-19

Note

Sterile kits are ready to use and do not require washing prior to first use.

Warning

Breast pump kits are a personal care item intended for use by a single user only and should not be shared between users.

5. Wash – Before first use and after each use*

Supplies needed:

- Mild dish soap
- Clean dish towel or soft brush
- Clean sink or bowl
- Drinking-quality water

Parts to wash:

- Breast shields
- Breast milk bottles
- Lids
- Connectors
- Yellow Valve Heads
- White Membranes

Note

*Sterile kits are ready to use and do not require washing prior to first use.

Warning

Breast pump kits are a personal care item intended for use by a single user only and should not be shared between users.

1

- Separate parts for washing

Caution

Separate and wash all parts that come in contact with breast and breast milk immediately after use to avoid dry up of breast milk residue and prevent growth of bacteria.

Note

- Wash breast pump kit parts after each use.
- Take care not to damage parts of the pump kit during cleaning.
- **It is not necessary to wash tubing, membrane cap or protective membrane prior to first use.**
- **Only wash tubing if milk is present or tubing is soiled.**

4

- Place parts on a clean surface and/or towel.
- Allow all parts to air dry.
- Store dry parts in a clean, cool place when not in use. Do NOT store wet or damp parts.

Note

Store the pump kit in a clean bag/container until next use.

Wash in sink

2

3

Rinse all separated parts that came in contact with breast and breast milk in cool water in order to remove breast milk residue.

- Soak all separated parts in warm, soapy water for 5 minutes.
- Wash each part with a clean dish-cloth or soft brush.
- Rinse all separated parts with clear water.

Caution

Only use drinking-quality tap or bottled water for cleaning.

You may also wash your kit parts by following the instructions on the Quick Clean™ Breast Milk Removal Soap bottle. (sold separately)

Alternatively, wash in dishwasher

1

- Wash all separated parts on top rack of dishwasher.
- Allow all breast pump parts to air dry in a clean area.
- Store dry parts in a clean, cool place when not in use.

Note

If the individual components of the pump kit are cleaned in a dishwasher, parts may be discolored by food pigments. This will not impact part function.

6. Sanitize – Before first use and once per day*

Supplies needed:

- Clean pot for boiling water
- Drinking-quality water
- Tongs

Parts to sanitize:

- Breast shields
- Breast milk bottles
- Lids
- Connectors
- Yellow Valve Heads
- White Membranes

Note

Sanitize breast pump kit parts once daily.

Warning

Breast pump kits are a personal care item intended for use by a single user only and should not be shared between users.

Medela pump kit parts can be sanitized by boiling or using the Medela Quick Clean™ Micro-Steam™ bag. (sold separately).

Note

*Sterile kits are ready to use and do not require washing or sanitizing prior to first use.

1

- Separate all parts that come in contact with breast milk.

4

Allow water to cool and gently remove parts from water with tongs.

2

- Wash hands thoroughly.

3

- Fill a pot with enough water to cover all parts.
- Bring water to a boil.
- Place parts in boiling water for 10 minutes.

5

- Place parts on a clean surface and/or towel.
- Allow all parts to air dry.
- Store dry parts in a clean, cool place when not in use. Do NOT store wet or damp parts.

Note

If you notice a white residue on your parts after boiling, you may have a high mineral content in your water. Remove residue by wiping parts with a clean washcloth and air dry. We recommend that you boil parts with distilled water to prevent substantial mineral build-up over time which may compromise your parts.

Note

Store the pump kit in a clean bag/container until next use.

If sanitizing your pump kit parts with the Quick Clean Micro-Steam bag, separate all parts that come in contact with your breast or breast milk and follow the instructions on the Quick Clean Micro-Steam bag.

7. Tubing & Breast Pump Care

Supplies needed:

- Mild dish soap
- Clean dish towel
- Clean sink or bowl
- Drinking quality water
- Damp (not wet) cloth

Parts needed:

- Tubing
- Symphony® Breast Pump
- Membrane cap
- Protective membrane

7A. Tubing care

- Inspect tubing after each pumping session.

Caution

- Do not store wet or damp parts as mold may develop.
- Inspect the tubing for condensation or milk after each use. Wash or replace the tubing immediately after use if soiled or milk is observed in the tubing.
- If tubing becomes moldy, discontinue use and replace tubing. Replacement parts can be purchased at www.ShopMedela.com.

Note

If there is condensation in the tubing: Shake water droplets out of the tubing and hang it up to dry or let the pump run with tubing attached for 1–2 minutes or until dry.

3

Remove, tubing, protective membrane, and membrane cap. Clean protective membrane and membrane cap by following steps in **Section 5**.

4

Remove tubing from breast shield.

!

Note

- Tubing should be washed if milk is present or tubing is soiled.
- Cleaning tubing is not necessary if condensation present is from previous washings or atmospheric conditions.

Warning

Breast pump kits are a personal care item intended for use by a single user only and should not be shared between users.

7

Use a damp (not wet) cloth to wipe the pump and membrane area.


Turn off breast pump.

Unplug breast pump from power source.

Rinse tubing in cool water to remove breast milk.

- Wash tubing in warm soapy water.
- Rinse tubing with clear water.
- Hang to air dry.
- Or shake out any drops of water from the tubing and let the pump run with tubing attached until dry.

Wipe off with a clean, damp (not wet) cloth. Water and basic dishwashing agent or non-abrasive detergent may be used.

Wipe off with a clean, damp (not wet) cloth. Alcohol based disinfecting agent may be used. Do not immerse the pump in water; do not run water over the pump.

Caution

Unplug breast pump before cleaning.

8. Powering Your Pump

Note

- Only use the pump with the appropriate program card (Standard, Preemie⁺ or PLUS).

Optional for pump with battery:

- Built-in battery must be charged for 12 hours on receipt of the device.
- The pump can be used during charging.
- If the pump is to remain unused, charge the battery every two months (for 12 hours each time).
- For optimum service life of the battery, charge the battery while breast pump is not in use and avoid fully discharging the battery on a regular basis.

Pumping time (fully-charged battery)	approx. 60 min.
Charging time (for full charge)	12 hours

8A. Card and Cord Protector

1

To remove, unscrew the card and cord protector and slide it out of the guide rail.

8B. Powering On

1

Insert the power cord into the power socket on the back of the pump.

2

Insert the cord into the card and cord protector in a loop around the pin.

- Make sure there is enough play to plug the cord in.

3

Push the card and cord protector over the guide rail on the back until it stops and tighten the screw.

8C. Changing the Program Card

2

Plug the power cord into the mains socket.

1

As soon as the card compartment is revealed, the card can be removed and replaced or changed.

Warning

- Only use the power cord that comes with the Symphony®.
- Make sure the voltage of the power cord is compatible with the power source.

8D. Traveling Outside the U.S.

While travelling internationally, we recommend purchasing a Universal Power Plug Adaptor (not included). Please consult with the country you are visiting to find out what type of adaptor will work best. Symphony power cord is two-pronged and does not require a grounded outlet.

For information about traveling with your breast pump and expressed breast milk, visit www.TSA.gov.

9. Assembly of Pump Kit

Tip

Carry out all the steps with care and assemble the pump kit correctly for optimum performance.

Parts Needed:

- Tubing
- Symphony® pump
- Breast shields
- Breast milk bottles
- Lids
- Connectors
- Yellow Valve Heads
- White Membranes
- Membrane cap
- Protective membrane

1

Caution

Wash hands thoroughly with soap and water before touching breast pump, kit and breasts and avoid touching the inside of containers or lids.

4


Note

- Use Authentic Medela Spare Parts only. See **Section 2** for details.
- Check pump kit components for wear or damage before use and replace if necessary.
- All components must be completely dry before use for proper performance and to prevent damage to the pump.

7

Place the protective membrane onto a flat surface and carefully push the membrane cap down until it engages.

2

Carefully push the white membrane onto the yellow valve head until it engages and is laying flat.

3

Push the breast shield onto the connector.

5

Screw the bottle into the connector.

6

Insert the tubing into the corresponding opening in the connector.


8

Insert the tubing adapter into the opening in the membrane cap.

9. Assembly of Pump Kit (cont.)

Open the lid of the breast pump.
→ Press the oval button on the top of the Symphony and lift the handle.

Position the membrane cap firmly on the breast pump membrane to form a seal. Center the tubing connector in the recess in the pump housing.
→ Make sure that the membrane cap engages and lies flush to the surface.

Note

If single pumping seat only one membrane cap.

Center the tubing connector in the recess in the pump housing.
→ Make sure that the membrane cap engages and lies flush to the surface.

Close the lid.

Note

- For step 12:
Always close the lid when pumping. The lid keeps the membrane cap pushed down to prevent a loss of vacuum that could occur if the membrane cap came off.

10. Breast Shield Sizing

PersonalFit™ Sizing

1. Medela's PersonalFit™ sizing

Pumping should not hurt. For maximum comfort and pumping efficiency, Medela offers five breast shield sizes.

21
mm

24
mm

27
mm




30
mm

36
mm

For additional information visit MedelaBreastshields.com.

2. Test your breast shield size

1. Start with the 24 mm that came with your pump, or the size determined by measuring.
2. Center nipple and gently hold breast shield against your breast.
3. Adjust for Maximum Comfort Vacuum™, to achieve optimum suction level.
4. Refer to images while pumping in expression.

3. Should you try a new size?

- Does your nipple rub sides of tunnel, to the point of causing discomfort?
- Do you see excessive areola being pulled into tunnel?
- Do you see any redness?
- Is your nipple or areola turning white?
- Do you feel unexpressed milk after pumping?

If you answered “**YES**” to any of these questions, consider trying a new size by following the measuring instructions above.

If you are still unsure if you selected the correct size, see a lactation consultant, breastfeeding specialist, or visit MedelaBreastshields.com for assistance in choosing the right size breast shield.

Caution

- While some discomfort may be felt when first using a breast pump, using a breast pump should not cause pain. If you are unsure about breast shield sizing, please contact a healthcare professional or breastfeeding specialist who can help you get a proper fit.

11. Determine the Program Card

There are three different program cards available for use with the Symphony®: Standard, Preemie+™, and PLUS™. The program cards are utilized to generate the pump's suction patterns.

The Symphony Standard program card is designed to build and maintain milk production with 2-Phase Expression® technology. The 2-Phase Expression technology mimics the baby's natural sucking rhythm. Short, stimulation cycles in the stimulation phase are followed by longer ones in the expression phase.^{1,2,3}

The Symphony Preemie+ program card was developed to help pump-dependent mothers of preterm and term infants initiate their milk production. Preemie+ program card contains two pump programs, Preemie+ and Standard. The Preemie+ program is designed to be used during the immediate postnatal period. The Standard program with 2-Phase Expression technology is utilized after initiation of milk production.

Symphony PLUS program card has the same technology as Preemie+, with INITIATE and MAINTAIN as the new interface descriptions being inclusive of the preterm and term population. Research has shown that Initiation Technology is beneficial for mothers of both preterm and term infants.^{11,12,13,14}

- The **INITIATE** program supports pump-dependent mothers to successfully initiate milk production. The program mimics the irregular sucking pattern, rapid sucking and pausing, of a term infant in the first days of lactation.¹¹
- The **MAINTAIN** program is designed to optimize milk output after secretory activation (milk coming in) has occurred, in order to build and maintain an adequate milk production. The program has 2-Phase Expression technology that mimics a term infant's sucking pattern during established lactation.^{1,2,3}

11A. Determining Program Card by LCD Screen


Standard Program Card:

The LCD screen will display "Symphony 2.0" when plugged in.

Preemie+ Program Card:

The LCD screen will display "Preemie/Standard" when plugged in.

PLUS Program Card:

The LCD screen will display "Press power" when plugged in.

11B. Determining Program Card by Card Image

Standard Card

Premie+ Card

PLUS Card

12. Pumping

Note

- If the pump has been running for 30 minutes without any user action (e.g. adjustment of the vacuum), the Symphony will automatically switch off.

Tip

- Do not fill bottle above maximum volume indicator line.

12A. Before pumping

Caution

Wash hands thoroughly with soap and water before touching breast pump, kit, and breast. Avoid touching the inside of containers or lids.

Switch on the breast pump with .

Place the second breast shield on the breast so that the nipple is properly **centered** in the tunnel.

12B. Preparing for pumping

Push membrane caps into place on the pump. A "click" will be heard once properly seated. (**Note:** if single pumping, only seat one membrane cap.)

Place the breast shield on your breast so that your nipple is properly **centered** in the tunnel. Hold the breast shield onto your breast with your thumb and index finger. Support your breast with the palm of your hand.

Caution

Pressing shields too firmly onto breasts may affect milk flow.

Set to Maximum Comfort
Vacuum™: Rotate the vacuum to the right and increase until there is slight discomfort. Rotate to the left to reduce the vacuum slightly.

12. Pumping (cont.)

12C. Finish Pumping

Switch off the breast pump with .

- Do not fill the bottle above maximum volume indicator line.
Always disconnect the Symphony® from the power source after expression.

Close the bottle with a lid.

- Follow the instructions in **Section 14** “Storing, preparing and feeding breast milk”.

Use the container stand to prevent the bottle from tipping over.

Caution

- Contact your healthcare professional or lactation professional if you can express only minimal or no milk or if expression is painful. See **Section 10 and 17** for more information.
- Do not try and express with vacuum that is too high and uncomfortable (painful). The pain, along with potential breast and nipple trauma, may decrease milk output.

If you are experiencing low or no suction, refer to Troubleshooting Section 15.

13. Program Cards

13A. Symphony® Standard Program Card

! Note

- Symphony Standard program card is designed to build and maintain milk production.
- Symphony Standard has 2-Phase Expression® technology, which mimics the baby's natural sucking rhythm.^{1,2,3}
- Symphony Standard has two pumping patterns, stimulation and expression.
 - Stimulation Phase:
Research-based fixed cycles per minute^{1,2,3}
 - Expression Phase:
supports milk production^{5,8,9}
- The combination of double pumping with 2-Phase Expression Technology, compared to single pumping, allows mothers the following benefits¹⁻¹⁰:
 - Time saving - more milk in less time
 - Maintenance of milk supply when not able to directly breastfeed
 - Milk with higher energy content (as measured by cream content)
 - comfort and efficiency


T Tip

- Press the "Let down" button  to revert from the expression phase to the stimulation phase if necessary. If milk has not begun to flow after two successive stimulation phases, take a break and try again in 15–30 minutes. Try breast massage (ask your lactation professional or health care practitioners) and then pump using the expression phase.
- The program stores the last vacuum setting in the stimulation phase.

Operation: Standard program

Press On/Off Button  to switch on the breast pump. The program will automatically start in Stimulation.
→ Stimulation phase will automatically change to the expression phase after two minutes.

If your milk starts to flow before the end of the two minutes, press "Let-down button  to change to the expression phase

→ Make sure that the milk flows into the bottle properly.

Set to Maximum Comfort Vacuum™ in both stimulation and expression phases. Refer to Section 12.

13. Program Cards (cont.)

13B. Symphony® Preemie+™ Program Card

Note

- Symphony Preemie+ program card helps initiate, build and maintain breast milk supply.¹¹
- The Preemie+ Program Card contains two pumping programs.
- **Preemie+:**
Program for the immediate postnatal period until initiation of milk production. Then use Standard program. The Preemie+ program runs uninterrupted for 15 minutes and comprises a sequence of various phases. “Done” appears in the display when the program is complete.
- **Standard:**
Use the standard program after secretory activation (milk has come in):
Once the mother has pumped at least 20 mL in total, in each of her last three pumping sessions (from both breasts combined),
OR
Beginning on day six, whichever comes first.

Operation: Preemie+ program

Press On/Off Button to switch on the breast pump. The display text shows “ for Preemie+”.

Operation: Standard program

Press On/Off Button to switch on the breast pump. The Standard program will start automatically with stimulation phase after 10 seconds.


2

Press the “Let-down” button  within 10 seconds. The display text will change to “start Preemie”.

It is important to complete the entire 15 minute Preemie⁺ program.

3

The Preemie⁺ program runs automatically for 15 minutes with stimulation, expression and pauses phases. An auditory signal indicates the end of the program. The display shows “Done”. The pump will automatically turn off.

2

If your milk starts to flow before the end of the two minutes, press “Let-down” button  to change to the expression phase

→ Make sure that the milk flows into the bottle properly.

Set to Maximum Comfort Vacuum™ in both stimulation and expression phases. Refer to Section 12.

13. Program Cards (cont.)

13C. Symphony® PLUS™ Program Card

Note

- Symphony PLUS program card has the same technology as Preemie⁺.
- Symphony PLUS program card helps initiate, build and maintain breast milk supply.¹¹
- Symphony PLUS program card contains two pumping programs, INITIATE and MAINTAIN.
- **INITIATE program:**
Program for the immediate postnatal period until initiation of milk production. Then use MAINTAIN program. The INITIATE program runs uninterrupted for 15 minutes and comprises a sequence of various phases. "Program complete" appears in the display when the program is complete.
- **MAINTAIN program:**
Use the MAINTAIN program after secretory activation (milk has come in):
 - Once the mother has pumped at least 20 mL in total, in each of her last three pumping sessions (from both breasts combined),
OR
 - Beginning on day six, whichever comes first.

Operation: INITIATE program


Press On/Off Button to switch on the breast pump. The display text shows "INITIATE press .

Operation: MAINTAIN program

Press On/Off Button to switch on the breast pump. The MAINTAIN program will start automatically with stimulation phase after 10 seconds.


Press the "Let-down" button within 10 seconds. The display text will change to "INITIATE running".

It is important to complete the entire 15 minute INITIATE program.

The INITIATE program runs automatically for 15 minutes with stimulation, expression and pauses phases. An auditory signal indicates the end of the program. The display shows "Program complete". The pump will automatically turn off.

If your milk starts to flow before the end of the two minutes, press "Let-down" button to change to the expression phase

- Make sure that the milk flows into the bottle properly.

Set to Maximum Comfort Vacuum™ in both stimulation and expression phases. Refer to Section 12.

14. Storing, Preparing & Feeding Breast Milk

14A. Storing breast milk

Freshly Expressed Breastmilk Storage Guidelines <i>(For Healthy Term Babies)</i>			
Room Temperature	Cooler with Ice Pack	Refrigerator	Freezer
4–6 hours at 66–78 °F (19–26 °C)	24 hours at 59 °F (15 °C)	3–8 days at 39 °F or lower (4 °C)	6–12 months 0–4 °F (-18–-20 °C)

References: www.BreastmilkGuidelines.com

14B. Preparing breast milk

Warning

- Do not thaw frozen breast milk in a microwave or in a pan of boiling water.
- Do not microwave breast milk. Microwaving can cause severe burns to baby's mouth from hot spots that develop in the milk during microwaving. (Microwaving can also change the composition of breast milk.)
- Thaw breast milk overnight in the refrigerator. Thawed breast milk is safe in the refrigerator for 24 hours. Do not refreeze thawed breast milk.
- Quickly thaw breast milk by holding the bottle under warm running water.
- Place the sealed bottle in a bowl of warm water for 20 minutes to bring it to body temperature.
- If you are adding expressed breast milk to a container of already frozen breast milk, make sure to add a lesser amount than the already frozen amount.

14C. Feeding breast milk

It is recommended that breastfeeding is well established prior to bottle feeding your baby.

- Always inspect the bottle, nipple and other components immediately before and after each use. If nipple appears cracked or torn, discontinue use immediately.
- To prevent possible choking hazard, test strength of nipple by pulling on bulb portion of the nipple.
- Do not attempt to enlarge the nipple hole.
- Infants must not be bottle fed without adult supervision.
- Nipple should not be used as a pacifier.

15. Troubleshooting

If the breast pump is not running	<p>Check whether electrical outlet power/battery is available.</p> <p>Check that the device is switched on.</p> <p>Check whether program card is inserted correctly in the card compartment on the back of the breast pump.</p>
If there is low or no suction	<p>Assemble the pump set in accordance with the instructions in Section 9.</p> <p>Check whether all connections of the pump set and breast pump are secure.</p> <p>Pump in accordance with Section 12. Make sure that the breast shield forms a complete seal around the breast.</p> <p>The white membrane must lie flat against the yellow valve head.</p> <p>Check that the yellow valve head and white membrane are clean and not damaged.</p> <p>See Section 9 "Assembly of Pump Kit".</p> <ul style="list-style-type: none">• The ends of the tubing must be firmly connected to the connector and the opening in the membrane cap.• The protective membrane and membrane cap should be correctly assembled and the membrane cap must be firmly fitted to the breast pump membrane to form a seal.• The protective membrane must be intact.• Do not kink or pinch the tubing during expression.• All components must be completely dry. <p>Yellow Valve Head and Membrane Check</p> <ul style="list-style-type: none">• Remove and separate yellow valve head from the white membrane.• Inspect the yellow valve head for chips or cracks. Inspect the membrane for holes or tears. If flawed or damaged, replace the yellow valve head and/or membrane prior to pumping.• Wash yellow valve head and membrane per Cleaning Instructions in Section 5.• Reassemble. Make sure the white membrane lays flat on the yellow valve head, and if it does not, replace the yellow valve head and/or membrane prior to pumping. <p>Breast Shield and Breast Shield Connector Check</p> <ul style="list-style-type: none">• Inspect breast shield and breast shield connector to make sure there are no cracks and they securely attach to each other.• Make sure the breast shield forms a complete seal around the breast. <p>Other Checks</p> <p>Inspect the following to check for damage and secure connections.</p> <ul style="list-style-type: none">• Secure tubing adapter to the breast shield connector.• Test vacuum. If suction is not improved, contact Medela Customer Service at 1-800-435-8316. <p>Tubing Check</p> <ul style="list-style-type: none">• If damaged replace prior to pumping.• Inspect tubing for damage and secure connections.• Check that the ends of the tubing are firmly connected to the connector and the opening in the membrane cap.• If single pumping, make sure only one membrane cap is seated. <p><i>Take care not to kink or pinch tubing.</i></p>

15. Troubleshooting (cont.)

If there is an overflow	<ol style="list-style-type: none">1. Place the pump set into the bottle holder.2. Turn off the pump and unplug the power cord from the electrical outlet.3. Use a damp (not wet) cloth to wipe the pump and membrane area.4. If there is milk in the tubing: Remove tubing, protective membrane and membrane cap and clean in accordance with Section 7. Shake out any drops of water from the tubing and hang up the tubing to dry or let the pump run with tubing attached until dry. If there is condensation in the tubing: Let the pump run with tubing attached until dry.5. When components are completely dry, reassemble breast pump and accessories set.
-------------------------	---

Condensation or Milk in Tubing	<ul style="list-style-type: none">• Turn off the pump and unplug from power source.• Use a damp (not wet) cloth to wipe the pump.• Refer to cleaning instructions "Tubing Cleaning" Section 7.• If necessary, clean the pump kit.• When components are completely dry, reassemble breast pump and kit.
No "Let-Down" or Breast milk Expression	<ul style="list-style-type: none">• Ensure your breast pump kit is assembled correctly and suction is available.• Relax and take a 10-15 minute break, if let down is not achieved after 2 consecutive pumping sessions.• Consult with a healthcare or lactation professional if expression does not occur.

If you have not resolved the problem with your breast pump or you have further questions, please contact Medela Customer Service at 1-800-435-8316 or fill out the contact us form at www.medelabreastfeedingus.com.

16. Error Codes

There are two ways in which a Symphony® breast pump will communicate that the unit has an issue. If the Symphony is experiencing an issue, an error code will be displayed on the LCD screen (see fig. 1). Previous problems are recorded in the Menu function under “Err:”. To retrieve a history of errors, hold down the “On” button while turning the vacuum control knob, then continue to move the knob until you reach “Err” (see fig. 2). If you continue to experience issues with your Symphony, please contact Customer Service at 1-800-435-8316.

LCD Display Messages

The following describes LCD messages that may be displayed indicating an issue with the unit, the potential causes, and action steps to correct the issue.

	Potential Issue	Actions
Invalid card	Program card isn't a Symphony card	Check if correct program card is used
	Program card is incorrectly inserted	Check that card is positioned correctly
	Program card is damaged and not working	Replace program card
No card	Program card isn't inserted	Insert program card
	Program card isn't completely inserted	Push card in all the way
Battery low	Battery needs charging soon, will beep every 20 seconds (15 minutes of pumping time remaining from when the indication started)	Be sure power cord is plugged into back of pump
Battery flat	Battery is completely discharged	Charge the batteries, but you may continue to pump
Motor overload	Motor is drawing too much current	Send Symphony pump to Medela service center if this happens repeatedly*
Configuration	Software configuration	Replace program card if experiencing error (See <i>Maintenance and Technical Guide 1907555</i>)
Error	Potential fault in electronics of gear control	Send Symphony pump to Medela service center if this happens repeatedly*
Power fault	Potential fault in electronics of power supply	Send Symphony pump to Medela service center if this happens repeatedly*

*All returns require prior authorization, including rental equipment. For the Return/Exchange/Repair form, please visit the Client Accounts page at www.medelabreastfeedingus.com – login, then click on Account Communication. Or, contact customer service at accountsupport@medela.com for the form. You will receive return instructions and a Return Authorization Number, which must be clearly written on the outside of the package in order to process your return.

17. Supplemental Information

Supplemental Information

Following are some common breastfeeding related conditions. If you experience any of these symptoms, contact a healthcare professional or lactation professional.

	Symptom	Potential Causes
Engorgement	Breasts are hard and uncomfortable, possibly with reddened areas.	Milk not draining properly from breast. Milk can drain back into the tissue causing areas to swell and become tender.
Blocked/Plugged Ducts	Area of breast will look reddened and may be tender to touch.	Milk not being drained from a specific duct. The area becomes "clogged" up and milk is then prevented from flowing.
Mastitis	It often follows engorgement but can just occur suddenly. The initial symptoms are similar to flu – fatigue, headache and muscle aches, fever and localized breast tenderness. It can be in one or even both breasts and needs immediate treatment.	A bacterial infection in the breast tissue which is often accompanied by cracked nipples.

A Healthcare Professional or Lactation Professional can give you guidance in breastfeeding your baby and using a breast pump. Their guidance is valuable in terms of the lifetime of health benefits for both you and your baby. A listing of lactation professionals are available at www.medelalocator.com or calling 1-800-TELL-YOU (1-800-835-5968).

18. Warranty

This product is warranted by Medela LLC to the original retail purchaser to be free from defects in material and workmanship for the period of three (3) years from the date of purchase. In the event of a defect or failure to conform to this warranty, Medela will, at Medela's option, repair or replace this product without charge for such replacement, on parts or labor, at the Medela office listed below. The Purchaser shall bear all responsibility and expense for returning this product to Medela, including risk of loss prior to receipt by Medela, and shipping, packing, and insurance costs.

Medela will have no obligation under this warranty to repair any malfunction or damage arising from any misuse, abuse, or alteration of this product. Without limiting the generality of the foregoing, bending, or dropping of this product or its components, visible cracking of the equipment housing will be presumed to be defects resulting from misuse or abuse. Accessories that are not expressly manufactured by Medela and subsequently attached to the equipment will immediately void this warranty.

ANY AND ALL IMPLIED WARRANTIES ARE LIMITED TO THE DURATION OF THREE YEARS FROM DATE OF PURCHASE. SOME STATES DO NOT ALLOW LIMITATIONS ON HOW LONG AN IMPLIED WARRANTY LASTS, SO THE ABOVE LIMITATIONS MAY NOT APPLY TO YOU.

THE LIABILITY OF MEDELA FOR BREACH OF ANY WRITTEN OR IMPLIED WARRANTY IS LIMITED TO REPAIR OR REPLACEMENT OF THIS PRODUCT, AND MEDELA WILL HAVE NO LIABILITY UNDER ANY CIRCUMSTANCES FOR ANY INCIDENTAL OR CONSEQUENTIAL DAMAGES. SOME STATES DO NOT ALLOW THE EXCLUSION OR LIMITATION OF INCIDENTAL OR CONSEQUENTIAL DAMAGES, SO THE ABOVE LIMITATIONS OR EXCLUSIONS MAY NOT APPLY TO YOU.

THIS WARRANTY GIVES YOU SPECIFIC LEGAL RIGHTS AND YOU MAY ALSO HAVE OTHER RIGHTS WHICH VARY FROM STATE TO STATE.

If you wish to make a claim under this warranty, you must return this product to Medela at the address below. You must contact Medela Customer Service and obtain a Return to Merchant Authorization (RMA) Number and reference the number on your return package. The original Invoice Number or Purchase Order Number must be referenced on the RMA.

Customer Service toll free number: 1-800-435-8316.

Medela LLC – Returns, Door 4501
1101 Corporate Dr.
McHenry, IL 60050
ATTENTION: RETURNS

Medela Canada Inc.
4160 Sladeview Cres., Unit #8
Mississauga, ON
Canada L5L 0A1

Call first for authorization number. Returns not accepted without an authorization number.

19. Meaning of Symbols

The warning symbol identifies all instructions that are important for safety. Failure to observe these instructions can lead to injury or damage to the breast pump. When used in conjunction with the following words, the warning symbol stands for:

Symbols on the packaging

This symbol indicates that the material is part of a recovery/recycling process.

This symbol indicates a carton package.

This symbol indicates to keep the device away from sunlight.

This symbol indicates to handle the fragile device with care.

This symbol indicates the temperature limitation for operation, transport and storage.

This symbol indicates the humidity limitation for operation, transport and storage.

This symbol indicates the atmospheric pressure limitation for operation, transport and storage.

This symbol indicates to keep the device dry.

This symbol indicates that the packaging contains products intended to come in contact with food according to regulation 1935/2004.

This symbol indicates unique GS1 Global Trade Item Numbers (GTIN).

Symbols on the device

-  This symbol indicates to follow instructions for use.
-  This symbol indicates the compliance with the essential requirements of the Council Directive 93/42/EEC of 14 June 1993 concerning medical devices.
-  This symbol indicates the manufacturer.
-  This symbol indicates the compliance with additional USA and Canada safety requirements for medical electrical equipment.
-  This symbol indicates a type BF applied part.
-  This symbol indicates manufacturer's serial number of the device.
-  This symbol indicates manufacturer's catalog number of the device.
-  This symbol indicates the protection against ingress of solid foreign objects and against harmful effects due to the ingress of water.
-  This symbol indicates the date of manufacture (four digits for the year and two digits for the month).
-  This symbol indicates that the Symphony is a class II device.
-  This symbol indicates fuse (to identify fuse boxes or their location).
-  This symbol indicates that the fuse box is safety tested.
-  This symbol indicates open the cover of the mains socket.
-  Let-down button
-  On/off button.
-  Vacuum adjustment knob.
-  Indication to press letdown button for Preemie+™ and PLUS program
-  Vacuum level indicator/Stimulation Phase - in Standard and Preemie+ program, and PLUS Maintain Program
-  Vacuum level indicator/Expression Phase - in Standard program and PLUS Maintain Program
-  Vacuum level indication/Expression Phase - in Preemie+ and throughout PLUS Initiate Program

20. EMC/Technical Description

The Symphony® breast pump needs special precautions regarding EMC (Electromagnetic compatibility) and needs to be installed and put into service according to the EMC information provided in these instructions for use.

Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations and walkie-talkies can affect the Symphony electric breast pump and should be kept at least a distance 1.0 m away from the equipment.

Electromagnetic Compatibility (EMC, IEC 60601-1-2:2007, Table 1)

Electromagnetic emissions

The Symphony electric breast pump is intended for use in the electromagnetic environment specified below. The customer or the user of the Symphony electric breast pump should assure that it is used in such an environment.

Emission tests	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Symphony electric breast pump uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	The Symphony electric breast pump is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Warning

Warning – The Symphony electric breast pump should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the Symphony electric breast pump should be observed to verify normal operation in the configuration in which it will be used.

Electromagnetic Compatibility (EMC, IEC 60601-1-2:2007, Table 2)

Electromagnetic immunity

The Symphony electric breastpump is intended for use in the electromagnetic environment specified below. The customer or the user of the Symphony electric breastpump should assure that it is used in such an environment.

The Symphony breastpump has no Essential Performance but was tested for immunity to electromagnetic disturbances and passed using the following criteria:

- Change of operation mode except reset into standby mode (pump is switched off)
- Change of pressure value more than 20 mmHg (the specified tolerance), when the EUT is operated as intended
- Change in programmable parameters, including mode and suction level
- Error codes, including the following:
 - Error indicated; no further restriction in device performance
 - Error indicated, the pump turns off but is restartable
 - No error indicated; the pump turns off and is not restartable
- Permanent hardware failure is induced in the EUT

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input / output lines	± 2 kV for power supply lines ± 1 kV for input / output lines	Mains Power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV differential mode ± 2 kV common mode	Mains Power quality should be that of a typical commercial or hospital environment.

20. EMC/Technical Description (cont.)

Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	$<5\% U_T$ $(>95\% \text{ dip in } U_T)$ for 0,5 cycle $40\% U_T$ $(60\% \text{ dip in } U_T)$ for 5 cycles $70\% U_T$ $(30\% \text{ dip in } U_T)$ for 25 cycles $<5\% U_T$ $(>95\% \text{ dip in } U_T)$ for 5 s	$<5\% U_T$ $(>95\% \text{ dip in } U_T)$ for 0,5 cycle $40\% U_T$ $(60\% \text{ dip in } U_T)$ for 5 cycles $70\% U_T$ $(30\% \text{ dip in } U_T)$ for 25 cycles $<5\% U_T$ $(>95\% \text{ dip in } U_T)$ for 5 s	<p>Mains Power quality should be that of a typical commercial or hospital environment.</p> <p>If the user of the Symphony electric breast pump requires continued operation during power interruptions, it is recommended that the Symphony electric breast pump is powered from an uninterruptible power supply or a battery.</p>
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE U_T is the a.c. voltage prior to application of the test level.

Electromagnetic Compatibility (EMC, IEC 60601-1-2:2007, Table 4)

Electromagnetic immunity			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	10 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the Symphony electric breast pump, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 0,35\sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz	3 V/m	$d = 1,2\sqrt{P}$ 80 MHz to 800 MHz $d = 2,3\sqrt{P}$ 800 MHz to 2,5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 

20. EMC/Technical Description (cont.)

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Symphony electric breast pump is used exceeds the applicable RF compliance level above, the Symphony electric breast pump should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the Symphony electric breast pump.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10 V/m.

Electromagnetic Compatibility (EMC, IEC 60601-1-2:2007, Table 6)

Recommended separation distances between portable and mobile RF communications equipment and the electric breast pump Symphony

The Symphony electric breast pump is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Symphony electric breast pump can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Symphony electric breast pump as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter M		
	150 kHz to 80 MHz $d = 0,35\sqrt{P}$	80 MHz to 800 MHz $d = 1,2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2,3\sqrt{P}$
0,01	0,04	0,12	0,23
0,1	0,11	0,37	0,74
1	0,35	1,2	2,3
10	1,11	3,7	7,4
100	3,50	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

21. Technical Specifications

vacuum (approx.)
-50 ... -250 mmHg
-7 ... -33 kPa
45 ... 120 cpm

100–240V~
50/60Hz
0.5A

T 0.8A, 250V
(slow blow)
5 x 20mm L = Low Breaking Capacity

Car adapter
T 2.5 A
(slow blow)

2x6V, 1200mAh
Pb (lead acid)
Yuasa NP 1.2-6

900 410 410 mm
↓ ↔ ↗
245 215 155 mm

3.5kg
2.9kg

Operation

Transport / Storage

Operation

Transport / Storage

Operation / Transport / Storage

22. Disposal

Disposal

The unit is made of various metal and plastics. Before disposal, the device is to be rendered unusable and it must not be disposed of as unsorted municipal waste in accordance with local regulations. Use your local return and collection system for waste electrical and electronic equipment (incl. batteries). Improper disposal may have harmful effects on the environment and on public health.

- (1) Wolff,P.H. The serial organization of sucking in the young infant. *Pediatrics* 42, 943-956 (1968).
- (2) Woolridge,M.W. The 'anatomy' of infant sucking. *Midwifery* 2, 164-171 (1986).
- (3) Kent,J.C., Ramsay,D.T., Doherty,D., Larsson,M., & Hartmann,P.E. Response of breasts to different stimulation patterns of an electric breast pump. *J Hum Lact* 19, 179-186 (2003).
- (4) Meier,P.P. et al. A comparison of the efficiency, efficacy, comfort, and convenience of two hospital-grade electric breast pumps for mothers of very low birthweight infants. *Breastfeed Med* 3, 141-150 (2008).
- (5) Kent,J.C. et al. Importance of vacuum for breastmilk expression. *Breastfeed Med* 3, 11-19 (2008).
- (6) Prime,D.K., Kent,J.C., Hepworth,A.R., Trengove,N.J., & Hartmann,P.E. Dynamics of Milk Removal During Simultaneous Breast Expression in Women. *Breastfeed Med* 7, 100-106 (2011).
- (7) Mitoulas,L., Lai,C.T., Gurrin,L.C., Larsson,M., & Hartmann,P.E. Effect of vacuum profile on breast milk expression using an electric breast pump. *J Hum Lact* 18, 353-360 (2002).
- (8) Burton, P., Kennedy, K., Ahluwalia, JS., Nicholl, R., Lucas, A., and Fewtrell, M.S. Randomized Trial Comparing the Effectiveness of 2 Electric Breast Pumps in the NICU. *J Hum Lact* published online 17 June 2013
- (9) Kent,J.C. et al. Volume and frequency of breastfeeds and fat content of breastmilk throughout the day. *Pediatrics* 117, e387-e395 (2006).
- (10) Prime,D.K., Garbin,C.P., Hartmann,P.E., & Kent,J.C. Simultaneous Breast Expression in Breastfeeding Women Is More Efficacious than Sequential Breast Expression. *Breastfeed Med* 7, 442-447 (2012).
- (11) Meier PP, Engstrom JL, Janes JE, Jegier BJ, Loera F. Breast pump suction patterns that mimic the human infant during breastfeeding: greater milk output in less time spent pumping for breast pump-dependent mothers with premature infants. *Journal of Perinatology* 2012;32 (2):103-10
- (12) Lawrence RA, Lawrence RM. *Breastfeeding: A guide for the medical profession*, 7th Edition. Elsevier Health Sciences 2010
- (13) Cregan MD, de Mello TR, Kershaw D, McDougall K, Hartmann PE. Initiation of lactation in women after preterm delivery. *Acta Obstetricia et Gynecologica Scandinavica* 2002; 81: 870-77
- (14) Post EDM, Stam, G. More milk with irregular breastpump suction (Abstract A14), *Early Human Development* 2013; 89 (S1): S79

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Medela wordmark and logo, Symphony, Lactina, Invented by Medela 2-Phase Expression and 2-Phase Expression are registered in the U.S. Patent and Trademark Office.
Plus, Initiation Technology, PersonalFit, Preemie+, Harmony, Quick Clean, Micro-Steam and Maximum Comfort Vacuum are trademarks of Medela.
Medela AG
Assembled by Medela LLC