INSTALLATION MANUAL • INSTALLATIEHANDLEIDING • INSTALLATIONSANLEITUNG  
MANUEL D'INSTALLATION • MANUAL DE INSTALACIÓN • MANUALE DI INSTALLAZIONE

# HANDICARE 2000 SIMPLICITY / STYLE

# Handicare 2000 Simplicity / Style installation manual

---

## Pre-installation check list

1 Check that all of the required components are available.

- |                        |                      |
|------------------------|----------------------|
| <b>a</b> Track Drawing | <b>f</b> Seat        |
| <b>b</b> Track         | <b>g</b> Fitting kit |
| <b>c</b> Stanchions    | <b>h</b> Footplate   |
| <b>d</b> Feet          | accessories kit      |
| <b>e</b> Power pack    |                      |

**Note:** Before assembly you must ensure that all of the required track is to the correct specification, the stanchions supplied are the correct length, the correct number have been supplied and the corresponding number of feet are available. In order to do this you need to cross reference the items supplied to the drawing supplied. Please refer to **Appendix 4** for further details.

## 2 Installation tools required

- Loading toggle
- Loading bars
- Shorting links
- 8mm Spanner
- 13mm Spanner
- 17mm Spanner
- No. 2 Pozidrive
- No. 3 Pozidrive
- Nylon electricians pull through – at least the length of the track
- Track Lubricant: HTEP grease (food safe grease)
- 2.5mm Allen key
- 3mm Allen key
- 4mm Allen key
- 5mm Allen key
- 6mm Allen key (ball ended)
- Superglue (inside track only)
- Torque wrench
- ½" in-hex socket driver – 8 x 120
- ½" Ratcheting socket wrench
- 2 x 250 mm socket extension bar
- 1 x 125mm socket extension bar
- T30 Torx driver
- 13mm podger spanner
- ¾" Sliding 'T' Bar

**Note:** The installation engineer must check with Handicare for the latest installation manual issue number.

Issue No of this manual is 0416-1  
ORIGINAL LANGUAGE ENGLISH

# Track assembly

## Items required:

- Track drawing
- Fitting kit
- Track
- Feet
- Stanchions
- Red cable
- Quantity of cable ties – 142 x 3.2mm and 300 x 4.8mm

## Tools required:

- 6mm Allen key
- 17mm Spanner
- Nylon electricians pull through
- 1/2" in-hex socket driver – 8 x 120
- 1/2" Ratcheting socket wrench
- 2 x 250 mm socket extension bar
- 1 x 125mm socket extension bar
- T30 Torx driver
- 13mm podger spanner
- 3/8" sliding 'T' bar

- 1 Study the drawing provided. This will be similar to the one shown below (**Figure 1**). See **Appendix 4** for a larger version of this example.


Figure 1

## Track assembly

### Internal curve tracks

An inside curve track is one that runs around the shortest side of the staircase such that the track bends away from you when you look at the stairs. An example of an inside curve track can be seen in **Appendix 5**.

The installation of an inside curve and outside curve track is identical until very near to completion of the installation. The additional steps that are required are featured in the corresponding part of the installation manual (**page 17**).

### Hinge track installations

A hinge track is a folding section of track that can be moved away from the bottom of the stairs so that it does not cause an obstruction or hazard.

**Note:** See **Appendix 3** for hinge track settings.

- 1** The hinge track section must be installed before the rest of the track.
  - a** Position the hinge track section in to the feet with the back of each foot 60 mm from the stringer or line of any obstruction, e.g. window sill (**Figure 2**).
  - b** Release the motor by pressing on the motor release to allow the arm of the hinge to be moved.
  - c** Lower the arm and position the hinge on the treads correctly so that the track reaches the floor.
  - d** The grey wire with a bullet connector must be connected to all of the communication contacts on the top track. In order to achieve this connection you should:
    - i** Pull some grey cable from the hinge to the first comms contact in the TOP tube.
    - ii** In the BOTTOM tube pull a red and a grey cable from the first +ve charge contact to the top of the stairs (**Figure 3**). The connection of these will be covered in the Track Electrical Connection section.

**Figure 2**


**Figure 3**

- e Follow the normal track installation procedure from step 2 to complete the installation.
- f Once the track installation is complete cut any excess threaded rod protruding from the top of the foot to prevent it stopping the lift early.

**Note:** If the top communications contacts are not connected to the hinge section the lift will not work.

2 In the drawing shown in **Appendix 1** the first section is actually a straight section with a gradient decrease section manufactured as a single piece of track. The straight section (2) has 152 teeth (152 S) and it is 1671mm long. The gradient decrease section (4) has 28 teeth (28 GD R665) and it is bent on a 665mm radius. The protective packaging that the track is delivered in will be marked to indicate the track section or sections that it contains and whether or not it is a top track or bottom track section (**Figure 4**).

**Note:** Do not remove the blue rack protection cover at this stage.

**Figure 4**

**Figure 5**

ED	REV	COMMENTS	BY	DATE											
	<b>A</b>	<b>1st Issue</b>													
General Table															
rtc-89955 TRACK DETAILS															
W/S	No.	SG	Joins	Teeth	TYPE	LENGTH	SB	No.	UF	RAIL	TOOTH	DISTANCE	LENGTH	TYPE	SETTING HT
	1				GROUND FLOOR			16	2	1	11	395	A	74	
A	2	B-W	28	S		308		17	2	24	264	470	A	149	
	3	W-W	48	HR17 R159			N	18	4	4		649	B	127	
	4	F W-X	37	GI R596				19	5	26	286	693	C	171	
	5	L X-C	44	S		484		20	6	5	55	697	C	175	
B	6	C-W	10	S		110		21	7	41		708	C	186	
	7	W-W	54	HR54 R223			N	22	9	16	176	649	B	127	
	8	F W-X	13	GD R665				23	9	73	803	661	C	139	
	9	L X-C	158	S		1737		24	9	129	1418	662	C	140	
C	10	F C-X	12	S		132		25	12	6	66	709	C	187	
	11	L X-W	13	GI R596				26	12	48	528	884	D	31	
	12	W-W	70	S		770		27	14	15		685	C	163	
	13	W-W	54	HR54 R223			N	28	15	14	154	620	B	98	
	14	F W-X	24	GD R665							TOTAL =	8482			
	15	L X-B	15	S		165									
						TOTAL =	6377								

**No.** – The foot number.

**Rail** – The section of rail that the stanchion fixes to.

**Tooth** – The corresponding tooth on the section of track that the stanchion should be positioned.

**Distance** – The distance of the track joint that the stanchion is positioned.

**Length** – The length of the stanchion.

**Type** – The type of stanchion required. There are four versions.

**Setting height** – using the size guide provided on each stanchion. This is the size that the stanchion needs to be set too.

### Stanchion information


- 1 The drawing shown in **Figure 5** has a table which shows the length of the stanchion that has to be inserted at each position. There are four different types of stanchion which will be adjusted to the correct height. To ease identification these have been labelled with a letter, A,B,C or D.
- 2 Each stanchion can be adjusted to the required height using the ‘setting height’ guide.
- 3 Grease the stanchion lugs as you go.

## Track assembly

### Installing the stanchions/track

**Note:** Hand ratchet the bolts – do not use power tools.

<b>Stanchion lock bolt</b>	M10 x 50 socket cap head screw	45Nm
<b>Stanchion/rail bolt</b>	M10 x 55 button head screw	15Nm
<b>Rail Joint bolt</b>	M8 x 55	25Nm
<b>Rail bolts at the end of the track – x 2 (with plastic cap)</b>	M10 x 55 button head screws	10Nm

### First section

- 1 Set the furthest bottom stanchion at the correct height. (Example:- Stanchion 16 – Type A – Setting height = 74mm)
- 2 Set the last section before the first joint to the correct height (**Figure 6**). Tighten the lock bolt to 45Nm (**Figure 7**).
- 3 Attach the top rail on to the two stanchions. Use stanchion/rail bolts supplied (**Figure 8**).

**Note:** See note for specified torques.

- 4 Attach the bottom rail.
- 5 Fasten the feet down to provide additional stability for the rest of track.

### Second section

- 1 Set the next stanchion at the correct height ready for the second section of rail. This needs to be placed just before the joint.
- 2 Add the next sections of rail. Top section first and then the bottom section. Use a podger spanner tool to help ease the two sections of rail together (**Figure 9**).
- 3 Set the next stanchion at the correct height two down from the previous one nearest the joint.


Figure 8


Figure 9


Figure 10

### Third section

- 1 Set the next stanchion one down from the last stanchion at the top of the stairs.
- 2 Add the top rail section.
- 3 Set the next stanchion from the one above.
- 4 Add the bottom rail section.
- 5 Install all of the remaining stanchions without adjusting the height.
- 6 Add all of the track bolts. These should be hand-tightened to a torque of 25Nm (**Figure 10**).
- 7 Check all stanchion heights and tighten the stanchions (**Figure 11**).

**Note:** For multi-flight stairs push cable through each separate section of rail.

- 8 Fasten some of the screws to the stanchion feet, typically two per metre, to hold the rail captive, allowing the power pack to be installed and checked.

**Note:** Care must be taken to ensure the rail is secure enough to carry out the installation and checks; add more screws if required.


Figure 11

# Track electrical connections

### Items required:

- Track drawing
- Fitting kit
- Track
- Feet
- Stanchions

### Tools required:

- 6mm Allen key
- 3mm Allen key
- 2.5mm Allen key
- 13mm Spanner

- 1** Use the pull through that you fitted inside the track to pull the red cable through the track **(Figure 12)**.
- 2** Decide if the power supply will be installed at the top or bottom of the staircase.
- 3** Remove the end caps from the bottom track at the top and bottom of the stairs **(Figure 13)**.
- 4** Connect the positive charge contact to the cable that you pulled through the track at the top and bottom of the stairs using the supplied terminal blocks.

## Hinge track installations

Each stanchion next to a charge comms contact should be drilled **(Figure 14)**. For the first occurrence of charge comms contacts you must:

- a** Connect the grey cable from the hinge track to the grey cable that you have pulled through the top track.
- b** Insulate this connection and push it through the hole in the front of the track.
- c** Feed the grey cable from behind the comms contact out of the track and into the stanchion and out of the top of the stanchion.
- d** Connect the comms contact to the grey wire in the top track using a connector block.
- e** Connect some more grey cable to this connector block and feed it into the stanchion and out of the top of the stanchion.
- f** Insulate the connector block behind the comms contact and then push it into the track.

- g** Connect the two grey cables that are protruding from the top of the stanchion with a connector block, insulate and then push in to the top of the stanchion.

### For every further intermediate comms contact:

- a** Pull the grey wire out of the bottom track behind the charge contact.
- b** Splice a length of grey cable from this and feed it in to the stanchion and out of the top of the stanchion.
- c** Feed the grey wire from the comms contact in to the stanchion and out of the top of the stanchion.
- d** Connect the two grey cables that are protruding from the top of the stanchion with a connector block, insulate and then push in to the top of the stanchion.

### For the top comms contact:

- a** Pull the grey wire out of the bottom track behind the charge contact.
  - b** Feed the grey wire into the stanchion and out of the top of the stanchion.
  - c** Feed the grey wire from the comms contact into the stanchion and out of the top of the stanchion.
  - d** Connect the two grey cables that are protruding from the top of the stanchion with a connector block, insulate and then push in to the top of the stanchion.
- 5** At the end of the track that you will be connecting the power supply to continue the installation as follows:
    - a** Remove top or bottom stanchion from rail.
    - b** Remove the paint from the stanchion and rail **(Figures 15 and 16)**.
    - c** Insert the negative charge connection ring terminal between the bottom stanchion lug and bottom rail section **(Figure 17)**. Insert stanchion bolt to hold in place.


Figure 12


Figure 15

Figure 13

Figure 16

Figure 14




Figure 17

## Track electrical connections

- d** Feed the grey negative charge connection loom into the large hole at the rear of the track section. Feed this wire through to the end of the rail section. **(Figure 18)**.
  - e** Feed the power supply cable through the same hole to the end of the rail section **(Figure 19)**.
  - f** Place a cable tie around the power supply cable **(Figure 20)**.
  - g** Pull the power supply cable tight at the rear of the rail section. Using a cable tie, attached the power supply cable to the stanchion lug **(Figure 21)**.
  - h** Using two cable ties, securely fasten the power supply cable to the back of the stanchion as shown in **(Figure 22)**
  - i** Connect the power supply cable and negative charge connection loom to the two-connector block **(Figure 23)**.
  - j** Wrap the two red charge contact wires together and feed this into the two-way connector block **(Figure 24)**.
  - k** Feed all of the wires into the end of the rail section **(Figure 25)**.
  - l** Replace the stanchion bolt **(Figure 26)**.
  - m** Replace the end cap at the **BOTTOM OF THE STAIRS ONLY** **(Figure 26)**.
- 6** Fit the remaining bolts in to the top and bottom track sections for the first and last stanchion.
- 7** Fit the power supply to the wall using the supplied bracket.
- 8** Terminate the power supply to the available AC supply.

**Note:** All mains wiring must comply with the local standards in force at the time of installation. Do not fix the power supply to the stringer.


Figure 21


Figure 24


Figure 22


Figure 25


Figure 23

Figure 26

# Loading the power pack and fitting the footplate

### Items required:

- Power pack
- Footplate accessories kit

### Tools required:

- Loading toggle
- Shorting links
- Loading bars
- 13mm socket
- No. 2 Pozi-drive
- Pliers
- T30 Torx driver

- 1 Fit the power pack loading bars to the top of the track (**Figure 27**).
- 2 Feed the power pack on to the loading bars (**Figure 28**).
- 3 Remove trunnion loading strip (**Figure 29**).
- 4 Take off the front cover from the power unit.
- 5 Assemble the footplate (footplate accessories kit)

## Manual (gas strut) footplate

- a Remove the footrest assembly from the footrest kit box. Loosen the screw and remove the shaft from this assembly.
- b Attach the footrest to the power pack (**Figure 30**) – make sure that the gas strut bracket remains in position. Tighten the screw to hold the shaft in place (**Figure 31**).
- c Loosen the top nut from the gas strut bracket (**Figure 32**).
- d Attach the upper part of the gas strut to the bracket (**Figure 33**).
- e Moving the footplate to the up position, push the bracket downwards (**Figure 34**).
- f Once the footplate is in the correct position, retighten the gas strut bracket (**Figure 35**).

Figure 27


Figure 28


Figure 29


Figure 30


Figure 33


Figure 31


Figure 34


Figure 32


Figure 35

## Power pack and footplate

### Powered (motor) footplate

- a Remove the footrest assembly from the footrest kit box. Loosen the screw and remove the shaft from this assembly.
- b Feed the footplate loom and spring around the back of the motor (**Figure 36**).
- c Attach the footrest to the power pack (**Figure 37**).

**Note:** Lubricate the three pivot points and shaft with HTEP grease (food safe grease).

- d Attach the spring to the main chassis bracket (**Figure 38**).
  - e Attach the connector (**Figure 39**).
- 6 Fit the batteries (footplate accessories kit) (**Figure 40**).
- 7 Connect the dummy seat toggle looms:
- a Driving toggle 8 way connector (Part No. 181001.50083) and conversion loom (SIM22900) – seat, direction, key switch and swivel (can only be fitted one way) (**Figure 41**).
  - b 2.5mm battery link (Part No. 181001.52036) across the flying battery leads – red link wire (**Figure 42**).

**Note:** Ensure that this connection is not earthed as it will blow a fuse.

- 8 Taking care not to damage the charge contacts, drive the power pack on to the track using the driving toggle.

**Note:** Ensure that the top and bottom trunnions are correctly engaged to the corresponding teeth on the top and bottom track (**Figure 43**).

**Note:** If the power pack is loaded without corresponding engagement of the top and bottom trunnions; and the power pack is operated it will cause permanent damage to the power pack. **Figure 44** shows how an incorrectly loaded power pack will appear.

**If this occurs a new power pack MUST be ordered as a matter of safety.**


Figure 36


Figure 37








Figure 38

## Power pack and footplate

---

- 9 When the power pack is fully and correctly engaged with the rack, remove the loading bars from the track (**Figure 45**).
- 10 Remove the dummy seat toggle looms.
- 11 Refit the end stop and end cap to the top and bottom tracks (**Figure 46**).

## Inside curve rail

**Items required:**

- Magnets

**Tools required:**

- Superglue

The procedure for the installation of an inside curve track is identical to that of an external curve track with the following additions.

1 Once the track has been assembled and fixed to the stairs it is necessary to fit magnets at the beginning and end of each bend in the track. This magnet will be detected by a reed switch on the power pack trunion. It will slow the lift down as it approaches a bend and then speed it up as it exits the bend.

a Position the magnets on the bottom track (Figures 47 and 48).

b Drive the lift from the bottom of the track up the stairs. The lift should slow down as it enters the bend (passes the first magnet) and speed up again when it exits the bend (passes the second magnet). If the lift does not behave as expected send it back to the bottom of the stairs, reverse the polarity of both magnets and retest.

c Once tested and aligned correctly superglue the magnets into place.

Figure 47


Figure 48

**Note:** The polarity of the magnets can be identified as shown in Figure 49. The polarity must be opposite at the entrance and exit of the bend.

**Note:** Ensure that the bolt is parallel to the track (Figure 49).

**NOTE: THE MAGNETS MUST BE PLACED IN THE CORRECT POSITIONS TO ENSURE THE SAFE TRAVEL OF THE STAIRLIFT.**


Figure 49

## Fitting the 2000 Simplicity seat (manual)

**Items required:**

- Seat

**Tools required:**

- 13mm Socket
- 4mm Allen key
- No. 2 Pozi-drive

**Note: For powered seat – see pages 28-33.**

- 1 Remove seat from the seat box.
- 2 Remove power pack carry handle (**Figure 50**) and fasten the seat post to the correct height (Torque 33Nm) (**Figure 51**).
- 3 Lay the seat on its back on the floor and remove the seat base bolt and washer (**Figure 52**).
- 4 Lift the seat and position it on to the seat stem (**Figure 53**).


**Note:** Be careful not to knock the swivel safety switch when lowering the seat onto the stem (**Figure 54**).

- 5 Gently rotate the seat backwards and forwards until the swivel locking peg locates.
- 6 Fully swivel the seat using the swivel levers.
- 7 Refit the seat base bolt and washer that you removed in step 3 (**Figure 55**).
- 8 Swivel the seat back to the drive position using the swivel levers.
- 9 Remove the drive toggle and link cables and make the seat cable connections by feeding through the seat stem/post (**Figure 56**).


**Note:** Ensure that no cables are pinched or trapped when fitting the power pack front cover.

**Figure 50**

**Figure 51**

**Figure 52**

Figure 53


Figure 55




Figure 54

Figure 56

**10** Fit the earth cable to the top of the seat plate  
**(Figure 57).**

**11** Fit the plastic cap to the seat frame  
**(Figure 58).**

## Electrical connections 2000 Simplicity seat (manual)

### Items required:

- None

### Tools required:

- None

#### 1 Make seat and footplate electrical connections:

- Connect the power switch (2 x red).
- Connect the diagnostic display (red, black and yellow).
- Connect the key switch loom from the power pack (orange and green) to the swivel switch loom connector (2 x white). **(Figure 59)**.
- Connect the seat to the swivel switch loom connector (white) **(Figure 60)**.
- Connect the toggle loom (orange, blue and green) **(Figure 61)**.
- Connect the infra-red receivers on the left and right sides of the pack (blue sleeved).

#### Powered footplate only

- Connect the powered footplate looms (twin yellow) **(Figure 62)**.
- Connect the motor cables (red and black) as per diagram on **page 52**.


Figure 60


Figure 61


Figure 59


Figure 62

## Electrical connections 2000 Simplicity seat (manual)

- 2 Connect the earthing terminal (black) to the pack and secure with the nut (**Figures 63 and 64**).
- 3 Connect the main seat loom (brown and pink) to the chassis loom (red and pink) – two black connectors. (**Figure 65**).
- 4 Test the function of the powered features, including the toggle (**see pages 39-41**).
- 5 Replace and refit the pack front cover (**Figure 66**).


Figure 65


Figure 63

Figure 64

Figure 66

## Fitting the Style seat (manual)

### Items required:

- Seat

### Tools required:

- 13mm Socket
- 4mm Allen key
- No. 2 Pozi-drive

- 1 Remove the seat stem from the packaging (**Figure 67**).
- 2 Attach the seat stem to the chassis (**Figure 68**).
- 3 When lowering the seat onto the stem, ensure the seat is in the ride position and the seat switch roller does not get damaged by the switch cam (**Figure 69**). There is a gap in the cam where the roller will fit (**Figure 70**).


Figure 68

Figure 69


Figure 70




Figure 67

## Fitting the Style seat (manual)


**Figure 71**

**Figure 72**


- 4 Drop the seat onto the stem **(Figure 71)**.
- 5 Gently rotate the seat backwards and forwards until the swivel locking peg locates.
- 6 Fully swivel the seat using the swivel levers.
- 7 Refit the seat base bolt and washer **(Figure 72)**.
- 8 Swivel the seat back to the drive position using the swivel levers.
- 9 Remove the drive toggle and link cables and make the seat cable connections by feeding through the seat stem/post **(Figure 73)**.

**Note:** Ensure that no cables are pinched or trapped when fitting the Power Pack front cover.

**Figure 73**

- 10 Fit the earth cable to the top of the seat plate **(Figure 74)**.
- 11 Fit the plastic cap to the seat frame **(Figure 75)**.
- 12 Tuck the seat cables into the slot of the plastic blanking plate **(Figure 76)**.

**Figure 74**


Figure 75

Figure 76

## Electrical connections Style seat (manual)

**Items required:**


- None

**Tools required:**


- None

**1** Make seat and footplate electrical connections:

- Connect the power switch (2 x red).
- Connect the diagnostic display (red, black and yellow).
- Connect the key switch loom from the power pack (orange and green) to the swivel switch loom connector (2 x white). **(Figure 77)**.
- Connect the seat to the swivel switch loom connector (white) **(Figure 78)**.
- Connect the toggle loom (orange, blue and green) **(Figure 79)**.
- Connect the infra-red receivers on the left and right sides of the pack (blue sleeved).


**Figure 77**

**Figure 78**

**Powered footplate only**

- Connect the powered footplate looms (twin yellow) **(Figure 80)**.
  - Connect the motor cables (red and black) as per diagram on **page 52**.
- 2** Connect the earthing terminal (black) to the pack and secure with the nut **(Figures 81 and 82)**.
- 3** Connect the main seat loom (brown and pink) to the chassis loom (red and pink) – two black connectors. **(Figure 83)**.
- 4** Test the function of the powered features, including the toggle **(see pages 39-41)**.
- 5** Replace and refit the pack front cover **(Figure 84)**.

**Figure 79**


Figure 80


Figure 83


Figure 81


Figure 82


Figure 84

# Fitting the Style seat (powered)

### Items required:

- Seat

### Tools required:

- 13mm Socket
- 4mm Allen key
- No. 2 Pozi-drive

- 1 Remove the seat stem from the packaging (**Figure 85**).
- 2 Attach the seat stem to the chassis (**Figure 86**).
- 3 When lowering the seat onto the stem, ensure the seat is in the ride position and the seat switch roller does not get damaged by the switch cam (**Figure 87**). There is a gap in the cam where the roller will fit (**Figure 88**).
- 4 Drop the seat onto the stem (**Figure 89**).
- 5 Feed the seat cables through the seat stem (**Figure 90**).
- 6 Fit the plastic cap to the seat frame (**Figure 91**).
- 7 Tuck the seat cables into the slot of the plastic blanking plate (**Figure 92**).


Figure 86

Figure 87


Figure 85






Figure 88

# Electrical connections Style seat (powered)




**Items required:**

- None

**Tools required:**

- None

- 1 You now have five seat cables to connect to the chassis (**Figure 93**).
- 2 Starting with the motor connector, connect to the in-line PCB. Make sure the flat surface of the connector is facing the outer edge of the PCB (**Figure 94**).
- 3 Gently push the motor loom connector into the in-line PCB (**Figure 95**).
- 4 Connect the motor loom earth to the chassis earth stud (**Figure 96**).
- 5 Connect the manual swivel over-ride loom to the motor loom (4 x Green 5) (**Figure 97**).
- 6 Connect the manual swivel over-ride loom earth to the chassis earth stud (**Figure 98**).
- 7 Seat Loom 1. Connect earth to chassis earth stud and secure with nut (**Figure 99**).
- 8 Make seat and footplate electrical connections:
  - a Connect the power switch (2 x red).
  - b Connect the diagnostic display (red, black and yellow).
  - c Connect the key switch loom from the power pack (orange and green) to the swivel switch loom connector (2 x white). (**Figure 100**).
  - d Connect the seat to the swivel switch loom connector (white) (**Figure 101**).
  - e Connect the toggle loom (orange, blue and green) (**Figure 102**).
  - f Connect the infra-red receivers on the left and right sides of the pack (blue sleeved).


Figure 96


Figure 99


Figure 97


Figure 100


Figure 98


Figure 101

## Electrical connections Style seat (powered)


Figure 102


Figure 103

### Powered footplate only

- g** Connect the powered footplate looms (twin yellow) (**Figure 103**).
- h** Connect the motor cables (red and black) as per diagram on **page 52**.
- 9** Connect the earthing terminal (black) to the pack and secure with the nut (**Figures 104 and 105**).
- 10** Connect the main seat loom (brown and pink) to the chassis loom (red and pink) – two black connectors. (**Figure 106**).
- 11** Test the function of the powered features, including the toggle (**see pages 39-41**).
- 12** Put the front chassis cover back on the pack and secure with the screws. Ensure the plastic cover plate is under the front cover (**Figure 107**).
- 13** Assemble the seat retaining washer. This is a safety feature and must be fitted (**Figure 108**).


Figure 104


Figure 105


Figure 106

Figure 107

Figure 108

# Remote control handsets

### Items required:

- Handsets

### Tools required:

- No. 2 Pozi-drive

The 2000 series stairlift may be supplied with infrared or radio frequency remote controls, dependant on the requirements of the installation. The infrared remote control has a red dipswitch cover on the rear, whilst the radio frequency remote control has a green dipswitch cover on the rear (**Figure 109**). If you are still unsure as to the type of remote control system you have fitted refer to your drawing or remove the pack underpan and look to see if a radio receiver board has been fitted (**Figure 110**).

## Programming

The infra-red remote control handsets supplied with the lift should work out of the box. If the lift does not travel up and down the stairs when the corresponding handset buttons are pressed, or if there are multiple lifts in one area, follow the procedure set out below.

- 1 Remove the cover from the power pack (**Figure 111**).
- 2 Bring both handsets to the lift and then programme the handset to the board in the following way.
  - a For single lift applications:
    - i Press and hold the red button on the PCB (**Figure 112**) – a yellow LED will illuminate on the PCB.
    - ii Press any button on the handset and the yellow light will go out.
    - iii Programming is complete.

**Note:** Low energy light bulbs can interfere with the infrared signal causing the lift to stop and start. Reprogramming the handsets with the lights on can resolve this problem.


Figure 109


Figure 110



Figure 111

**Note:** If the red button on the control board is depressed for more than 10 seconds then the handset will not be programmed. You must repeat the programming again.

**b** For multiple lift applications:

- i** Remove the red dipswitch protection cap with a small screwdriver (**Figure 109**).
  - ii** Set the dipswitches on the PCB on both handsets to the same settings – note that the next pair of handsets will require the dipswitches to be set to a different combination (**Figure 113**).
  - iii** Press and hold the red button on the PCB (**Figure 112**) – a yellow LED will illuminate on the PCB.
  - iv** Press any button on one handset and the yellow light will go out.
  - v** Programming is complete.
- 3** Test that the handsets have programmed correctly by pressing the up and down button on each and check that the lift moves in the corresponding direction.
  - 4** Finally fit the power unit front cover.

**Figure 112**

**Figure 113**

### Radio frequency remote control

If the stairlift is fitted with radio frequency remote handsets they will automatically program themselves to the lift when first powered on. Some large multi-storey installations may require RF repeaters to be installed in order to ensure a strong and constant control signal reaches the lift at all points along its travel. Handicare's drawing team will normally be able to determine if RF and repeaters will be required based on the staircase length and configuration and therefore an appropriate number of repeaters will be supplied.

RF repeaters can also be requested as a spare to fit retrospectively to a site if required.

If necessary the RF handsets can be programmed as follows:

- 1 Turn on the stairlift power – the receiver lights blink then should go solid green. (By default the dip switches should be OFF-OFF. The dipswitches on the receiver PCB and handset must match).
- 2 Hold the remote control handset near to the lift.
- 3 Remove the cover from the power pack.
- 4 Press and hold the red button on the main PCB – a yellow LED will illuminate on the PCB.
- 5 Press the up direction on the handset and the yellow light will go out.
- 6 Programming is complete.
- 7 Test that the handsets have programmed correctly by pressing the up and down button on each and checking that the lift moves in the corresponding direction.

**Note:** If a second stairlift using RF handsets is fitted in close proximity to the first the dipswitches on the handsets and receiver board for that lift must be set to ON-ON.

### Positioning the RF repeaters

The RF handsets have a built in signal strength indicator; when the signal is strong the green led lights, the orange led indicates an acceptable signal strength, whilst the red led shows that the lift is out of range.

- 1 Drive the lift up from the bottom of the stairs until the red led shows on the handset.
- 2 Position a RF repeater in the nearest appropriate point approximately one metre back down the stairs from where the lift has stopped.
- 3 The repeaters are powered via the red cable in the bottom track rail. Find the nearest pre-drilled hole in the track and splice into the cable.
- 4 The positive connection comes from the red cable in the track. The negative must be connected back to the rail in the same way as the standard negative connection is made (**see page 10**).
- 5 Run the two core cable to the repeater and connect to either pair of positive and negative terminals.
- 6 The RF repeater does not need to be in sight of either the stairlift or remote handset so the RF repeaters can be positioned discreetly, for example behind the track.
- 7 The repeaters will automatically register with the RF board in the lift.
- 8 Repeat steps 1 to 6 for any further repeaters required until the lift can travel the full distance of the rail.

## Fitting RF repeaters

If the lift will not drive to (or cannot be called from) the furthest end of the track when using the RF remotes then a RF repeater(s) can be fitted which will extend the range of the remotes and enable it to do so.

- 1 Whilst standing by the lift use the remote handset to drive the lift toward the other end of the track i.e. from the bottom to the top of the track.
- 2 When the red light shows on the handset stop the lift.
- 3 A RF repeater should be fitted approximately 1m back down the track from where the lift has stopped. The repeater does not require line of sight so can be placed discreetly, for example behind the track.
- 4 After turning off the power to the system, wire the RF repeater into the 33V lift charging supply from the bottom rail.
- 5 Should the stairlift still not reach the end of the track once the RF repeater has been fitted, repeat the above process to fit further RF repeaters. Each additional RF repeater should have different dipswitch settings from the last. Start by standing next to the lift at the point at which you have fitted the last repeater.


Figure 114

Figure 115

Figure 116

## Fitting the handset holder

**Note:** The handset **MUST** be fixed into the holder by screwing through the holder into the handset (**Figure 114**).

- 1 Fit the hand control holder to the wall in the required position using the fixings supplied (**Figures 115 and 116**).

# Fix the track to the stairs

### Items required:

- Screws

### Tools required:

- 6mm Allen key
- No. 3 pozi-drive

- 1 Check that the height of each stanchion is correct.
- 2 Drive the power pack up and down the stairs to ensure that it will not encounter any unforeseen obstacles.

**Note:** Look for any other potential obstacles such as window ledges or exposed pipework that could obstruct the movement of the lift at this stage. If there are any other obstacles then move the track away from the wall until the rear edge of the power pack can travel freely past them.

- 3 When all the stanchions are vertically aligned tighten the screws fully to fasten the foot to the floor or riser (**Figure 117**).

**Figure 117**

Figure 118

Figure 119

## Testing

### Items required:

- Handsets

### Tools required:

- None

### All variants

- 1 Sit on the lift and drive up and down the track to ensure that the footplate clears the nose of the risers along the complete length of the track.
- 2 Check that the lift charges correctly:
  - a Drive the lift to the top charge contacts and ensure that '-' shows in the diagnostic display.
  - b Drive the lift to the bottom charge contacts and ensure that '-' shows in the diagnostic display.
- 3 Check the key switch operation:
  - a Turn the key to the off position and ensure that the lift will not drive - '0' should show in the diagnostic display when you try to drive the lift with the key switch in the off position.
  - b Turn the key switch back on.
- 4 Ensure that the under chassis safety edge on the power pack is functioning (**Figure 118**):
  - a Drive the lift in the down direction.
  - b Press the under chassis safety edge – the lift should stop.
- 5 Ensure that all of the safety edges on the trunnions are functioning (**Figure 119**):
  - a Drive the lift in the up direction.
  - b Press the upside top trunnion safety edge – the lift should stop.
  - c Drive the lift in the up direction.
  - d Press the upside bottom trunnion safety edge – the lift should stop.
  - e Drive the lift in the down direction.

## Testing

- f** Press the downside top trunnion safety edge – the lift should stop. Drive the lift in the down direction.
  - g** Drive the lift in the down direction.
  - h** Press the downside bottom trunnion safety edge – the lift should stop.
- 6** Check the safety edges on the footplate (**Figure 120**):
- a** Drive the lift in the up direction.
  - b** Press the upstairs edge of the footplate – the lift should stop.
  - c** Drive the lift in the down direction.
  - d** Press the downstairs edge of the footplate – the lift should stop.
  - e** Drive the unit in the down direction.
  - f** Press the underside of the footplate – the lift should stop.
- 7** Check the swivel interlock (**Figure 121**):
- a** Drive the lift in the up direction.
  - b** Swivel the seat – the lift should stop.
  - c** Drive the lift in the down direction.
  - d** Swivel the seat – the lift should stop.


- 8** Check the stopping distance of the lift:
  - a** Using the handset, drive the lift in the down direction.
  - b** After the lift has travelled for about 2 seconds, release the down button.
  - c** Check that the lift stopped within 20mm of the point where the button was released.

### Power footplate only

- 1** For seat operated versions:
  - a** Lift the front seat squab.
  - b** The footplate should lift.
  - c** Lower the front seat squab.
  - d** The footplate should lower.

### Manual swivel only

- 1** Check the swivel interlock:
  - a** Drive the lift in the up direction.
  - b** Swivel the seat – the lift should stop.
  - c** Drive the lift in the down direction.
  - d** Swivel the seat – the lift should stop.

## Power swivel only

- a** Drive the lift up to the top stop/charge point.
- b** Keep up-direction toggle held on.
- c** The powered swivel will now swivel out after a 2-3 second delay.
- d** Now move and hold the toggle to the down direction.
- e** The seat will now rotate to the 'ride' position.
- f** Check the manual override switch works by holding down, manually rotating the seat out and then back to 'ride' position.

## Power hinge track only

### 1 Using the handset:

- a** Drive the lift to the top of the stairs.
- b** Depress the down button and keep it depressed.
- c** Observe the function of the lift.  
The lift should:
  - i** Drive down the stairs.
  - ii** Stop at the park position.
  - iii** Beep whilst lowering the hinge.
  - iv** After a short pause continue to the bottom of the track.
- d** Depress the up button and keep it depressed.
- e** Observe the function of the lift.  
The lift should:
  - i** Drive up the stairs.
  - ii** Stop at the park position.
  - iii** Beep whilst raising the hinge.
  - iv** After a short pause continue to the top of the track.

### 2 Sit on the stairlift:

- a** Drive the lift to the top of the stairs.
  - b** Push the toggle switch in the down direction and hold it in that position.
  - c** Observe the function of the lift.  
The lift should:
    - i** Drive down the stairs.
    - ii** Stop at the park position.
    - iii** Beep whilst lowering the hinge.
    - iv** After a short pause continue to the bottom of the track.
  - d** Push the toggle switch in the up direction and hold it in that position.
  - e** Observe the function of the lift.  
The lift should:
    - i** Drive up the stairs.
    - ii** Stop at the park position.
    - iii** Beep whilst raising the hinge.
    - iv** After a short pause continue to the top of the track.
- 3** Turn off the power supply to the charge contacts.
  - 4** Repeat steps 1 and 2 above.
  - 5** Turn the power supply back on.

# Handover

**Items required:**

- None

**Items required:**

- None

**Note:** Ensure that the site is cleaned and tidied before demonstrating the lift to the customer.

Demonstrate the stairlift to the customer, carer and any users or potential users before leaving the installation site.

Please use the following checklist to ensure that all of the items that should be demonstrated are covered:

Feature	Explanation	Done?
<b>Keyswitch</b>	Used to disable the stairlift against unauthorised use – especially useful to prevent children from playing with the stairlift.	
	The stairlift will still charge with the keyswitch disabled.	
<b>On/off switch</b>	DO NOT turn the stairlift off using this switch unless you will not be using the stairlift for a prolonged period such as a holiday.	
<b>Seat-belt</b>	The seat-belt should be used every time the stairlift is used.	
<b>Operating lever</b>	Which way is up / which way is down.	
	How the toggle can be used, e.g. with the fingers, palm of the hand, etc.	
	The delay in pressing the lever before the stairlift will move.	
<b>Operating the stairlift</b>	Always keep your FEET ON THE FOOTPLATE whilst the stairlift is in motion, and try to avoid your feet hanging over the edges of the footplate.	
	Always sit fully back in the chair when the stairlift is in motion.	
	Demonstrate the ‘normal’ noise that a stairlift will make in operation.	

Feature	Explanation	Done?
<b>Remote control</b>	How to call and send the stairlift. The Handicare 2000 Simplicity can be parked on any additional charge contact.	
<b>Folding the stairlift</b>	How to fold and unfold the stairlift.	
	The stairlift should be folded when not in use.	
<b>Operating the swivel</b>	How to operate the swivel.	
	Never swivel the seat whilst the stairlift is in motion.	
	Never remove the seat belt until the chair is swivelled.	
	Never dismount the chair unless the seat is in a locked position.	
	Never dismount the chair whilst the stairlift is in motion.	
<b>Power Footplate (if fitted)</b>	How to operate the power footplate.	
<b>Power Hinge (if fitted)</b>	How to operate the hinge section using the lift toggle control.	
	How to operate the hinge section using the handsets.	
	The lift beeps before and whilst the hinge operates to warn the client that it is moving.	
	How to use the park function on the handset to park the lift and automatically raise the hinge section.	
	The importance of parking the lift and lifting the track when not in use.	
<b>Diagnostic codes</b>	Show the section in the user manual about self help, using the fault codes and the location of the display on the lift.	
<b>Reporting a fault</b>	What number should the customer call and what information do they need to have available when they call.	
<b>Emergency hand winding</b>	When should this be done (overrun, with keyswitch off and power on).	
	How to report a stairlift that has overrun repeatedly.	
<b>Trapped articles</b>	Reverse the stairlift away from the trapped article and remove the item before use.	

## Handover

---

Feature	Explanation	Done?
<b>Other warnings</b>	Never allow more than one person to use the stairlift at any one time. The maximum carrying capacity is 115kg (18 stone/252lb).	
	The stairlift is designed for carrying people only.	
	NEVER allow children to play on or with the stairlift.	
	NEVER allow water to come into contact with the components in the stairlift. If you have to transport liquids DO SO WITH CARE.	
	NEVER place objects in or on the track, or leave objects on the stairs, where they could come into contact with the stairlift in operation. Your stairlift is fitted with sensitive side edges and undertray on the footplate, which will automatically stop the stairlift if it detects any obstructions.	
	NEVER use the stairlift in a standing position.	
	Maintenance and repairs should only be undertaken by a qualified engineer to maintain the validity of the warranty.	
Under no circumstances attempt to repair or re-site the stairlift yourself.		
<b>Servicing</b>	Recommend that the stairlift is serviced by a qualified engineer after 12 months and every 12 months after that.	
<b>Cleaning</b>	Turn the lift off using the key switch and clean with a damp, not wet, lint free cloth and a small quantity of washing up liquid.	
	Do not use abrasive cleaners, bleach or solvent based cleaners as they can damage the stairlift.	

## Maintenance

### Items required:

- Track lubricant:  
HTEP grease  
(food safe grease)
- Cleaning materials

### Tools required:

- None

- 1 Lubricate the rack with a small quantity of lubricant:
  - a Apply a small amount of petroleum jelly at 4 equally spaced points along the rack of both the top and bottom rails (**Figure 122**).
  - b Run the lift up and down the track several times to distribute the lubricant.

**Note:** Do not use an oil based lubricant such as WD40. It will damage the lift.

- 2 Lubricate the swivel mechanism with a small quantity of petroleum jelly:
  - a Remove the seat by unscrewing the retaining bolt (**Figure 123**).
  - b Apply lubricant to the nylon washer.
  - c Refit the seat and refit the retaining bolt (**Figure 123**).
- 3 Check all of the safety features on the lift as described in the Testing section of the installation manual.
- 4 Check that the lift still stops as per the original design.
- 5 Sit on the lift and ride it up and down the stairs several times:
  - a Listen for any unusual noises.
  - b Check for poor ride quality, especially across the spliced joints.
- 6 Clean the track.
- 7 Clean the seat.
- 8 Clean the power pack.

If anything does not function correctly further investigation and corrective action should be undertaken.

Figure 122


Figure 123

# Diagnostic Codes

**Code:**  
None

**Meaning:**  
No display

**Description:**

**Telephone fault finding action:**

Ask client to press toggle and see if display comes on. Ask if main power switch on the lift is turned on. If does not correct fault – send engineer.

**On site fault finding action:**

Check batteries. Check display board is working correctly. Check comms circuit for short circuit.

**Code:**  
Okay

**Meaning:**  
Charging

**Description:**

**Telephone fault finding action:**

If no other code is displayed and lift does not drive – are the arms fully down? Ask client if 8 or 9 show in the display when the toggle switch is activated – if no send an engineer.

**On site fault finding action:**

Send engineer to check toggle and arm circuit.

**Code:**  
0

**Meaning:**  
Final limit activated

**Description:**

Keyswitch, seat swivel, OSG or overrun activated, E-stop, handwind switch.

**Telephone fault finding action:**

Ask client to ensure keyswitch is on and seat is fully swivelled to drive position. If this does not correct the fault – send engineer.

**On site fault finding action:**

Test each item in turn to determine fault.

**Code:**  
1

**Meaning:**  
Requires charge

**Description:**

**Telephone fault finding action:**

Ask client to drive lift to charge contacts – if 'Okay' code does not show send engineer.

**On site fault finding action:**

Drive lift to charge contacts. Check charging circuit if 'Okay' code not displayed.

**Code:**  
2

**Meaning:**  
Off charge

**Description:**

**Telephone fault finding action:**

Ask client to drive lift to charge contacts – if 'Okay' code does not show send engineer.

**On site fault finding action:**

Drive lift to charge contacts. Check charging circuit if 'Okay' code not displayed.

**Code:**  
3

**Meaning:**  
Top stop – Right hand  
Bottom stop – Left hand

**Description:**

End limits activated.

**Telephone fault finding action:**

Ask client to tap trunnion guards in case they are stuck. If this does not correct the fault – send engineer.

**On site fault finding action:**

Check end limit circuit and trunnion guards.

**Code:**  
4

**Meaning:**  
Top pressure – Right hand  
Bottom pressure – Left hand

**Description:**

Safety edge (up/down direction) activated – footplate.

**Telephone fault finding action:**

Ask client to tap footplate edges. If this does not correct the fault – send engineer.

**On site fault finding action:**

Check footplate leading edge.

**5** **Code:** 5 **Meaning:**  
Bottom stop – Right hand  
Top stop – Left hand

**Description:**  
End limits activated.

**Telephone fault finding action:**  
Ask client to tap trunnion guards in case they are stuck. If this does not correct the fault – send engineer.

**On site fault finding action:**  
Check end limit circuit and trunnion guards.

**6** **Code:** 6 **Meaning:**  
Bottom pressure – Right hand  
Top pressure – Left hand

**Description:**  
Safety edge (down/up direction) activated – footplate, underpan.

**Telephone fault finding action:**  
Ask client to tap footplate edges. If this does not correct the fault – send engineer.

**On site fault finding action:**  
Check footplate bottom and downside and chassis underpan.

**7** **Code:** 7 **Meaning:**  
Low voltage

**Description:**  
Batteries approaching critical level.

**Telephone fault finding action:**  
Ask client to drive lift to charge contacts – if 'Okay' code does not show send engineer. If the lift will not drive – send engineer.

**On site fault finding action:**  
Drive lift to charge contacts. Check charging circuit if 'Okay' code not displayed. If lift will not drive replace batteries and check charging circuit.

**8** **Code:** 8 **Meaning:**  
Up direction – Right hand  
Down direction – Left hand

**Description:**  
Shows when toggle switch is activated in up/down direction.

**Telephone fault finding action:**  
None.

**On site fault finding action:**  
None.

**9** **Code:** 9 **Meaning:**  
Up direction – Left hand  
Down direction – Right hand

**Description:**  
Shows when toggle switch is activated in down/up direction.

**Telephone fault finding action:**  
None.

**On site fault finding action:**  
None.

**A** **Code:** A **Meaning:**  
Hinge open

**Description:**

**Telephone fault finding action:**  
None.

**On site fault finding action:**  
None.

**b** **Code:** b **Meaning:**  
Switch active

**Description:**  
A switch is active at power on i.e. Toggle Up, Toggle down, IR Up IR down or learn.

**Telephone fault finding action:**  
Send engineer.

**On site fault finding action:**  
A short circuit has occurred on toggle circuit – check circuit.

**C** **Code:** C **Meaning:**  
IR Address fail

**Description:**  
The IR Address (Dip switch) does not match.

**Telephone fault finding action:**  
Send engineer.

**On site fault finding action:**  
Re-program handsets.

## Appendix 2

---

**Code:**  
d

**Meaning:**  
Relay not open (stopped)

**Description:**

The main power relay is welded closed.

**Telephone fault finding action:**

Send engineer.

**On site fault finding action:**

Replace main board.

**Code:**  
E

**Meaning:**  
Relay not closed

**Description:**

The main power relay did not close.

**Telephone fault finding action:**

Send engineer.

**On site fault finding action:**

Check battery voltage and replace if necessary, otherwise replace main board.

**Code:**  
F

**Meaning:**  
Brake semi-conductor failed

**Description:**

**Telephone fault finding action:**

Send engineer.

**On site fault finding action:**

Replace main board.

**Code:**  
g

**Meaning:**  
Brake not connected

**Description:**

**Telephone fault finding action:**

Send engineer.

**On site fault finding action:**

Check brake.

**Code:**  
H

**Meaning:**  
Relay not open (Pre delay)

**Description:**

The main power relay is closed early.

**Telephone fault finding action:**

Send engineer.

**On site fault finding action:**

Replace main board.

**Code:**  
J

**Meaning:**  
Hinge interlock switch error

**Description:**

**Telephone fault finding action:**

Send engineer.

**On site fault finding action:**

Replace hinge interlock roller switch.

**Code:**  
L

**Meaning:**  
Current limit exceeded

**Description:**

**Telephone fault finding action:**

Check client understands loading limits.  
Send engineer if persistent fault.

**On site fault finding action:**

Check motor and track for obstructions.

**Code:**  
n

**Meaning:**  
Half speed

**Description:**

**Telephone fault finding action:**

Send engineer if permanently showing.

**On site fault finding action:**

Check reed switches and dipswitch settings.  
Check battery voltage.

**Code:**  
o**Meaning:**  
Default Eeprom**Description:**

The Eeprom has been reset to its default value (Flash Failure).

**Telephone fault finding action:**

Send engineer.

**On site fault finding action:**

Replace main board.

**Code:**  
P**Meaning:**  
PS no reply**Description:**

Powered swivel did not respond.

**Telephone fault finding action:**

NOT YET ACTIVE.

**On site fault finding action:**

NOT YET ACTIVE.

**Code:**  
r**Meaning:**  
No charge current**Description:**

Activated.

**Telephone fault finding action:**

Check that power supply is switched on from mains supply. If yes send engineer.

**On site fault finding action:**

Check power supply.

**Code:**  
U**Meaning:**  
No float indication V1.19 software and later only**Description:**

Faulty power supply.

**Telephone fault finding action:**

Send engineer.

**On site fault finding action:**

Replace power supply.

**Code:**  
y**Meaning:**  
Software error**Description:**

The main control board has a software fault.

**Telephone fault finding action:**

Send engineer.

**On site fault finding action:**

Replace main board.

**Code:**  
N/A**Meaning:**  
Over current**Description:**

N/A.

**Telephone fault finding action:**

Send engineer with a service exchange power pack.

**On site fault finding action:**

Replace power pack.

**Code:**  
N/A**Meaning:**  
Lift is in "sleep mode" to conserve battery power**Description:****Telephone fault finding action:**

None.


**On site fault finding action:**

Press toggle or handset button to wake lift.


**Note: If the display is showing anything but the diagnostic codes listed it requires resetting.**

Drive the lift off the charge contacts.  
Turn the main switch on the unit off for 10 seconds and then back on.


2000 Simplicity / Style block schematic – LH stairlift – HAN125


# 2000 Simplicity / Style block schematic – RH stairlift – HAN125


## Power footplate PCB connections


## HAN126 Style powered swivel PCB connections

# Style power swivel overview – brushless

## HAN125 Internal RTC reed switch connections

(INTERNAL CURVE CONNS)

WHEN LOOM 181001. 50062 IS PLUGGED INTO THE 4 WAY CONNECTOR, 2 ORANGE + 2 RED WIRES WILL BE AVAILABLE TO CONNECT INTO.


**CONNECTION WITH ONE SWITCH:**

CONNECT ORANGE (1) TO ONE SPADE CONNECTOR  
 CONNECT RED (4) TO SECOND SPADE CONNECTOR  
 SET JUMPER ON PCB TO POSITION A  
 SAME CONNECTION REGARDLESS OF LIFT HAND

**CONNECTION WITH TWO SWITCHES:**

LEFT HAND  
 ORG (1 & 2) = UPPER REED, RED (3 & 4) = LOWER REED  
 RIGHT HAND  
 ORG (1 & 2) = LOWER REED, RED (3 & 4) = UPPER REED  
 SET JUMPER ON PCB TO POSITION B  
 O/C ON CIRCUIT = HALF SPEED  
 S/C ON CIRCUIT = FULL SPEED


## Standard track wiring – LH and RH


NOTE: ON A 1000 LIFT, POWER SUPPLY CAN BE WIRED FROM BOTTOM TRACK SECTION IF REQUIRED. SIMPLY FIT END CAP PRESSING (WITH EARTH CABLE) ON TRACK SECTION DOWNSTAIRS.

ON A 2000 LIFT, POWER SUPPLY CAN BE WIRED FROM TOP, BOTTOM OR MIDDLE OF TRACK SECTION IF REQUIRED.


# Hinge PCB connections – LH and RH


# Hinge schematic – LH and RH

## HAN125 2000 DIP settings **LH** standard


### FUNCTION OF SWITCHES


**IF ANY DIP SWITCHES ARE CHANGED, PLEASE TURN OFF LIFT, WAIT 10 SECONDS, THEN POWER THE LIFT BACK ON.**

## HAN125 2000 DIP settings **RH** standard


### FUNCTION OF SWITCHES


**IF ANY DIP SWITCHES ARE CHANGED, PLEASE TURN OFF LIFT, WAIT 10 SECONDS, THEN POWER THE LIFT BACK ON.**

# HAN125 2000 DIP settings **LH** internal bend


## FUNCTION OF SWITCHES


**IF ANY DIP SWITCHES ARE CHANGED, PLEASE TURN OFF LIFT, WAIT 10 SECONDS, THEN POWER THE LIFT BACK ON.**

# HAN125 2000 DIP settings **RH** internal bend


## FUNCTION OF SWITCHES


**IF ANY DIP SWITCHES ARE CHANGED, PLEASE TURN OFF LIFT, WAIT 10 SECONDS, THEN POWER THE LIFT BACK ON.**

# HAN125 2000 DIP settings **LH** auto hinge HAN126


**IF ANY DIP SWITCHES ARE CHANGED, PLEASE TURN OFF LIFT, WAIT 10 SECONDS, THEN POWER THE LIFT BACK ON.**

## DIP SWITCH SETTING HAN125


## DIP SWITCH SETTING HAN126


THIS SETTING WILL ENABLE 2000 HINGE MODE AND ALSO SET THE CURRENT LIMIT FOR AN RTC HINGE

# HAN125 2000 DIP settings **RH** auto hinge HAN126


**IF ANY DIP SWITCHES ARE CHANGED, PLEASE TURN OFF LIFT, WAIT 10 SECONDS, THEN POWER THE LIFT BACK ON.**

## DIP SWITCH SETTING HAN125

## DIP SWITCH SETTING HAN126

THIS SETTING WILL ENABLE 2000 HINGE MODE AND ALSO SET THE CURRENT LIMIT FOR AN RTC HINGE

# External track details

1 FITTING KIT DETAILS

2 TRACK PARTS REQUIRED

3 SEAT TYPE

4 SIDE OF STAIRS FOR INSTALLATION LOOKING FROM THE BOTTOM

5 TRACK COLOUR

6 APPLICATIONS ENGINEER

7 CUSTOMER REFERENCES

8 DRAWING NUMBER

9 STANCHION DETAILS

10 TRACK DETAILS

11 TRACK LAYOUT

**DO NOT SCALE ALL DIMENSIONS IN MM UNLESS OTHERWISE STATED**

DRAWING NUMBER: rtc-91602

Product type: RTC-2000

**Fitting Kit**  
 Hand Sbs: **RF**  
 Number needed: **2**  
 Stanchion Feet Standard: **11**  
 Lower: **1**  
 Located on: **9**  
 Fan: **0**  
 Located on:

Wood Screws for feet: **36**  
 Stanchion bolts: **24**  
 Wall Brackets: **0**  
 Located on:

Heavy Duty Kit: **NO**  
 Wire needed (m): **7**  
 End Stops: **MK2**  
 Magnets: **0**  
 Brace Bar: **0**

**Parts for Build**  
 Track Balls: **22**  
 Joining plugs: **12**  
 Charge Points: **Standard**  
 Additional Charge Points (section/ by stanchion)

Swivel Bracket: **NO**  
 UnderFloor Lower: **0**  
 UnderFloor Standard: **0**

Seat Type: **2000 STYLE**

Hand of Rail: **RIGHT**

Rail Colour: **Light Grey -**  
 Build party / track format checked  
 By:

WS No.	ISG	Loans	Teeth	TYPE	LENGTH	SB No.	UF	RAIL	TOOTH	DISTANCE	LENGTH	TYPE	SETTING HT
1				GROUND FLOOR									
A	2	B-B	152	S	1671	10	2	37	407	615	B	114	
	3	F	B-X	140	S	1539	11	2	97	1067	659	B	158
B	4	L	X-W	28	GD R665		12	2	151	1660	658	B	157
	5	W-B	31	S	341		13	3	52	572	654	B	153
C	6	B-B	65	HL14.25 R447		N	14	3	105	1155	648	B	147
	7	F	B-X	10	GI R596		15	4	16		615	B	114
D	8	L	X-B	44	S	484		5	30	330	541	B	13
				TOTAL =	5168		17	6	30		619	B	118
							18	7	1		718	B	217
							19	8	21	231	648	B	147
							20	8	43	473	556	B	55
										TOTAL =	7281		

MARKED FOR WELD BY:	Marked + Checked M/C SHOP BY:	Date:
Track Assembled BY:	Track Checked BY:	

DRAWN: Kamaldeep.Bilku	CONTRACT / CUSTOMER: Companion Stairlifts	DATE: 18/03/15
REFERENCE: SMITH54033	DRAWING NUMBER: rtc-91602	PG NUMBER: 8

# External track details

**1** FITTING KIT DETAILS

**2** TRACK PARTS REQUIRED

**3** SEAT TYPE

**4** SIDE OF STAIRS FOR INSTALLATION LOOKING FROM THE BOTTOM

**5** TRACK COLOUR

**6** APPLICATIONS ENGINEER

**7** CUSTOMER REFERENCES

**8** DRAWING NUMBER

**9** TRACK SIDE VIEW

**10** CLIENT DETAILS

**11** TRACK ISOMETRIC VIEW

**12** TRACK PLAN VIEW

**13** SPECIAL INSTRUCTIONS

**14** TRACK FRONT VIEW

**15** RISER HEIGHT DETAILS

DO NOT SCALE ALL DIMENSIONS IN MM UNLESS OTHERWISE STATED		REV	COMMENTS	BY	DATE
No	Riser Height Left				
1	199				198
2	202				199
3	200				202
4	196				197
5	199				197
6	194				198
7	199				198
8	199				197
9	197				198
10	201				196
11	200				197
12	195				208
13	197				192
14	200				200

**12**

WARNING!  
Less than 50mm knee clearance. Modifications may be required.

**14**

WARNING!  
Less than 100mm bulkhead clearance. Modifications may be required.

**11**

**13**

**Style Seat Applied Due To Knee Clearance**

USER DETAILS	
Standing Ht:	1650 mm
Seated Ht:	840 mm
Back To Knee:	950 mm
Seat Height:	485 mm

**10**

	DRAWN: Kamaldeep.Bilku	DATE: 18/03/15
	CONTRACT / CUSTOMER: Companion Stairlifts	PG NUMBER: 7
	REFERENCE: SMITH54033	DRAWING NUMBER: rtc-916028

**1** Fitting Kit

Hand Sets: **RF**

Number needed at Stanchion Feet: **2**

Standard: **11**

Lower: **1**

Located on: **9**

Fan: **0**

Located on:

Wood Screws for feet: **36**

Stanchion bolts: **24**

Wall Brackets: **0**

Located on:

Heavy Duty Kit: **NO**

Wire needed (m): **7**

End Stops: **MK2**

Magnets: **0**

Brace Bar: **0**

Parts for Build

Track Bolts: **22**

Joining plugs: **12**

Charge Points: **Standard**

Additional Charge Points (section/ by stanchion)

Swivel Bracket: **NO**

UnderFloor Lower: **0**

UnderFloor Standard: **0**

Seat Type: **2000 STYLE**

Hand of Rail: **RIGHT**

Rail Colour: **Light Grey -**

Build parts/ track format checked

By:

This drawing and any information or descriptive matter set out herein are the confidential and copyright property of HANDICARE LTD and must not be disclosed, copied or used for manufacturing, marketing or any other purpose without their written permission.

# Inside track details

**DRAWING NUMBER:** rtc-89955  
**Product type:** RTC-2000

**DO NOT SCALE ALL DIMENSIONS IN MM UNLESS OTHERWISE STATED**

**REV**      **COMMENTS**      **BY**      **DATE**  
**A**      **1st Issue**

**Fitting Kit**  
 Hand Sets: **N/A**  
 Number needed: **0**  
 Stanchion Feet  
 Standard: **7**  
 Lower: **1**  
 Located on: **16**  
 Fan: **5**  
 Located on: **19, 20, 21**  
**27, 28**  
 Wood Screws for feet: **39**  
 Stanchion bolts: **26**  
 Wall Brackets: **0**  
 Located on:

**Heavy Duty Kit: NO**  
 Wire needed (m): **8**  
 End Stops: **MK2**  
 Magnets: **6**  
 Brace Bar: **1**

**Parts for Build**  
 Track Bolts: **44**  
 Joining plugs:  
 Charge Points:  
**Standard**  
 Additional Charge Points (section/ by stanchion)  
 Swivel Bracket: **NO**  
 Underfloor Lower: **0**  
 Underfloor Standard: **0**

Seal Type:  
**Track Only**

Hand of Rail:  
**RIGHT**

Rail Colour:  
**Light Grey**

Build parts/ track format checked  
 By:

rtc-89955 TRACK DETAILS																	
W	S	No	SG	Joins	Teeth	TYPE	LENGTH	SB	No	UF	RAIL	TOOTH	DISTANCE	LENGTH	TYPE	SETTING	
GROUND FLOOR																	
		1							16	2	1	11	395		A	74	
		2	B-W	28	S		308		17	2	24	264	470		A	149	
		3	W-W	48	HR17 R159				18	4	4			649	B	127	
<b>A</b>		4	F	W-X	37	GI R596			19	5	26	286	693		C	171	
		5	L	X-C	44	S	484		20	6	5	55	697		C	175	
		6	C-W	10	S		110		21	7	41		708		C	188	
		7	W-W	54	HR54 R223				22	9	16	176	649		B	127	
<b>B</b>		8	F	W-X	13	GD R665			23	9	73	803	661		C	139	
		9	L	X-C	158	S	1737		24	9	129	1418	662		C	140	
		10	F	C-X	12	S	132		25	12	6	66	709		C	187	
		11	L	X-W	13	GI R596			26	12	48	528	884		D	31	
		12	W-W	70	S		770		27	14	15		685		C	163	
<b>C</b>		13	W-W	54	HR54 R223				N	28	15	14	154	620	B	98	
		14	F	W-X	24	GD R665								TOTAL =		8482	
		15	L	X-B	15	S	165										
																TOTAL =	6377

MARKED FOR WELD BY: \_\_\_\_\_ Date: \_\_\_\_\_  
 Marked + Checked M/C SHOP BY: \_\_\_\_\_  
 Track Assembled BY: \_\_\_\_\_  
 Track Checked BY: \_\_\_\_\_

DRAWN: mark.staniland      DATE: 26/01/15  
 CONTRACT / CUSTOMER: Handicare Stairlifts Ltd      PG NUMBER: 30122      DRAWING NUMBER: rtc-89955  
 REFERENCE: R&D Position 5

This drawing and any information or descriptive matter set out hereon are the confidential and copyright property of HANDICARE LTD. and must not be disclosed, loaned, copied or used for manufacturing, forwarding or any other purpose without their written permission.

**Product type:** RTC-2000

**Fitting Kit**  
 Hand Sets: **N/A**  
 Number needed: **0**  
 Stanchion Feet  
 Standard: **7**  
 Lower: **1**  
 Located on: **16**  
 Fan: **5**  
 Located on: **19, 20, 21**  
**27, 28**  
 Wood Screws for feet: **39**  
 Stanchion bolts: **26**  
 Wall Brackets: **0**  
 Located on:

**Heavy Duty Kit: NO**  
 Wire needed (m): **8**  
 End Stops: **MK2**  
 Magnets: **6**  
 Brace Bar: **1**

**Parts for Build**  
 Track Bolts: **44**  
 Joining plugs:  
 Charge Points:  
**Standard**  
 Additional Charge Points (section/ by stanchion)  
 Swivel Bracket: **NO**  
 Underfloor Lower: **0**  
 Underfloor Standard: **0**

Seal Type:  
**Track Only**

Hand of Rail:  
**RIGHT**

Rail Colour:  
**Light Grey**

Build parts/ track format checked  
 By:

BY	DATE	
No	Riser Height	Riser Height
1	21	16
2	199	206
3	197	193
4	197	200
5	195	199
6	201	200
7	198	200
8	204	201
9	203	199
10	201	200
11	199	200
12	200	203
13	202	201
14	200	200
15	198	199
16	201	202
17	199	196
18	200	200
19	200	202
20	200	200
21	203	201
22	200	200

**USER DETAILS**  
 Standing Ht: 1850 mm  
 Seated Ht: 850 mm  
 Back To Knee: 600 mm  
 Seat Height: 485 mm

DRAWN: mark.staniland      DATE: 26/01/15  
 CONTRACT / CUSTOMER: Handicare Stairlifts Ltd      PG NUMBER: 30122      DRAWING NUMBER: rtc-89955  
 REFERENCE: R&D Position 5

This drawing and any information or descriptive matter set out hereon are the confidential and copyright property of HANDICARE LTD. and must not be disclosed, loaned, copied or used for manufacturing, forwarding or any other purpose without their written permission.

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Deze tekening en eventuele informatie of beschrijvingen zamen dienen tevens te worden gezien zijn verantwoordelijk en het copyright is in handen van HANDICARE LTD. Het is niet toegestaan deze tekening of informatie te kopiëren, verspreiden of anderszins openbaar te maken. Het is niet toegestaan deze tekening of informatie te kopiëren, verspreiden of anderszins openbaar te maken.
Handicare Treppenlifte GmbH, Hausberger Straße 16