Protocole CAO 3shape

Sans SCAN POST

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The screenshot displays the 3shape software interface for dental crown design. The main 3D view shows a tooth with a crown preparation in yellow. A 2D cross-section window, titled 'Coupe 2D', shows the internal profile of the preparation with a 0.102 mm gap indicated. The 'Paramètres avancés' (Advanced Parameters) panel is open, with the 'Espace ciment suppl.' (Additional Cement Space) parameter highlighted in red, set to 0.100 mm. Other parameters include 'Espace ciment' (0.010 mm), 'Dist. à la ligne cervic.' (1.20 mm), and 'Distance de lissage' (0.35 mm). The interface includes a top menu bar, a toolbar, and a left sidebar with navigation and design tools.

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Commande: Couronne 27, Inlay-Core 27
Cabinet:
Patient: Test work flow complet. TRI05 Basic (TRI05 Classic)

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Lisser les artefacts de surface

Commande: Couronne 27, Inlay-Core 27
 Cabinet:
 Patient: Test work How complet TRIOS Basic (TRIOS Classic)

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



Commande: Couronne 27, Inlay-Core 27
Cabinet:
Patient: Test work flow complet. TRIOS Basic (TRIOS Classic)

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
The screenshot displays a software interface with a 'Commandes' window. On the left, a sidebar lists various filters such as 'Par date' (today, two days ago, last week, two weeks ago) and 'Par état' (Created, Digitized, Conceived, Sent, Prepared, In production, Approved, Failed, Closed). The main window title is 'Résultats pour 'inlay-core' (cliquer pour renommer)'. Below the title is a search icon and a close button. The central area contains a table with three columns: 'Files', 'Status', and 'Progress'. One row is highlighted in blue, showing a file name starting with '82799_20190211_1025_TRIOS_Basic_(TRIOS_...', a status of 'In Progress', and a progress bar at 2%. At the bottom of the table area, there are two buttons: 'Stop Current' and 'Stop All'.

Files	Status	Progress
82799_20190211_1025_TRIOS_Basic_(TRIOS_...	In Progress	2%

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