

Wireless bedwetting alarm system

Product information

User manual

Warranty



Thanks, you've selected the Rodger bedwetting alarm

Thank you for your confidence in the Rodger bedwetting alarm. We assume that, if you use the bedwetting alarm in a good way, you will have success. Please, read all instructions before using the bedwetting alarm.

- The Rodger bedwetting system is suitable for children over the age of 5 years.
- When beginning the therapy, be motivated, relaxed and confident.
- Perform the alarm test before the start of the treatment. You and your child will get familiar with the bedwetting alarm. Repeat this before your child goes to sleep.
- Make no contest of the therapy. This will cause stress and is counterproductive.
- Place the bedwetting alarm in a certain distance from the bed, so your child has to get out of bed to switch off the alarm.
- Give your child an occasional reward. A small present is a nice positive reinforcement that your child is doing well.
- Use the Rodger bedwetting system every night until there have been no "wet beds" for 14 consecutive nights.
- Successful treatment with the Rodger bedwetting system takes from a few weeks to a few months. Persistence is the key to success.
- Do your best to use the alarm every night without interruption.
- If your child's bedwetting condition does not improve within 4 months, visit your doctor to rule out any medical conditions.
- Do not start the therapy if important events will take place soon, i.e. moving the family or the birth of a child.
- Do not wake your child at night to go to the toilet.
- Never punish your child; your child is not responsible for occasional accidents.
- Do not prohibit your child from drinking before going to bed. Drinking less does not decrease the likelihood of bedwetting. However, it is preferable that your child does not drink coffee, tea, or beverages containing carbonation or chocolate before going to bed.

- In the early days of the treatment, a child who is a heavy sleeper might not wake up. You could choose another sound of the alarm. Another good solution could be to sleep next to your child and help your child to wake up as soon as the alarm sounds, but make sure that your child stops the alarm and goes to the toilet. In time, your child will “tune in” and do this more independently.
- If relapse occurs, another short course will usually succeed. Store all parts together for future use. Remove battery before storing.

Transmitter & receiver

Range:

The range of the receiver and the transmitter is maximum 25 meters (80 feet) in open air. Inside buildings, the range could be less than 25 meters (80 feet).

Interference:

Other transmitters nearby (e.g.: wireless baby-monitor, garage door openers, weather station, headset, etc.) can disturb the bedwetting alarm system. Change the frequencies of the other devices or remove them.

Power supply:

The receiver works either plugged in or with batteries (never use in combination as this may damage the circuitry!).

Usage:

- Always remove the transmitter immediately from a wet pant. This extends the life of the battery and prevents for damage to the pants.
- When a red light is blinking inside the **transmitter**, please replace the battery.
- Make sure to remove the transmitter from the pants before laundering. The warranty is voided if the transmitter has been in the washing-machine.
- The transmitter and receiver are programmed to be used together. If either part is defective, both must be replaced.
- Remove the batteries if the bedwetting alarm system will not be used for a longer time.

Alarm test:

- Put the receiver into an electrical outlet (or use batteries).
- Adjust the volume **(D)** to maximum.
- The green light should shine continuously.
- Connect the two push-buttons of the transmitter with the metal strip **(F)** on the receiver. This will cause the alarm to sound.
- The alarm can be turned off by pushing the stop button **(A)**. If the receiver does not sound, then check the battery in the transmitter.

Pants

Laundering:

The pants need to be washed before use. They become more absorbent. The pants can be machine washed and tumble dried. Do certainly not wring! Avoid fabric softener. Over time, residue may build up on the sensor contacts inside the pant. You simply need to soak the sensor pants for about 30 minutes in vinegar solution (1 cup vinegar and 2 cups water).

Pantstest:

- Take a cup of lukewarm saltwater (1 teaspoon salt in 1 cup water).
- Click the transmitter on the pant.
- Put the receiver into an electrical outlet (or use batteries).
- Adjust the volume **(D)** to maximum.
- The green light should shine continuously.
- Pour two teaspoons of this salty solution on the front of the pant. This should cause the alarm to sound.
- The alarm can be turned off by pushing the stop button **(A)**.

Setting sounds

See inside cover for the reference of the letters for the operation.

Option 1: For every alarm to be the same sound > set switch to 1;

Select the sound of your choice by pressing the “select” button **(C)**. The last selected sound will appear every time when the alarm goes off.

To test the sound, connect the two push-buttons of the transmitter with the metal strip **(F)** on the front of the receiver. This will cause the alarm to sound.

Option 2: For every alarm to be a different sound > set switch to 8;

With every alarm 1 of the 8 available sounds will sound.

Manual instruction pictograms




































See inside cover for an overview of the instruction pictograms belonging to the explanations below.

1. Plug receiver into the wall. Then turn it on and adjust the volume.
2. Pull on the underpants.
3. Attach the transmitter to the underpants.
4. Sleep well.
5. The receiver will sound an alarm when urination occurs.
6. Get up, turn off the alarm and go to the toilet.
7. Repeat steps 2 and 3, sleep well.
8. Wake up dry!!!

Online success calendar

The online success calendar is an interactive version of the standard success calendar and is very motivating during the treatment. By using the online success calendar, you receive informative messages about the progress of the training. After a dry night, your child may choose a game of its choice as a reward. Further instructions: www.rodger.nl.




Explanation success calendar

	1	2		4			
							
							
							
							

Dry: Slept through

Dry: Visited toilet

Wet: Woke by alarm

Wet: Slept through alarm



Child: Dry-night; child sticks on a sticker.

Parents: Tick where appropriate.

Vibrating cushion

For children who are deaf or hard of hearing, a vibrating cushion can be ordered to attach to the alarm system **(E)**. This vibrating attachment can also be useful for children who need extra stimulus to respond to the alarm system. The vibrating cushion can be placed under the pillow to alert your child that the alarm is sounding.

Mattress covers

Mattress covers protect your bedlinen and keep your mattress dry if the child wets the bed. Further information: www.rodger.nl.

Warning

- Keep sufficient ventilation while applying AC power source.
- To prevent fire or shock hazard, do not expose this product to rain or any type of moisture.
- Please ensure yourselves of your local laws, by throwing away the batteries.

Help desk

For questions about the bedwetting treatment or the use of the bedwetting alarm, please get in contact with your supplier.

Warranty

Rodger carries a one-year warranty, after receipt of the device by the purchaser, against manufacturer's defects and against breakage or failure under normal use - at our discretion. This warranty does not cover batteries, battery lids, misuse or immersing the transmitter in water. We shall not be liable for more or other than the aforesaid warranty obligations. This warranty is valid only if you are the first buyer and have the receipt. The warranty is not valid if the alarm is not used in accordance with the instructions or if repairs to the device are made by third parties. Rodger reserves the right to change specifications and warranty without prior notice.