



One Step Multi-Drug Oral Fluid Test

Catalogue No. See Box Label

For in vitro diagnostic use.

PrimeScreen® One Step Multi-Drug Oral Fluid Test offers qualitative detection of the following drugs of abuse and their principal metabolites in human oral fluid at specified cut-off levels for use in employment and insurance testing: Amphetamine (AMP), Barbiturates (BAR), Cocaine (COC), Methylenedioxyamphetamine (MDMA), Methamphetamine (MET), Methadone (MTD), Opiate (OPI), Oxycodone (OXY), Phencyclidine (PCP), and Marijuana (THC).

INTENDED USE

PrimeScreen® One Step Multi-Drug Oral Fluid Test is a rapid oral fluid screening test. The test is a lateral flow, one-step immunoassay for the qualitative detection of specific drugs and their metabolites in human oral fluid at the following cut off concentrations for use in employment and insurance testing.

Test	Calibrator	Cut off (ng/mL)
Amphetamine (AMP)	D-Amphetamine	50
Barbiturates (BAR)	Secobarbital	20
Cocaine (COC)	Cocaine	20
Methylenedioxyamphetamine (MDMA)	3,4-Methylenedioxyamphetaminel	50
Methamphetamine (MET)	D-Methamphetamine	50
Methadone (MTD)	Methadone	30
Opiate (OPI)	Morphine	40
Oxycodone (OXY)	Oxycodone	20
Phencyclidine (PCP)	Phencyclidine	10
Marijuana (THC)	Δ 9-THC	40

This test will detect other related compounds, please refer to the Analytical Specificity table in this package insert.

The assay provides a qualitative, preliminary test result. A more specific analytical method must be used in order to obtain a confirmed result. Gas Chromatography/Mass Spectrometry (GC/MS) or Liquid Chromatography/Tandem Mass Spectrometry (LC/MS-MS) are preferred confirmatory methods. Professional judgment should be applied to any drug test result, particularly when preliminary results are positive.

Amphetamine (AMP): Amphetamine is a sympathomimetic amine with therapeutic indications. The drug is often self-administered by nasal inhalation or oral ingestion.

Barbiturates (BAR): Barbiturates are a class of central nervous system depressants. Abuse of barbiturates can lead not only to impaired motor coordination and mental disorder, but also to respiratory collapse, coma and even death. Barbiturates are taken orally, rectally, or by intravenous and intramuscular injections.

Cocaine (COC): Cocaine derived from leaves of coca plant, is a potent central nervous system stimulant and a local anesthetic. Among the psychological effects induced by using. Cocaine are euphoria, confidence and a sense of increased energy, accompanied by increased heart rate, dilation of the pupils, fever, tremors and sweating.

Methylenedioxyamphetamine (MDMA): Methylenedioxyamphetamine (ecstasy) is a designer drug first synthesized in 1914 by a German drug company for the

treatment of obesity. Those who take the drug frequently report adverse effects, such as increased muscle tension and sweating. MDMA is not clearly a stimulant, although it has, in common with amphetamine drugs, a capacity to increase blood pressure and heart rate. MDMA does produce some perceptual changes in the form of increased sensitivity to light, difficulty in focusing, and blurred vision in some users.

Methamphetamine (MET): Methamphetamine is a potent stimulant chemically related to amphetamine but with greater CNS stimulation properties. The drug is often self-administered by nasal inhalation, smoking or oral ingestion.

Methadone (MTD): Methadone is a synthetic analgesic drug that is originally used in the treatment of narcotic addict. The drug is often administered orally or intravenously and is metabolized in the liver and excreted in urine.

Opiates (OPI): The opiates such as heroin, morphine, and codeine are derived from the resin of opium poppy. The principal metabolites of opiates are morphine, morphine-3-glucuronide, normorphine and codeine with a half-life of about 3 hours. Heroin is quickly metabolized to morphine. Thus, morphine and morphine glucuronide might both be found in the saliva of a person who has taken only heroin. The body also changes codeine to morphine. Thus, the presence of morphine (or the metabolite, morphine glucuronide) in the saliva indicates heroin, morphine and/or codeine use. The window of detection varies for different opiates. Codeine can be detected within one hour and up to 7-21 hours after a single oral dose. Morphine is detectable for several days after a dose.

Oxycodone (OXY): Oxycodone is known as Oxycontin, Roxicodone and is an ingredient of Percodan, Percocet, Roxicet and Tylox. Oxycodone is a semi-synthetic opiates derived from opium. Like other opiates, oxycodone is characterized by its analgesic properties, and the tendency for users to form a physical dependency and develop tolerance with extended use. Oxycodone is usually administered in combination with non-opiate analgesics such as acetaminophen and salicylates for the relief of moderate to severe pain. Oxycodone is a central nervous system depressant that may cause drowsiness, dizziness, lethargy, weakness and confusion. Toxicity in an overdose of oxycodone can lead to stupor, coma, muscle flaccidity, severe respiratory depression, hypotension, and cardiac arrest.

Phencyclidine (PCP): Phencyclidine is an arylcyclohexylamine that was originally used as an anesthetic agent and a veterinary tranquilizer. Phencyclidine can produce hallucinations, lethargy, disorientation, loss of coordination, trance-like ecstatic states, a sense of euphoria and visual distortions. It has many street names, such as "angel dust" and "crystal cyclone," etc. Phencyclidine can be administered orally, by nasal ingestion, smoking, or by intravenous injection.

Marijuana (THC): Tetrahydrocannabinol, the active ingredient in the marijuana plant (cannabis sativa), is detectable in oral fluid shortly after use. The detection of the drug is thought to be primarily due to the direct exposure of the drug to the mouth (oral and smoking administrations) and the subsequent sequestering of the drug in the buccal cavity.

PRINCIPLE

PrimeScreen® One Step Multi-Drug Oral Fluid Test is a competitive immunoassay that is used to screen for the presence of drugs in oral fluid. It is a chromatographic absorbent device in which drugs or drug metabolites in a sample competitively combine to a limited number of antibody-dye conjugate binding sites.

When the sponge end of the collector is immersed into the oral fluid sample, the sample is absorbed into the device by capillary action, mixes with the antibody-dye conjugate, and flows across the pre-coated membrane. When sample drug levels are zero or below the target cutoff (the detection sensitivity of the test), antibody-dye conjugate binds to the drug/protein conjugate immobilized in the Test Region (T) of the device. This produces a colored band that, regardless of its intensity, indicates a negative result.

When sample drug levels are at or above the target cutoff, the free drug in the sample binds to the antibody-dye conjugate preventing the antibody-dye conjugate from binding to the drug-protein conjugate immobilized in the Test Region (T) of the device. This prevents the development of a distinct colored band in the Test Region (T), indicating a potentially positive result.

To serve as a procedure control, a colored band will appear at the Control Region (C), if the test has been performed properly.

PRECAUTIONS

- Not to be used for clinical diagnosis.
- Do not swallow.
- Discard after first use. The test cannot be used more than once.
- Do not use the test kit beyond expiration date.
- Do not use the test if the pouch is punctured or not sealed.
- Keep out of the reach of children.
- Do not read results after 5 minutes.
- The used collector and cube should be discarded according to local regulations.

MATERIAL

Materials Provided

- 25 PrimeScreen® One Step Multi-Drug Oral Fluid Test
- 5 Additional Sponge Collectors
- 25 Sponge Collectors
- One (1) Package Insert
- One (1) Procedure Card

Material Required but Not Provided

- Timer

STORAGE AND STABILITY

- Store at 4°C-30°C (39°F-86°F) in the sealed pouch up to the expiration date
- Keep away from direct sunlight, moisture and heat.
- DO NOT FREEZE.
- Preferably open the pouch only shortly before collection and testing.

SPECIMEN COLLECTION AND PREPARATION

Collect the oral fluid sample using the sponge collector provided. Instruct the donor not to place anything in the mouth including food, drink, gum, or tobacco products for at least 10 minutes prior to collection. No other collection devices should be used with this assay. Oral fluid collected at any time of the day may be used.

TEST PROCEDURE


Allow the kit and specimen to come to room temperature (65°F-86°F/18°C-30°C) prior to testing. AVOID PLACING ANYTHING IN THE MOUTH 10 MINUTES PRIOR TO TESTING.

- Remove the test cube and the sponge collector from the foil pouch by tearing at the notch. Place the test cube upright on a level surface.
- Put the sponge end of the collector on your tongue or near cheek to collect oral fluid for about 3 minutes until color on saturation indicator strip appears RED in the indicator window. If color on saturation indicator has not turned red at 7 minutes, repeat the collection using one additional sponge collector provided, beginning with Step 1.
- Open the test cube and place the fully saturated sponge collector inside the test cube. Press the sponge collector down firmly until it reaches the bottom of the test cube, then

close the cube lid tightly while compressing the collector. Keep test cube upright on flat surface and follow Step 4.

Note: Make sure the sponge collector is inserted vertically and the handle of collector is put into the clamp.

- Interpreting Drug Test Results:
Read results at 5 minutes. Do not read after 5 minutes.



INTERPRETATION OF RESULTS

Negative (-)


A colored band is visible in the Control Region (C) and the appropriate Test Region (T). It indicates that the concentration of the corresponding drug of that specific test zone is zero or below the detection limit of the test.

Preliminary Positive (+)

A colored band is visible in the Control Region (C). No colored band appears in the appropriate Test Region (T). It indicates a positive result for the corresponding drug of that specific Test Region (T).

Invalid

If a colored band is not visible in the Control Region (C), the test is invalid. Another test should be run to re-evaluate the specimen. If test still fails, please contact the distributor with the lot number.



Note: There is no meaning attributed to line color intensity or width

QUALITY CONTROL

Though there is an internal procedural control line in the test device of Control Region (C),

the use of external controls is strongly recommended as good laboratory testing practice to confirm the test procedure and to verify proper test performance. Positive and negative control should give the expected results. When testing the positive and negative control, the same assay procedure should be adopted.

LIMITATIONS OF PROCEDURE

- The test provides only a qualitative, preliminary result. A secondary analytical method must be used to obtain a confirmed result. Gas Chromatography/Mass Spectrometry (GC/MS) or Liquid Chromatography/Tandem Mass Spectrometry (LC/MS-MS) are preferred confirmatory methods.
- A positive test result does not indicate the concentration of drug in the specimen or the route of administration.
- A negative result may not necessarily indicate a drug-free specimen. Drug may be present in the specimen below the cutoff level of the assay.

PERFORMANCE CHARACTERISTICS

A. Analytical Sensitivity

Standard drugs were spiked into negative PBS pool to the concentration of 0% Cut-off, -50% Cut-off, -25% Cut-off, Cut-off, +25% Cut-off and +50% Cut-off. The results were summarized below.

Drug Conc. (Cut-off Range)	N	AMP		BAR		COC		MDMA		MET	
		-	+	-	+	-	+	-	+	-	+
0% Cut-off	30	30	0	30	0	30	0	30	0	30	0
-50% Cut-off	30	30	0	30	0	30	0	30	0	30	0
-25% Cut-off	30	28	2	25	5	25	5	25	5	28	2
Cut-off	30	12	18	10	20	10	20	10	20	10	20
+25% Cut-off	30	8	22	6	24	6	24	6	24	8	22
+50% Cut-off	30	0	30	0	30	0	30	0	30	0	30

Drug Conc. (Cut-off Range)	N	MTD		OPI		OXY		PCP		THC	
		-	+	-	+	-	+	-	+	-	+
0% Cut-off	30	30	0	30	0	30	0	30	0	30	0
-50% Cut-off	30	30	0	30	0	30	0	30	0	30	0
-25% Cut-off	30	25	5	14	16	14	16	26	4	14	16
Cut-off	30	12	18	10	20	14	16	14	16	14	16
+25% Cut-off	30	6	24	5	25	5	25	5	25	5	25
+50% Cut-off	30	0	30	0	30	0	30	0	30	0	30

B. Analytical Specificity

The following table lists the concentration of compounds (ng/mL) above which PrimeScreen® One Step Multi-Drug Oral Fluid Test identified positive results at the read time of 5 minutes.

Amphetamine (AMP)		Methadone (MTD)	
D-Amphetamine	50	Methadone	30
D,L-Amphetamine	125	Doxylamine	5,000
β-Phenylethylamine	4,000		
Tryptamine	1,500	Opiate (OPI)	
p-Hydroxyamphetamine	800	Morphine	40
(+)-3,4-Methylenedioxyamphetamine (MDA)	2,500	Codeine	100
Methamphetamine	11,000	Ethyl morphine	100
3,4-	100,000	Hydromorphone	1,000

Methylenedioxyamphetamine			
Dopamine hydrochloride	8,000	Hydrocodone	2,000
		Levorphanol	400
Barbiturates (BAR)		Morphine 3-β-D-Glucuronide	50
Secobarbital	20	Norcodeine	1,500
Amobarbital	30	Normorphine	12,500
Alphenol	15	Nalorphine	10,000
Aprobarbital	20	Oxycodone	>300,000
Butabarbital	10	Oxymorphone	25,000
Butathal	10	Thebaine	1,500
Butalbital	250		
Cyclopentobarbital	60	Oxycodone (OXY)	
Pentobarbital	30	Oxycodone	20
Phenobarbital	10	Dihydrocodeine	4,000
		Codeine	10,000
Cocaine (COC)		Hydromorphone	300,000
Cocaine	20	Morphine	11,000
Benzoylcegonine	100	Acetylmorphine	>100,000
Cocaethylene	25	Buprenorphine	>100,000
Ecgonine	40,000	Ethyl morphine	>100,000
Ecgonine methylester	12,500		

Phencyclidine (PCP)			
		Phencyclidine	10
3,4-Methylenedioxyamphetamine	50	4-Hydroxyphencyclidine	12,500
3,4-Methylenedioxyamphetamine HCl	300		
3,4-Methylenedioxyethylamphetamine	60	Marijuana (THC)	
		11-nor-Δ 9-THC-9-COOH	25
Methamphetamine (MET)		11-nor-Δ 8-THC-9-COOH	60
D-Methamphetamine	50	11-hydroxy-Δ 9-THC	2,500
Fenfluramine	10,000	Δ 8-THC	7,500
p-Hydroxymethamphetamine	400	Δ 9-THC	40
Methoxyphenamine	25,000	Cannabinol	1,000
3,4-Methylenedioxyamphetamine	500	Cannabidiol	10,000
L-Phenylephrine	4,000		
Procaine	2,000		
(1R,2S) - (-) Ephedrine	400		

C. Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds spiked into drug-free PBS stock. The following components show no cross-reactivity when tested with PrimeScreen® One Step Multi-Drug Oral Fluid Test at a concentration up to 100 µg/mL.

Acetaminophen	Ketoprofen
Acetophenetidin	Loperamide
N-Acetylprocainamide	Maprotiline

Acetylsalicylic Acid	Meprobamate
Aminopyrine	Labetalol
Amoxicillin	Meperidine
Ampicillin	Meprobamate
Ascorbic Acid	Methylphenidate
Apomorphine	Nalidixic Acid
Aspartame	Naloxone
Atropine	Naltrexone
Benzilic Acid	Naproxen
Benzoic Acid	Niacinamide
Benzphetamine	Nifedipine
D,L-Brompheniramine	Norethindrone
Caffeine	D-Norpropoxyphene
Chloralhydrate	Noscapine
Chloramphenicol	D,L-Octopamine
Chlorothiazide	Oxalic Acid
(±) Chlorpheniramine	Oxolinic Acid
Chlorpromazine	Oxymetazoline
Chloroquine	Papaverine
Cholesterol	Penicillin-G
Clonidine	Pentazocine
Cortisone	Perphenazine
(-) Cotinine	Phenelzine
Creatinine	D,L-Propranolol
Deoxycorticosterone	D-Propoxyphene
Dextromethorphan	D-Pseudoephedrine
Diclofenac	Quinidine
Diflunisal	Quinine
Digoxin	Ranitidine
Diphenhydramine	Salicylic acid
(-) Ephedrine	Serotonin (5-Hydroxytyramine)
β-Estradiol	Sulfamethazine
Ethyl-p-aminobenzoate	Sulindac
Fenoprofen	Tetracycline
Furosemide	Tetrahydrocortisone, 3 Acetate
Gentisic Acid	Thiamine
Hemoglobin	Thioridazine
Hydralazine	D, L-Tyrosine
Hydrochlorothiazide	Tolbutamide
Hydrocortisone	Triamterene
O-Hydroxyhippuric Acid	Trifluoperazine
p-Hydroxytyramine	Trimethoprim
Ibuprofen	D, L-Tryptophan
Iproniazid	Tyramine
Isoproterenol	Uric Acid
Isosuprine	Verapamil
Ketamine	Zomepirac

ASSISTANCE

If you have any question regarding to the use of this product, please call our Toll Free Number 1-888-444-3657 (9:30 a.m. to 5:00 p.m. CDT M-F).

BIBLIOGRAPHY OF SUGGESTED READING

- Moolchan, E., et al, "Saliva and Plasma Testing for Drugs of Abuse: Comparison of the Disposition and Pharmacological Effects of Cocaine", Addiction Research Center, IRP, NIDA, NIH, Baltimore, MD. As presented at the SOFT-TIAFT meeting October 1998.

- Kim, I, et al, "Plasma and oral fluid pharmacokinetics and pharmacodynamics after oral codeine administration", Clin Chem, 2002 Sept.; 48 (9), pp 1486-96.
- Schramm, W. et al, "Drugs of Abuse in Saliva: A Review," J Anal Tox, 1992 Jan-Feb; 16 (1), pp 1-9.
- McCarron, MM, et al, "Detection of Phencyclidine Usage by Radioimmunoassay of Saliva," J Anal Tox. 1984 Sep-Oct.; 8 (5), pp 197-201.

ADDITIONAL INFORMATION AND RESOURCES

The following list of organizations may be helpful to you for counseling support and resources. These groups also have an Internet address which can be accessed for additional information. National Clearinghouse for Alcohol and Drug Information www.health.org 1-800-729-6686 Center for Substance Abuse Treatment www.health.org 1-800-662-HELP The National Council on Alcoholism and Drug Dependence www.ncadd.org 1-800-NCA-CALL American Council for Drug Education (ACDE) www.acde.org 1-800-488-DRUG

INDEX OF SYMBOLS



Keep away from sunlight



Store between 4°C - 30°C (39°F - 86°F)



Keep dry



Do not re-use

Distributed by Wondfo USA Co., Ltd.

545 Willowbrook Centre Parkway, Unit B
Willowbrook, IL 60527

www.wondfousa.com

Made in China

Rev. A1
Rel.: 2021/11/02