


Butterfly Transfer Board

DRAWING

SKU: ABEZ012



Specially designed transfer board that is used to transfer patients between bed and wheelchair or from wheelchair to car. Its special curve design can work well around wheelchair arms. It is lightweight and easy to carry.

PRODUCT SPECIFICATIONS

Size: 75.5 x 34.5 x 1cm ($\pm 3\%$)

Weight Capacity: 130Kg

MDA Product Reg. No.: GA10776223-139094

SAFETY INFORMATION

- Before using the transfer board, please carefully bend the board to check its strength and durability. Visually check to see whether there is any crack shown. Do not use the board if there is any damage occurred.
- Before transferring the patient, be sure both bed and wheelchair are securely locked and two surfaces are levelled.
- The max weight is tested for 130Kg when the gap < 30cm, 150Kg when the gap < 20cm.
- Do not sterilize EZ-510 at high temperatures.

OPERATING INSTRUCTION

Place the wheelchair as close to the bed as possible. Butterfly Transfer Board's special curve design works well around the wheelchair's armrest, so it can be used for those wheelchairs where their armrests cannot be removed. Butterfly Transfer Board works as a bridge between the bed and wheelchairs. Place the board underneath the patient's hip area and the other side on the other surface. To position the Butterfly Transfer Board under the patient, lean the patient slightly in the opposite direction. Remove the board after transfer. The same application can be used for transferring from wheelchair to car seat.

