

■ ■ wissner-  
■ ■ bosserhoff



## sentida

Universal low nursing beds





# Safety and comfort in perfection

sentida – universal low nursing beds

Falling is one of the greatest risks for the physical health and psychological well-being of elderly people. For this reason, systematic fall prevention is of crucial importance in care for the elderly, making it an important quality feature of a nursing home. The consequences of falls are often very serious for elderly people and, in the worst case, even irreversible. The fear of further falls is one of the most common psychological consequences, leading to limited physical activity, a loss of self-confidence and social isolation.

With the sentida range of low nursing beds, wissner-bosserhoff has achieved a consequent and above all uncompromising implementation of the requirements of modern nursing care: High safety standards hand in hand with a cozy design.

Innovative solutions such as the 3-stop strategy and the patented SafeFree® side guard concept reduce falls as well as the consequences associated with them and set standards in terms of safety, mobility, height adjustability and coziness.

5

million  
falls of elderly  
people per year:

Almost one third of people aged 65 and above, as well as 50% of those aged 80 and above fall at least once a year.

In Germany alone, there are about 5 million cases of elderly people falling each year.

(Source: "Gesundheit und Krankheit im Alter" (= health and illness in old age).  
Publisher: Karin Böhm, Federal Statistical Office of Germany / Clemens Tesch-Römer, German Centre of Gerontology / Thomas Ziese, Robert Koch Institute, Berlin 2009)



Prevention  
with sentida  
low nursing beds


for use in institutions with high standards of universal care close to the ground: Nursing homes and homes for the elderly, care resorts, palliative wards, dementia and geriatric wards and rehabilitation institutions.



## Living concept Feathers

sentida | Design O10 | Wood dekor Okapi Walnut

Nature as an example for the choice of materials and colors. As harmoniously coordinated as a bird's plumage – this is how a comfortably furnished room was created which provides residents with security and protection. From the window seat you can hear the morning birdsong. This concept is just a suggestion – we will happily give you individual advice.



# sentida Portfolio

The entry model



The mobile one



The flexible one



The cozy one



The extra mobile one



Model	sentida 1	sentida 3	sentida 4	sentida 5	sentida 6
Movability	In low position	At every bed height	At every bed height	At every bed height	At every bed height
Castors	Two fixed castors at the foot end Two 100 mm castors at the head end	Four 50 or 75 mm castors at the foot end double castors Two 100 or 160 mm castors at the head end fixed castors	50 or 75 mm double castors	50 mm double castors, completely covered by cozy undercarriage cladding	100, 125 or 150 mm easy-roll castors or 125 mm double or design castors
Braking capacity	Electric brake with handset	Electric braking with handset foot-end brake bar	Head-end and foot-end brake bar	One pedal brake for each axle with 2-castor central braking system	Foot-end brake bar with 4-castor central braking system
Safe working load	225 kg	225 kg	225 kg	270 kg	270 kg
Height adjustment	25-80 cm	25-80 cm	25-80 cm	27-80 cm	27-80 cm with 100 mm easy-roll castors 29,5-82,5 cm with 125 mm easy-roll castors 28-81 cm with 125 mm double or design castors 32-85 cm with 150 mm easy-roll castors
Patient surface	87/100 x 200/210/220 cm	87/100 x 200/210/220 cm	87/100 x 200/210/220 cm	87 x 200/210/220 cm 100 x 200/210/220 cm (sentida 5-xl) 120 x 200 cm (sentida 5-xxl)	87/100 x 200/210/220 cm (sentida 6 & 6-d*) 100 x 200/210/220 cm (sentida 6-xl) 79 cm x 200 cm (sentida 6-s) 79 cm x 190 cm (sentida 6-xs)

Also available as dismantlable version, see p. 23





## Benefits at a glance



Protection and mobility without FDM* SafeFree® – the patented side guard concept	12 - 13
Side guards – individual according to requirements SafeFree® FLex	14 - 15
Safe nursing care, mobilization & sleep 3-stop strategy	16
Flexible & time-saving 3 bed lengths that can be adjusted without tools	17
Bedsore prevention Patient surfaces – ergonomics & modularity	18
<b>NEW:</b> Manual mobilisation assistant MobiStick 2	19
<b>NEW:</b> Clearly arranged. Ergonomic. Intuitive. SafeControl®   Handset	20
Improving and maintaining mobility SafeLift   Revolutionary bed operation for residents	21
Safety by day and night SafeSense® – the intelligent bed exit sensor system	22
<b>NEW:</b> SafeSense® 3 From bed exit sensor technology to digitisation platform	23
Energy efficiency & stability Thanks to innovative drive concept	26
LED-Designer reading light soluna & lymera	27
The one easy to disassemble sentida 6-d	28
We think ahead memoriana   Room concept for dementia patients	29
Additional options Design   Colors and decors   Features and accessories	30



## Living concept **Urban**

sentida | Design HD-3 | Wood decor Havana cherry | Fabric Stone | Padding Deep buttoned

More and more people choosing to forgo retiring to a quiet cottage in the country and are instead choosing to spend their time in the city surrounded by life and culture. Retiring in the city is not about relaxing in tranquillity after years of hard work but about making the most of new found freedom, experiencing new things, forging new friendships and re-discovering yourself in the heart of the ever changing urban landscape. This is just a suggestion – we will happily give you individual advice.



# Protection and mobility without FDM\*

SafeFree® the patented side guard concept

The SafeFree® divided side guard system is patented and voluntarily certified. The technology has been tested in practice more than 1,000,000 times. You can lower the space-saving side guards quickly in only 2 seconds and individually adjust them to four heights, depending on the care situation and size of the resident. Obtaining freedom instead of taking it away, in agreement with the Werdenfelser Weg (Werdenfelser care pathway)\*\*.

**1**  
No protection needed



Lowered, the divided side guards look homely and provide unhindered access to the resident for daily routine nursing tasks. No need to take off and store accessories that are in the way, such as a middle post or fixed side guards.

**2**  
Soft protection (no FDM\*)



For safety at night, it is often sufficient to combine the low position of the bed with the first height level of the head-end side guard. This avoids unnecessary barriers and effectively reduces the risk of falls.

**3**  
¾ protection (no FDM\*)



Creates a deliberate bed egress gap at the foot end to reduce freedom-depriving measures (FDM) while at the same time maintaining a high safety standard.

**4**  
Full protection (FDM\*)  
for mattresses up to 18 cm



The highest position provides comprehensive protection of the resident. The distance between the two parts of the side guards is minimal, making them just as safe as a continuous side guard and exceeding the requirements of the standard IEC 60601-2-52. It is not necessary to stock up on accessories to close a central gap when using an abdominal belt for restraint.



If you use a thicker mattress, you can raise the side guards by an additional 6 cm.


Learn more:



Videos on the side guard concept




More safety due to smaller central gap\* < 60 mm




- Targeted compliance with the standard IEC 60601-2-52
- Safety in working position, i.e. also with raised backrest or leg rest
- Safety without using an additional central gap protector as accessory for restraint with an abdominal belt in accordance with the BfArM\*\* (German Federal Institute for Drugs and Medical Devices)

All-round access for nursing tasks = with 4 divided side guards



- Flexible use of the bed for standard care, even in limited spaces
- Flexible access from both sides of the bed in case of increased need for care
- Flexible use of the bed for bed egress on the left or right side to mobilize the resident
- Avoiding follow-up costs for upgrading divided side guards

Safe mobilization regardless of body size



Divided side guards are very useful for effective mobilization, as they are easy to grab and clasp. They can be used as a safe standing-up aid and verifiably support the resident's mobility.

2 Automatic positions



Watching television, reading or eating breakfast in bed with automatic comfort seating



Bedsores prevention and keeping heels free with automatic positioning

# SafeFree® Flex

Side guards as individual as your needs

In addition to permanently mounted side elements, individual and future-oriented SafeFree® Flex offers as much freedom as possible and as much protection as necessary. At the same time SafeFree® Flex retains the homely design which distinguishes all wissner-bosserhoff nursing home beds.

0

No protection needed



Comfort bed, without side guards

1|2

Soft protection (no FDM\*)



Mobilization aid only the head end elements

3

¾ protection (no FDM\*)



Wall protection with foot-end exit

4

Full protection (FDM\*) for mattresses up to 18 cm



4 fully adjustable side guards

## Safe

Complete protection



- Secure fastening of elements
- Prevents an accidental loosening or removal of the side elements

Seamless safety



- Gap closure between the side panel and the patient surface
- Protection against accidents through getting trapped or wedged

Maximum mobilization



- Optimum length and individually adjustable height of the head-end side guard
- Perfect hip protection when sleeping as well as comfortable and ergonomic mobilization of the resident



## Free

Fewer FDM\*




- Conscious "omission" of side guards
- Investment security for the operator Legal certainty regarding FDM\* for the nursing staff

## Flex

More economic efficiency

**Side elements can be adjusted without tools**  
Your budget can be used as flexibly as the side elements can be attached, because you only pay for the elements that your residents really need.

**No technician needed**  
Quick modification for the changing resident situation – with gap closure without (1) and with a handrail (2) or fully adjustable side guards (3)



Benefit from a homely design with a continuous side panel, independent of the required mobilization or resident protection situation:

Learn more:



Video






## Safe nursing care, mobilization and sleep

### 3-stop strategy

- Comfortable adjustment range from 27 to 80 cm
- First automatic stop at 41 cm at chair height for bed egress
- Collision protection when adjusting the bed from the low position and when using a bed extension


### 80 cm

Optimal nursing position




### 41 cm\*

Safe mobilization at chair height



### 27 cm\*

Fall prevention during sleeping phases




## Flexible and time-saving

### 3 bed lengths that can be adjusted without tools

- All standard wissner-bosserhoff nursing home beds are equipped with an integrated bed extension.
- Adjusted flexibly and without tools to different bed lengths.
- The requirements of the standard IEC 60601-2-52 are complied with in each configuration.


### 200 cm

Bed length of 200 cm:  
The divided side guards offer protection without a central gap > 318 mm\* and without a protector




### 210 cm

Alternative bed length of 210 cm:  
Protection with longer side guards at the foot end without accessories and without a gap > 318 mm\*\*



### 220 cm

Bed extension to 220 cm:  
is temporarily adjustable with divided side guards for sentida - complete protection through foot-end protector (accessory) is guaranteed



# Bedsore prevention

## Patient surfaces – ergonomics & modularity

The so-called double retraction has been a component of the wissner-bosserhoff patient surfaces for years. Here, not only the backrest but also the leg rest is moved away from the seat section, so that the surface area is increased. Results of pressure measurements clearly show that the pressure load with the double retraction is considerably lower than with a standard lying surface.

approx. **37.5 %**  
of all residents in German old people's or nursing homes are subject to an increased **risk of decubitus**


Source: Own survey, wissner-bosserhoff 2020

approx. **2.7 m.**  
People suffer in Germany from chronic wounds like **decubitus**


Source: Study BVMed 2012

## Our solutions

WIBO patient surface with notable relief




Double retraction = decubitus prevention




Thermal image of a BiNetic patient surface

Standard patient surface without relief




Single retraction = risk of decubitus



Thermal image of a standard patient surface


Patient surface modularity  
3 different surfaces



**wire mesh patient surface Aero (Standard)**  
→ good aeration



**Synthetic mats EasyClean**  
→ easy to clean



**Comfort patient surface**  
with its 55 spring elements  
→ bedding comfort and simple pressure relief

↓  
Double retraction =  
**17cm**  
more space in bed

# Manual mobilisation assistant

## MobiStick 2

The MobiStick 2 enables clever bed operation for residents. Its ergonomic shape makes it easy and intuitive to use anywhere. Residents can effectively support and move themselves.

**NEW!**



Due to its ergonomic shape MobiStick 2 is easy and intuitive to use and universally applicable. The hand knob is rotatable and has a simple locking.



When not in use, it can simply be folded away and is therefore not in the way, but quickly ready for use when needed.

↓  
**Without tools**  
MobiStick 2 can be attached without tools on both sides of the bed, is suitable for every resident size and can be retrofitted for wissner-bosserhoff beds from 2004 onwards.


Learn more:





**Advantages at night**

- SafeControl® handset is activated by movement and displays all selectable functions
- Unnecessary lighting at night is avoided
- Respect for the night's rest of the residents



## Clearly arranged. Ergonomic. Intuitive.

SafeControl® | Handset



In addition to its clarity and ergonomic design, the SafeControl® handset stands out due to its intuitive operation. Thanks to motion sensors and controlled backlighting, all selectable functions are displayed to the operator, ensuring easy operation. This ensures easy operation, especially at night.



### 3 safety levels

Key combinations can be used to different handset modes can be set:

**Care mode:**

All selection buttons are unlocked

**Locking mode:**

All functions of the handset are locked

**Resident mode:**

The light grey selection keys are available to the resident

**Thanks to the integrated locking function an additional magnetic key is not needed.**

### Comfort & Safety

- Clear and ergonomic
- First error-proof thanks to preselection
- Intuitive operation
- Motion sensors
- Backlighting
- No additional key/accessories necessary
- Flexible use due to long cable

## Simple bed operation & safe mobilization

SafeLift | Product benefits at a glance



SafeLift offers the simplest operation so that residents can adjust their bed independently and set it to the best mobilization position for them. At the same time, they can support themselves when standing up on the SafeLift controller. This promotes residents' self-mobilization, and reduces staff calls for bed adjustment.



- Adjustment buttons for bed height and backrest
- Large buttons
- Firm position on a SafeFree® side guard
- GO button
- 360° rotation
- Ergonomic design
- Strong stability



\*Exceptions: carisma 300, sentida 1-p / 4-p, sentida 7-i

# Safety by day and night

## SafeSense® – the intelligent bed exit sensor system

wissner-bosserhoff, your reliable partner for the reduction of freedom-depriving measures (FDM): With SafeSense®, wissner-bosserhoff presents an intelligent assistance sensor system, which reliably helps reducing freedom-depriving measures day and night. It is equipped with an automatic night light, which is activated when the resident leaves the bed and is deactivated when the resident is back in bed again. SafeSense® is compatible with standard nurse call systems and can be retrofitted to wissner-bosserhoff nursing home beds manufactured in 2002 or later.



### Automatic night light

- Night light is automatically switched on when the resident leaves the bed
- Night light is switched off again as soon as the resident is back in bed

### Individual nurse call

- Freely adjustable immediate call from 0 seconds to 30 minutes
- Adaptation to the resident's mobility

### Integrated or available as retrofit kit\*

- Accessory option for current WIBO nursing home beds
- Retrofit kit for electrical WIBO nursing home beds manufactured in 2002 or later
- Available with cable or wireless



Video

## More safety – less freedom-depriving measures (FDM)

- Fast response time in the event of falls
- Safety during night and weekend shifts
- Fewer inspections
- Ideal in combination with sentida low nursing beds and SafeFree® side guards

# SafeSense® 3

## From bed exit sensor technology to digitisation platform



After more than 6 years of experience with integrated bed exit sensor technology, the time is right for a real digital care assistant. Since 2014, wissner-bosserhoff has been able to gain valuable experience with digital products with the bed exit sensor system SafeSense® and the intelligent nursing bed sentida 7-i. During this time, SafeSense® has matured into a care assistant with further supporting functionalities and is recommended as a future-proof digitalisation platform thanks to its open architecture. This can be easily retrofitted with further functions that additionally relieve the burden of care.



### Satisfied residents through higher quality of care

- 24 h movement monitoring provides valuable information about lack of movement, nocturnal restlessness and lack of sleep
- Support in the prevention of risk of pressure sores and falls
- Preservation of independence and mobility through bed-exit monitoring
- Adaptation of the safety level to the individual resident risk



### Technology to „touch“-the intelligent care bed

- Assistance system for continuous **resident monitoring in real time**
- Contactless **data acquisition and visualisation**
- Differentiated **call management** for faster **prioritisation**
- **Bed presence** - Resident in bed YES or NO

## Intelligently networked - digital care in real time

- Pressure ulcer prevention through movement monitoring
- Fall prevention through bed-exit monitoring
- Future-proof sensor platform through open system architecture
- Mobile care monitoring through connectivity and mobile devices



# Living concept Lilac

sentida | Design G10 | Wood decor Salisbury Elm natural

Timeless and classic – this choice of colors and materials gives the room a special atmosphere. The trend color of lilac is reminiscent of spring – everything is blossoming and fresh colors replace the winter gray. This concept is just a suggestion – we will happily give you individual advice.



# Energy efficiency and stability

Thanks to innovative drive concept

With the new motor, the sentida 5 and 6 nursing home beds are even more user-friendly, stable and energy-efficient. The integrated LED technology, energy-saving features and low-maintenance design are vital requirements for the cost-effective operation of a nursing home.



## Even more user-friendly

The new sentida 5|6 is equipped with an LED bedside light connection integrated in the motor. You only need one single cable from the bed to the power supply. The LED light can also be switched on and off using the handset. You can now also watch video tutorials on how to operate the bed by scanning the QR code on the bed.



Video

## Stable and safe

The new motor is especially powerful. Its performance is especially impressive under high load. It has a safe working load of up to 270 kg (sentida 5|6) and is suitable for a maximum resident weight of up to 235/255 kg (depending on the accessories used).



## Economical and energy-efficient

With continuously increasing electricity prices, energy efficiency is an increasingly important factor. With an extremely low stand-by consumption of only 0.5 W, the new motor is especially economical. A clever cable management reduces wear and tear. The power cable can be exchanged separately, causing less storage costs.

The following calculation example demonstrates how the stand-by consumption of nursing home beds can influence the costs:



Product	Stand-by operation	Electricity costs over the course of 15 years*
The NEW sentida	0.5 W	€1,950
Bed A	3.76 W	€15,345
Bed B	2.5 W	€10,200
Bed C	1.7 W	€6,900

\*Example: Institution with 100 beds

# LED lights

soluna & lymera

The undulating shape of the lampshade, with its partially translucent character, is prominent as a special design element of soluna lights. Equipped with a swivel arm, the light allows the light focus to be precisely positioned close to the object to be illuminated, regardless of whether the light is attached to the bed or the bedside cabinet.



Comfort LED reading light lymera

## Advantages of LED

- Long service life
- Low energy consumption
- No heat build-up
- Safe operation



Design LED light soluna with wavy lampshade



# The one easy to disassemble

sentida 6-d

At first glance, the sentida 6-d cannot be distinguished from a standard sentida bed. The technical parameters of both beds also are identical. You can only see the difference once you take a closer look – the sentida 6-d can be completely disassembled and reassembled by one person – quickly and without tools. None of the individual components weighs more than 18 kg, complying with all occupational health and safety requirements.

The sentida 6-d is used in home care or nursing homes with narrow doors – without compromising design, quality or stability.



sentida 6-d dismantled  
No assembly weighs more than 18 kg.

# We think ahead

memoriana | Room concept for dementia patients

Adequately taking care of people with dementia is a particular challenge for nursing staff and relatives. In cooperation with dementia experts from Great Britain and Germany, wissner-bosserhoff has developed a room concept especially for dementia patients. Specifically adapted to the needs of the patient, the new concept offers a considerable reduction in the strain on those providing care. The focus of the development includes three targets for the care of residents with dementia:



## Improving orientation

The resident is offered reference points for orientation in time and space, as well as colored orientation aids. The targeted integration of colors, forms and objects helps in terms of spatial orientation – worries and frustrations are reduced.



## Increasing independence

The more independently the resident can move in his surroundings, the less strain on the nursing staff. Functional furnishing elements such as handrails and open storage areas promote the dementia patients' mobility and independence. Other areas of the room are, however, only accessible to those providing care.



## Maximizing safety

Resident safety is of utmost importance. Low nursing beds with side guards, handrails on furniture and walls as well as visual orientation aids provide the person diagnosed with dementia a maximum level of safety without feeling restrained.



# Colors | Wood décors\*



NEW: Crystal white (U11026)\*2



NEW: Magnolia (U11509)\*2



NEW: Norway maple (R27044)\*2



Lindberg Oak (R20021)



NEW: Nordic Teak (R50094)\*2



NEW: Italian Oak (R20286)\*2



NEW: Salisbury Elm natural (R37016)\*2



Havana Cherry (R42006)




NEW: Altmühl beech (R24006)\*2



Steamed Beech (R24034)



NEW: Okapi Walnut (R30135)\*2



Color combinations for the memoria room concept\*3

# Certified security

Voluntarily, wissner-bosserhoff has its manufacturing processes and management regularly audited by TÜV Süd. The company is certified in accordance with ISO 9001 for quality management systems, DIN EN ISO 13485 for the manufacture of medical products and ISO 14001 for environmental management.



ISO 9001:2015 | ISO 13485:2016 | ISO 14001:2015 | CB certificate for product safety, TÜV Süd | Gutachten zu Seitensicherungen TÜV Süd



\*1 = Some decors available at a surcharge | \*2 = Available from 07.2021 | \*3 = Matching furniture for this can be found in the „memoriana“ catalogue  
K1 = Nordic Teak / Italian Oak | K2 = Pointed Maple / Elm Salisbury Natur | K3 = Cherry Havana / Okapi Walnut

# Synthetic leather & fabric covers


New combination options  
We will be happy to advise you!



With our new covers made from high-quality materials by renowned partners we offer an excellent surface material with the following properties:

- Resistant to moisture and dirt
- Breathable
- Water-repellent
- Blood and urine resistant
- Disinfectant resistant
- Extremely hard-wearing
- Environmentally friendly and non-toxic
- Durable
- Flame retardant (in accordance with the ignitability criteria in accordance with DIN EN 1021 Part 1 + 2 and BS 5852 Crib 5\*)

## Synthetic leather



Silver gray



Pearl



Birch



Pebble



Cashmere



Smoke



Cherry



Corn




Lime




Royal

## Fabrics



Green, floral (Aurora 6520)




Blue, striped (Agon 5520)




Gray, striped (Agon 8520)



Gray, polka dot (Flint 8520)



Meadow (Mailo 6520)




Sea (Mailo 5520)



Sand (Mailo 1521)



Stone (Mailo 8520)



Rose (Mailo 4520)

\* Flame retardancy depends on the foam material used. | Item numbers: 50-0818, fabric covers for head section Df | 50-0819, fabric covers for foot section Df



# Design | head / foot sections

## Head section



Design A  
Head-/Foot section



Design D



Design Df



Design G



Design K



Design O



Design T



Design Ts

## Foot section



Design Ak  
Head-/Foot section



Design D



Design Df



Design G



Design K



Design O



Design T



Design Ts

## Low foot section



Design Ak  
Folded down



Design Dx



Design Dfx



Design Gx



Design Kx



Design Ox




Design Tx




Design Tsx

# Accessories | SafeSense® & SafeLift

## SafeSense® sensor mat 2.1



**02-000876**  
SafeSense® sensor mat 2.1  
80 x 20 cm




**10-0812**  
Digital bed exit box 2  
cable or wireless version

Also available  
as a set:

Sensor mat 2.1  
& bed exit box 2.1  
**10-0871**  
Cable version  
**10-0872**  
Wireless version

## SafeSense® connection kit



**04-1508**  
Adapter cable with on/off switch and  
suitable plug for connection to  
the nurse call system (customer-specific)  
for the cable and wireless version




**04-1608**  
Adapter cable without on/off switch;  
with suitable plug for connection to  
the nurse call system (customer-specific)  
for the cable and wireless version

## SafeLift bed operation



**10-0811**  
SafeLift Controller




**10-0808**  
SafeLiftBox


As a set:

**10-0799**  
SafeLift set

## Retaining plate SafeSense® and SafeLift




**01-011193**  
Retaining plate for DU09 motor




**01-011194**  
Retaining plate for Q7 motor

# Furnishings and accessories

## Dimensions




sentida with divided (10) side guards, without gap  
Outer dimensions: L = 208 cm, W = 103/108 cm  
Mattress dimensions: L = 200 cm, W = 87 cm  
H1 25/27 cm (low position), H2 80 cm (highest position)




sentida with continuous 2-part  
09 side guard with trigger release


## Patient surfaces



Comfort




EasyClean




Aero (Standard)


## Bed extensions




Temporary bed extension divided  
10 side guard (standard)




Permanent bed extension divided  
10 side guard (optional)



Permanent bed extension  
09 side guard (optional)




Bed extension with express release  
(optional)



10-0570/10-0580  
Protector for bed extension  
by 20 cm (optional)


## Options




Vertical wall stopper castors



3D horizontal and vertical wall stopper castors




CPR release on the backrest




Plug & play use of the handset on one or both sides of the bed

01-010812, 01-010813  
Underbed light switchable (not shown)


## Accessories




10-0419  
Handset holder for divided side guards  
10-0421  
Handset holder, short SD-Handset  
10-0790  
Handset holder, short SafeControl




10-0422  
Accessories and handset holder for continuous side guards  
10-0421  
Handset holder, short SD-Handset  
10-0790  
Handset holder, short SafeControl




50-0345  
Shortened urine bottle basket  
50-0164  
Urine bag holder (without illustration)




10-0817  
Gap closure base (Image similar)




10-0818  
Gap closure design




10-0819  
Side element Flex




10-0634  
Mobilizations-Assistant-MobiStick 2




10-0715  
Standing-up aid for 09 side guards as well as 3/4 side guards




50-0550 Design LED light soluna without plug power supply




50-0677 Comfort LED reading light lymera (with & without power supply) Transparent lampshade




50-0771  
Fall prophylaxis mat with anti-slip W 90 x L 200 cm, height approx. 4.5 cm




50-0174  
Cushioning for divided 10 side guards




50-0338  
Cushioning for 09 two-part side guards




10-0563/10-0533/10-0521  
Segufix bracket for thigh-backrest- for the seat part




10-0617  
Integrated design linen holder Compatible with all beds (except 09 3-part side guards and A-kk design)



10-0719  
Linen holder



Head-end wall spacer, telescopic, only for sentida 5 and 6  
10-0801 short, 1 x 20 mm  
10-0802 medium, 4 x 20 mm  
10-0803 llong, 8 x 20 mm  
Head-end wall spacer, only for sentida 5 and 6 (no image)



Side-mounted wall spacer, telescopic, only for sentida 5 and 6  
10-0804 short, 1 x 20 mm  
10-0805 medium, 4 x 20 mm  
10-0806 long, 8 x 20 mm  
Side-mounted wall spacer, only for sentida 5 and 6 (no image)

**wissner-bosserhoff GmbH**

Hauptstr. 4-6 | 58739 Wickede (Ruhr) | Germany

Phone: +49 2377.784-0 | Fax: -163 | [info@wi-bo.de](mailto:info@wi-bo.de) | [wi-bo.com](http://wi-bo.com)

