

T-POD™

Instruction and Use

Simple, Easy and Fast 5-Step Application Complete, safe stabilization and circumferential compression of the pelvic region is quickly and easily achieved in 5 steps:

1 Sliding belt into position

Medical personnel slides T-POD belt under the patient.



2 Making proper size adjustment

Trim the belt to leave a 15-20 cm gap in the center of the abdomen.

This will provide for simultaneous circumferential compression.



3 Applying pulley system

The velcro-backed mechanical advantage pulley system is easily attached.



4 Applying Compression


Medical personnel effortlessly draws the pull-tab on the mechanical advantage pulley system, creating simultaneous circumferential compression around the pelvic region.

This results in the reduction and stabilization of the open-book fracture and attenuation of blood loss.



5 Securing Compression

The velcro-backed pull-tab is easily attached to the belt, securing compression.



Re-application

Application time and date are recorded on the writing space provided.

