STATIVWAAGE
STAND SCALE
BALANCE A COLONNE

7730, 7731

Stativwaage 7730

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STAND SCALE

Thank you for choosing this Soehnle Professional product. This product is equipped with all features of the very latest technology and is optimised to make operation as simple as possible. If you have any questions, or if there is a problem with the equipment that cannot be resolved with this user manual, please contact your dealer or service partner, or visit our website at www.soehnle-professional.com.

1. MODELS OF SCALES

7730.01.001 Stand scale without height rod, approved

7731.01.001 Stand scale with height rod, approved

2. FUNCTIONS

2.1 Mounting the scale and start-up

Set up the scale on place of use and ensure that it is placed horizontal on a solid, level ground. According to the bubble adjust the feet and turn the balance to a horizontal position.

Connect the mains adapter to the scale.

Please pay attention to the free accessibility of the power socket.

Before the initial start-up of the scale, the battery must first be charged for 5 hours!

Calling up the approval count:

Display with mains operation:

After connecting the mains plug with the power supply the approval count is shown for two seconds.

Display with rechargeable battery operation:

After pressing the ON/TARA button the approval count is shown for two seconds.

2.2 Switching on

Display with mains operation:

By stepping on the platform the scale turns on automatic.

Display with rechargeable battery operation:


The scale turns on through pressing the ON/TARA button.

2.3 Weighing

Stand on the scale. The weight appears automatically after the scale is loaded.

2.4 Taring

Through pressing the function key on the back of the display the scale can be tared out. Place towels or other coatings on the scale and press the ON/TARA button. The display goes back to »0.0 kg« and the scale is ready for use. The taring can be done over the entire weight range.

2.5 Height rod with BMI-display

Fold up and slide up the height rod carefully. Then place the head stopper on the head of the patient. After stabilisation of the weight, the values of the patient - weight, height and BMI are shown rotational.

2.6 Weighing small children


An adult steps on the scale.

> The display shows »57.6 kg« for example

Press the ON/TARA button.

> The display changes to »0.0 kg«

The adult steps off the scale ...

> The display shows »-57.6 kg« or shows underload stripes 

... then steps on the scale again holding the child.






> The display now indicates the child's weight, for example »5.1 kg«

2.7 Switching off / End the operation safely

Display with mains operation: The scale automatically switches off, when the display has remained at »0.0 kg« for about 60 seconds. To switch off the display unit completely, disconnect the plug from the power socket. By pulling out the plug from the power socket, the device is all-pole disconnected in all from the power supply. So the operation is safely completed.

Display with rechargeable battery operation: The scale switches off automatically, when the display has remained at »0.0 kg« for about 60 seconds. Otherwise you can press the ON/TARA button, if the display shows »0.0 kg«. The display can be reactivated by pressing the ON/TARA button on the back.

3. FAULTS – CAUSES AND REMEDIES

DISPLAY	CAUSES	REMEDY
	The scale is not reset to »0.0 kg« when switched on.	Unload the scale and restart the display unit. If the error remains, please contact your service partner.
	Underload: only dashes appear at the bottom of the display field.	Switch the display off and on again, or pull out the mains plug. If the error occurs again, please contact your service partner.
	Overload: only dashes appear at the top of the display field. Maximum weighing range is exceeded.	Unload the scale, or switch them off and on again. Check the maximum weighing range.
	The weighing value does not stabilise.	Check the surface under the scale, and check that the weighing platform is level. If the error persists after readjustment, please contact your service partner.
	The weighing value does not stabilise.	Stand still or stabilise the weight on the scale or restart the scale.
	Height rod is not in reference position during power-up	Fold up the height rod during power-up.
	The rechargeable battery is empty.	Charge the battery.

If the fault cannot be cleared, or other error messages appear, please contact your service partner.

4. GENERAL INFORMATION

4.1 Technical data

The weighing range and scale interval are indicated on the data plate on the scale.

Operating temperature:	+ 10 °C to + 40 °C
Functioning temperature:	- 10 °C to + 50 °C
Storage temperature:	- 20 °C to + 65 °C
Moisture:	20 % to 85 % relative moisture non-condensing
Air humidity:	950 ... 1050 hPa
Rechargeable battery operation:	7,2 Volt
To use power adapter:	Power adapter 618.020.059 (plug-on)
Mains operation:	100 ... 240 V / 50-60 Hz, 250 mA

The technical specifications of the power supply are printed on the power adapter.

4.2 Intended use

The device is used to weigh people and may only be used for this purpose. Any other use is prohibited. This scale is designed for calibrated applications. It meets the applicable requirements of the EC guidelines 2014/30/EU, 2014/31/EU, 2014/35/EU and 93/42/EEC.

Messages about errors that can endanger the patient, as well as errors that lead to false measurement results are to be notified to the manufacturer's medical advisors in writing. This user manual is part of the device. Strict adherence to these instructions is a precondition for the proper use and correct operation of the device.

Soehnle Professional only accepts the responsibility for the safety of the device, if the instructions are followed and the device will be operated in accordance with the user manual. The device is a medical device and should only be used by persons who possess knowledge based on their training, which guarantees the correct handling. The user has to check the safety function and the condition of the appliance before using the device. The operator must be familiar with the operation of the device. The device is not intended for potentially explosive areas of medical rooms. Explosive danger areas may result from the use of flammable anesthetics, skin cleansing and skin disinfectant.

4.3 Classification

Medical product of class I with measuring function.

Electrical Safety Class II (safety isolated, no protective conductor connection).

Water protection according to EN60529: IPX2 for the whole device.

Application part

The medical product is a Type B application part and intended for direct contact with the patient. The leakage currents correspond to the classification of type B application parts.

4.4 Safety information

This device is designed for professional use. Please read carefully the information listed in the user manual before using the device. It contains important information regarding installation, intended use and maintenance of the device. The manufacturer is not liable if the following instructions are not followed: When using electrical components under increased safety requirements, the relevant provisions are complied with. Never work on the device, if it is energized.

There are no special requirements for the staff, the scale can be used anywhere. The use in critical environments is not provided (CT, defibrillator). The product is not for use in explosion-prone environments. The scale with mains operation must be set up, so that the separation of power supply and scale can be easily done (ensure accessibility of the socket). The display and the scale may not get wet from below.

Raise the scale in order to move it; do not drag it along the ground. Afterwards check the levelling and adjust if necessary.

In accordance with EC guideline 2014/30/EC, this device is equipped with interference suppression. However, the display reading can be affected by extreme electrostatic or electromagnetic interference, for example if a radio or mobile phone is used immediately beside it. The display should return to normal once the source of interference is removed. Possibly a restart of the scale is necessary. If there is permanent electrostatic interference in the area, please contact your responsible service partner.

This equipment is a measuring instrument. Air draught, vibrations, rapid temperature changes and direct sunlight may affect the weight results. The scale conforms to protection class IPX2. Avoid high levels of air humidity, steam, aggressive fluids and heavy soiling.

Improper installation will invalidate the warranty. The electrical connection conditions must be consistent with the values printed on the power supply unit. The equipment is designed to be used indoors. Note the permitted ambient temperatures for operation (see the Technical data). This equipment meets the requirements for electromagnetic compatibility. Do not exceed the maximum values specified in the standards. If the device is used properly, there is no risk for other equipment being affected by this product. If you have any problems, please contact your service partner.

Warning: This equipment must not be modified without the manufacturer's permission. If the medical device will be modified, suitable tests and inspections must be performed to ensure that it can still be used safely.

4.5 Cleaning

Before cleaning, disconnect the equipment from the power supply by pulling out the mains plug. All you need to clean the equipment is a damp cloth and a commercially available cleaning agent. Do not use abrasive cleaning agents, which can damage the glass surface.

The following disinfectants are allowed:

Methylated spirits; Isopropanol; 2% Kohrsolin; 1% aqueous Sokrena solution; 5% Sagrotan; 5% Gigasept. Do not spray the equipment or the connecting plug.

4.6 Maintenance and service

The Metrological inspection was carried out by the initial approval during manufacturing. More regular metrological controls (Subsequent approvals) are carried out according to national rules of the competent approval authorities. Repairs must only be performed by a service station authorized from Soehnle Professional using original spare parts.

If the device does not work as intended, there is a suspicion of damage. In this case the device must be supplied necessarily to a service station authorized from Soehnle Professional. For a repair by an authorized service station only original spare parts must be used. These original parts are described in the service documentation with different order numbers.

4.7 Warranty / Liability

If a defect attributable to the manufacturer is found in the supplied goods, the manufacturer will be entitled to choose either to rectify the defect or to provide a replacement. Replaced parts become property of the manufacturer. If the defect cannot be rectified or the supply of a replacement is unsuccessful, the statutory regulations will apply.

The warranty covers a period of **2 years** starting from the date of purchase. **Please keep your invoice as proof of purchase.** If your equipment needs servicing, please contact your dealer or the manufacturer's Customer Service.

In particular, liability will not be accepted for any damage resulting from any of the following: Unsuitable, improper storage or use, faulty assembly or start-up by the customer or a third party, natural wear and tear, modifications or tampering, faulty or careless handling, particularly excessive load, chemical/electrochemical/electrical influences or dampness, if these cannot be attributed to the manufacturer. The warranty overall correct function of the equipment will be reserved, if operating, climatic or other influences cause a substantial change in the conditions or the state of the material.

Wear parts (e.g. batteries) are covered by a 6 month warranty.

Please keep the original packaging in case a return shipment is necessary.

4.8.1 Disposal of batteries

Batteries that contain hazardous substances are identified by a **symbol showing a crossed out waste bin and should not be disposed of as normal household waste.** As a consumer, you are required by law to return used batteries. Old batteries can be handed in as special waste at the public recycling points in your borough, or in shops in which batteries of similar kinds are sold. These symbols are found on batteries containing hazardous substances: Pb = battery contains lead, Cd = battery contains cadmium, Hg = battery contains mercury.

4.8.2 Disposal of scale

According to the current state of knowledge, this equipment does not contain any substances that are particularly harmful for the environment. It must not be treated as normal waste, but must be taken to a collection point where electric and electronic devices are recycled. You can obtain further information from your municipal authorities, local waste disposal services, or the company from which this product was purchased.


4.9 Service interface

Applicable for medical licensed devices according to EN 60601-1 for service purposes.

Optional: USB-Interface

5. APPROVAL

5.1 EC-Mark of conformity

The product carries the CE mark according to the following guidelines:

EMC Directive: 2014/30/EC
 Low Voltage Directive: 2014/395/EC
 Scale Directive: 2014/31/EC
 Medical Directive: 93/42/EEC

Explanation of the symbols:

CE 0483 EC mark of conformity

 Safety class of the product

 Product number

 Manufacturer of the product

0483 official number of the „Notified Body“ for the Medical Directive

GB

Device of protection class II according to standard 60601-1 Application part type B

This device has a USB port. This USB port is only intended for service by trained service personnel.

Use of the measuring rod

Do not load the balance on the side

Switch on the scale, tare, reset

Do not sit on the scale.

Pay attention to manual

5.2 Explanation of the approval label and symbols

	↔1	↔2	↔3		SOEHNLE PROFESSIONAL 773x
Max	150kg	250kg			
Min	1kg	2kg			
e =	50g	100g		-10°C/+40°C	
(III)	TCM 128/10-4760	SN	7730.0116xxxx	CE M16	0122

Explanation of the symbols:

↔1	Active weighing range	7730	Product number of the manufacturer
Max	Maximum capacity of the weighing range	0122	e.g. 0122, official number of the „Notified Body“ (this authority has conducted the initial approval)
Min	Minimum capacity of the weighing range		
e =	Approval value (resolution)	001	The approval count indicates how often a scale has been approved. The stored count must be secured with the approval count (see sticker mark) to match the approval connector.
-10°C /+40°C	Working temperature		
CE M16	EC mark of conformity with Year of labelling and in same time the beginning of the approval validity	SN	Serial number of the scale (scale model, final digit of manufacturing year, counting number)
(III)	Accuracy class		

5.3 Explanation of the symbols on the packaging

Attention breakable

Respect the position direction during the transport

Protect from humidity and moisture

Comply with storage temperature

6. WARNINGS

- › In the event of a malfunction, please contact your dealer or the manufacturer. Unauthorised modifications or repairs can damage your scale and will cause the manufacturer's warranty to become void.
- › Display units should be protected from moisture. Liquids (e.g. water) can damage the display unit. Use a dry cloth, e.g. a towel, to wipe the display.
- › Disconnect the power supply before installation, cleaning or maintenance of the equipment. Otherwise, the equipment could be damaged.
- › If you are not using the display for long periods, disconnect the power supply, or switch off the display if it is operating on battery power. Recharge the batteries after 4 weeks at the latest, to compensate for self-discharge.
- › Avoid stacking any materials on top of the display, or placing any weights on it. This can damage the display.
- › Place the scale on a solid, stable and level surface: this will ensure that measurements are accurate. Accuracy is affected, if the scale is tilted or placed on a carpeted floor.
- › Do not connect the scale to an unstable power source.
- › Only use the original equipment. The use of other brands can cause damage to the display.

7.GUIDANCE AND MANUFACTURER'S DECLARATION – ELECTROMAGNETIC EMISSIONS

Guidance and manufacturer's declaration—electromagnetic emissions		
The Model 773X is intended for use in the electromagnetic environment specified below. The customer or the user of the Model 773X should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment—guidance
RF emissions CISPR 11	Group 1	The Model 773X uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The Model 773X is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

The Scales are defeated by special precautions concerning the electromagnetic compatibility and must be installed as mentioned in the installation instructions.

Portable and mobile HF-Communication devices might be interfere at short distance

Guidance and manufacturer's declaration—electromagnetic immunity			
The Model 773X is intended for use in the electromagnetic environment specified below. The customer or the user of the Model 773X should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment—guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines Not applicable	±2 kV for power supply lines ±1 kV for input/output lines Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.

<p>Voltage dips, short interruptions and voltage variations on power supply input lines</p> <p>IEC 61000-4-11</p>	<p><5 % <i>UT</i> (>95 % dip in <i>UT</i>) for 0.5 cycle</p> <p>40 % <i>UT</i> (60 % dip in <i>UT</i>) for 5 cycles</p> <p>70 % <i>UT</i> (30 % dip in <i>UT</i>) for 25 cycles</p> <p><5 % <i>UT</i> (>95 % dip in <i>UT</i>) for 5 s</p>	<p><5 % <i>UT</i> (>95 % dip in <i>UT</i>) for 0.5 cycle</p> <p>40 % <i>UT</i> (60 % dip in <i>UT</i>) for 5 cycles</p> <p>70 % <i>UT</i> (30 % dip in <i>UT</i>) for 25 cycles</p> <p><5 % <i>UT</i> (>95 % dip in <i>UT</i>) for 5 s</p>	<p>a typical commercial or hospital environment. If the user of the Model 004 Image Intensifier requires continued operation during power mains interruptions, it is recommended that the Model 004 Image Intensifier be powered from an uninterruptible power supply or a battery.</p>
<p>Power frequency (50/60 Hz) magnetic field</p> <p>IEC 61000-4-8</p>	<p>3 A/m</p>	<p>0,3 A/m</p>	<p>If image distortion occurs, it may be necessary to position the Model 004 Image Intensifier further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.</p>
<p>NOTE—<i>UT</i> is the a.c. mains voltage prior to application of the test level.</p>			

Guidance and manufacturer's declaration—electromagnetic immunity			
The Model 773X is intended for use in the electromagnetic environment specified below. The customer or the user of the Model 773X should assure that it is used in such an electromagnetic environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment—guidance
Conducted RF IEC 61000-4-6	3 V _{eff} 150 kHz bis 80 MHz	10 V _{eff} 150 kHz bis 80MHz	Portable and mobile RF communications equipment should be used no closer to any part of the Model 005, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 0,4 \sqrt{P}$ $d = 0,4 \sqrt{P}$ für 80 MHz bis 800 MHz $d = 0,7 \sqrt{P}$ für 800 MHz bis 2,5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^c should be less than the compliance level in each frequency range. ^d Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz bis 2,5 GHz	10V/m 26MHz bis 2,7 GHz	
NOTE 1—At 26 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2—These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.			
a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Model 773X is used exceeds the applicable RF compliance level above, the Model 773X should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the scale.			
b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Recommended separation distances between portable and mobile RF communications equipment and the Model 005

The Model 773X is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled.

The customer or the user of the Model 773X can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Model 773X as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz bis 80 MHz $d = 0,4 \sqrt{P}$	80 MHz bis 800 MHz $d = 0,4 \sqrt{P}$	800 MHz bis 2,5 GHz $d = 0,7 \sqrt{P}$
0,01	0,04	0,04	0,07
0,1	0,13	0,13	0,22
1	0,40	0,40	0,70
10	1,30	1,30	2,20
100	4,00	4,00	7,00

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1—An additional factor of 10/3 is used in calculating the recommended separation distance for transmitters in the ISM frequency range 80 MHz to 2.5 GHz to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas.

NOTE 2—These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

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Soehnle Industrial Solutions GmbH