

SLING INSTRUCTION SHEET

DISPOSABLE FULL STANDING SLING


GENERAL GUIDELINES

Under no circumstances should any Waverley Glen Lift System be put in the control of persons who have not received proper training in its operation. If operated improperly, this lift system can cause serious injury.

The safety of both the attendant and the individual depends upon strict adherence to the following guidelines:

1. All instructions pertaining to the operation of the lift system should be read prior to the system being put into use.
2. Do not, under any circumstances, exceed the system's standard working load.
3. Prior to using the sling, make a visual inspection for signs of wear. Look for fraying, cuts, or tears in the material of both the sling's body and straps.
4. The sling must be positioned properly under and around the individual, to ensure a balanced lift.
5. The lift system is not a toy; do not use it for unsafe practices.

In facilities when more than one staff member is responsible for operating the lift system, it is imperative that all such members be trained in its proper use. A training and orientation program should be established by the facility to acquaint new caregivers with the system.



Safe Working Load: 625 lbs

LIFTING FROM CHAIR



1. Take hold of the sling at base of back section ensuring that the quilted side will fit against the person's body. The handle loops at the back will then be visible when the sling is fitted.



2. Introduce the sling by sliding it down between the chair backrest and the person's back so that it touches the chair seat.



3. Draw the complete leg and aperture sections under the person's upper thighs and buttocks. This is made easier if the person can be turned to one side.



4. The leg sections should then be brought up between the thighs. The outer flaps should wrap around the legs to cover the inner flaps. Connect the four buckles and tension the leg straps to give a firm snug fit by pulling on the loose straps ends.


Note: Disposable Full standing Sling is compatible with full range of Lift equipment and carry bars designed for loops options of Prism medical Ltd. Make sure to check Working Load of sling and use compatible equipment for lifting. Should you have any questions or concerns regarding the compatibility of the sling with Lifting equipment, Please contact your Ergosafe representative should you have any further questions.



SLING INSTRUCTION SHEET

DISPOSABLE FULL STANDING SLING

LIFTING FROM A CHAIR CONTINUED




5. Bring the chest sections around the body and overlap the short chest flap section with the long chest flap. Connect the two buckles and tension the chest straps to give a firm snug fit.



6a. Lower the carry bar of the ceiling lift to approximately the seated person's eye level height. Attach all four straps to the carry bar using loops of the same colour. Block the person's knees with a bed or frame. (Instructions also apply to fig. 6b). The sling is not suitable with all ceiling lifts. Please check with your supplier.



6b. When using a ceiling lift, the chair must be positioned at right angles to the track. (The sling is not suitable for use with all ceiling lifts. Please check with your supplier.)



7. Raise the carry bar just enough to tension the straps and ensure that they are still securely attached. The person can now be raised to the standing position. To remove the sling, carry out the fitting procedure in reverse.

SPECIAL INSTRUCTIONS

Prior to proceeding with the actual use of this sling, the following information and instructions should be reviewed carefully to ensure a proper understanding of the use of this sling is obtained.

This sling may not be practical for use by all individuals.

Please consult your authorized representative should you require this sling for individuals who are unable to assist with this process.

The standing sling should be checked on a regular basis to ensure cleanliness and that no fraying has occurred in the body of the sling or straps. Please contact your authorized representative if you have any concerns about the sling.

APPLICATION

This sling is only designed for use as a rehabilitation aid to assist a person to be raised from a seated position to a standing position. It must be used in conjunction with a knee block such as a bed against a wall or with a standing frame. The person must be assessed by a healthcare professional who should be fully trained in the application and fitting of the sling.

WASHING INSTRUCTIONS



SINGLE PATIENT USE ONLY
WARNING
DO NOT WASH

“ As with all patient handling slings, Ergosafe recommends that the slings be inspected on a regular basis for signs of abnormal wear and tear. Please refer to Ergosafe's Sling Inspection Checklist. Disposable or Patient Specific slings are designed such that laundry is not required or allowed. Therefore the sling should be disposed of anytime there is visible soiling of the sling, and replaced with a new disposable/patient specific sling. Should you have any questions or concerns regarding the integrity of a sling you should err on the side of caution and no longer use the sling in question. Please contact your Ergosafe representative should you have any further questions.

Note: Client-specific slings last for one to three months – If not soiled or laundered. Disposable full standing Support Sling is made from non-woven polyester material that conforms to body contour and improves comfort and safety. And because these slings are intended for only one client they help prevent infections from spreading, as well as eliminating laundry costs.



SLING INSTRUCTION SHEET

DISPOSABLE FULL STANDING SLING


Below is the detailed description of the “DO NOT USE” LABEL
SLING SHOWN BELOW IS A SAMPLE REPRESENTATION

THIS PRODUCT IS DESIGNED FOR SINGLE-PERSON USE.

IT IS NOT DESIGNED TO BE LAUNDERED


In the event of a sling being accidentally laundered, a prominent 'DO NOT USE' label will show and the sling should be discarded. Light soiling may be carefully sponged and the sling 'air dried'.

DO NOT TUMBLE DRY



When the 'DO NOT USE IF WASHED' label is visible, the sling is safe to use following a visual inspection as specified at the beginning of this manual.

When the 'DO NOT USE' label is visible, the sling must be withdrawn from use and discarded immediately.



The sling can be used for bathing, but the 'DO NOT USE IF WASHED' label must remain dry. The Ink used for 'IF WASHED' is water soluble and will dissolve should it become immersed. Soiled slings should be discarded.