
User Manual

Blood Pressure Monitor

Model: TMB-1872-B

Arm Type

- Thank you very much for selecting RENPHO Blood Pressure Monitor TMB-1872-B.
- Please read the user manual carefully and thoroughly to ensure the safe usage of this product. Keep the manual well for further reference in case you have problems.

Guangdong Transtek Medical Electronics Co., Ltd.
 Zone A, No.105 ,Dongli Road, Torch Development District,
 Zhongshan, 528437, Guangdong, China


EC REP MDSS-Medical Device Safety Service GmbH
 Schiffgraben 41, 30175 Hannover, Germany

UK REP MDSS-UK RP Ltd.
 6 Wilmslow Road, Rusholme, Manchester M14 5TP UNITED KINGDOM

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COMPLIED STANDARDS LIST

♥ Complied Standards List

| | |
|---------------------------------|---|
| Risk management | EN ISO 14971:2012 / ISO 14971:2007 Medical devices - Application of risk management to medical devices |
| Labeling | EN ISO 15223-1:2016 / ISO 15223-1:2016 Medical devices. Symbols to be used with medical device labels, labelling and information to be supplied. Part 1: General requirements |
| User manual | EN 1041:2008 +A1:2013 Information supplied by the manufacturer of medical devices |
| General Requirements for Safety | EN 60601-1:2006+A1:2013+A12:2014/ IEC 60601- 1:2005+A1:2012 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance EN 60601-1-11:2015/ IEC 60601-1-11:2015 Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment |
| Electromagnetic compatibility | EN 60601-1-2:2015/ IEC 60601-1-2:2014 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic disturbances - Requirements and tests |
| Performance requirements | EN ISO 81060-1:2012 Non-invasive sphygmomanometers - Part 1: Requirements and test methods for non-automated measurement type IEC 80601-2-30:2018 Medical electrical equipment- Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers |

| | |
|-------------------------------|--|
| Clinical investigation | EN 1060-4:2004 Non-invasive sphygmomanometers - Part 4: Test procedures to determine the overall system accuracy of automated non-invasive sphygmomanometers ISO 81060-2:2013 Non-invasive sphygmomanometers - Part 2: Clinical validation of automated measurement type |
| Usability | EN 60601-1-6:2010+A1:2015/IEC 60601-1-6:2010+A1:2013 Medical electrical equipment - Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability IEC 62366-1:2015 Medical devices - Part 1: Application of usability engineering to medical devices |
| Software life-cycle processes | EN 62304:2006 + A1:2015/IEC 62304: 2006+A1:2015 Medical device software - Software life-cycle processes |
| Bio-compatibility | ISO 10993-1:2018 Biological evaluation of medical devices- Part 1: Evaluation and testing within a risk management process ISO 10993-5:2009 Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity ISO 10993-10:2010 Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization |

♥ EMC Guidance

The ME EQUIPMENT or ME SYSTEM is suitable for home healthcare environments

Warning: Don't near active HF surgical equipment and the RF shielded room of an ME system for magnetic resonance imaging, where the intensity of EM disturbances is high.

Warning: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Warning: Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the equipment TMB-1872-T, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Technical description:

1. all necessary instructions for maintaining BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to electromagnetic disturbances for the excepted service life.
2. Guidance and manufacturer's declaration -electromagnetic emissions and Immunity

Table 1

| Guidance and manufacturer's declaration - electromagnetic emissions | |
|---|-------------|
| Emissions test | Compliance |
| RF emissions CISPR 11 | Group 1 |
| RF emissions CISPR 11 | Class [B] |
| Harmonic emissions IEC 61000-3-2 | Class A |
| Voltage fluctuations/ flicker emissions IEC 61000-3-3 | Comply |

Table 2

| Guidance and manufacturer's declaration – electromagnetic Immunity | | |
|--|---|---|
| Immunity Test | IEC 60601-1-2 Test level | Compliance level |
| Electrostatic discharge (ESD) IEC 61000-4-2 | ±8 kV contact ±2 kV, ±4kV, ±8 kV, ±15 kV air | ±8 kV contact ±2 kV, ±4kV, ±8 kV, ±15 kV air |
| Electrical fast transient/burst IEC 61000-4-4 | ±2 kV for power supply lines ±1 kV signal input/output 100 kHz repetition frequency | ±2 kV for power supply lines ±1 kV signal input/output 100 kHz repetition frequency |
| Surge IEC61000-4-5 | ±0.5 kV, ±1 kV differential mode ±0.5 kV, ±1 kV,±2 kV common mode | ±0.5 kV, ±1 kV differential mode ±0.5 kV, ±1 kV,±2 kV common mode |

| Immunity Test | IEC 60601-1-2 Test level | Compliance level |
|--|--|--|
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11 | 0 % UT; 0,5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0 % UT; 1 cycle and 70 % UT; 25/30 cycles; Single phase: at 0°. 0 % UT; 250/300 cycle | 0 % UT; 0,5 cycle. At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0 % UT; 1 cycle and 70 % UT; 25/30 cycles; Single phase: at 0°. 0 % UT; 250/300 cycle |
| Power frequency magnetic field IEC 61000-4-8 | 30 A/m 50Hz/60Hz | 30 A/m 50Hz/60Hz |
| Conducted RF IEC61000-4-6 | 3 V 0,15 MHz – 80 MHz 6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz | 3 V 0,15 MHz – 80 MHz 6 V in ISM and amateur radio bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz |
| Radiated RF IEC61000-4-3 | 10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz | 10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz |
| NOTE U T is the a.c. mains voltage prior to application of the test level. | | |

Table 3

| Guidance and manufacturer's declaration - electromagnetic Immunity | | | | | | | |
|--|----------------------|------------|--|----------------------------------|----------------|--------------|---------------------------|
| | Test Frequency (MHz) | Band (MHz) | Service | Modulation | Modulation (W) | Distance (m) | IMMUNITY TEST LEVEL (V/m) |
| Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment) | 385 | 380-390 | TETRA 400 | Pulse modulation b) 18Hz | 1.8 | 0.3 | 27 |
| | 450 | 430-470 | GMRS 460, FRS 460 | FM c) ± 5kHz deviation 1kHz sine | 2 | 0.3 | 28 |
| | 710 | 704-787 | LTE Band 13, 17 | Pulse modulation b) 217Hz | 0.2 | 0.3 | 9 |
| | 745 | | | | | | |
| | 780 | | | | | | |
| | 810 | 800-960 | GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5 | Pulse modulation b) 18Hz | 2 | 0.3 | 28 |
| | 870 | | | | | | |
| | 930 | | | | | | |

| | | | | | | | |
|--|----------------------|------------|---|---------------------------|----------------|--------------|---------------------------|
| Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment) | Test Frequency (MHz) | Band (MHz) | Service | Modulation | Modulation (W) | Distance (m) | IMMUNITY TEST LEVEL (V/m) |
| | 1720 | 1700-1990 | GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS | Pulse modulation b) 217Hz | 2 | 0.3 | 28 |
| | 1845 | | | | | | |
| | 1970 | | | | | | |
| | 2450 | 2400-2570 | Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7 | Pulse modulation 217 Hz | 2 | 0.3 | 28 |
| | 5240 | 5100-5800 | WLAN 802.11 a/n | Pulse modulation 217 Hz | 0.2 | 0.3 | 9 |
| | 5500 | | | | | | |
| | 5785 | | | | | | |