



Bioenno Tech, LLC

12630 Westminster Ave., Suite H


Santa Ana, California 92706, USA

www.bioenno.com

Email: contact@bioenno.com

Phone: +1 714 234-7363, +1 949 310-9899

Proven results: The sliceGolgi Kit was used to perform Golgi-staining alone as well as combined Golgi-staining and immunohistochemistry/immunocytochemistry (ICC) on free-floating sections.



A-C: The **sliceGolgi Kit** was used to impregnate and stain neurons in 100 – 200 micron thick, **free-floating sections** of postnatal day 12 (P12) C57BL mouse. **(A)** a stained neuron in the lateral hypothalamic area with numerous filopodia-like protrusions (arrowheads); **(B)** a neuron in the frontal cortex with mature spines (arrows) and filopodia (arrowheads); and **(C)** the main (stem) and oblique dendrites of a pyramidal neuron in the motor area of the frontoparietal cortex.

D-E: The combination of **Golgi-staining** and **ICC** allows simultaneous visualization of dendritic spines and immuno-reactive products. Golgi-staining was first performed on 50 – 100 micron thick, free-floating sections, followed by ICC. The impregnation time was 2 – 4 days at $22 \pm 1^\circ\text{C}$. Images were taken from the subiculum of a 4-month old C57 mouse **(D)** and from the frontoparietal cortex (somatosensory area) of a 2-month old C57BL mouse **(E)**. Boxed areas in **E1** were magnified (63x) **(E2)** to highlight the dendritic spines (arrows) and immuno-labeled neurons (arrowheads).