

Treatment Preparation

AquaBoost Back Facial



PRODUCTS REQUIRED FOR SERVICE:

SKU	PRODUCT NAME & PRODUCT AMOUNT
PSN130 OR PSN131	Serum 1: Cleanse & Exfoliate Serum 2: Deep Cleanse
PSN132	Serum 3: Soothe & Hydrate
AQB206	Treatment Tip
NB525	Gloves

CONTRAINDICATIONS
<ul style="list-style-type: none"> • Open lesions on the skin. • Pregnant & nursing women should consult with a physician. • Accutane or other skin-thinning medications or topical gels/creams. <p><i>**Always follow your company guidelines on who should and should not receive a treatment**</i></p>

SERVICE TIPS
<ul style="list-style-type: none"> • Guests may find their sides to be sensitive and prefer not to have their sides included in the treatment. • The AquaBoost Back Facial is suitable for the entire back however the 60-minute protocol includes the Mid Back, Upper Back, & Shoulders only in regards to product quantity. • It is appropriate to perform the back facial by spot treating versus the entire back if necessary.

HEMOCARE RECOMMENDATIONS
<p>It is recommended guests go home with or replenish the following to maintain results</p> <ul style="list-style-type: none"> • Purifying and/or Exfoliating Cleanser • Spot treatment for acne challenges

POST TREATMENT INSTRUCTIONS
<ul style="list-style-type: none"> • Wear Sunscreen Daily • Resume your normal skin care next day • Avoid exfoliation for the next 7 days. • Schedule next visit.

Detailed Protocol

AquaBoost Back Facial



STEPS	PRODUCT NAME	INSTRUCTIONS
Cleanse	Select Cleanser of Your Choice	Thoroughly cleanse the back to prepare the skin for the AquaBoost treatment.
Exfoliate	Serum 1 Cleanse & Exfoliate	<ul style="list-style-type: none"> • Work from the inside to the outside of the back. (ie from the spine to the side) • Hold the skin taut with one hand while gliding the treatment wand with the other. This will help avoid skipping and rapid release of the wand. • Work in sections from inside to outside completing: <ul style="list-style-type: none"> Mid Back Upper Back & Shoulders <p>**The full back can receive the AquaBoost service, however will require additional serum than accounted for in the protocol. It is also feasible to spot-treat specific areas such as upper back only, only acne challenged locations, etc.**</p>
Extractions	Serum 2 Deep Cleanse	<ul style="list-style-type: none"> • If your guest requires extractions or excessive sebum removal, setup Serum 2. • Perform passes from the inside to the outside of the back in the areas needing extractions and deep cleansing. • Use the twist method on stubborn blemishes to help with extraction. <p>**Apply Desincrustation Solution to areas with more congestion or stubborn extraction sites to help prep the skin. Apply solution prior to setting up Serum 2 to allow time to take action.**</p> <p>**Perform manual extractions if necessary at this time.**</p>
Hydrate	Serum 3 Soothe & Hydrate	<ul style="list-style-type: none"> • Work from the inside to the outside of the back. (ie from the spine to the side) • Hold the skin taut with one hand while gliding the treatment wand with the other. This will help avoid skipping and rapid release of the wand. • Work in sections from inside to outside completing: <ol style="list-style-type: none"> 1. Mid Back 2. Upper Back & Shoulders • Press any serum into the skin that is not immediately absorbed to avoid dripping.
CryoWand	Select Gel Mask or Heavy Gel Serum	<p>Perform 1 pass with CryoWand</p> <ul style="list-style-type: none"> • Turn on the CryoWand • Alert the guest they will experience a cold feeling from the CryoWand • Apply Gel Mask / Heavy Gel Serum • Work in sections from inside to outside completing: <ol style="list-style-type: none"> a. Mid Back b. Upper Back & Shoulders • Move in circular motions upward and outward. • Keep the CryoWand moving at all times. • Remove excess product with a warm towel.
Moisturize	Select Moisturizer of your choice	<ul style="list-style-type: none"> • Apply additional moisturizer if needed.