

Treatment Preparation

AquaBoost Fanny Facial



PRODUCTS REQUIRED FOR SERVICE:

SKU	PRODUCT NAME & PRODUCT AMOUNT
PSN130 OR PSN131	Serum 1: Cleanse & Exfoliate Serum 2: Deep Cleanse
PSN132	Serum 3: Soothe & Hydrate
AQB206	Treatment Tip
NB525	Gloves
PSN131 (Optional)	Serum 2: Deep Cleanse

CONTRAINDICATIONS
<ul style="list-style-type: none"> • Open lesions on the skin. • Active Rosacea flare-up. • Pregnant & nursing women should consult with a physician. • Accutane or other skin-thinning medications or topical. <p><i>**Always follow your company guidelines on who should and should not receive a treatment**</i></p>

SERVICE TIPS
<ul style="list-style-type: none"> • Please follow local and state guidelines on draping and treatment locations policies. • It is strongly recommended that you wear gloves throughout the entire treatment.

HEMOCARE RECOMMENDATIONS
<p>It is recommended guests go home with or replenish the following to maintain results</p> <ul style="list-style-type: none"> • Cleanser • Purifying and/or exfoliating product • Moisturizer

POST TREATMENT INSTRUCTIONS
<ul style="list-style-type: none"> • Wear Sunscreen Daily • Resume your normal skin care next day • Avoid exfoliation for the next 7 days.

Treatment Preparation

AquaBoost Fanny Facial



STEPS	PRODUCT NAME	INSTRUCTIONS
CLEANSE	Select Cleanser of Your Choice	<ul style="list-style-type: none"> Thoroughly cleanse the fanny to prepare the skin for the AquaBoost treatment. Steam is optional during cleansing. <p>**Disposable panties are recommended for modest guests</p> <p>**For states that allow fully uncovered fannies, a rolled towel vertically along the center of the glutes is recommended, allowing both glutes to be fully exposed for treatment.**</p>
SERUM 1	Cleanse & Exfoliate	<ul style="list-style-type: none"> Hold the skin taut with one hand while gliding the treatment wand with the other. This will help avoid skipping and rapid release of the wand. Perform one pass: <ul style="list-style-type: none"> Begin on the left. Pull the treatment wand vertically from the bottom to the top of the glute. Repeat one pass on the right glute. <p>For areas of congestion, repeat a pass horizontally in that section.</p>
SERUM 2 OPTIONAL	Deep Cleanse	<ul style="list-style-type: none"> Spot treat acne and areas of congestion using Serum 2 Use small vertical and horizontal strokes Use the twist method on stubborn blemishes to help with extraction. <p><i>**Apply Desincrustation Solution to areas with more congestion or stubborn extraction sites to help prep the skin. Apply solution prior to setting up Serum 2 to allow time to take action.**</i></p>
SERUM 3	Soothe & Hydrate	<ul style="list-style-type: none"> Hold the skin taut with one hand while gliding the treatment wand with the other. This will help avoid skipping and rapid release of the wand. Perform one pass: <ul style="list-style-type: none"> Begin on the left. Pull the treatment wand vertically from the bottom to the top of the glute. Repeat one pass on the right glute. Press any serum into the skin that is not immediately absorbed to avoid dripping.
CRYOWAND	Select Gel Mask or Heavy Gel Serum	<p>Perform 1 pass with CryoWand</p> <ul style="list-style-type: none"> Turn on the CryoWand. Alert the guest they will experience a cold feeling from the CryoWand. Apply Gel Mask / Heavy Gel Serum. Begin on the left glute rotating in circular motions, up the glute. Keep the CryoWand moving at all times. Remove excess product with a warm towel.
MOISTURIZE	Select Moisturizer of your choice	<ul style="list-style-type: none"> Apply additional moisturizer if needed.