



Compass Chair

It's the perfect inclusion chair. The rock-solid frame is built to withstand years of demanding use. The gracefully contoured seat and sturdy armrests provide reassuring boundaries, giving mildly involved students the support they need to feel secure and comfortable. At the same time, its simple and elegant construction lets it blend easily into any classroom environment. It looks almost like a regular classroom chair – it's designed to.

The five available sizes accommodate everyone from small children to adults, and the adjustable legs allow fine-tuning for optimal positioning.



Find important details and informative facts about the Compass Chair at:
www.rifton.com/compasschair

"We have had trouble getting one of our students to sit even for a few moments, which, of course, makes learning very difficult. We've tried different seating options but nothing seemed to work. Imagine our surprise when he saw the Compass Chair for the first time and, entirely on his own, went and sat in it until we had to ask him to get up!"

*Kathleen Marshall
Anderson Center for Autism
Staatsburg, NY*

Compass Chair features



Unique built-in support
The contoured sides of the seat, together with the inside surfaces of the armrests, provide subtle lateral support that can be particularly helpful for students with Autism spectrum disorder, giving security and sensory boundaries.



The Compass Chair comes in five sizes.



Stability feet attach to the back legs to accommodate more active, tip-prone students.



Adjusting the rear legs shorter than the front legs gives the chair a slight tilt, which can help some students relax.

Components



Seat belt

The optional seat belt clips into the seat when needed for safety or postural support.



Stability feet


The stability feet enlarge the chair's base of support. They are adjustable like the standard chair legs.



Adjustable legs

The legs provide 5 cm of height adjustment to ensure optimal seat height and allow for growth. Both chairs in this picture are size 3.

Compass Chair dimensions

User dimensions (cm)	R310 size 1	R320 size 2	R330 size 3	R340 size 4	R350 size 5
Height	91 – 109	109 – 122	122 – 140	140 – 160	160+
Knee to heel	23 – 28	27 – 32	32 – 37	37 – 42	42 – 47
 Key user dimension: height The user's overall height is a general guide to help you select the appropriate chair. Choose the model that allows for growth. Important: Make sure that seat width, depth and height are adequate for each user, and that the user's weight does not exceed the maximum weight recommended.					
Item dimensions (cm)	R310	R320	R330	R340	R350
Seat height	23,25,28	27,29,32	32,34,37	37,39,42	42,44,47
Seat depth	24	28	32.5	37.5	42
Seat width	27	31	36	41.5	46.5
Backrest height	24	28	32	36.5	41
Armrest height above the seat	14	15	18	21	22.5
Armrest length	17	19	22	24.5	28
Max. working load (kg)	45	55	80	100	113

2021 order form

Use dimension chart to select appropriate size.

	R310 size 1	R320 size 2	R330 size 3	R340 size 4	R350 size 5
* Compass chair	R311 <input type="checkbox"/>	R321 <input type="checkbox"/>	R331 <input type="checkbox"/>	R341 <input type="checkbox"/>	R351 <input type="checkbox"/>
Seat belt	R312 <input type="checkbox"/>	R312 <input type="checkbox"/>	R332 <input type="checkbox"/>	R342 <input type="checkbox"/>	R342 <input type="checkbox"/>
Stability feet	R316 <input type="checkbox"/>	R316 <input type="checkbox"/>	R316 <input type="checkbox"/>	R316 <input type="checkbox"/>	R316 <input type="checkbox"/>
Total	R310	R320	R330	R340	R350

* Required