

SWEDE-O® PLANTAR F3™ FOOT ROLLER

- Developed by a Certified Athletic Trainer
- Helps relieve plantar fasciitis and heel/foot pain
- Use hot or cold while stretching plantar fascia
- Contoured design targets painful areas
- Latex free Sorbothane® material is firm, yet flexible
- Easy to clean and reusable
- Use at work, at home, or while traveling





Available From:



DEVELOPED BY A CERTIFIED ATHLETIC TRAINER

The Swede-O® Plantar F3™ Foot Roller was developed by a Certified Athletic Trainer. The rounded, segmented shape provides optimal stretch and targets painful areas.

The most common symptoms of Plantar Fasciitis are stabbing, burning foot and/or heel pain. These symptoms are usually most intense when you first step out of bed or after long periods of rest.

PLANTAR FASCIITIS RELIEF - 3 WAYS

1. STRETCH: Roll back and forth under the foot to gently stretch the plantar fascia (the fibrous tissue along the bottom of the foot that connects the heel bone (calcaneus) to toes). Stretching the plantar fascia prior to placing weight on it helps prevent damage and micro tears to the tissue.

2. COLD THERAPY: Freezer safe roller provides cold therapy to help relieve pain and reduce inflammation. Ideal following activities that irritate the plantar fascia.

3. HEAT THERAPY: Use roller warm to provide heat therapy to help increase elasticity and blood flow to promote healing of the plantar fascia. Ideal prior to walking or standing.

The roller is made of latex free Sorbothane® material which is firm, yet flexible and comfortable against the skin. When warmed or chilled, the roller maintains its density and stays flexible.

SPECIFICATIONS

Product Selection: Swede-O® Plantar F3™ Foot Roller, OMN-3120

Product Weight: .7 lbs

Dimensions: 6" x 2"

Materials: Sorbothane

Patent Number: 6,793,636

Case Quantity: 6

Case Dimensions: 12" x 10" x 8"

Case Weight: 5 lbs

Latex: This product is not made with natural rubber latex

Country of Origin: Made in USA

California Residents Only:

▲WARNING: Cancer and Reproductive Harm — www.P65Warnings.ca.gov