

HPLC column

*SunArmor*

サンアーマー **SunArmor**



**Innovations United**



# The evolution of further surface modification \*




**Development of a novel silyl-reagent which bonded with multifunctional end-capping reagents**

**Final TMS treatment**



**★ C18 phase can be used at pH range from 2 to 12 as well as hybrid C18s.**

**☆ An excellent peak for acidic, basic and metal chelating compounds without effect of residual silanol groups**



Characteristics of SunArmor


Stationary Phase	Particle size (μm)	Pore diameter (nm)	Specific surface area (m <sup>2</sup> /g)	Carbon loading (%)	Ligand	pH range for usage	USP Categor
SunArmor C18	3 and 5	12	340	17	C18	2 - 12	L1



# Stability under basic pH condition

## SunArmor C18

Almost same stability to compare with the hybrid C18s.




Durable test condition  
Column size: 50 x 2.1 mm  
Mobile phase:  
CH<sub>3</sub>OH/10mM Ammonium bicarbonate  
(pH10.5)=30/70  
Flow rate: 0.8 mL/min  
Temperature: 60 °C

Measurement condition  
Column size: 50 x 2.1 mm  
Mobile phase: CH<sub>3</sub>CN/H<sub>2</sub>O=60/40  
Flow rate: 0.2 mL/min  
Temperature: 40 °C  
Sample: 1 = Butylbenzene



# Stability under neutral pH condition at 80 °C

## SunArmor C18



Durable test condition  
Column size: 50 x 2.1 mm  
Mobile phase:  
CH<sub>3</sub>OH/10mM Ammonium acetate (pH6.8)=30/70  
Flow rate: 0.2 mL/min  
Temperature: 80 °C

Measurement condition  
Column size: 50 x 2.1 mm  
Mobile phase: CH<sub>3</sub>CN/H<sub>2</sub>O=60/40  
Flow rate: 0.2 mL/min  
Temperature: 40 °C  
Sample: 1 = Butylbenzene

The result of the comparison data in this catalog is not the representative example of all application.


# Parameter using standard sample



Evaluation of hydrogen bonding, hydrophobicity and steric selectivity

Reproducibility in retention

Separation of standard sample of SunArmor C18



Condition

Column: SunArmor C18 5  $\mu$ m, 150 x 4.6 mm

Mobile phase: CH<sub>3</sub>OH/H<sub>2</sub>O=75/25

Flow rate: 1.0 mL/min

Temperature: 40 °C

Sample:

1 = Uracil




2 = Caffeine



3 = Phenol




4 = Butylbenzene




5 = o-Terphenyl



6 = Amylbenzene




7 = Triphenylene



	Hydrogen bonding (Caffeine/Phenol)	Hydrophobicity (Amylbenzene/Butylbenzene)	Steric selectivity (Triphenylene/o-Terphenyl)
SunArmor C18	0.40	1.54	1.35

Reproducibility in retention of SunArmor C18



Retention time of amitriptyline

Batch	Retention time
101	7.69 min
102	7.97 min
103	8.12 min
104	7.85 min
105	7.93 min
Average (Av)	7.91 min
Standard deviation ( $\sigma$ )	0.14 min

Condition (amitriptyline)

Column dimension: 150 x 4.6 mm


Mobile phase: Acetonitrile/20mM phosphate buffer pH7.0=(60:40)

Flow rate: 1.0 mL/min


Temp.: 40°C

Sample: 1=Uracil


2=Propranolol



3= Nortriptyline




4=Amitriptyline






# Peak Shape of basic compound

## Comparison of amitriptyline (4) using phosphate buffer



<tr: retention time, TF: tailing factor>

## Comparison of amitriptyline (4) using ammonium acetate buffer



<tr: retention time, TF: tailing factor>

### Condition (amitriptyline)

Column dimension: 150 x 4.6 mm

Mobile phase:

- A) Acetonitrile/**20mM phosphate buffer pH7.0**=(60:40)
- B) Acetonitrile/**10mM ammonium acetate pH6.8**=(40:60)

Flow rate: 1.0 mL/min


Temp.: 40°C

Sample: 1=Uracil


2=Propranolol

3= Nortriptyline

4=Amitriptyline



## Comparison of pyridine (2) using methanol/water mobile phase



### Condition (pyridine)

Column dimension: 150 x 4.6 mm


Mobile phase: CH<sub>3</sub>OH/H<sub>2</sub>O=30/70

Flow rate: 1.0 mL/min

Temperature: 40 °C

Detection: UV@250nm


Sample: 1 = Uracil, 2 = Pyridine, 3 = Phenol



The result of the comparison data in this catalog is not the representative example of all application.



# Peak shape of acidic and metal chelating compounds




< $t_r$ : retention time, TF: tailing factor>

### Condition (formic acid)

Column dimension: 150 x 4.6 mm  
 Mobile phase: CH<sub>3</sub>CN/0.1% H<sub>3</sub>PO<sub>4</sub>=2/98  
 Flow rate: 1.0 mL/min  
 Temperature: 40 °C  
 Detection: UV@210nm  
 Sample: 1 = Formic acid  
 2 = Acetic acid  
 3 = Propionic Acid


\*Hybrid C18s showed a very poor peak shape for formic acid. It is doubted that some amines as a by-product remained on the surface of packing materials.




< $t_r$ : retention time, TF: tailing factor>

### Condition (oxine)

Column dimension: 150 x 4.6 mm  
 Mobile phase: CH<sub>3</sub>CN/20mM H<sub>3</sub>PO<sub>4</sub>=10/90  
 Flow rate: 1.0 mL/min  
 Temperature: 40 °C  
 Detection: UV@250nm  
 Sample: 1 = 8-Quinololinol (Oxine)  
 2 = Caffeine





\*Distributor network

**Innovations United**

Innovations United  
 300 East 57th Street Suite 11J  
 New York, NY 10022  
 Tel: 212-204-0075  
 Fax: 484-313-6368  
 Email: [info@innovationsunited.com](mailto:info@innovationsunited.com)  
 Web: [www.innovationsunited.com](http://www.innovationsunited.com)

**Ordering information of SunArmor**

Packings	Inner diameter (mm)	2.0	3.0	4.6	10	20	USP category
	Length (mm)	Catalog No.	Catalog No.	Catalog No.	Catalog No.	Catalog No.	
SunArmor C18, 3 μm	50	HB2241	HB2341	HB2441	-----	-----	L1
	75	HB2251	-----	-----	-----	-----	
	100	HB2261	HB2361	HB2461	-----	-----	
	150	HB2271	HB2371	HB2471	-----	-----	
	250	-----	HB2381	HB2481	-----	-----	
SunArmor C18, 5 μm	50	HB3241	HB3341	HB3441	-----	-----	L1
	100	HB3261	HB3361	HB3461	-----	-----	
	150	HB3271	HB3371	HB3471	-----	-----	
	250	HB3281	HB3381	HB3481	HB3781	HB3881	

Guard Cartridge column of SunArmor

product	Particle size	Catalog No.
SunArmor C18, 5 μm Guard cartridge column (1-pak + Holder) 4 x 10mm	5 μm	HB3A1H
SunArmor C18, 5 μm Guard cartridge (4-pak) 4 x 10mm	5 μm	HB3A1C
SunArmor Guard cartridge holder	---	HOLA1C

**Manufacturer**

**ChromaNik Technologies Inc.**

1512

6-3-1 Namiyoke, Minato-ku, Osaka, 552-0001 Japan

TEL: +81-6-6581-0885 FAX: +81-6-6581-0890

E-mail: [info@chromanik.co.jp](mailto:info@chromanik.co.jp)

URL: <http://chromanik.co.jp>