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FR

Patient instructions for use

EN

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DE

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CS

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Patient instructions for use

1P700998 2024-01

Read before use, in the presence of your prosthetic orthopedic prosthodontist
Corresponding certified prosthodontist instructions: 1P700999 2024-01

ENGLISH

Table of contents

1. DESCRIPTION OF THE SYMBOLS	3	8. OTHER USES OF THE SYNSYS SMARTPHONE APP	12
2. PACKAGE.....	4	9. BATTERY.....	13
3. DESCRIPTION, PROPERTIES AND MECHANISM OF ACTION .	4	A. BATTERY CHARGE	13
A. DESCRIPTION	4	B. BATTERY DRAINED	14
B. PROPERTIES	4	C. SAVING THE BATTERY LIFE ..	14
C. MECHANISM OF ACTION.....	4		
4. INTENDED USE / INDICATIONS	5	10. DETECTION OF MALFUNCTIONS.....	15
5. CLINICAL BENEFITS	6	11. WARNING, CONTRAINDICATIONS, SIDE EFFECTS	15
6. OPERATION / INSTRUCTIONS FOR USE	6	A. WARNING	15
A. SETTING THE HEEL HEIGHT	7	B. CONTRAINDICATIONS.....	17
B. STANCE PHASE RESISTANCE ...	8	C. SIDE EFFECTS.....	17
C. SWING PHASE RESISTANCE	8	D. ELECTRO-MAGNETIC COMPATIBILITY	18
D. SWING PHASE INITIATION	9	<i>Generalities</i>	18
7. USING SPECIAL MODES ..	10	<i>Emissions.....</i>	18
A. WALKING MODE 2	10	<i>Immunity.....</i>	18
B. LOCKING MODE	10		
C. CYCLING MODE.....	11		
D. FREE MODE	11		
		12. CLEANING.....	20
		13. CARE, MAINTENANCE, AND SERVICE LIFE	20
		14. CONDITIONS OF USE AND STORAGE	21
		15. REGULATORY INFORMATION.....	21

1. Description of the symbols

	On / Off		Manufacturer
	Serial number		Consult the instructions
	Danger		Caution
	BF applied part = socket (not part of the device). The buttons and the upper part of the prothesis are regarded as applied parts	IP54	Protection against dust and water splashes in all directions
	Lithium - ion battery		Charger insulation class: class II
	CE marking and year of 1 st declaration		Single patient, Multiple use
	Fragile, handle with care		Up (Indicates the correct vertical position)
	Keep in a dry area		Temperature limits
			Contains electronic components – Separate and recycle

2. Package

Designation	Reference	Included / Sold separately
SYNSYS	1P700LEG	Included: Knee, ankle cover, ankle, foot, foot shell
On/Off tool	1P700132-S	Included
Charger	1P700234-S	Included
Maximum flexion stops	1P700204-[105 to 120]	Included
USB stick	M6X111	Included
Protection cover (option)	1P700218-[color]-S	Not included Sold separately

3. Description, properties and mechanism of action

A. Description

SYNSYS is a microprocessor-controlled knee – ankle – foot assembly.

B. Properties

Reference	1P700LEG
Weight	3.2 kg
Maximum bend	125°
Maximum patient weight (Including carried load)	125 kg / 275 lbs.

C. Mechanism of action

The link between the knee and ankle enables hip/knee/ankle triple-flexion movements, making it safe in addition to walking, for going downstairs, upright/seated transitions, picking up objects from the ground, or doing up your laces, for example.

This knee/ankle link also makes it possible to lift the forefoot during the swing phase (dorsiflexion), to reduce the risk of tripping or walking problems associated with compensation movements. This applies even on upward slopes.

The mobility of the ankle enables the foot to be placed flat on the ground as soon as the heel strikes (plantar flexion).

The onboard sensor in the prosthesis automatically detects the gait phases and manages real-time changes to the prosthesis' resistance.

The knee resistance will be adapted:

- For going downstairs or down a slope
- For sitting or bending down
- In case of a trip (knee flexion locked)
- Other uses besides walking, thanks to Biking, Free and Locking modes (knee locked)

In case of alteration in performance, the knee or ankle may not behave as expected; for example, the knee will remain stiff, in extension.

A certified prosthetist will train you in how to use SYNSYS correctly.

4. Intended use / Indications

This medical device is supplied to healthcare professionals (orthopedic prosthetists), who will train the patient in its use. The prescription is drawn up by a doctor, together with the orthopedic prosthetist, who assess the patient's ability to use it.

This device is for a SINGLE PATIENT. It must not be reused on another patient.

This microprocessor knee-ankle-foot is to be used for lower limb exoprosthetic fittings for transfemoral or knee disarticulated amputees, having a knee-to-ground height between 43 and 55 cm and a foot size between 23 and 30 cm. It was designed for an active person (**IFC d4602** and **d4608**, or **Medicare K3**, and K2 ambulators showing a potential to become a **K3**), as it allows them to carry out their everyday activities on all types of terrain as well as negotiate slopes and stairs.

This medical device intended use is everyday activities, like walking and cycling.

IP54: SYNSYS is splash resistant (fresh water) and can be used in rain without risk of malfunction.

Maximum weight (including carried load): 125 kg / 275 lbs

Not suitable for bilateral femoral amputees or people with hip disarticulation

Not suitable for children

The knee's maximum flexion angle is 125°. It may however be limited by the socket size and shape.

5. Clinical benefits

The device enables the patient to walk, and also go downstairs naturally by allowing free positioning of the prosthetic foot on the stair step, thereby enhancing safety. It provides the same safety when going downhill.

The ankle joint enables the patient to wear shoes with a heel of up to 5 cm, or to walk barefoot.

Typically, thanks to SYNSYS, at an ambient temperature of 24°C, it will be possible to take the stairs down 4 floors (60 alternate steps), walk 1 km on flat ground, go down a 40 m slope, and continue on flat ground for another 1 km.

Once a smartphone has been paired with the prosthesis, the SYNSYS Smartphone app provides access to activity information, and enables the patient to select footwear or change modes (§6 and §7). This app can be downloaded from the store.

6. Operation / instructions for use

Some features of this prosthesis can be viewed / modified via the Bluetooth link on a smartphone and the SYNSYS Smartphone app. You must always be very careful when adjusting via this app. The connection is secure: pairing authorization + confirmation, via the buttons on the prosthesis. In case of a connection problem, check that the prosthesis is not already connected by Bluetooth to another device.

The prosthesis can communicate by Bluetooth up to 10 m with a tablet or smartphone, for the purposes of adjusting the settings, viewing the parameters and activity history; the distance will vary according to any obstacles between the prosthesis and the smartphone.

It will take 2 minutes after disconnecting the Bluetooth for the static standing stabilization to be active.

The smartphone used for Bluetooth communication with the prosthesis must comply with the IEC or ISO security standards applicable to these devices (IEC60950-1 or IEC62368-1).

The App. SYNSY is intended for Android and Apple iOS devices.

SYNSYS App has been designed to connect to iPhone 13, iPhone 12, iPhone SE (2nd generation), iPhone SE (1st generation).

Apple, iPad, iPad Pro, and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries.

Minimum configurations: Android ≥ 6.

To connect:

- The first time: activate Bluetooth on the smartphone and start the SYNSYS app; press and hold the right-hand button for 3 seconds, with the knee bent. Once the prosthesis has been detected, confirm the pairing by briefly pressing the right-hand button on the prosthesis.
- The following times: activate Bluetooth on the smartphone and start the SYNSYS app; activate Bluetooth on the prosthesis by briefly pressing the right-hand button.

For greater security, it is recommended to protect access to the app by the locking options in the phone's parameters.

In case of difficulties using the SYNSYS Smartphone app, contact your certified prosthetist.

A. Setting the heel height

The ankle angle can be adjusted to the height of the shoe heel.


To make this setting, put on the prosthesis and the desired shoes, and then adjust the ankle angle using the buttons.

Setting range: 0 to 5 cm. Precision: 1 long press ≈ 1 to 1.5 mm.

To reduce the heel height:

RIGHT-HAND button


Repeat until you obtain a satisfactory angle.

To increase the heel height:


LEFT-HAND button

Take a few steps or press on the heel so that the ankle position adjusts properly. Repeat until you obtain a satisfactory angle.

Each button press is acknowledged by a single vibration. In case of doubt, you can restore the initial setting with your reference shoes.

On the SYNSYS app, you will be able to:

Associate the shoe photo with the set heel height, rename it, and select the associated walking mode


B. Stance phase resistance

The knee provides stance resistance, e.g. for walking, sitting down, going downhill or downstairs.

- If the knee does not bend fast enough, **decrease** this resistance.

- + If the knee bends too fast, **increase** the resistance.

Setting range: +/-2.5% around the value set by the orthopedic prosthodontist.

 The knee stance resistance is effective from **0° to 95°**flexion. Take care during standing-seated transitions and going down tall stair steps.

C. Swing phase resistance

The swing phase resistance defines the swing flexion damping, especially at slow walking speeds.

- If the foot does not push off backwards properly, **decrease** this value.
- + If the heel rise is too high even at slow walking speed, leaving you "waiting for" your prosthesis, **increase** this value.

Setting range: +/-5% around the value set by the certified prosthodontist.

D. Swing phase initiation

The knee is in safe mode in stance phase; it is freed before toeing off.

- If the knee is not released quickly enough, **reduce** this value
- + If the knee is released too quickly, **increase** this value.

Possible malfunction on soft ground, snow, sand...

Setting range: +/-5% around the value set by the certified prosthetist.

7. Using special modes

Remember to return to Walking mode when this special mode is no longer necessary, except for Biking mode, where the return to Walking mode is automatic.

In Walking mode, the knee joint by default applies a flexion resistance level defined by the certified prosthetist's settings.

A. Walking mode 2


This mode can be set by the orthopedic prosthodontist to adjust the electronic settings to a particular situation (e.g.: wearing heavy shoes, carrying loads...). Select this mode when this sort of situation is encountered.

B. Locking mode


The knee can be locked:

Manual lock: hydraulic locking of the knee joint in full extension:

And, if made accessible to the user by the certified prosthodontist:

- **Active flexion locking ("automatic locking") - between 0° and 45° max.** - useful to maintain a fixed position, with the knee bent


Before walking again, check that the knee is presenting the right flexion resistance.

C. Cycling mode

This mode ensures that the knee joint is completely free (no resistance) while the knee is bent. As soon as the leg is straightened, the prosthesis will automatically return to walking mode, in order to secure the stance for getting off the bicycle.

Knee flexion may have been limited to less than 125° by the certified prosthetist adding a stop, to prevent socket contact: this may make biking more difficult.

D. Free mode

This mode ensures that the knee joint is completely free (no resistance). It is not accessible via the buttons, but can be selected in the SYNSYS app.

It may be useful, for example, in a car (seated in a confined space).


Never use this mode for walking.

To exit it, select another mode in the SYNSYS app, or press & hold the two buttons simultaneously for 1 second.

8. Other uses of the SYNSYS smartphone app

On the SYNSYS app, you will also be able to:


- Look up the battery charge

- Select the alert types: either by vibrations or in "silent" mode, only reporting changes in the prothesis' behavior (see Indicators table)
- "Silent" mode is temporary. It is deactivated automatically for safety reasons.
- Set the buzzer volume.
- Test the operation of the vibrator and buzzer and play with the alerts to memorize them.

- Check your activity data: total number of steps, distance travelled, walking speed
- Check the alerts / notifications.

In case of abnormal use (impact, excessive vibration, or electromagnetic discharge), **the prosthesis will alert you via a series of 10 vibrations, alternating between short and long, accompanied by an audible signal.**

Err on the side of caution, and switch to safety mode. Then check the operation of the device: if the knee can bend, continue to use it without fear. If the knee remains stiff, **contact your orthopedic prosthetist to repair the default**, and use the spare prosthesis.

9. Battery

A. Battery charge

Do not wear the prosthesis while charging the battery and respect the usual electrical safety instructions. Have the charger accessible during charging, to be able to quickly disconnect it if necessary.

Battery charging temperature: +5°C to +40°C; air humidity: < 93%.

Only charge the battery with the 1P700234-S charger supplied by PROTEOR.

An unapproved charger may cause prosthesis failures, or an electrocution risk.

In case of any doubt over this charger (breakage...), do not use it, order a new one from your certified prosthetist.

This charger is not protected from water splashes (IP21).

Remove the cap from the connector, and then connect the charger, cable downward. Plug the power supply block into the mains (conventional wall outlet). A vibration will indicate that the charger is correctly connected.

Disconnecting the charger from the prosthesis or mains will be confirmed by a vibration. Reclose the cap properly, since there is a risk of magnetization of metal particles, which would disrupt the connector's connection.

The charge level will be automatically indicated upon disconnecting the charger, after a few seconds.
It can also be viewed on the SYNSYS Smartphone app, or by **briefly pressing the left-hand button**:


	x 4 =		$\geq 75\%$
	x 3 =		< 75% and > 50%
	x 2 =		$\leq 50\%$ and > 25%
	x 1 =		$\leq 25\%$

The prosthesis has an energy management system, which gives it a battery life of several days. The battery life will vary with use. However, it is recommended to recharge the battery at least once a week.

The charging time is 6 hours, for a full charge. The battery is guaranteed for 500 charging cycles.

B. Battery drained

There are three alert thresholds, indicating:

When the battery charge no longer allows normal operation of the prosthesis, it will be automatically placed in a mode defined with the orthopedic prosthodontist (free or resistant knee).

Stay still for 10 seconds while the mode switches and recharge the battery as soon as possible. Check the knee's flexion resistance before walking again.

Without the electronic control, the prosthesis has:

- Constant knee resistance (stance = swing), the level of which can be set as the patient wishes: "free" for higher mobility (no brake) or "braked" for higher safety (no swing phase)
- End of swing phase damping
- Free ankle movement, especially placing the foot flat on the ground.


It no longer enables:

- Adaptation to the knee resistance, under load or free
- Adaptation to the walking speed
- Adaptation to the shoe heel height
- Switching mode (cycling, lock...)

C. Saving the battery life

In case of a prolonged period out of use (several weeks), switch off the prosthesis:

1. Bend the knee then press & hold the two buttons simultaneously for 10 seconds.
2. OR via the ON/OFF button using the tool 1P700132-S supplied: press & hold for 10 seconds.

Think about switching the prosthesis back on and checking its charge, before reusing. The knee will be ready to be worn with the reference shoes.

If the prosthesis is shut down or the battery heavily drained, the heel height reference may be lost. When reconnected, the prosthesis will be automatically reset

➔ Repeat the procedure for setting a suitable ankle angle for walking (see §6.A).

10. Detection of malfunctions

 If you notice any abnormal behavior or feel any changes in the characteristics of the device, or in case of a severe impact, or exposure to temperatures beyond defined limits, consult your Certified Prosthetist.

11. Warning, contraindications, side effects

A. Warning

 Any modification of this device is prohibited.

 Do not use non-PROTEOR accessories – risk of malfunction and electrocution.

 **Young children:** Risk of inhalation or ingestion of small parts, which could detach due to wear. Risk of strangulation with the charger cables. Keep the prosthesis and its accessories out of the reach of children.

The heel height can be adjusted. An inappropriate setting may cause loss of comfort or balance. The settings also enable mode switches, to modify the knee's behavior in stance phase. Using the wrong mode may cause a fall.

Pay attention to the vibration and audible signals.

In case of doubt, stop walking and check that the prosthesis is working properly.

Look up the alert indicators on the SYNSYS Smartphone app, if applicable.

To reduce the risk of falling when going downstairs, it is advisable to hold on to the handrail.

If the prosthesis emits an alert, stop immediately and check that the knee has the right resistance level.

You should be particularly careful when carrying a heavy load, especially a child: the behavior of the prosthesis may be altered (knee resistance / freed).

Despite all the safety mechanisms provided by this prosthetic solution, there are still risks involved in walking on loose, slippery or rough ground (sand, wet soil...)

Do not try to open or alter the prosthesis.

Only qualified personnel approved by PROTEOR are authorized to carry out maintenance operations.

Do not change the battery.

Do not use a battery other than the original.

Do not plug in anything else to the magnetic connector.

Only certified prosthetists are approved to modify the settings.

Any kind of overloading, excessive strain or excessive vibrations (e.g.: transport) should be avoided. This may cause a mechanical, hydraulic or electronic failure (see Alerts Table on the last page).

Driving:

Have your fitness to drive tested by an approved authority.

You need to ensure that it is completely safe to drive when the prosthesis is not operational.

Air travel:

Do not use a smartphone to connect to the prosthesis. SYNSYS contains a 9.4Wh lithium battery.

According to current regulations, the prosthesis can be carried in the hold when it has been switched **OFF** (see §9.C). NB: the ON/OFF tool is regarded as a sharp tool, i.e. not permitted in airplane cabins

Hydraulic unit overheating:

In case of excessive activity (e.g. going more than 4 floors downstairs), the knee resistance may be altered.

There is a risk of burns from touching overheating components (internal mechanics). It is recommended to wear trousers, to avoid burns.

This risk is increased in case of exposure to direct sunlight (UV).

Pay attention to the vibration and audible signals which indicate a danger of overheating. In this case, stop any activity, to let the prosthesis cool down for a few moments.

The same alert may indicate that the temperature is too low (< -5°C): do around ten knee bends before starting walking again.

Risk of falling. In case of doubt over the alert, look up the SYNSYS smartphone app.

Maximum measured prosthesis (pyramid) connector temperature: 45°C.

Prevent dirt and moisture from entering:

In case of contact with a liquid, wipe down the prosthesis with a clean cloth, and let it dry.

In case of accidental contact with salt water, wipe down the prosthesis with a cloth soaked in fresh water, and let it dry.

In case of any doubt, after unadvised use or a fall, in case of noise, visible signs of wear (impact mark, crack, breakage...) or abnormal behavior of your prosthesis, consult your certified prosthetist to check the functions.

Safety not enabled: in case of a prosthesis failure causing low knee resistance, make sure to initiate the stance phase with the knee extended. Contact your orthopedic prosthetist as soon as possible.

Risk of pinching between moving parts, especially in the knee bending area: take care to prevent any part of the body or clothing from getting caught. Watch out for other people (children, care staff...)

The ankle cover protects the device without disrupting ankle movement and prevents any pinching risk to the fingers.

Handle the prosthesis with care, to prevent any injury risk (cuts...).

Avoid wearing clothing that could impair the movement of the knee and ankle joints. Take care with heavy shoes, which might affect the movements. Avoid wearing high shoes (boots, ankle boots).

Mountaineering: max. altitude = 3000 m / 9800 ft.

Battery auto-drain: if the prosthesis is out of use for prolonged periods (several weeks), the battery will drain as normal.

Switch off the prosthesis to save the battery (see §9.C).

By way of precaution, check the prosthesis charge prior to use (see §9.A).

Risk of irritation due to skin contact with oil: wash your hands in case of contact. If a greasy substance is detected on the prosthesis, contact your orthopedic prosthetist.

In case of failure, the knee's behavior may change, in particular locking in extension.

B. Contraindications

This device is **not** designed for sports use (risks of impacts, bumps, falls...).

The prosthesis can be splashed by water and used in heavy rain without any risk of malfunctioning.

However, it must **not** be immersed. In case of exposure, wipe down the prosthesis with a clean cloth, and let it dry.

C. Side effects

There are no side effects directly associated with the device.

Any serious incident that occurs relating to the device must be reported to the manufacturer and to the competent authority.

D. Electro-magnetic compatibility

Generalities

Environment for use: SYNSYS is designed for use at home and places open to the general public. Do not stand near powerful magnetic or electrical sources (high-voltage lines, transmitters, transformers, MRIs, scanners, metal detectors, welding stations), ATEX zones (=explosive atmospheres, oxygen-rich) or inflammable anesthetic mixtures. The device is incompatible with defibrillation shocks.

⚠ RF portable communication devices (including peripherals such as antenna cables and external antennas) should not be used closer than 30 cm (12 inches) from any part of SYNSYS, including the charger cable. Otherwise, SYNSYS' performance could be affected. Risk of interference causing altered prosthesis behavior (e.g.: knee resistance), which could lead to a fall. The prosthesis may also disrupt the operation of other nearby electronic devices.

⚠ This device should not be used next to other devices or stacked on them, because this might cause a malfunction. If this use is required, check that this device and others are working properly.

⚠ Using accessories, transducers, or cables others than those specified or supplied by the manufacturer of this device may increase the electromagnetic emissions from this device, or reduce its immunity, and lead to improper operation.

List of cables and accessories: charger 1P700234-S (manufacturer: PROTEOR) 155 cm.

Emissions

The medical device is designed for use in an electromagnetic environment as described in the table below. So the user and certified prosthetist must make sure that the medical device is used in the environment described below.

Emissions test	Compliance	Electromagnetic environment – notes
Electromagnetic radiation disturbance (radiated emissions) (CISPR 11)	Group 1	This medical device uses RF energy for its internal / external operation.
Disturbance voltage on the supply terminals (Conducted emissions) (CISPR 11)	Class B	
Harmonic current emissions (IEC61000-3-2)	Class A	
Voltage changes, voltage fluctuations and flicker (IEC61000-3-3)	Compliant	Home health care environment, and environment of a professional healthcare establishment.

Immunity

This medical device is designed for use in a home health care environment, and the environment of a professional healthcare establishment. The user and certified prosthetist must ensure the electromagnetic compliance of the environment.

Immunity test	Test level as per IEC60601-1-2 Ed4		Compliance level	
Electrostatic discharges (DES) (IEC61000-4-2)	± 8 kV in contact $\pm 2/4/8/15$ kV in air		± 8 kV in contact $\pm 2/4/8/15$ kV in air	
Radiated radiofrequency electromagnetic fields. (IEC61000-4-3)	10 V/m 80 MHz - 2.7 GHz 80 % AM at 1 kHz		12 V/m 80 MHz – 2.7 GHz 80 % AM at 1 kHz	
Proximity fields emitted by RF wireless communication devices. (IEC 61000-4-3)	Frequency (MHz)	Modulation	Required level (V/m)	Compliance level (V/m)
	385	Pulsed modulation: 18 Hz	27	27
	450	Pulsed modulation: 18 Hz	28	28
	710 – 745 - 780	Pulsed modulation: 217 Hz	9	9
	810 – 870 - 930	Pulsed modulation: 18 Hz	28	28
	1720 – 1845 - 1970	Pulsed modulation: 217 Hz	28	28
	2450	Pulsed modulation: 217 Hz	28	28
Electrical fast burst transients (IEC61000-4-4)	Power supply: ± 2 kV Input/output lines: ± 1 kV Repetition frequency: 100 kHz		Power supply: ± 2 kV Input/output lines: ± 1 kV Repetition frequency: 100 kHz	
	Between phases: ± 0.5 kV, ± 1 kV Between earth and phases: ± 0.5 kV, ± 1 kV, ± 2 kV		Between phases: ± 0.5 kV, ± 1 kV Between earth and phases: ± 0.5 kV, ± 1 kV, ± 2 kV	
	3 V 0.15 MHz – 80 MHz 6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz		3 V 0.15 MHz – 80 MHz 6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	
Magnetic field at network frequency (IEC61000-4-8)	30 A/m		30 A/m	
Voltage dips and interruptions (IEC61000-4-11)	0 % UT; 0.5 cycles At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle at 0° 70 % UT; 25/30 cycles at 0° 0 % UT; 250/300 cycles		0 % UT; 0.5 cycles At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle at 0° 70 % UT; 25/30 cycles at 0° 0 % UT; 250/300 cycles	

12. Cleaning

Clean the prosthesis using a wet soft cloth and mild soap, following exposure to dust or moisture. In particular, keep the articulated area of the knee clean (magnet). Take care with the small mobile covers (access to socket for charging). Never use any solvents, which could damage the device. Check that there are no metal particles on the magnetic connector and the articulated area.

In case of high perspiration in the socket, and use of vacuum suspension, make sure the prosthesis does not get excessively wet.

13. Care, maintenance, and service life

Electrical network compatible with charger 1P700234-S (manufacturer: PROTEOR):

Input current: 0.3A max. at 100VAC	Output current: 2A max.
Input voltage: AC 100-240 V	Output voltage: 5V DC
Low Power Supply (LPS) compliant	Frequency: 50-60 Hz
Degree of protection: IP21	

Compatibility with standard Bluetooth v2.1 + EDR – class 2 (2.5 mW)
[2400 – 2483.5] MHz

EMC emission class:

Class B – use in residential environment as per standard IEC 60601-1-2. Immunity test level = 12v/m (as per standard EN22523).

Degree of protection as per standard IEC 60529:

IP54: Protected against dust and other microscopic residue.
Protected against water splashes from any directions.

On request, your orthopedic prosthodontist will tell you the prosthesis software reference and version.

The device must undergo a check by your certified prosthodontist at least once a year.

Service life: 3 to 6 years, depending on activity level.

The device needs to be replaced after 6 years' use or after 3 million steps.

If the predicted service life is exceeded, the risk of device failure is increased. This device is prescribed for one user only; any resale is prohibited.

This device is made of components made from various materials: elastomer, plastic, composite (carbon fiber), aluminum, titanium, steel, brass. It also contains oil, a Lithium-ion battery and electronic components, which must not be disposed of with household waste. Avoid any contact with oil. At the end of the device's life, return it to your certified prosthodontist.

14. Conditions of use and storage

Temperature for use: SYNSYS is designed for use in places where the ambient temperature is between -10 and 40°C. In case of exposure to excessively low or high temperatures, your device's operation may be altered.

Atmospheric pressure: 700 – 1060 hPa.

Air humidity (RH): 15% - 90%. No condensation

Storage and transport temperature: RH 15% - 90% with no condensation

- Long term (more than 1 week): between 0°C and 30°C
- Short term (less than one week): between -20°C and +60°C

Following storage at -20°C, wait for 16 minutes for the knee to become operational ($\geq -10^{\circ}\text{C}$).

Following storage at 60°C, wait for 35 minutes for the knee to become operational ($\leq 40^{\circ}\text{C}$).

For transport, switch off the SYNSYS and pack it properly. Avoid leaving the prosthesis in places with excessive temperatures.

Preferably store the prosthesis in the vertical position.

Following prolonged storage, leave the prosthesis vertical for 2 hours, recharge the prosthesis, test the brake and check that the vibrator is operating properly; dust off.

15. Regulatory information


This product is tested and certified compliant with Regulation (EU) 2017/745, and with Class B under standards IEC 60601-1, 60601-1-2, 60601-1-4, 60601-1-11, ETSI EN 300-328, EN 62311 and ISO 10328 - P6 - 125 kg.


The device complies with part 15 of the FCC regulations, and bears the CE marking

RECAP OF BUTTON ACTIONS


Battery charge level





Activate
Bluetooth
(Active time: 2
mins)

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PROTEOR SAS