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SYNSYS

1P700LEG

Certified Prosthetist instructions for use

1P700999 2024-01

Read before use.

Inform end user about all safety precautions before they start using the device.

ENGLISH

Table of contents

1.	Package	4	L.	Limited heel rise	23
2.	Description, properties, and mechanism of action.....	4	M.	Slope adaptation	23
	A. Description	4	8.	Special modes	24
	B. Properties	5	A.	Lock mode	24
	C. Mechanism of action	6	B.	Cycling mode	25
3.	Intended use / Indications.....	7	C.	Free mode	25
4.	Clinical benefits	7	D.	Setting the heel height	26
5.	Accessories and compatibility	8	E.	Walking mode 2.....	27
6.	Preparation	8	F.	Knee resistance in battery critical state	27
	A. Battery charge	8	G.	Battery flat.....	28
	B. Prerequisites.....	9	H.	End of settings checks	28
	C. Preparation.....	9	9.	Use via smartphone Bluetooth.....	28
	D. Prosthesis length adjustment	10	10.	Alerts.....	31
	E. Rod and tibial tube replacement procedure.....	11	11.	Detection of malfunctions.....	31
	F. Foot selection	15	12.	Warning, contraindications, side effects	32
	G. Replacing foot 1A70020.....	15	A.	Warning.....	32
	H. Replacing the ankle cover	16	B.	Contraindications	34
	I. Bench alignment.....	16	C.	Side effects	34
	J. Ankle neutral position (plantar / dorsal flexion).....	17	D.	Electro-magnetic compatibility...	34
7.	Settings and use	17	13.	Maintenance, storage, disposal, and service life.....	36
	A. Static alignment.....	17	A.	Maintenance / cleaning.....	36
	B. Dynamic alignment	17	B.	Environment for use:.....	38
	C. Static standing stabilization	18	C.	Storage	38
	D. Foot internal / external rotation.	18	D.	Disposal	39
	E. Stance phase resistance.....	19	E.	Service life	39
	F. Transition to swing phase threshold	20	F.	Additional information	39
	G. Swing flexion resistance	20	14.	Regulatory information	39
	H. Setting extension return	21			
	I. Terminal impact damping.....	21			
	J. Small Steps Option.....	22			
	K. Plantar flexion resistance	22			

Description of the symbols

	On / Off		Manufacturer
	Serial number		Consult the instructions
	Danger		Caution
	BF applied part = socket (not part of the device). The buttons and the upper part of the prosthesis are regarded as applied parts	IP54	Protection against dust and water splashes in all directions
	Lithium - ion battery		Charger insulation class: class II
	CE marking and year of 1 st declaration		Single patient, Multiple use
	Fragile, handle with care		Up (indicates the correct vertical position)
	Keep in a dry area		Temperature limits
			Contains electronic components – Separate and recycle

1. Package

Designation	Reference	Included / Sold separately
SYNSYS	1P700LEG	Included: Knee, ankle cover, ankle, foot, foot shell
On/Off tool	1P700132-S	Included
Charger	1P700234-S	Included
Maximum flexion stops	1P700204-[105 to 120]	Included
USB stick	M6X111	Included
Protection cover (option)	1P700218-[color]-S	Not included Sold separately

2. Description, properties, and mechanism of action

A. Description

SYNSYS is an all-in-one solution incorporating a knee joint and an ankle joint working in synergy, as well as a prosthetic foot, providing electronic control of the stance phase. Its behavior adapts to the walking speed and situations encountered (going downstairs, going up or down slopes...). Furthermore, the ankle joint enables the patient to wear shoes with a heel of up to 5 cm, or to walk barefoot.


This device will be fitted by a certified prosthetist who has received specific certifying training. The certified prosthetist will in turn train their patient how to use the device correctly, in particular by reading the instructions together.



SYNSYS comprises a proximal part equipped with a male pyramid proximal connection, a tibial tube and rod, available in different lengths, and a foot suited to the user's size and activity level, along with its sock and foot shell.

The SYNSYS SOFT configuration software will enable the certified prosthetist to make the settings on the prosthesis, while guiding them through the sequence of these settings; and view their patient's activity data.

The patient dedicated SYNSYS app will enable the patient to adjust certain settings on their prosthesis and view their activity data.

B. Properties

Reference	1P700LEG
Weight	3.2 kg
Maximum bend	125°
Maximum patient weight (Including carried load)	125 kg
Foot size	23 to 30 cm
A: Knee-ground height	430 to 550 mm
B: Build heights = BH	453 to 573 mm
C: Total heights = TH	467 to 587 mm
	<p>Knee axis – connector base = 21 mm</p> <p>In the seated position: knee axis – front of patella = 64 mm</p>

C. Mechanism of action

This knee-ankle assembly provides control of the stance phase and swing phase, thanks to a microprocessor system.

The knee resistance will be suitable for walking, as well as:

- For going downstairs or downhill
- For sitting or bending down
- Upon a trip
- To its use, with the cycling, free and locking modes.

The ankle makes it possible to:

- Place the foot flat during the heel strike phase (plantar flexion)
- Automatically adapt the ankle when going up a slope, to facilitate step completion (dorsiflexion)
- Adjust the ankle angle to the height of the shoe heel.

The hydraulic link between the knee and ankle enables hip/knee/ankle triple-flexion movements, which help the patient with going downstairs, upright/seated transitions, picking up objects from the ground, or doing up their laces.

This knee/ankle synergy also makes it possible to lift the toe during the swing phase, to reduce the risk of tripping or walking problems associated with compensation movements.

The onboard sensors are able to automatically recognize the walking cycle phases and manage real-time adaptation of the prosthesis' resistance.

SYNSYS provides optimal safety: the knee joint braking mechanism is continuously active except upon automatic recognition of forward swing phase. Combined with the knee/ankle hydraulic link, this enables the patient to go downstairs naturally, with the freedom to place their prosthetic foot on the step, thereby enhancing safety. It provides the same safety when going down a slope.

Safety, in case of a stumble: knee flexion is locked as soon as the knee returns from extension, enabling the patient to recover by catching their weight on the prosthesis.

In case of alteration in performance, the knee or ankle may not behave as expected; for example, the knee will remain stiff, in extension.

3. Intended use / Indications

This prosthetic component is a medical device, supplied to healthcare professionals (certified prosthetist), and patients trained by their CP to use this medical device properly. The prescription is drawn up by a doctor, together with the certified prosthetist, who assess the patient's ability to use it. In case of doubt over the patient's physical or mental capacities, the certified prosthetist must inform the prescribing physician.

 This device is for a SINGLE PATIENT. It must not be reused on another patient.

This microprocessor knee-ankle-foot is to be used for lower limb exoprosthetic fittings for transfemoral or knee disarticulated amputees, having a knee-to-ground height between 43 and 55 cm and a foot size between 23 and 30 cm. It was designed for an active person (IFC d4602 and d4608, or Medicare K3, and K2 ambulators showing a potential to become a K3), as it allows them to carry out their everyday activities on all types of terrain as well as negotiate slopes and stairs.

This medical device intended use is everyday activities, like walking and cycling.

IP54: SYNSYS is splash resistant (fresh water) and can be used in rain without risk of malfunction.

Maximum weight (including carried load): 125 kg

 **Not suitable for bilateral femoral amputees or people with hip disarticulation**

 **Not suitable for children**

 The knee's maximum flexion angle is 125°. It may however be limited by the socket size and shape.

It can withstand water splashes and be used in heavy rain without any risk of malfunctioning.

4. Clinical benefits

The device enables the patient to go downstairs naturally by allowing free positioning of the prosthetic foot on the stair step, thereby enhancing safety. It provides the same safety when going downhill.

The ankle joint enables the patient to wear shoes with a heel of up to 5 cm, or to walk barefoot.

For example, thanks to SYNSYS at an ambient temperature of 24°C / 75°F, it will be possible to take the stairs down 4 floors (60 alternate steps), walk 1 km on flat ground, go down a 40 m slope, and continue on flat ground for another 1 km.

Once a smartphone has been paired with the prosthesis, the SYNSYS Smartphone app provides access to activity information, and enables the patient to select footwear or change modes (§9 and §10). This app can be downloaded from the store.

5. Accessories and compatibility

 **Any modification to this device, other than as described in these instructions, is prohibited.**

 **Do not use non-PROTEOR accessories: risk of malfunction.**

Please refer to the PROTEOR catalogue. (References: §14.A)

Male pyramid proximal link, which can be connected directly to the socket fitted with a female pyramid, or 1D41 or 1D41-HD tube holder. Please refer to the connector's instructions for use.

6. Preparation

A. Battery charge

 **Before using for the first time, charge the knee battery for 6 hours and check the charge level.**

Do not wear the prosthesis while charging the battery.

Comply with the usual electrical safety advice.

Keep the charger accessible during charging, so as to be able to easily disconnect it in case of any concerns. This charger is not protected from water splashes (IP21).

Remove the cap from the connector, and then connect the charger, cable downward. Plug the power supply block into the mains (conventional wall outlet). A vibration will indicate that the charger is correctly connected.

Disconnecting the charger from the prosthesis or mains will be confirmed by a vibration. Reclose the cap properly (risk of magnetization of metal particles, which would subsequently disrupt the connector's connection).

Battery charging temperature: +5°C to +40°C; air humidity: < 93%.

The charge will be automatically indicated upon disconnecting the charger, after a few seconds; or by **briefly pressing the left-hand button**:

	x 4 =		≥ 75%
	x 3 =		< 75% et > 50%
	x 2 =		≤ 50% et > 25%
	x 1 =		≤ 25%

SYNSYS has an energy management system, which provides more freedom to the user and requires charging only every few days. The battery life will vary with use. It is recommended to recharge the battery at least once a week.

The charging time is 6 hours, for a full charge. The battery is guaranteed for 500 full charging cycles.

Only charge the battery with the 1P700234-S charger supplied by PROTEOR. An unapproved charger may cause prosthesis failures, and lead to an electrocution risk. Connecting this device to other devices not described in the instructions for use may be dangerous. In case of any doubts over this charger, order a new one.

If the prosthesis is shut down or the battery heavily drained, the heel height reference may be lost. When the prosthesis is reconnected, it will automatically be reset to its factory position.

➔ Return the ankle to an appropriate angle for walking (see §7.1.).

B. Prerequisites

SYNSYS SOFT is intended for Windows and Apple iOS devices.

SYNSYS Soft has been designed to connect to iPad Pro 12.9-inch (5th generation).

Apple, iPad, iPad Pro, and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries.

Minimum configuration for installing the SYNSYS SOFT is Windows 10. The software can be directly downloaded on Proteor EBOX.

In case of any difficulty installing the software, retry with the firewall and antivirus deactivated. To ensure that SYNSYS SOFT operates correctly, and safeguard patient data, you are requested to adhere to good practice and a high-level security policy on your computer / tablet (anti-virus, firewall...).


You must update the software and firmware as soon as you receive the message. The firmware reference and version are indicated on the SYNSYS SOFT Patient Folder page, after connecting to the prosthesis.

For better security, it is recommended to protect access to the app. and software by the locking options in the smartphone's parameters, or by Windows session automatic locking on a PC.

C. Preparation

The knee can bend to 125°; however, this angle may be limited by the socket volume. Pay particular attention:

- To the socket design: **avoid any contact between socket – covers / hydraulic unit / frame.**
- If using a rotator, make sure that the button does not come into contact with the frame.
- Use appropriate flexion stops (available every 5°, from 100° (1P700204-100) to 120° (1P700204-120) e.g., to ensure kneeling is completely safe.

Fit the stop 1P700204-~ on the knee back part. Screw together the assembly using a screw WMCX3008Z (Allen key 2.5mm), to **1 Nm**. Repeat the operation for the 2nd stop.


When first using the prosthesis, switch it on by plugging it in. The prosthesis has no wake-up delay, it is usable immediately after storage.

D. Prosthesis length adjustment

The 1P700LEG prosthesis is delivered assembled. If the length is not suitable, you will need to exchange the carbon tube and aluminum rod to obtain the desired **knee-to-ground height (A)**, every 5mm between **43 and 55 cm**.


Selecting the right kit:

**Knee-ground height = center of knee axis – under the heel with the foot shell*

E. Rod and tibial tube replacement procedure

Before disassembly, check that the *Heel Height* is set to a “factory” value of between 10 and 20, equivalent to a 10 mm shoe heel. Perform a plantar flexion movement.

 Use the appropriate torque wrench for the points below.

1) Removing the outer structure:


- ① Unscrew the tibia cover screw (**from KTG height 460**), and then remove the tibia cover.
- ② Gently remove the foot shell + ankle cover.
- ③ Unscrew the ankle rear axis (Allen key 4 mm).
- ④ Remove the screw guard, and then unscrew the tibial tube retaining screws (Allen key 3 mm).

CAUTION do not lose the set washer and O-rings.

2) Replace the tibial tube

- ① Unscrew the tibial tube screws (Allen key 3 and 4 mm).
- ② Push out the axis with the pin punch. Pull out the ankle parts.

Insert the new carbon tube, insert the axis with the center pin guide and tighten:

- The 4 tibial tube screws: **6Nm** for the 2 lower screws, and **2Nm** for the 2 upper screws
- For **KTG heights 430 to 440**: Just tighten both screws to **6Nm**.

3) Change the rod

CAUTION Rod tightening tool to be secured in a vice during fitting.

1) Fasten the rod tightening tool on a vice.

2) Mount SYNSYS on the tool, and if necessary, adjust the sensor positioning so that the frame reaches the end stop.

3) Tighten both tibial tube screws to **6 Nm. (3x each screw alternately)**

Unscrew the rod locknut and then the rod, using a **13 mm open-ended spanner**.

DO NOT USE ANY HEAT SOURCE, RISK OF DAMAGING THE DEVICE.

Screw the new rod into the knee.

Tighten the rod to **20 Nm** using the torque wrench, and the 13 mm open-ended bit.

Then screw in the lower ball pivot, without tightening the nut. Position the tool slider at the required **[knee-to-ground]** mark.

Adjust the nut and tighten until there is slight contact with the gauge pin.

Lock the ball pivot by folding down the tool.

Screw on the locknut to **16 Nm**.

Open the locking system, check the ball pivot is well aligned, unscrew the two frame screws, and remove the tool assembly.

4) Assemble the tibial tube with the knee:

① Fit the tibia/ankle/foot blade assembly onto the bottom stop of the frame, and alternately screw in **(3 times)** the 2 clamping screws to **6N.m**.

② Refit the ankle rear axis with the green set seals (one on each side) + the washer (threaded side) and washer chamfer facing the ball pivot, and then tighten the rear ankle axis to **7N.m**.

Once fitting is complete, make sure that, at maximum plantar flexion, the carbon tibial tube comes into contact with the aluminum ankle part.

Refit the foot shell, ankle cover and then the tibia cover.

For KTG heights ≥ 460 mm, screw in the tibia cover to contact. For KTG heights < 460 mm, they are fitted by clipping.

F. Foot selection

Table specifying the foot categories according to patient weight and activity:

Patient weight \ Activity	45-59 kg 99-130 lbs.	60-74 kg 131-163 lbs.	75-89 kg 164-196 lbs.	90-104 kg 197-229 lbs.	105-125 kg 230-275 lbs.
Normal	1	2	3	4	5
High	2	3	4	5	N/A

Category	ISO 10328 tested load level
1	ISO 10328 - P3 - 60 kg ^{*)}
2	ISO 10328 - P4 - 80 kg ^{*)}
3	ISO 10328 - P5 - 100 kg ^{*)}
4	ISO 10328 - P5 - 110 kg ^{*)}
5	ISO 10328 - P6 - 125 kg ^{*)}

^{*)} Mandatory maximum body weight limit

Table of module availability by foot size:

	23	24	25	26	27	28	29	30
1								
2								
3								
4								
5								

G. Replacing foot 1A70020

- ① Remove the carbon tibia cover, foot shell and ankle protection
- ② Unscrew a small side screw on the ankle by one turn (Allen key 2 mm).
- ③ Unscrew the lower screw on the foot blade (Allen key 8 mm). Remove the foot, retrieve the screw and washer, then replace the foot with the desired model. Screw the side screw back in to contact.

Use suitable torque wrenches.

- Retighten the lower screw (Allen key 8 mm) to a torque of **35 Nm**.
- Retighten the side screws (2 Allen key) to a torque of **2 Nm**.
- Put the sock and foot shell back on.
- Refit the tibia cover.


H. Replacing the ankle cover

Remove the tibia cover.

Gently detach the ankle cover from the tibial ring (**pin by pin**).

Remove the assembled Foot Shell and Ankle cover, together.

Follow the procedure in reverse to refit the ankle cover.


The SYN SYS SOFT configuration software guides you through the settings, step by step, so you can be sure of covering all the parameters.

⚠ Check the device to ensure it is operating correctly after any length modifications or foot replacement.

I. Bench alignment

⚠ Pay attention to socket contact at maximum flexion: the contact should only occur with the internal flexion stop. Do not remove the protective covers, this could damage the prosthesis and will void the warranty. Use the provided flexion stops (1P700204-100, -105, -110, -115, -120) if needed.

⚠ The socket must not come into contact with the elastomeric front upper cover (battery cover underneath). The warning label may be removed for delivering the prosthesis to the patient.

The ankle axis is close to the natural joint center, which differs from other prosthetic feet.

The prosthesis is factory set for shoes with a 10 mm heel. Depending on the shoes, adjust the ankle neutral position (§7.I.) so that the tibial tube is vertical (0/+1°).

On a vertical transfer fixture, with shoe, align the prosthesis so that the vertical construction line passes through the greater trochanter, the knee axis, and the ankle axis.

Align via the upper pyramid (tightening to the torque specified on the connector or in its instructions), and the ankle setting buttons or the SYN SYS SOFT configuration software (this will affect the heel height).

Note: knee stability is affected by the position of the knee AND of the ankle relative to the load line. If the ankle is too far back, it will make the knee unstable. If the ankle is too far forward, it may make bending the knee difficult when going down a slope or stairs.

J. Ankle neutral position (plantar / dorsal flexion)


To set ankle flexion: proceed via the SYNSYS SOFT configuration software, or **with the knee in extension, press and hold**.

- + The left-hand button (plantar flexion) = heel higher
- The right-hand button (dorsal flexion) = heel lower

Setting range: 0 → 5 cm.

Precision: 1 vibration ≈ 1 to 1.5 mm (according to foot size)

This setting may also be made via the SYNSYS SOFT configuration software. There is an option for disabling this setting via the buttons: risk of de-alignment of the prosthesis if misused by the patient.


7. Settings and use

A. Static alignment

The load line, passing vertically through the greater trochanter, should pass 5 to 10 mm anterior to the knee center, and between the rear third and halfway point of the foot.

For a more dynamic alignment (easier knee flexion), move the knee forward so that the load line passes over the knee center.

For more stability, move the knee backwards so that the load line passes up to 25 mm anterior to the knee axis.

B. Dynamic alignment


On its own, the **foot** inversion/eversion is not settable.

Have the patient walk between parallel bars. If the knee gives out too quickly at the end of the stance, increase the stability = move the prosthesis backward: the load line will pass at most 25

mm in front of the knee. For more dynamic walking (easier knee flexion), move the knee forward until the load line passes over the knee axis; in this configuration, it is possible to walk with stance flex.

The prosthesis can communicate using Bluetooth up to 10 m away, with a tablet or computer, to adjust the settings, view the parameters and activity; this distance may decrease according to any obstacles between the prosthesis and the tablet. This is a secured connection: pairing authorization + confirmation, via the prosthesis buttons. In case of a connection problem, check that the prosthesis is not already BT connected to another device.

C. Static standing stabilization

When activated, this option provides support in a prolonged standing position. The sensors detect that the user is standing without moving, the prosthesis helps to maintain this position without effort. The knee will be released by the movement.

It will take 2 minutes after disconnecting the Bluetooth for the static standing stabilization to be active.

Notify this info to patients using the App. SYN SYS and having this function active.

D. Foot internal / external rotation

To rotate the foot:


- 1- Remove the shell and sock,
- 2- Unscrew the lower screw ② (8 Allen key),
- 3- Adjust the internal/external rotation via the side screws ③ (Allen key 2 mm).

Firmly tighten the lower screw ② to a torque of 35 Nm, before refitting the sock and foot shell.

Setting range: +/- 9°


E. Stance phase resistance

Between parallel bars, have the user try this triple flexion movement, with knee resistance. Determine whether the resistance level is suitable, for sitting down on a chair for instance. If the knee does not bend quickly enough, **reduce** this resistance. If the knee bends too quickly, **increase** the resistance.


Setting range: 40-80%, recommended range: 55-70%, Precision: 1%

This same setting must be readjusted for the downstairs and slope phases, depending on the patient's activity level and the environment.

The knee stance resistance is **effective from 0° to 95°** flexion. Take care when going down high stair steps.

F. Transition to swing phase threshold


The knee is in safe mode in stance phase; it is released by pressing down with the forefoot. Have the patient start by standing between the bars.

- If the knee is not released when walking (no swing phase), **reduce** this value.
- + If the knee is released too quickly, **increase** this value.

Possible malfunctions on soft ground, snow, sand... Train the patient to walk on these different surfaces.

Setting range: 30-100N.m, recommended range: 50-80Nm Precision: 1


G. Swing flexion resistance

This defines the swing flexion damping, especially at slow walking speeds.

- If you observe practically no swing flexion, reduce this value.

- + If the heel rise is excessive, even at slow walking speed, increase this value.


Setting range: 0-40%, recommended range: 10-25%, Precision: 1%

H. Setting extension return

- 1- Open the cover by unscrewing the outer screw.
- 2- Use the pump provided.
 - If the swing extension is too fast, even at slow walking speeds, **decrease the air pressure (min.: 8 bar)**.
 - + If the knee extension is too slow, **increase the air pressure (max. 10 bar)**.

Setting range: [8-10] bar, varying in 1-bar steps.

- 3- Firmly reclose the cover, tightening to 0.5 Nm


I. Terminal impact damping

- Open the cover by unscrewing the visible screw.
- Use the Torx T10 key.
 - If the knee does not return to full extension, **turn it anti-clockwise**.
 - + If the user feels an impact at the end of swing extension, **screw in**.

Setting over 2 turns, Precision 1/4-turn

- Firmly reclose the cover - 0.5 Nm

SYNSYS-1863639
v. 012a 21 N 20°C

Softortho Demo

Alignment Walking mode 1 Special modes Walking mode 2 Alerts Checking

Swing phase initiation 70 Nm

Swing resistance 15 %

Terminal impact damping

Option small steps Disabled Enabled

Threshold speed 30 °/s

Heel rise damping 50 %

Stance resistance 70 %

Heel height 18

Plantar flexion resistance

Dampening of Terminal Impact


Use a T10 torx wrench.
Open the rear top cover.

Dampening of Terminal Impact

If the user feels an impact at the end of the swing extension, turn clockwise.
If the knee doesn't return fast enough to full extension, turn counter clockwise.
Sensitivity: operate by 1/8th of a turn.

J. Small Steps Option

If the option is deactivated, the knee will stay extended when trampling. Activating this option will make the knee more mobile and increase safety: at slow or comfortable walking speed, the switch to swing phase will happen when the toes are unloaded.

Speed threshold: This speed transition (normal vs fast) is adjustable.

+ If, at comfortable walking speed, the knee releases too early, increase this speed threshold.

- If, at comfortable walking speed, the knee does not release fast enough ("the foot is stuck to the ground"), decrease this threshold.


K. Plantar flexion resistance

This mechanical setting (using the tool 1P700139-S) is used to adjust the speed of when foot flat occurs, at the beginning of stance phase (at heel strike).

! Do not scratch the piston rod with the tool.

- Turn **left** to **reduce** the resistance (for a quicker foot midstance phase); (min.: -3)
0 Default setting.
- Turn **right** to **increase** the resistance (more damped foot midstance phase) (max.: +3)

Setting range: 1 turn [-3;+3]; precision: 1/14th of a turn


L. Limited heel rise

The heel rise can be limited proportionally to the walking speed.

- If the patient is having trouble walking at different speeds ("waiting for" their prosthesis), **decrease** this value.
- If the heel rise is too high in swing flexion, **increase** this value.

Setting range: 40-80; Recommended range: 40-60%, Precision: 1

M. Slope adaptation

This option enables the ankle to adapt to the slope angle: this facilitates rollover when going uphill.


The initial setting is made with the reference shoes, on flat ground, in a straight line, over five steps at a comfortable speed.

- To adapt to shallow slopes, decrease the value.
- + To adapt to steep slopes only, increase the value.

8. Special modes

This device enables the user to walk safely but can also adapt to other activities thanks to the Special Modes: Biking, Free or Lock modes. The access to some modes can be limited by the certified prosthetist depending on the user's needs, via the SYNSYS SOFT configuration software.

Remember to return to normal mode.

The patient can view / modify some prosthesis features, via the Bluetooth connection on their smartphone, in the patient dedicated SYNSYS app.


A. Lock mode

In Walking mode, the knee joint has flexion resistance at a level defined during tests with the certified prosthetist.


The knee can be flexion locked:

Manual lock: hydraulic locking of the knee joint in full extension:

And, if the certified prosthetist enables it for the user:


Active flexion locking ("automatic locking") - between 0 and 45° - useful to maintain a fixed position, with the knee bent.

Setting ranges:

Static force: 10-50N.m, recommended range: 20-40N.m; Precision: 1

Static duration: 0.5-2 seconds, recommended range: 0.8-1.5 seconds; Precision: 0.1 seconds


Before walking again, please check the knee resistance.

B. Cycling mode

This mode enables the knee joint totally free flexion (no resistance) unless the knee is extended.

C. Free mode

This mode ensures that the knee joint is completely free (no resistance). For example, it can be used when sitting in a car. **It is not accessible via the buttons.**

- Unlocking: select this mode in the patient app.
- Confirmation: short vibration.
- Exit: select another mode in the SYN SYS patient-dedicated app or press & hold the two buttons simultaneously for 1 second.

D. Setting the heel height

The ankle angle can be adjusted to the shoe heel height, or for walking barefoot; 10 different shoes can be saved. The settings are accessible via the buttons on the prosthesis, the SYNYSYS SOFT configuration software and the patient dedicated SYNYSYS app.


Risk of slipping when walking barefoot

To make this setting, ask the patient to put on the desired shoes, and then adjust the ankle angle using the buttons or the SYNYSYS SOFT configuration software.

Setting range: [0-5] cm, Precision: 1 to 1.5 mm

Each button press is acknowledged by a single vibration. Only modifications made via the SYNYSYS SOFT configuration software or the patient dedicated SYNYSYS app will be saved.

To reduce the heel height:

- Knee extended.
 - Press and hold the **right-hand button** on the prosthesis.
- Repeat until you obtain a satisfactory angle.

To increase the heel height:


- Knee extended.
 - Press and hold the **left-hand button** on the prosthesis.
 - Take a few steps or press on the heel so that the ankle position adjusts properly.
- Repeat until you obtain a satisfactory angle.

E. Walking mode 2

SYNSYS offers the opportunity to adjust all the walking parameters settable via SYNSYS Soft to a particular situation encountered by the patient: carrying loads or wearing heavy shoes, in a professional context for instance.

By default, the parameters are the factory settings. Choose the shoes associated with Walking mode 2 and adjust the parameters. Check this step carefully: slow walking, quick walking, stairs, and slope. The mechanical settings cannot be specific to this mode.

F. Knee resistance in battery critical state


When battery is out of charge, the knee resistance is no longer managed electronically; one constant resistance level will be used for both swing and stance phases (standing, going downstairs, downhill...).

To test the behavior of the prosthesis without battery, activate the "Low battery" mode.


Depending on the user, select the desired resistance level within the range 0-100%: 0% = prosthesis completely free, and 100% = prosthesis fully locked.

Setting range: 0-100%;
precision: 1%

G. Battery flat

There are three alert thresholds, indicating:

When the battery charge no longer allows normal operation of the prosthesis, it will be automatically placed in “battery critical” mode.

Tell the patient to stay still for 10 seconds while the mode switches, and to recharge the battery as soon as possible. Check the knee’s flexion resistance before walking again.

Without the electronic control, the prosthesis has:

- + Constant knee resistance (stance = swing), the level of which can be set as the patient wishes: “free” for higher mobility (no brake) or “braked” for higher safety (no swing phase)
- + End of swing phase damping
- + Free ankle movement, especially placing the foot flat on the ground.

It no longer enables:

- o Adaptation to the knee resistance, under load or free
- o Adaptation of swing phase to walking speed
- o Adaptation to the shoe heel height
- o Switching mode (cycling, lock...)
- o Locking upon a trip

H. End of settings checks

To check that this prosthesis has been correctly installed, and limit any risk, check:

- + That the stance resistance level is correct
- + That it re-extends correctly, with the right terminal impact damping.
- + That the transition to swing phase threshold is suitable
- + If the Static Standing Stabilization is activated, two minutes after BT disconnection, check the expected function.


9. Use via smartphone Bluetooth.

Inform the patient: two profiles are available: **normal profile** (buzzer and/or vibrator) or **silent** (only alerts entailing a change in the prosthesis’ behavior are retained). They can switch profile and set the buzzer volume via the SYNSYS patient-dedicated app.

Note: Silent mode lasts 2 hours, SYNSYS will switch back to normal mode after this time.

Check the battery level.

Change the language or units:

Select the walking mode (1 or 2)


Engage the prosthesis lock, in extension.

Engage cycling mode,


with the knee bent.

Switch to free mode,

and then exit.

Select the shoe to wear or set the ankle angle to the heel height of a new shoe.


Associate the picture of your shoe to the adjusted heel height


Select the user profile, play with the alerts to memorize them, adjust the volume, and test the operation of the vibrator and buzzer.

Look up the prosthesis' internal counters: total number of steps, distance travelled, number of steps sorted by walking speed...

Adjust the stance resistance, swing resistance and transition to swing phase threshold, by a few percent around the value set by the orthopaedic prosthetist.

 Any change to the settings other than those above are prohibited for the patient.

The smartphone, touchscreen tablet or PC used for Bluetooth communication with the prosthesis must comply with the IEC or ISO security standards which are applicable to these devices (IEC60950-1 or IEC62368-1).

The App. SYNSY is intended for Android and Apple iOS devices.

SYNSYS App has been designed to connect to iPhone 13, iPhone 12, iPhone SE (2nd generation), iPhone SE (1st generation).

Apple, iPad, iPad Pro, and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries.

Minimum configurations: Android \geq 6.

10. Alerts

The prosthesis alerts the user in case of duplication errors on the SYNSYS SOFT configuration software.

If the SYNSYS SOFT configuration software accidentally closes (application crash, PC power outage...) or the prosthesis is accidentally disconnected, make sure the patient is safe, and restart the settings. Do not let them leave without verifying all the settings.

In case of abnormal use (impact, excessive vibration, or electromagnetic discharge), **the prosthesis will alert you via a series of 10 vibrations, alternating between short and long, accompanied by an audible signal.**

Err on the side of caution, and switch to safety mode. Then check the operation of the device: if the knee can bend, continue to use it without fear. If the knee remains stiff, **contact your orthopedic prosthetist to repair the default**, and use the spare prosthesis.

Always disconnect the SYNSYS SOFT configuration software from the prosthesis once the settings are completed.

To remove the warning indicator, press & hold (>3 secs) the left-hand button, with the knee bent.

11. Detection of malfunctions

If the patient notices any abnormal behavior or feel any changes in the characteristics of the device, in case of a severe impact or exposure to temperatures above the permitted thresholds, the patient must consult their certified prosthetist.

12. Warning, contraindications, side effects

Only qualified personnel approved by PROTEOR are authorized to carry out maintenance operations. **Do not change the battery.** Do not use a battery other than the original.

A. Warning

Risk of falls on stairs: always use the banister or handrail. In case of an alert, stop the user immediately, and have them check that they have the right braking level.

You should be particularly careful when carrying a heavy load, especially a child: the behavior of the prosthesis may be altered (knee resistance / freed). Reminder: maximum load capacity of the prosthesis: 125 kg.

Despite all the safety mechanisms provided by the prosthesis solution, there are still risks involved in walking on loose, slippery, or rough ground (sand, wet soil...).

+ Any excessive handling, mechanical overloads, excessive vibrations (e.g.: some forms of transport) must be avoided. This may cause a mechanical, hydraulic, or electronic failure.

Driving: remind the user to undergo a fitness to drive test by an approved authority. You need to ensure that it is completely safe to drive when the prosthesis is not operational.

Air travel: Do not use a smartphone to connect to the prosthesis. SYNSYS contains a 9.4Wh lithium battery. So according to current regulations, the prosthesis can be carried in the hold when it has been switched **OFF**.

Mountaineering: max. altitude = 3000 m / 9800 ft.

Hydraulic unit overheating: in case of sustained activity (e.g., going more than 4 floors downstairs), the braking performance may be altered. There is also a risk of burns by touching overheating components. It is recommended to wear trousers, to avoid burns.

If a socket is used in direct contact with the skin, we advise using a silicone pad.

This risk is increased in case of exposure to direct sunlight (UV). Pay attention to the vibration and audible signals which indicate a danger of overheating. In this case, stop any activity, to let the prosthesis cool down for a few moments: **risk of falling**. If the temperature is too low (< - 5°C), the same alert is triggered: the patient must do around ten knee bends before starting walking again.

Maximum measured prosthesis (pyramid) connector temperature: 45°C.

Low resistance: make sure to initiate the stance phase with the knee extended. If this behavior is abnormal, and not due to the battery being drained, contact After-Sales.

Safety measures

Always check that the user is safe before making any setting: non-slippery ground, handrail / support points / chair.

Always warn the user orally before altering a setting.

The prosthesis enables heel height adjustment (see §9.D). An inappropriate setting may cause loss of comfort or balance. Take care with high shoes (e.g.: boots), the rigidity may affect the freedom of movement of the ankle.

The settings also enable mode switches, to modify the knee's response in stance phase. Using the wrong mode may cause a fall.

Test that the vibrator and buzzer are working properly on the Alerts page of the SYNSYS SOFT configuration software.

Pay attention to the vibration and audible signals emitted by the prosthesis.

In case of doubt, stop walking and check that the prosthesis is working properly. Look up the alert indicators in the SYNSYS SOFT configuration software, if applicable.

Prevent dirt and moisture from entering. In case of contact with a liquid, wipe down the prosthesis with a clean cloth, and let it dry. In case of accidental contact with salt water, wipe down the prosthesis with a cloth soaked in fresh water, and let it dry.

Risk of pinching in the knee bending area: take care to prevent any part of the body or clothing from getting caught in this area. Remind the patient to watch out for other people (children, care staff...). Handle the prosthesis with care, to avoid injuries (cuts...)

The ankle cover 1P700256-S / 1P700257-S must be worn to protect the device without disrupting ankle movement, and to prevent any pinching risk to the fingers.

Check the settings if a trim is worn, especially for full re-extension of the knee. Similarly, avoid wearing clothing that could impair the movement of the knee and ankle joints.

Battery auto-drain: if the prosthesis is out of use for prolonged periods (several weeks), the battery will drain as normal. Switch off the prosthesis to save the battery. By way of precaution, check the prosthesis charge prior to use (§7.A).

There is a risk of irritation due to skin contact with oil: wash your hands in case of contact. If the patient observes a greasy substance on the prosthesis, they must contact their certified prosthetist, who will contact the After-Sales service.

B. Contraindications

This device is **not** designed for sports use (risks of impacts, bumps, falls...). The prosthesis can be splashed by water and used in rain without any risk of malfunctioning.

However, it must **not** be immersed. In case of exposure, wipe down the prosthesis with a clean cloth, and let it dry.

C. Side effects

There are no side effects directly associated with the device. Any serious incident that occurs relating to the device must be reported to the manufacturer and to the competent authority.

D. Electro-magnetic compatibility

General

Environment for use: SYNSYS is designed for use at home and places open to the public. Do not stand near powerful magnetic or electrical sources (high-voltage lines, transmitters, transformers, MRIs, scanners, metal detectors, welding stations), ATEX zones (=explosive atmospheres, oxygen-rich) or inflammable an aesthetic mixture. **The device is incompatible with defibrillation shocks.**

RF portable communication devices (including peripherals such as antenna cables and external antennas) should not be used closer than 30 cm (12 inches) from any part of SYNSYS, including the charger cable. Otherwise, SYNSYS' performance could be affected. Risk of interference causing altered prosthesis behavior (e.g.: knee resistance), which could lead to a fall. The prosthesis may also disrupt the operation of other nearby electronic devices.

This device should not be used next to other devices or stacked on them, because this might cause a malfunction. If this use is required, check this device and others are working properly.

Using accessories, transducers, or cables others than those specified or supplied by the manufacturer of this device may increase the electromagnetic emissions from this device, or reduce its immunity, and lead to improper operation.

Charger reference: 1P700234-S (manufacturer: PROTEOR) 155 cm.

Emissions

The medical device is designed for use in an electromagnetic environment as described in the table below. So the user and certified prosthetist must make sure that the medical device is used in the environment described below.

Emissions test	Compliance	Electromagnetic environment – notes
Electromagnetic radiation disturbance (radiated emissions) (CISPR 11)	Group 1	This medical device uses RF energy for its internal / external operation.
Disturbance voltage on the supply terminals (Conducted emissions) (CISPR 11)	Class B	Home health care environment, and environment of a professional healthcare establishment.
Harmonic current emissions (IEC61000-3-2)	Class A	
Voltage changes, voltage fluctuations and flicker (IEC61000-3-3)	Compliant	

Immunity

This medical device is designed for use in a home health care environment, and the environment of a professional healthcare establishment. The user and installer must ensure the electromagnetic compliance of the environment.

Immunity test	Test level as per IEC60601-1-2 Ed4		Compliance level	
Electrostatic discharges (DES) (IEC61000-4-2)	± 8 kV in contact ± 2/4/8/15 kV in air		± 8 kV in contact ± 2/4/8/15 kV in air	
Radiated radiofrequency electromagnetic fields. (IEC61000-4-3)	10 V/m 80 MHz - 2.7 GHz 80 % AM at 1 kHz		12 V/m 80 MHz – 2.7 GHz 80 % AM at 1 kHz	
Proximity fields emitted by 450 MHz wireless communication devices (IEC 61000-4-3)	Frequency (MHz)	Modulation	Required level. (V/m)	Compliance level (V/m)
	385	Pulsed modulation: 18 Hz	27	27
	450	Pulsed modulation: 18 Hz	28	28
	745 - 780	Pulsed modulation: 217 Hz	9	9
	810 – 870 - 930	Pulsed modulation: 18 Hz	28	28
	1720 – 1845 - 1970	Pulsed modulation: 217 Hz	28	28
	2450 – 5240 – 5500 - 5785	Pulsed modulation: 217 Hz	28 9	28 9
Electrical fast burst transients (IEC61000-4-4)	Power supply: ± 2 kV Input/output lines: ± 1 kV Repetition frequency: 100 kHz		Power supply: ± 2 kV Input/output lines: ± 1 kV Repetition frequency: 100 kHz	
Shock waves (IEC61000-4-5)	Between phases: ± 0.5 kV, ± 1 kV Between earth and phases: ± 0.5 kV, ± 1 kV, ± 2 kV		Between phases: ± 0.5 kV, ± 1 kV Between earth and phases: ± 0.5 kV, ± 1 kV, ± 2 kV	
Conducted disturbances, induced by RF fields (IEC610004-6)	3 V 0.15 MHz – 80 MHz 6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz		3 V 0.15 MHz – 80 MHz 6 V in ISM and amateur radio bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz	
Magnetic field at network frequency (IEC61000-4-8)	30 A/m		30 A/m	

Immunity test	Test level as per IEC60601-1-2 Ed4	Compliance level
Voltage dips and interruptions (IEC61000-4-11)	0 % UT; 0.5 cycles At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle at 0° 70 % UT; 25/30 cycles at 0° 0 % UT; 250/300 cycles	0 % UT; 0.5 cycles At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % UT; 1 cycle at 0° 70 % UT; 25/30 cycles at 0° 0 % UT; 250/300 cycles

13. Maintenance, storage, disposal, and service life

A. Maintenance / cleaning

For the patient's safety, adhere to the prescribed overhaul and inspection intervals:

- Orthopedic prosthetist checking interval: 1 year.
- Overhaul interval (PROTEOR After-Sales): 24 months maximum, or after 1 million steps (activity data visible in the SYNSYS SOFT configuration software).

Before sending the prosthetic component to After-Sales, remove the foot (see §7.B.), and refit it on the Loaner Unit. SYNSYS will be sent back fitted with a new foot.

The foot shell may be changed every year.

The parameters of a replacement knee or After-Sales return have been preset, but you should check their configuration upon receiving and installing any new knee.

The certified prosthetist should be the first contact in case of failure or for servicing. However, the patient may approach any PROTEOR certified prosthetist.

On request, PROTEOR will provide the diagrams, component lists, descriptions, calibration instructions or any other information useful for the certified prosthetist to replace the spare parts or perform annual maintenance operations.

To contact us (assistance, if necessary, with assembly, use or maintenance; or to report an unforeseen event): manufacturer's e-mail or website: www.proteor.com

Spare parts / accessories

Use only PROTEOR spare parts and accessories. Using accessories or spare parts not listed in the instructions may be dangerous.

Instructions in §7.F. :

Foot shell kit (including sock and ankle cover):


X= G or D; YY = [23-30]

Caucasian	1A70020-XYX-S
Black	1A70020-XYYN-S
Brown	1A70020-XYXB-S

 Make sure that the ankle has good mobility after replacing these components.

- Cover: 1P700218-XX-S (XX = colors = 01, 02, 03...)

Instructions:

- Charger 1P700234-S (manufacturer: PROTEOR) – 155 cm.
- Plantar flexion resistance setting tool kit: 1P700139-S
- Pump kit: 1P700230-S
- ON/OFF tool kit 1P700132-S

Spare parts availability period: 6 years.

Ongoing maintenance

Make sure that the device is not operating, before performing any servicing or maintenance.

Cleaning

Clean the prosthesis using a wet soft cloth and mild soap, following exposure to dust or moisture. In particular, keep articulated area of the knee clean (magnet). Take care with the small mobile covers (access to socket for charging). Never use any solvents, which could damage the device. Check that there are no metal particles on the magnetic connector and the articulated area.

In case of high perspiration in the socket, and use of vacuum suspension, make sure the prosthesis does not get excessively wet.

Repressurizing the hydraulic unit

The hydraulic circuit is kept under pressure to ensure smooth operation of the prosthesis. As in any pneumatic circuit, the air pressure will decrease over time. The rate of pressure loss will depend on the use and environment of use of the prosthesis (e.g.: frequent temperature variations). We recommend checking the pressure or having it checked once a year. Use the pump supplied during certification (1P700230-S). Do not use an electric pump.

DO NOT EXCEED A PRESSURE OF 10 BARS, nor a flowrate of more than 1L/min.

Procedure: knee extended: 1- unscrew the cover, 2- unscrew the cap, 3- screw down the pump, making sure that the connection is fully tightened. Adjust the pressure.

 Min. max. pressure range: 8 – 10 bar.

Once the correct pressure has been set, unscrew the pump, and screw the cap back on when in contact with the seal, tighten it by another $\frac{3}{4}$ of a turn. Then firmly reclose the cover. If this cover 1P700162 is missing or damaged, replace it.

B. Environment for use:

Temperature for use: -10°C to $+40^{\circ}\text{C}$.

Atmospheric pressure: 700 – 1060 hPa.

Air humidity (RH): 15% - 90%. No condensation

C. Storage

Storage and transport temperature: RH 15% - 90% with no condensation

- Long term (more than 1 week): between 0°C and 30°C
- Short term (less than one week): between -20°C and $+60^{\circ}\text{C}$

Following storage at -20°C , wait for 16 minutes for the knee to become operational ($\geq -10^{\circ}\text{C}$).

Following storage at 60°C (140°F), wait for 35 minutes for the knee to become operational ($\leq 40^{\circ}\text{C}$ / 104°F).


For transport, switch off the SYNSYS and pack it properly. Avoid leaving the prosthesis in places with excessive temperatures.

Preferably store the prosthesis in the vertical position.

After prolonged storage, leave the prosthesis vertical for 2 hours, charge the prosthesis, test the knee resistance and check that the vibrator is operating properly; dust off.

To switch the SYNSYS off:

1. Bend the knee then press & hold the two buttons simultaneously for 10 seconds. Check if SYNSYS is effectively off by pushing buttons. There should be no reaction from the device.

2. **OR** via the ON/OFF button using the tool 1P700132-S supplied: press & hold for 10 seconds.

D. Disposal

This device is made of components made from various materials: elastomer, plastic, composite (carbon fiber), aluminum, titanium, steel, brass. It also contains oil, a Lithium-ion battery, and electronic components, which must not be disposed of with household waste. Avoid any contact with oil. The device and packaging must be disposed of in accordance with current local or national regulations.

E. Service life

This device is tested as per standard NF EN ISO 10328, based on 3 million cycles. According to the user's activity level, this corresponds to 3 to 6 years. Refer to the prosthesis' activity data. If the predicted service life is exceeded, the risk of device failure is increased. This device is prescribed for one user only; any resale is prohibited.

Battery lifetime: 500 complete charge cycles (approximately 2 years).

F. Additional information

Electrical network compatible with charger 1P700234-S (manufacturer: PROTEOR):

Input current: 0.3A max. at 100VAC	Output current: 2A max.
Input voltage: AC 100-240 V	Output voltage: 5V DC
Low Power Supply (LPS) compliant	Frequency: 50-60 Hz
Degree of protection: IP21	

Compatibility with standard Bluetooth v2.1 + EDR – class 2 (2.5 mW)
[2400 – 2483.5] MHz

EMC emission class:

Class B – use in residential environment as per standard IEC 60601-1-2. Immunity test level = 12v/m (as per standard EN22523).

Degree of protection as per standard IEC 60529:

IP54: Protected against dust and other microscopic residue.
Protected against water splashes from any directions.

Materials or substances in contact: carbon fiber frame, thermoplastic covers and foot shells, oiled piston rod, socket as chosen by the certified prosthetist.

14. Regulatory information


This product is tested and certified compliant with Regulation EU MDR 2017/745, and with Class B under standards IEC 60601-1, 60601-1-2, 60601-1-4, 60601-1-11, ETSI EN 300-328, EN 62311 and ISO 10328 - P6 - 125 kg.


The device complies with part 15 of the FCC regulations, and bears the CE marking

RECAP OF BUTTON ACTIONS


Battery
charge
level


Activate
Bluetooth
(Active
time: 2
mins)


Adjust
heel
height


Switch to
cycling
mode

Activate /
deactivate
locking
mode

Switching
OFF

ALERTS

Low battery


Diagram illustrating low battery alerts:

- A yellow triangular warning sign with a battery icon and a lightning bolt through it.
- Three battery icons in a row, each with a small black dot in the center, followed by "x 3 = [battery icon] < 10%".
- Five battery icons in a row, each with a small black dot in the center, followed by "x 5 = [battery icon] < 5%".
- Ten battery icons in a row, each with a small black dot in the center, followed by "x 10 = [battery icon] → [moon icon] [warning sign icon]".

Overheating


Diagram illustrating overheating alerts:

- A thermometer icon with a red liquid level and a yellow sunburst above it, indicating high temperature.
- A yellow triangular warning sign with three wavy lines above it.
- A battery icon with a small black dot in the center.
- An arrow pointing to a grey circular button with two vertical bars (pause symbol).
- A second yellow triangular warning sign with three wavy lines above it.
- A battery icon with a small black dot in the center, followed by a double battery icon with a "5''" below it.
- An arrow pointing to a red octagonal STOP sign.

Fault


Diagram illustrating fault alerts:

- A cracked glass icon with a red question mark overlaid on it.
- A yellow triangular warning sign with an exclamation mark inside.
- A battery icon with a small black dot in the center, followed by a double battery icon.
- An arrow pointing to a grey circular button with two vertical bars (pause symbol).
- A second yellow triangular warning sign with an exclamation mark inside.
- A battery icon with a small black dot in the center, followed by a double battery icon, followed by "x10".
- An arrow pointing to a red octagonal STOP sign.

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The screenshot shows the SYNSYS software interface for 'Gehmodus 2'. The main settings area includes:

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On the right, there is an 'Anleitung' (Instructions) panel with the following text:

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The screenshot shows the SYNSYS software interface for 'Alarme abgeben' (Give alarms). The main settings area includes:

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On the right, there is an 'Alarme abgeben' panel with the following text:

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Below the panel, there is a list of alarm levels:

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PROTEOR SAS