# **RUSH** FOOT<sup>®</sup>

by  PROTEOR

**RUSH Foot<sup>®</sup>**  
Instructions for

07 2021

## ***Table of Contents / Tabla de contenido / Sommaire / Inhaltsangabe Indice / Obsah***

**General Product Statement / Información del producto /  
Information produit / Produktinformationen /  
Informazioni sul prodotto / Informace o výrobku**

Page 3

**English**  
Pages 4 - 7

**Spanish / Español**  
Pages 8—11

**French / Français**  
Pages 12 - 15

**German / Deutsch**  
Pages 16 - 19

**Italian / Italiano**  
Pages 20 - 23

**Czech / Český**  
Pages 24 - 27

**Japanese / 日本語**  
Pages 28 - 31

## General Product Statement

**ENGLISH** : The RUSH Foot<sup>®</sup>'s superior engineering delivers the most realistic and responsive foot and ankle motion available. It's unique glass composite material eliminates prevalent weak points and is exceedingly more flexible than carbon fiber. The RUSH Foot<sup>®</sup> delivers an unmatched range of motion from the roll-through design which allows for a smooth transition from heel strike through toe-off with no "dead spots" in between. The innovative design helps maintain a natural, smooth gait, with superior energy return, even in rugged or uneven terrain. The radically simplistic design demonstrates that stress points won't fail when subjected to pressure, water, or impact. Virtually indestructible even in the most extreme conditions.

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### H2O Collection

#### *Instructions for Use*

**This instructional booklet is intended for the Prosthetist and the patient.  
Read it carefully and completely before any use.**

**Please transmit to the patient the safety instructions and all the information  
necessary for the proper use of the product (Sections 8 & 9).**

### 1—Description and Functions

The fiberglass composite material of the RUSH HiPro H2O and Rampage H2O feet provide great flexibility, resistance and durability.

Rampage LP H2O feet offer the same great flexibility, resistance and durability as the HiPro H2O and RAMPAGE H2O in a lower build height configuration.

Each H2O RUSH foot<sup>®</sup> is comprised of the following components:

- An upper blade that allows energy storage and return.
- A continuous lower blade with a patented design that ensures progressive flattening of the foot and smooth step rollover.
- A unique structural bond connecting the two blades that eliminates fasteners.
- A high energy return rubber heel.
- A Vibram Sole and a protective rubber toe cap.

### 2—Contents at Arrival

1. H2O RUSH Foot<sup>®</sup>

### 3—Indications for Use and Maximum Weight 365lbs/166kg

The HiPro H2O, RAMPAGE H2O, and RAMPAGE LP H2O feet are designed exclusively for lower limb prostheses.

RUSH H2O feet are intended for barefoot activities only. Acceptable activities include:

showering, pool activities, surfing, kayaking, stand up paddle boarding, indoor yoga, karate, gymnastics, other miscellaneous barefoot activities.

They are appropriate for patients with low, moderate, or high impact activity levels whose weight (carried load included) complies with the table below.

Advice: Patients of similar weight and size often require different category selections due to the activities they perform in their day to day life. Patients who perform high impact activities such as sprinting and basketball benefit from a slightly stiffer and more responsive foot compared to a patient who only performs low activities such as walking and playing golf. Proteor USA has created four different activity levels with examples to aid Prosthetists in selecting the correct activity level for their patient.

The four activity levels and their examples are the following:

- Low: Walking, Showering
- Moderate: Paddle Boarding, Yoga
- High: Water Skiing, Wakeboarding
- High Impact: Gymnastics, Karate.

The authorized maximum load is 365 lbs / 166 kg.

Each H2O RUSH foot® is resistant to fresh, sea, and chlorinated water provided that the recommendations of Sections 8 and 9 are respected.

		HiPro, & RAMPAGE / Rampage LP (x)								
Patient Weight	lbs	0 to 105	106 to 140	141 to 175	176 to 210	211 to 245	246 to 280	281 to 315	316 to 350	351 to 365
	kg	0 to 48	49 to 64	65 to 79	80 to 95	96 to 111	112 to 127	128 to 143	144 to 159	160 to 166
Low		1	2	3	4	5	6	7	8	9
Moderate		1	2	3	4	5	6	7	8	9
High		1 (2)	2 (3)	3 (4)	4 (5)	5 (6)	6 (7)	7 (8)	8 (9)	9 (-)
High Impact		2 (3)	3 (4)	4 (5)	5 (6)	6 (7)	7 (8)	8 (9)	9 (-)	- (-)

They include a male pyramid that allows the use of adult pyramid receivers.

To select the correct category for the patient, the Prosthetist must know the

### 4—Contraindications for Use


RUSH H2O Feet are intended for bare feet activities only. Usage during activities where shoes would typically be worn will prematurely wear the Vibram sole and void the warranty of the foot. Activities that are not permitted include:

trail hiking, mountain biking, dirt biking, running, skateboarding, and use in extreme/aggressive terrain.

### 5—Bench Alignment

#### Sagittal Plane:

After determining the appropriate socket flexion and heel height, position the weight line so that it falls through the center of the foot adapter. As the unique glass composite material provides substantially more flexibility than the other prosthetic feet, this represents a good starting point for a bench alignment.

#### Coronal Plane:

After determining the appropriate socket adduction/abduction, position the weight line so that it falls through the midline of the foot in a neutral M-L position.

### 6—Alignment (With the Patient)

#### Static Alignment:

Due to the shape of the rocker sole (bottom blade), patients may experience a new sensation when seeking the midpoint of the foot. The rocker sole allows patients to find their own comfortable static or standing position. The rocker bottom shape of the sole serves 2 primary functions :

- Provide a continuous and progressive point of contact throughout the entire step.
- Eliminate any “flat” or “dead” spot.

Adjustment of the A/P set screws at the proximal adapter is the more appropriate place to make changes regarding plantar flexion or dorsiflexion.

#### Dynamic Alignment:

The use of a slide adapter with the H2O RUSH Foot® is highly encouraged for dynamic alignment purposes, as it provides the best solution to troubleshooting the most common alignment issues listed below:

- Hard or soft heel
- Hard or soft toe
- Varus or Valgus movements during Stance Phase.

Once the optimal relative socket/foot alignment is determined, H2O RUSH Foot® encourages the use of plantar flexion / dorsiflexion at the proximal foot adapter to optimize heel-to-toe comfort and energy return.

The patient should wear the H2O RUSH Foot® for at least one week to fully experience the responsiveness of a non-carbon fiber device.

## 7—Product Life

This component has been tested cyclically, which corresponds to at least 3 years of use depending on the activity of the patient.

An annual check by the Prosthetist is recommended.

## 8—Maintenance

No maintenance operation such as lubrication, work on the screws, or other parts is required for this foot.

The patient is encouraged to periodically clean the foot to remove any dirt or debris which could alter the foot behavior or damage it.

## 9—Advice for Use and Safety

 Even though the foot is meant for activities in water, it does not guarantee that the patient will not slip and fall while using it (example: slipping in the bathtub)

 Use temperature : - 20°C / + 60°C (- 4°F / 140°F) .

 Use of the foot that does not adhere to your Prosthetists recommendations can damage foot components (carrying heavy loads for example).

 In case of an abnormal behavior of your foot, or if you feel changes in your foot features, or if the foot sustained a hard impact, contact your Prosthetist.

 After use in water or in a very dirty environment, rinse the foot and the other parts thoroughly with clean water to eliminate debris such as dirt, sand, rock, etc.

**Failure to follow these guidelines may result in premature failure and void the foot warranty.**

## 10—Recycling

The different items of the foot are special wastes and must be handled according to local laws.

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Proudly handcrafted in the U.S.A.
**PROTEOR SAS**