## RUSH HiPro®

## RUSH RAMPAGE®

## RUSH RAMPAGE LP®

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EN Instruction for use

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NO Bruerveiledning

SV Bruksanvisning

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SK Návod na používanie

HR Upute za uporabu

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Rev. B

2023-10

**Pass on § 11, 12, 13, 14, 15, 17, 18 and 19 of these instructions to the patient.**
**1. INCLUDED ITEMS**

RUSH HiPro, RAMPAGE, RAMPAGE LP		
Part description	Part number	Included / Sold separately
RUSH HiPro	HIP-xx-x-xx*	Included
RUSH RAMPAGE	RAM-xx-x-x*	Included
RUSH RAMPAGE LP	RAML-xx-x-xx*	Included
Spectra sock	S0-NPS-200xx-00*	Suitable sock included
Foot shell	FS-xx-xx* FTC-4R-1-xxxxx-xx*	Sold separately
RUSH HiPro, RAMPAGE, RAMPAGE LP EVAQ8		
Part description	Part number	Included / Sold separately
RUSH HiPro EVAQ8	EVQH-xx-x-xx*	Included
RUSH RAMPAGE EVAQ8	EVRAM-xx-x-xx*	Included
RUSH RAMPAGE LP EVAQ8	EVQRAML-xx-x-xx*	Included
Spectra sock	S0-NPS-200xx-00*	Suitable sock included
Foot shell	FS-xx-xx* FTC-4R-1-xxxxx-xx*	Sold separately
RUSH EVAQ8 Rebuild Kit	EVRB	Sold separately
RUSH EVAQ8 Release Valve	EV RV	Sold separately
RUSH HiPro, RAMPAGE, RAMPAGE LP H2O		
Part description	Part number	Included / Sold separately
RUSH HiPro H2O	H2H-xx-x-xx*	Included
RUSH RAMPAGE H2O	H2RAM-xx-x-xx*	Included
RUSH RAMPAGE LP H2O	H2RAML-xx-x-xx*	Included

*\* see catalog*
**2. DESCRIPTION**

RUSH HiPro, RAMPAGE and RAMPAGE LP are energy-return prosthetic feet consisting of:

- A fiberglass composite keel
- A fiberglass composite sole plate
- A rubber heel shock bumper
- A male pyramid connection
- A Spectra sock

RUSH HiPro EVAQ8, RAMPAGE EVAQ8 and RAMPAGE LP EVAQ8 incorporate an elevated vacuum system inside the rubber heel shock bumper. RUSH HiPro H2O, RAMPAGE H2O and RAMPAGE LP H2O integrate a Vibram<sup>®</sup> sole and custom molded rubber toe-wrap.

**3. PROPERTIES**

Model		RUSH HiPro	RUSH RAMPAGE	RUSH RAMPAGE LP
Weight*		920 g / 2.03 lb	850 g / 1.87 lb	805 g / 1.77 lb
Build height*	22 to 24 cm	165 mm / 6.50"	168 mm / 6.63"	114 mm / 4.50"
	25 to 27 cm	175 mm / 6.88"	181 mm / 7.13"	121 mm / 4.75"
	28 to 29 cm	181 mm / 7.13"	184 mm / 7.25"	121 mm / 4.75"
Heel height		10 mm / 3/8"		
Model		RUSH HiPro EVAQ8	RUSH RAMPAGE EVAQ8	RUSH RAMPAGE LP EVAQ8
Weight*		930 g / 2.05 lb	860 g / 1.90 lb	815 g / 1.80 lb
Build height*	22 to 24 cm	165 mm / 6.50"	168 mm / 6.63"	114 mm / 4.50"
	25 to 27 cm	175 mm / 6.88"	181 mm / 7.13"	121 mm / 4.75"

	28 to 29 cm	181 mm / 7.13"	184 mm / 7.25"	121 mm / 4.75"
<b>Heel height</b>	10 mm / 3/8"			
<b>Model</b>	<b>RUSH HiPro H2O</b>	<b>RUSH RAMPAGE H2O</b>	<b>RUSH RAMPAGE LP H2O</b>	
<b>Weight*</b>	765 g / 1.69 lb	695 g / 1.53 lb	650 g / 1.43 lb	
<b>Build height*</b>	22 to 24 cm	152 mm / 6.00"	159 mm / 6.25"	105 mm / 4.13"
	25 to 27 cm	159 mm / 6.25"	165 mm / 6.50"	105 mm / 4.13"
	28 to 29 cm	175 mm / 6.88"	178 mm / 7.00"	111 mm / 4.38"
<b>Heel height</b>	10 mm / 3/8"			

\*Weight based on size 26, cat. 4, with foot shell and Spectra sock

Build height based on sizes 23, 26 or 29, cat. 4, with foot shell, Spectra sock, and 10 mm heel height

H2O weight and build height does not include foot shell or Spectra sock

These devices have been tested according to ISO 10328 for a maximum patient weight up to 166 kg for 2 million cycles.

Selection of foot category based on patient's weight and impact level – RUSH HiPro & RUSH RAMPAGE										
<b>Weight*)</b>	<b>lb</b>	0-105	106-140	141-175	176-210	211-245	246-280	281-315	316-350	351-365
	<b>kg</b>	0-48	49-64	65-79	80-95	96-111	112-127	128-143	144-159	160-166
<b>Impact level</b>	<b>Low</b>	1	2	3	4	5	6	7	8	9
	<b>Moderate</b>	1	2	3	4	5	6	7	8	9
	<b>High</b>	2	3	4	5	6	7	8	9	-

\*) Body mass limit not to be exceeded (ISO 10328)

Selection of foot category based on patient's weight and impact level – RUSH RAMPAGE LP										
<b>Weight*)</b>	<b>lb</b>	0-105	106-140	141-175	176-210	211-245	246-280	281-315	316-350	351-365
	<b>kg</b>	0-48	49-64	65-79	80-95	96-111	112-127	128-143	144-159	160-166
<b>Impact level</b>	<b>Low</b>	1	2	3	4	5	6	7	8	9
	<b>Moderate</b>	2	3	4	5	6	7	8	9	-
	<b>High</b>	3	4	5	6	7	8	9	9	-

\*) Body mass limit not to be exceeded (ISO 10328)


#### 4. MECHANISM OF ACTION

The half rocker sole plate provides constant contact to the ground while eliminating dead spots, providing a smooth rollover, and natural energy return. The heel shock bumper stores energy during initial contact and loading response which is then released during midstance. The instep absorbs energy during terminal stance and releases the stored energy during pre-swing.

#### 5. INTENDED USE/INDICATIONS

These medical devices are supplied to healthcare professionals (prosthetists), who will train the patient in their use. The prescription is drawn up by a doctor who assesses the patient's ability to use the device.

⚠️ These devices are for multiple use on a **SINGLE PATIENT**. They must not be used on another patient.

These devices are intended to be integrated in a custom-made external lower limb prosthesis to ensure the function of the foot in patients with unilateral or bilateral lower limb amputation or deficiencies (transtibial/transfemoral amputation, knee/hip disarticulation, congenital limb deficiencies). They are intended for patients who would benefit from the dynamic toe response and vertical compression.

These devices are indicated for patients with moderate to high activity level (K3 to K4) for walking and physical activities without excessive overload.

Maximum weight (load carrying included): 166 kg / 365 lb (See table §3)

#### 6. CLINICAL BENEFITS

- Walking comfort
- Possible ambulation on uneven ground
- Shock absorption and reduction of socket forces

#### 7. ACCESSORIES AND COMPATIBILITY

For RUSH HiPro, RAMPAGE, RAMPAGE LP and EVAQ8 models, an appropriate foot shell must be mounted on the foot (refer to our catalog).

The foot includes a male pyramid connection designed to be compatible with standard female pyramid connectors (refer to our catalog).

## 8. ALIGNMENTS

### Bench alignment


**Recommendation:** Add a 10 mm ( $\frac{3}{8}$ " ) wedge under back of the foot during bench alignment for all but the H2O models, as shown below.

Sagittal plane:

- Define the appropriate socket flexion and heel height.
- Position the load line so that it falls through the center of the foot pyramid connection.

Frontal plane:

- Define the appropriate socket adduction/abduction.
- Position the load line so that it falls through the midline of the foot in a neutral mediolateral position.

### Static alignment

Due to the shape of the rocker sole (bottom blade), patients may experience a new sensation when seeking the midpoint of the foot. The rocker sole allows patients to find their own comfortable static or standing position.

Adjustment of the anteroposterior set screws at the proximal connector is the more appropriate place to make changes regarding plantarflexion or dorsiflexion.

⚠ The heel wedges must not be introduced at this stage of alignment.

### Dynamic alignment

The use of a slide adapter is recommended for dynamic alignment purposes, as it provides the best solution to troubleshooting the most common alignment issues listed below:

- Hard or soft heel
- Hard or soft toe
- Varus or valgus movements during stance phase

Once the optimal relative socket/foot alignment is determined, use plantarflexion or dorsiflexion at the proximal foot adapter to optimize heel-to-toe comfort and energy return.

## 9. ASSEMBLY

After dynamic alignment, tighten the pyramid adjustment screws according to the specifications of the connector manufacturer. Secure pyramid adjustment screws with a thread locking adhesive (e.g., Loctite 242).

### EVAQ8 models

The straight barb, exhaust filter, tubing, inline filter, socket right angle barb and housing, tube securing hook and loop tape are included with the foot and will require assembly prior to use.

To connect the vacuum system to the socket:

- Connect a short piece of tubing to the socket right angle barb.
- Install the inline filter into the distal end of the tubing (the inline filter may be placed anywhere in the tubing between the socket and vacuum pump heel).
- Connect a second piece of tubing to the distal end of the inline filter and route the tube to the medial side of the pylon or wrap the tubing around the pylon (to prevent damage to the tube or snagging while walking).
- Cut the tubing to desired length and connect to the straight barb located in the recessed area of the vacuum pump heel.
- Secure the tubing to the pylon using the included hook and loop tape or other appropriate tape.

### Spectra sock

A Spectra sock is included to protect the foot shell and fiberglass components and minimize noise. It must be placed over the keel and sole plates before mounting the foot shell.

**NOTE:** The H2O models are not intended to be worn with a foot shell. Therefore, they do not include a Spectra sock.

### Foot shell

To install and remove the foot shell, use a shoehorn to prevent damage to the foot module.

- ⚠ Never remove the foot from the foot shell by pulling manually. Never use a screwdriver or any other inappropriate instrument to remove it. This could damage the foot.
- ⚠ The H2O models are not intended to be worn with a foot shell. Doing so may damage the foot module.

## 10. ADJUSTMENTS

If the patient still requests additional heel stiffness, use the heel wedges provided. These heel wedges use a double-sided sticker so the lower face of the sole plate should be degreased before using. See figures below for correct wedge placement.

**NOTE:** The H2O models do not include heel wedges.

CORRECT HEEL WEDGE PLACEMENT	INCORRECT HEEL WEDGE PLACEMENT
	
<p>Install Heel Wedge on plantar surface of the sole plate approximately 3.2 mm (1/8") forward from posterior end.</p>	<p>Placement between components will adversely affect foot performance.</p>

## 11. TROUBLESHOOTING

- ⚠ If the patient notices any abnormal behavior or feel any changes in the characteristics of the device (noise, play, excessive wear...), or if the device has received a severe impact, they should stop using the device and consult their prosthetist.

## 12. WARNINGS

- ⚠ In case of damaged packaging, check the integrity of the device.
- ⚠ Never use the foot without a foot shell and a Spectra sock (except for H2O models).
- ⚠ Never loosen the pyramid fastening screws.
- ⚠ The patient must inform their prosthetist if they gain or lose significant weight.
- ⚠ Make sure that the foot and inside of the foot shell are free of impurities (e.g., sand). The presence of impurities causes the fiberglass parts and the foot shell to wear out. Clean the foot according to the instructions (see §17).
- ⚠ After swimming or using in water, the foot including its foot shell must be cleaned (see §17).
- ⚠ Failure to follow the instructions for use is dangerous and will void the warranty.

## 13. CONTRAINDICATIONS

- ⚠ Do not use for a patient whose maximum weight (load carrying included) may exceed 166 kg / 365 lbs.
- ⚠ Do not use for a patient who does not meet the requirements of an activity level of K3 or higher.
- ⚠ Do not use for activities associated with a risk of significant impact or excessive overloading.

## 14. SIDE EFFECTS

There are no known side effects directly associated with the devices.

**Any serious incident that has occurred in relation to the devices should be reported to the manufacturer and to the competent authority of the Member State in which the user is established.**

**15. MAINTENANCE AND CONTROL**

No maintenance operation such as lubrication, work on the screws or other parts is required.




The foot module must be inspected by the prosthetist at least every six months. Inspections at shorter intervals are required if the user is more active.

The Spectra sock and the foot shell must be replaced by the prosthetist at regular intervals, depending on the patient’s level of activity. If these parts are damaged, it can lead to premature foot wear.

The **EVAQ8 models’** components (tubing, inline filter, one-way valves housed inside the vacuum heel) may need periodic cleaning or replacement during the life cycle of the system and are not replaceable under the warranty as it is considered normal wear.

**16. PERIODIC INSPECTION OF THE EVAQ8 SYSTEM**

- Visually inspect the tubing for kinks, cracks, or wear that may leak air into the system. Replace tubing if any of these conditions exist.
- Remove the inline filter from tubing and look through it. If light can be seen, the filter is clean. If light is blocked, blast air from a syringe through the inline filter from distal to proximal end (reverse of normal flow) to attempt to clear the blockage. If blockage persists, the filter needs to be replaced.
- The one-way valves contained in the vacuum heel may need to be cleaned and flushed with distilled water or isopropyl alcohol to ensure proper function. This procedure should be done only by a qualified professional.
- To flush the one-way valves and vacuum heel:

					
1. Exhaust Assembly 2. Exhaust Filter 3. Large O-Ring 4. Valve Body Adapter 5. Small O-Ring	2. Exhaust Filter 6. Duckbill Valve	7. Straight Barb	8. Duckbill Valve	9. Paperclip (foot on side)	9. Paperclip (foot upright)

- Remove the vacuum hose from the user’s socket keeping it attached to the EVAQ8 foot.
- Remove the EVAQ8 foot from the user’s socket.
- Remove the vacuum hose from the EVAQ8 foot.
- Using a 5/16” socket, remove the (1) Exhaust Assembly from the heel, the (4) Valve Body Adapter will most likely stay attached to the (2) Exhaust Filter.  
**NOTE:** If (4) Valve Body Adapter does not remove with (2) Exhaust Filter, use a flat blade screwdriver to remove it.
- Place the (4) Valve Body Adapter into a soft jaw vice or grip with soft face pliers and remove the (2) Exhaust Filter with 5/16” socket.  
**NOTE:** (6) Duckbill Valve will be slightly pressed into the bottom of the exhaust filter.
- Using a 1/4” deep socket, remove the (7) Straight Barb from the other side of the heel.
- Inside the heel beneath where the (7) Straight Barb was located is another (8) Duckbill Valve. Remove the (8) Duckbill Valve by tapping the foot against your hand or by straightening out a (9) Paperclip and inserting it into the other side of the heel to push out the (8) Duckbill Valve.
- Inspect both (3 & 5) O-rings on the (4) Valve Body Adapter. There is one at the base of the threads, and one in a groove on the body. Replace both if they show any wear.
- Clean the female threads in both sides of the heel with a cotton swab and isopropyl alcohol or distilled water.
- If you are reusing the (1) Exhaust Filter, (7) Straight Barb, in-line filter, and (8) Duckbill Valve, clean them with isopropyl alcohol or distilled water. Take extra care and make sure the (8) Duckbill Valve are clean and free from debris (a magnifier is helpful for inspection). Flush the in-line filter from both directions to ensure it is clean. **Allow valves to air-dry. DO NOT use a towel or cloth.**
- Once the parts are dry, or if you are using new parts from a rebuild kit, set all the parts out on a clean surface.
- Insert the turquoise (8) Duckbill Valve into the threaded opening of the (1) Exhaust Filter so that the (8) Duckbill Valve rim is flush against the opening of the (1) Exhaust Filter and the tip of the (8) Duckbill Valve is inside the (1) Exhaust Filter. Thread the (1) Exhaust Filter into the (4) Valve Body Adapter **BY HAND** until snug.
- **BY HAND**, thread the assembled (1) Exhaust Assembly into side of the heel that will **NOT** be connected to the prosthetic socket by vacuum hose.
- Once the (1) Exhaust Assembly is threaded in snugly **BY HAND**, torque to 15 in-lbs. Do not over torque. Over torquing will snap the threads and will not be covered under warranty.

**NOTE:** If you do not have a torque wrench, thread the (1) Exhaust Assembly until you feel a hard stop and then turn  $\frac{1}{16}$  of a revolution more.

- Insert the white (8) Duckbill Valve into the recessed side of the heel so the tip of the (8) Duckbill Valve points into the heel. Using a small screwdriver or a straightened-out (9) paperclip, make sure the (8) Duckbill Valve is seated all the way into the recess.
- **BY HAND**, thread the (7) Straight Barb into the side of the heel that is recessed.
- Once the (7) Straight Barb is threaded in snugly **BY HAND**, torque it to 15 in-lbf. This is a very low torque value, and over torquing will snap the threads on the (7) Straight Barb and will not be covered under warranty.
- If you do not have a torque wrench, thread the (7) Straight Barb until you feel a hard stop and then turn  $\frac{1}{16}$  of a revolution more.
- Reattach vacuum hose to EVAQ8 foot by sliding the vacuum hose over the (7) Straight Barb.
- Put the Spectra Sock and foot shell over the EVAQ8 foot.
- Reattach the EVAQ8 foot to the user’s socket.
- Reattach the other end of the vacuum hose to the user’s socket. The vacuum hose may be routed to the prosthetist’s preference.

**17. CLEANING**

Remove the foot shell and the Spectra sock, rinse the foot in clear water, clean it with neutral soap and dry carefully. The foot shell can be cleaned with a damp cloth or sponge. It must be dry before reuse.

The devices are not resistant to solvents. Exposure to solvents may cause damage.

**18. ENVIRONMENTAL CONDITIONS**

Temperature range for use and storage: -20°C to 60°C [-4 to 140° F]

Relative air humidity: no restrictions

**Waterproof: the devices are resistant to fresh, sea, and chlorinated water.**

The foot shell is not resistant to ultraviolet light (UV). Do not store it under direct sunlight.

**19. DISPOSAL**

The devices are made of fiberglass composite material, rubber, Vibram (H2O), metal and plastic (EVAQ8). The devices and their packaging must be disposed of in accordance with local or national environmental regulations.

**20. DESCRIPTION OF THE SYMBOLS**

	Manufacturer		Identified risk		CE marking and year of 1 <sup>st</sup> declaration
	Authorized representative in the European Union		Single patient, multiple use		

**21. REGULATORY INFORMATION**

These products are CE-marked medical devices and are certified as conforming with Regulation (EU) 2017/745.

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### **PROTEOR SAS**