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
**I - USE**

These components are designed for orthopedic fitting of patients :

- Transtibial amputees weighing (carried load included) no more than 150 Kg for 1K204 & 1K203
- Transtibial amputees weighing (carried load included) no more than 125 Kg for 1K202
- Transfemoral amputees weighing (carried load included) no more than 125 Kg for 1K173

**II - PACK CONTENTS**

- 1 lamination connector (4 short prongs for 1K202, 1K204, 1K203 and 3 long prongs for 1K173),
- 1 spacer for lamination,
- This instruction sheet.

**III - SHAPING THE LAMINATION CONNECTOR**

It is highly recommended to adapt the shape of the connector prongs according to the positive.

**IV- RECOMMENDATIONS FOR LAMINATION**

A) Cover positive (1) with one or more layers of glass or carbon fabric (2) (according to the desired weight and strength).

B) For reinforcement purpose, install a glass or carbon disc (3) where the lamination connector will be placed.

C) In order to reinforce the assembly, make braids (5) of glass or carbon fabrics pass through the holes of lamination connector (4) .

D) Position lamination connector (4) onto positive (1) and check perpendicularity.

E) Install second disc (6) of glass or carbon fabric. The pyramid must necessarily be placed through this disc.

F) Clip spacer (7) on lamination connector (4) .

G) Finish covering with one or more layers of glass or carbon fabric (8). Enclose the assembly in PVA bag (9) and inject resin. Before resin becomes hard, put string (10) around groove (11) located at the base of the lamination connector prongs. This will ensure a better strength.

H) When the resin is hard enough, cut off material excess and remove the spacer.

**V - CONNECTION**

1K173, 1K202, 1K203 and 1K204 can be used with any standard female pyramidal (e.g. 1K172, 1D41-HD , etc...)

**VI - RECYCLING**