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


A) PACK CONTENTS (Fig.1)

- **1K15001:**
 - . A threaded extension piece in light alloy (1) with 4 headless hexagonal socket screws M8x16 (2)
 - . 2 headless hexagonal socket screws M8 x 20 (not shown) [See use §F => Fig. 5a and 5b]
 - . This leaflet of instructions for use.
- **1K155:**
 - . A threaded extension piece in steel (1) with 4 headless hexagonal socket screws M8x16 (2)
 - . This leaflet of instructions for use.
- **1K183:**
 - . A rotating pyramid receiver M36x1.5 in steel (1) with 4 headless hexagonal socket screws M8x16 (2)
 - . This leaflet of instructions for use.

B) USE


These components are designed exclusively for lower limb prosthesis. They are appropriate for patients weighing (carried load included) less than :

- 100 kg for **1K15001**
- 150 kg for **1K155** and **1K183**

C) CONNECTIONS (Fig.2)


The **1K15001** and **1K155** adjustable threaded extension pieces and the **1K183** female pyramid connector can be assembled :

- at their threaded end **1a**, to any tapped component (A) M36x1.5
- at their end **1b**, to any male pyramid component (B).

D) WHICH LENGTH SHOULD BE CUT WHEN USING 1K15002 OR 1K184 (Fig.3) (**1K183 : not applicable**)

- Determine the build height of the product according to its location in the prosthesis (this dimension will be called Y).
- Screw the threaded extension piece **1K15001** or **1K155** into the tapped sleeve **1K15002** or **1K184** until you reach the threaded root.
- Measure dimension X.
- The length to be cut will be calculated as follows : **Length to be cut = X-Y.**


E) CUTTING TO LENGTH THE THREADED EXTENSION PIECE (Fig.4) (**1K183 : not applicable**)

Using the 5C413 cutting tool is highly recommended to cut the threaded extension piece to length


Caution : The cutting tool is 12 mm thick. This thickness is the minimal length of the re-cut threaded piece, NEVER re-cut the threaded piece shorter than 12 mm because the component would get fragile.

- Screw the threaded extension piece (1) into the **5C413** cutting tool so that the extension piece protrudes from side **F** by the length that should be cut.
- Install the assembled threaded extension piece and cutting tool into a vice and firmly tighten the jaws to lock the threaded extension piece.
- Use a hacksaw to cut the threaded extension piece (1) along side **F** of the **5C413** cutting tool.
- Loosen the jaws of the vice and remove the assembled threaded extension piece and cutting tool from the vice.
- Unscrew the threaded extension piece and remove it from the cutting tool.

F) ANGULATION ADJUSTMENT (Fig.5)

- a) Adjust angulation with the 4 screws (2) and a 4 mm hex wrench.


If one of the screws is too short, replace it by one of the 20 mm long screws included in the pack (Fig. 5a and 5b)

- b) Secure screwing by applying Proteor XC047 adhesive on these screws.
- c) Tighten the screws up to 15Nm.

G) RECYCLING

- **1K15001:** The threaded extension piece is in aluminium alloy and the screws in galvanized steel. Each of these components should be recycled according to the laws in force.
- **1K155** and **1K183:** The threaded extension piece and the rotating pyramid receiver are in stainless steel and the screws in galvanized steel. Each of these components should be recycled according to the laws in force.

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