

Class II Pre-Operative View

Rubber dam isolation is always preferred. Groom contacts with Bioclear TruContact sander before rubber dam to expedite placement and avoid tearing the rubber dam.



Cavity Preparation and Beveling

The pre-wedge helps to protect the dam, soft tissue, and neighboring tooth.



Blast Tooth Inside and Out

The blaster has 3 purposes here, first to remove biofilm from the infinity edge area, second to remove the pellicle from enamel past the bevel, and third to cleanse the inside of the cavity.



Measure Cavity Depth with Go/No Go Probe

Most heated Filtek One Restoratives allow 5 mm depth of cure in the interproximal if the matrix is clear. The top of the yellow stripe is 5 mmand should be compared to the marginal ridge area of the matrix and neighboring tooth. A No Go for cavities deeper than 5mm require 2 layers.




Air Thin the Adhesive

The Scotchbond Universal Adhesive should be air thinned until it stops moving.



Disclose the Tooth to Identify Biofilm

Place Bioclear Dual Color Disclosing Solution with a disposable brush onto the dried tooth.



Finished Cavity Preparation

For most cavities, no mechanical retention via box shape is needed when the proper Clark Class II preparation is



Place Biofit Matrix and Largest Wedge

The Biofit HD matrix is rugged and will withstand heavy seating pressure. Trap the "tab" against the neighboring marginal ridge during insertion to keep matrix fully seated and stable. Tab should be centered as shown.



Place Acid Etch

Make sure to inject acid etch all the way to the edge of the matrix which is the infinity edge area. Etch for 20 seconds, rinse and dry.



diameter, go inside, over, and outside the corners of the matrix.



Pre-Wedge with Diamond Wedge

Wedge can be placed from buccal or lingual direction. If possible the wider end near the handle should be placed into the widest, deepest embrasure.



Run Sanding Strip Through Contact

The contact area must be lightly broken. This step is much safer using a 3M equivalent lightening strip, not a bur, to protect the neighboring tooth.



Place TwinRing

The TwinRing has four "hands" that will straddle the wedge. Once the TwinRing is placed, gently rock the TwinRing apically, like a rubber dam clamp, to assure full seating.



Place and Massage Adhesive Into Dentin

Massage Scotchbond Universal Adhesive into dentin tubules according to manufacturer's recommendations.




Place Spot Welds of Flowable Composite

Two Spot Welds of heated Filtek Flowable composite, 1-2 mm in



BIOCLEARMATRIX.COM

Refer to manufacturer instructions for use (IFU) for complete product information.



Place Push/Pull, Quick Light Cure

The Push/Pull instrument is placed inside of the cavity prep, with pressure mid coronal-apical, to approximate the matrix for a broad, tight contact. The assistant light cures while the doctor holds the Push/Pull Instrument.



Air Thin the Adhesive

Gently air thin the Scotchbond Universal Adhesive to drive off solvents and to fully wet the tooth. - Do Not Light Cure Yet



Quickly Shape the Composite

A central groove should be placed and composite massaged toward the cusp tips. Do not over-manipulate or marginate the composite. Leave excess.




Remove Wedge and Matrix

Remove wedge, then release matrix from tooth with explorer. Then spin the matrix out with hemostat.



Occlusal Margination

A rubber abrasive point such as a "brownie" is used to marginate the occlusal infinity edge Tooth Restoration Interface (TRI).




Cure Again After Removing Push/Pull

Do a full light cure of cavity area once the Push/Pull Instrument has been removed, to ensure unrestricted light penetration.



Inject Heated Flowable Composite

The heated Filtek Flowable composite is injected slowly, at 90° to the matrix tooth interface, displacing the adhesive/surfactant. - Do Not Light Cure Yet



Light Cure from the Occlusal

This is the first step in 3-point curing, generally for 10 seconds.



Sculpt Restoration with Disc

A Sof-Lex XT Coarse Disc is used to round the marginal ridge and to marginate. Keep disc out of deep interproximal areas, maintaining the smooth mylar finish.



Magic Mix Pre-Polish


Bioclear Magic Mix is applied in a disposable cup for further margination and pre-polish.



2nd Application of Adhesive

The cavity area and infinity edge areas are wetted with Scotchbond Universal Adhesive which acts as a surfactant allowing good flow of composite in the next steps.

- Do Not Light Cure Yet



Inject Heated Bulk Fill Paste Composite

The heavier paste will displace most of the heated Filtek Flowable composite leaving a strong, seamless restoration.
- Do Not Light Cure Yet



Remove TwinRing, Light Cure from Buccal, then Lingual

Steps 2 and 3 of 3 point curing are to cure through the clear matrix from buccal and then from lingual, generally for 10 sec each



Grind Simple Anatomy

A diamond or carbide bur is used to create simple fossa near contact area and central groove.



Rock Star Polish

The Bioclear RS Polisher is used with light pressure in the dry phase, and then heavy pressure in the wet phase. Assistant blows gentle air stream during dry phase, and constant full air-water during the wet phase as the doctor applies heavy pressure to the cup.