





# COMPOSITE POLISHING GUIDE

Cosmedent gives you everything you need to polish composites. With our conveniently color-coded easy to use polishers, you can create a smooth and esthetic polish that mimics enamel long-term.

Cosmedent Polishing Tool	Speed Range	Use Notes
<b>Coarse FlexiDisc® Mini</b> 	3,000 to 15,000 RPM for 3/8" disc 5,000 to 15,000 RPM for all other sizes	<ul style="list-style-type: none"> <li>• Use light pressure</li> <li>• Gross reduction</li> <li>• Initial contouring</li> <li>• Buffing out imperfections</li> </ul>
<b>Medium FlexiDisc® Mini</b> 	3,000 to 10,000 RPM around gingival margin with 3/8" disc 10,000 to 15,000 RPM for all other sizes	<ul style="list-style-type: none"> <li>• Completion of finish</li> <li>• Removal of scratches</li> <li>• Use flex and pressure</li> </ul>
<b>Fine FlexiDisc® Mini</b> 	10,000 to 15,000 RPM for all sizes	<ul style="list-style-type: none"> <li>• Beginning final polish</li> <li>• Use flex and pressure</li> </ul>
<b>Superfine FlexiDisc® Mini</b> 	10,000 to 15,000 RPM for all sizes	<ul style="list-style-type: none"> <li>• Enhancing final polish</li> <li>• Use flex and pressure</li> </ul>
<b>Medium FlexiCups® &amp; FlexiPoints®</b> 	2,000 to 10,000 RPM MAX with Light Pressure for all cups and points	<ul style="list-style-type: none"> <li>• Use with water</li> <li>• Achieve labial characterization</li> </ul>
<b>Superfine FlexiCups® &amp; FlexiPoints®</b> 	2,000 to 10,000 RPM MAX with Light Pressure for all cups and points	<ul style="list-style-type: none"> <li>• Use with water</li> <li>• Polish gingival margin and labial surfaces</li> </ul>
<b>FlexiBuff® Mini</b> 	15,000 RPM for all sizes	<ul style="list-style-type: none"> <li>• Add Enamelize® to tooth first and use light pressure</li> <li>• Available sizes: 3/8", 1/2", &amp; 5/8"</li> </ul>
<b>Enamelize™ Aluminum Oxide Polishing Paste</b> 	15,000 RPM with Light Pressure	<ul style="list-style-type: none"> <li>• Creates a beautiful final polish on composites</li> <li>• Retains surface texture</li> <li>• Can also be used on porcelain, metal and natural tooth structure</li> <li>• Available sizes: 36 g Tube &amp; 3 Syringe</li> </ul>

5880\_4\_2022.09.16