EN

**COSMEDENT®** www.cosmedent.com | 800.621.6729  
**Renamel® Microfill** | LIGHT-CURED MICROFILL COMPOSITE

## INSTRUCTIONS FOR USE

## GENERAL CHARACTERISTICS

Renamel® Microfill is the only composite that truly mimics the enamel surface. Where enamel is to be replaced, Renamel® Microfill is the composite resin of choice. Renamel® Microfill has been and continues to be the number one rated microfill for more than two decades. When applied, sculpted, finished and polished accurately, no other composite resin can match the tooth's natural appearance. In essence, it is the artistic dentist's dream material.

## DESCRIPTION

Renamel® Microfill is a light-cured, microfilled, highly polishable dental filling material, according to EN ISO 4049 (Type I/Class 2/Group 1). Renamel® Microfill cures with an LED or Halogen light with a light intensity more than 500mW/cm<sup>2</sup> at wave- length range of 400 – 500 nm. Renamel® Microfill has a compressive strength of >320 MPa, flexural strength of > 80 MPa and 2 mm depth of cure.

## CREATING AN ENAMEL-LIKE SURFACE

Renamel® Microfill is a pure microfill, and as such can be polished as highly as glass. Using the appropriate instruments, a dentist can shape and texture the surface to exactly duplicate the surrounding dentition. This step should not be taken lightly, for this is the part that gives you the ultimate in esthetic dentistry. Renamel® Microfill has the best refractive and reflective index of any product in the marketplace. When finished and polished properly, Renamel® Microfill will create results that are beautiful and long-lasting, that both the dentist and the patient will be astounded.

## COMPOSITION

Renamel® Microfill contains UDMA and 1,4-Butanediol dimethacrylate (28 – 30%). The fillers contain silicon dioxide and prepolymerized composite (70% by weight and 60% by volume). Initiators, stabilizers and pigments are additional ingredients (< 1%). The particle size of the fillers range from 0.04 – 0.2 microm.

Renamel® Microfill has been in the marketplace for more than twenty-five years, and its handling, color stability, and long term polishability is unsurpassed.

## LONG-TERM COLOR STABILITY

Renamel® Microfill offers color stability that is unmatched in the marketplace. Renamel® Microfill holds its shade and polish so that you can rest assured that your patient will be happy for many years.

## INDICATIONS FOR USE

- Direct anterior veneers
- Top layer of Class I and II restorations
- Anterior Class III and IV restorations
- Class V restorations (cervical caries, root erosion, wedge-shaped defects)
- Closing of diastemas
- Repair of composite and ceramic veneers
- Reshape, lengthen and restoratively reposition anterior teeth

## CONTRAINDICATIONS

Renamel® Microfill is contraindicated if patient is known to be allergic to any of the ingredients of Renamel® Microfill.

## SIDE EFFECTS

Can cause irritation in sensitive patients.

**INTERACTIONS**  
Do not use Renamel® Microfill in combination with eugenol containing materials, as these can impede the polymerization process.

**NOTICE:**  
During time-consuming procedures, avoid overexposure from the operator's light to protect against premature polymerization of the material.

The light curing unit must be maintained and calibrated according to the manufacturer's specifications. Serious incidents must be reported to DeltaMed GmbH, the responsible importer and the responsible authority.

**WARNING**  
Renamel® Microfill contains 1,4-Butanediol dimethacrylate and Urethane dimethacrylate.

May cause an allergic skin reaction. Wear protective gloves. In case of contact with the eyes, rinse immediately with plenty of water and seek medical attention.

## DISPENSING COMPOSITE MATERIAL

- Extrude composite onto a mixing pad by turning the handle slowly in a clockwise manner. Slice an increment of proper size with a sterile instrument onto a mixing pad.
- Immediately replace syringe cap to prevent polymerization. If the extruded material is not used right away, cover the composite with a light protecting box such as our ResinKeeper™.

**Unit dose delivery**  
- Insert unit dose into dispensing gun and dispense composite material. Dispense of unit dose after use.

## PLACEMENT DIRECTIONS

Prior to the application of Renamel® Microfill, the use of opaques or microhybrids, or a combination of both may be necessary to achieve excessive color block. If replacing dentin, place a micro- hybrid or nanofill composite first for strength and opacity. Polymerize before the addition of your microfill layer.

## SHADE SELECTION

When selecting a shade, one should remember that teeth are usually polychromatic, having three zones: a translucent darker incisal, a lighter central body portion and a darker cervical area.

Assessing shade on a clean, moist tooth using magnification will help you to easily identify these areas.

- ① Prepare tooth, etch and dry according to manufacturer's instructions.
- ② Place dental adhesive of choice, such as Cosmedent's Complete™, G-3 or Rite-Lite 2™.
- ③ Renamel® Microfill should be placed in a light-protecting box, such as the ResinKeeper, to prevent prepolymerization by ambient light. Slice an increment of proper size and roll the composite into a ball, a cylinder, or an elliptical shape (convenience form) as needed.

- ④ Place with a clean metal instrument onto the area to be restored and shaped. Cosmedent's 8A and BAL (long and short bladed) and G-3 are excellent instruments for placement.
- ⑤ If resin still sticks to the instrument, the instrument may be scratched or dirty. If your room is excessively warm and some stickiness occurs, simply wipe the instrument with a light alcohol swab and this will reduce stickiness.

- ⑥ Take your time and sculpt carefully. Slightly overcontour every restoration, except in the interproximal areas. Sculpt exactly to contour in these areas and it will greatly aid in developing the final contour of your restoration and shorten your finishing and polishing time. Cosmedent's IPC instruments are excellent for interproximal sculpting.
- ⑦ Make sure your patient removes all makeup prior to treatment.
- ⑧ When applying composite, place convenience form in one piece, then sculpt off excess. This is much easier than placing too little and then adding more.

- ⑨ The use of Cosmedent's Multi-Purpose instrument is very effective for incisal groove and lobe placement, which allows enough room for addition of the incisal material.
- ⑩ Use a Caliper instrument for measurements.
- ⑪ To avoid contamination, wash and dry your gloves thoroughly prior to placement. Do not use powdered gloves.
- ⑫ Make sure your patient removes all makeup prior to treatment.

- ⑬ If placing an incisal shade, the use of a titanium coated instrument, such as Cosmedent's Multi-Purpose instrument, will help you create incisal lobes and grooves in your body composite to make room for the application of your incisal shade.

## HYGIENE / DISINFECTION

- Do not reuse a unit dose capsule.

- Disinfection of the syringe and the unit dose dispenser can be carried out using a commercially available disinfecting solution.

**Note:** Dentists tend to under-cure composites. You will get best results if you slightly over-cure.

- ⑭ Finishing should be done with a combination of polishing materials such as diamonds, carbide burs and aluminum oxide discs in coarse and medium grits (Mini FlexiDiscs® from Cosmedent).

- ⑮ High-gloss polish can be achieved using aluminum oxide impregnated instruments, such as fine and superfine flexible discs (Mini FlexiDiscs) and (FlexiCups), followed by Enamelize® polishing paste applied with Mini FlexiBuff®s, all from Cosmedent. For more detailed polishing techniques, please visit Cosmedent's Learning Center at [Cosmedent.com/Learning\\_Center](http://Cosmedent.com/Learning_Center) or call a Cosmedent Customer Care Representative at 800.621.6729.

## COLOR MATCHING

Is there anything more difficult than matching a tooth shade against an existing dentition? The Renamel® Microfill System™ is designed to take the guesswork out of shade selection. It recognizes that there are different zones of color and translucency. It understands that color selection should be about achieving the correct enamel shade with predictability.

Cosmedent's Renamel® Microfill System will help you easily achieve consistently reliable results when selecting a shade.

**Color matching is a difficult science, and many people have trouble with color selection.**

**What are the problems?**  
- Many materials on the market do not match their shade guides.  
- It is often difficult for an individual to see shade differences.  
- The color, opacity, translucency and polishability of many materials do not accurately simulate the tooth surface.

- Materials often have a refractive and reflective index that is dissimilar to that of the enamel surface.  
- Many materials do not have long-term color stability which can result in large color changes over time.

Cosmedent has consistently developed reliable and color stable materials that accurately match each other and the Classic VITA® Shade Guide.

Renamel® Microfill has been matched to VITA shades more accurately than any other material in the marketplace. The 16 different shades in the VITA shade guide encompass four shade ranges.

Cosmedent has added several additional useful shades to give the clinician still more options to develop the correct color. These shades include some darker geriatric shades for cervical restorations (A5, A6, and C5) and SuperBrite® shades (SB1, SB2, and SB3) to match intensely bleached teeth, or to create the Hollywood smile. This system is also complemented by unique half shades of A1, A2, A3, and A3.5, which are commonly missing in other systems. The Renamel® Microfill System also contains three vividly translucent incisal shades: light, medium and dark.

## SHADE SELECTION

When selecting a shade, one should remember that teeth are usually polychromatic, having three zones: a translucent darker incisal, a lighter central body portion and a darker cervical area.

Assessing shade on a clean, moist tooth using magnification will help you to easily identify these areas.

The body shades in The Renamel Restorative System correspond to the central body zone of the tooth surface. Primary shade selection should be taken against the body area of the tooth. A good suggestion is to remove the incisal and gingival portion on each shade tab which may help you achieve a more accurate match to the body area. You can also make your own shade guide out of composite material.

## Some additional hints on shade selection are:

- Understand how light affects the overall aspect of your shade selection. Assess shade in fluorescent light (where we view one another most of the time) and augment with other types of light, i.e., shade matching light such as Rite-Lite 2™ (from Cosmedent) or viewing in natural light.
- Your materials should potentially match an existing shade guide. If not, make your own shade guide.

- You can include the patient in shade changing, but never include them in shade matching.

## USEFUL TIPS FOR PLACEMENT

- ① Use of a clean, dry, gloved finger is very effective in rapidly forming the desired surface prior to your sculpting techniques.

- ② Place with a clean metal instrument onto the area to be restored and shaped. Cosmedent's 8A and BAL (long and short bladed) and G-3 are excellent instruments for placement.

- ③ It is convenient to place Renamel® Microfill in a light-protecting box, such as the ResinKeeper, to prevent prepolymerization by ambient light. Slice an increment of proper size and roll the composite into a ball, a cylinder, or an elliptical shape (convenience form) as needed.

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