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INSTRUCTIONS FOR USE

DESCRIPTION

Renamel® Posterior is a light cured, hybrid filled, radiopaque, dental filling material according to EN ISO 4049 (Type 1/Class 2/ Group 1). The Renamel Posterior System consists of opaque Body shades. It cures with a halogen light with a light intensity more than 500mW/cm² at wavelength range of 400 – 500 nm.

COMPOSITION

Renamel Posterior contains UDMA, Bis-GMA and 1,4-Butanediol dimethacrylate (22 %). The fillers contain strontium glass and silicon dioxide (77% by weight and 60% by volume). Initiators, stabilizers and pigments are additional ingredients (< 1%). The particle size of the fillers range from 0.04 – 5 micron. Renamel Posterior has a flexural strength of 120 MPa, radiopacity of > 1.0 mm Al and 1.5 mm depth of cure.

INDICATIONS FOR USE

(Areas of application)
- Restorations in the posterior region (Class I and II)

PRETREATMENT

Before treatment, tooth color should be taken using a standard color guide. Keep the teeth moist when taking the shade.

DISPENSING COMPOSITE MATERIAL

Syringe delivery:
- Extrude composite onto a mixing pad by turning the handle slowly in a clockwise manner. Slice an increment of proper size with a sterile instrument onto a mixing pad.
- Immediately replace syringe cap to prevent polymerization. If the extruded material is not used right away, cover the composite with a light protecting box such as our ResinKeeper™.

Unit dose delivery:
- Insert unit dose into dispensing gun and dispense composite material. Dispose of unit dose after use.

POSTERIOR PLACEMENT: CLASS I

- 1 Prepare the cavity according to general standards in composite filling therapy.
2 Recommendation for deep cavities: Cover very deep areas close to the pulp with a calcium hydroxide material and subsequently use a pressure-resistant cement (e.g. glass ionomer cement). Do not cover other cavity walls, since they can be used to support the bond with an enamel / dentin adhesive.
3 After preparation, cleanse, etch, wash, dry and place bonding adhesive of choice, such as Cosmedent's Complete™ for excellent bond strengths and marginal seal.
4 Place Renamel Posterior, adding enough composite to cover the entire base of the preparation to a depth about 1 mm and not exceeding 2 mm and polymerize each surface for 20 – 40 seconds (Halogen or LED Light).
5 Continue to buildup in incremental layers and slightly over-contour the occlusal. The use of Cosmedent's condensing and contouring instruments and round ball burnishers will work well to gently contour and move the material to where you want it. Do not over-manipulate material because this could cause voids and later pitting of the finished restoration.
6 Final cure each restoration on all surfaces for 60 seconds.
7 Finishing and polishing should be done with a combination of polishing materials such as diamonds, carbide burs, aluminum oxide discs (Mini FlexiDiscs® from Cosmedent), finishing cups and points, (FlexiCups®) and (FlexiPoints® from Cosmedent) and polishing paste (Enamelize® from Cosmedent).

POSTERIOR PLACEMENT: CLASS II

- 1 Prepare the cavity according to general standards in composite filling therapy.

- 2 Recommendation for deep cavities: Cover very deep areas close to the pulp with a calcium hydroxide material and subsequently use a pressure-resistant cement (e.g. glass ionomer cement). Do not cover other cavity walls, since they can be used to support the bond with an enamel / dentin adhesive.
3 Following cavity preparation, cleanse, etch, wash and dry and place a matrix band. You need an ultra-thin, yet highly resilient, matrix band such as Cosmedent's Ultra-Thin Precontoured Matrix Bands™, or a matrix system of your choice.
4 After placement, wedge thoroughly and make sure the gingival margin is completely sealed. Often, it is necessary to wedge from both the buccal and lingual to obtain complete closure.

- 5 Burnish the band using a round ball burnisher and the convex side of a spoon excavator in order to create a contact area and not just a contact point. Use heavy pressure and observe the burnished area closely. You should easily be able to see the contact area you have created.

- 6 Rinse the preparation and dry thoroughly and apply base if so desired. Then, etch and apply the dental adhesive of your choice. Cosmedent's Complete dental adhesive system is an excellent choice.

- 7 Apply a flowable cavity liner such as Cosmedent's hyper radiopaque Renamel De-Mark™ at the gingival seat of posterior restorations to differentiate your restoration from tooth structure. Using a nonstick composite placement instrument apply a small amount of material and gently tap to place with a small round ball burnisher and cover the gingival seat only. Light-cure for 40 – 60 seconds. Note: Cosmedent's posterior titanium-coated instruments are excellent for placing and developing contour in posterior restorations.

- 8 Place Renamel Posterior in desired shade. Remember, do not try to condense. This is not amalgam. Make sure all aspects of the gingival seat are covered and all margins are completely sealed. You can use an IPC Carver to contour the proximal areas to reduce flashing. Once the first layer has been placed, polymerize for 20 – 40 seconds (Halogen or LED Light).

- 9 Apply and sculpt further material as needed, and cure. The rule of thumb is not to exceed 2 mm of material at a time. You should polymerize each incremental layer for 20 – 40 seconds. This will give optimal strength and hardness to the final restoration. Polymerize the final layer for 60 seconds. Note: To prevent the "Cuspal-flexion" in large restorations, it might be desirable to develop one cusp at a time.

- 10 Finishing and polishing should be done with a combination of polishing materials such as diamonds, carbide burs, aluminum oxide discs (Mini FlexiDiscs from Cosmedent), finishing cups and points, (FlexiCups from Cosmedent) and (FlexiPoints from Cosmedent) and polishing paste (Enamelize® from Cosmedent).

HYGIENE / DISINFECTION

- Do not reuse a unit dose compule.
- Disinfection of the syringe and the unit dose dispenser can be carried out using a commercially available disinfecting solution.
- All instruments used should be thoroughly washed and autoclaved.

STORAGE

- Store at 39 – 75°F (4 – 24°C).
- Close syringe immediately after usage. Exposure to light causes premature polymerization.
- Use at room temperature.
- Do not use after expiration date.
- Keep away from sunlight.
- Keep material out of reach of children.

IMPORTANT PRODUCT INFORMATION

Federal (U.S.) law restricts this device to sale by or on the order of a dentist.

Note: This product was developed for use in the dental sector and should be used in accordance with the instructions for use. The manufacturer does not accept liability for damage caused by its use for any other purpose. The user is also obliged to accept responsibility for checking before use that the material is suitable and can be used for the intended purpose, particularly if this is not listed in the instructions for use.

CONTRAINDICATIONS

Placement of Renamel Posterior restorations is contraindicated if a patient is known to be allergic to any of the ingredients of Renamel Posterior.

SIDE EFFECTS

May cause irritation in sensitive patients.

INTERACTIONS

Do not use in combination with eugenol containing materials as these may impede the polymerization process.

NOTICE

During time-consuming procedures, avoid overexposure from the operatory light to protect against premature polymerization of the material. The light curing unit must be maintained and calibrated according to the manufacturer's specifications. Serious incidents must be reported to DeltaMed GmbH, the responsible importer and the responsible authority.

WARNING

Renamel Posterior contains 1,4 Butanediol dimethacrylate and Urethane dimethacrylate. May cause an allergic skin reaction. Wear protective gloves. In case of contact with the eyes, rinse immediately with plenty of water and seek medical attention.

DISPOSAL INSTRUCTIONS

Dispose of contents/ container in accordance with local regulation.

PACKAGING

Renamel Posterior is available in 4.5 gram syringes or 10 unit dose compules of 0.25 g, each.

SHADES

Renamel Posterior is available in the following shades: Body shades: A1, A2, A3, B-Zero.

Canada Site Licence Holder: CLINICAL RESEARCH DENTAL 167 Central Ave, Suite 200 London, Ontario, Canada N6A 1M6 Establishment License No. 666

Renamel Posterior is available in the following shades: Body shades: A1, A2, A3, B-Zero.

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Cosmedent® Inc.