CLEARFIL™ PHOTO CORE

ENGLISH INSTRUCTIONS FOR USE

I. INTRODUCTION

CLEARFIL PHOTO CORE is a light-cure, radiopaque core build-up composite material.

II. INDICATION

Core build-up of vital or non-vital tooth.

III. CONTRAINDICATION

Patients with a history of hypersensitivity to methacrylate monomers

IV. INCOMPATIBILITY

Do not use eugenol containing materials for pulp protection or temporary sealing, since the eugenol could retard the bonding system curing process.

V. PRECAUTIONS

1. Safety precautions

- Avoid use of the product for patients with a history of hypersensitivity to methacrylate monomers.
- If any hypersensitivity occurs, such as dermatitis, discontinue use of the product and consult a physician.
- Wear gloves or take other appropriate protective measures to prevent the occurrence of hypersensitivity that may result from contact with methacrylate monomers.
- Use caution to prevent the product from coming in contact with the skin or getting into the eye. Before using the product, cover the patient's eyes with a towel to protect the patient's eyes from splashing material.
- If the product comes in contact with human body tissues, take the following actions;
 - <If the product gets in the eye> Immediately wash the eye with copious amounts of water and consult a physician.
 - <If the product comes in contact with the skin> Immediately wipe it off with a cotton pledget moistened with alcohol or gauze and wash with copious amounts of water.
- Use caution to prevent the patient from accidentally swallowing the product.
- Avoid looking directly at the dental curing unit when curing the product.

2. Handling and manipulation precautions

- Do not use the product in conjunction with the other composite resins. Mixing materials may cause a change in physical properties, possibly a decrease, from the properties expected.
- When light curing the product, note the depth of cure in this instructions for use.
- The product polymerize under an operating light or natural light (sunlight from windows). Avoid operating light or natural light, and use within 5 minutes from dispensing.
- The emitting tip of the dental curing unit should be held as near and vertical to the resin surface as possible. If a large resin surface is to be light cured, it is advisable to divide the area into several sections and light-cure each section separately.
- Low intensity of light causes poor adhesion. Check the lamp for service life and the dental curing unit guide tip for contamination. It is advisable to check the dental curing light intensity using an appropriate light evaluating device at appropriate intervals.
- The use of the product is restricted to a licensed dental professional.

3. Storage precautions

- Do not use after the expiration date. Note the expiration date on the outside of package.
- Keep away from extreme heat or direct sunlight.
- The product must be stored at 2-25°C/36-77°F when not in use.
- Replace the cap as soon as possible after the composite has been dispensed from the syringe.
- The product must be stored in proper places where only dental practitioners can access it.

VI. COMPONENTS

1) Components

Please see the outside of the package for contents and quantity.

2) Ingredients

- Principal ingredients:
- Silanated silica filler
 - Silanated barium glass filler
 - Bisphenol A diglycidylmethacrylate (Bis-GMA)
 - Triethyleneglycol dimethacrylate
 - dl-Camphorquinone

The total amount of inorganic filler is approx. 68 vol%.
The particle size of inorganic fillers ranges from 0.49 to 75 µm.

VII. CLINICAL PROCEDURES

A. CORE BUILD-UP OF VITAL TOOTH

A-1. Cleaning tooth structure

Be sure the cavity is adequately cleaned. An adequately cleaned cavity assures maximum adhesive performance.

A-2. Moisture control

Avoid contamination of the treatment area from saliva or blood to produce optimal results. A rubber dam is recommended to keep the tooth clean and dry.

A-3. Cavity preparations

Remove any infected dentin and prepare the cavity in the usual manner.

A-4. Pulp protection

Any actual or near pulp exposure could be covered with a hard setting calcium hydroxide material. There is no need for cement lining or basing. Do not use eugenol materials for pulp protection.

A-5. Applying a matrix strip and wedges

Apply a matrix strip and wedges as needed.

A-6. Etching

Apply etching agent (e.g. K-ETCHANT GEL) into the entire cavity surface with a sponge or disposable brush. After about 30 to 60 seconds, wash the etching agent away completely with water and dry with an air flow.
Take care to prevent saliva from coming in contact with the etched tooth surface. Any contaminated tooth surface must be re-etched.

A-7. Applying the bonding agent

Mix one drop each of CLEARFIL NEW BOND Catalyst liquid and Universal liquid, and apply the mixture to the entire cavity surface with a sponge or disposable brush. Evaporate the ethanol of the mixture with a gentle air flow. (For detailed instructions, see the CLEARFIL NEW BOND instructions for use.)

A-8. Preparation of CLEARFIL PHOTO CORE

Dispense the necessary amount of the composite onto the paper pad by rotating the paste syringe plunger. After dispensing the composite, turn the plunger counter-clockwise a half-turn to prevent excess composite from oozing out. Replace the syringe cap immediately to prevent premature setting of the composite.

A-9. Placing CLEARFIL PHOTO CORE into the cavity

Place the composite into the cavity preparation using hand instruments. It is recommended that transparent matrix strips be used to help form the core.

A-10. Curing

Cure the composite with a dental curing unit. Hold the emitting tip as close to the composite as possible. The light curing time and depth of cure for each dental curing unit is used is shown below.

Light curing time and depth of cure

depth of cure	Dental curing unit	Light curing time			
		5 sec.	10 sec.	20 sec.	40 sec.
	Conventional halogen*	—	—	3.5mm	4.0mm
	Fast halogen*	4.0mm	4.5mm	—	—
	Plasma arc*	5.0mm	6.0mm	—	—
	LED*	—	—	4.0mm	4.5mm

*Dental curing unit

Type	Light source	Wavelength range and light intensity
Conventional halogen	Halogen lamp	Light intensity ¹⁾ of 300 - 550 mW/cm ² in wavelength range from 400 - 515 nm
Fast halogen	Halogen lamp	Light intensity ²⁾ of more than 550 mW/cm ² in wavelength range from 400 - 515 nm
Plasma arc	Xenon lamp	Light intensity ²⁾ of more than 2000 mW/cm ² in wavelength range from 400 - 515 nm, and light intensity ²⁾ of more than 450 mW/cm ² in wavelength range from 400 - 430 nm.
LED	Blue LED ¹⁾	Light intensity ²⁾ of more than 300 mW/cm ² in wavelength range from 400 - 515 nm

¹⁾ Peak of emission spectrum: 450 - 480 nm

²⁾ Evaluated according to ISO 10890-1.

³⁾ Wavelength distribution and light intensity values measured with a spectro-radiometer calibrated using an IEC or the NIST (National Institute of Standards and Technology) standard lamp

A-11. Preparation of the abutment tooth

After curing the composite resin, prepare the abutment tooth.

B. CORE BUILD-UP OF NON-VITAL TOOTH

B-1. Cleaning tooth structure

Be sure the cavity is adequately cleaned. An adequately cleaned cavity assures maximum adhesive performance.

B-2. Moisture control

Avoid contamination of the treatment area from saliva or blood to produce optimal results. A rubber dam is recommended to keep the tooth clean and dry.

B-3. Cleaning cavity and preparing root canal

Remove the temporary sealing material from the cavity and filling material from the root canal. Using a Pizo reamer, prepare and clean the root canal opening.

B-4. Try-in of the post

A post of the proper diameter is fitted to the prepared root canal and the length of the post adjusted. Sand blast the post as needed.

B-5. Either B-5a or B-5b

B-5a. When PANAVIA F 2.0 is used.

Cement the post into the root canal according to the instructions for Use of PANAVIA F 2.0.

[NOTE]

Do not use CLEARFIL SE BOND when PANAVIA F 2.0 is used. Since PANAVIA F 2.0 plays the role as adhesive bonding agent, apply excess PANAVIA F 2.0 paste to the remaining tooth structure on which ED PRIMER II, a component of PANAVIA F 2.0 is applied and the surface of the post.

B-5b. When another resin luting cement is used.

Cement the post into the root canal according to Instructions for Use of the resin luting cement.

B-6. Preparation of CLEARFIL PHOTO CORE

Dispense the necessary amount of the composite onto the paper pad by rotating the paste syringe plunger. After dispensing the composite, turn the plunger counter-clockwise a half-turn to prevent excess composite from oozing out. Replace the syringe cap immediately to prevent premature setting of the composite.

B-7. Build-up of CLEARFIL PHOTO CORE

Using appropriate instruments, build-up the composite to the appropriate anatomy being careful to avoid incorporation of air bubbles. When building-up an anterior tooth core, after the matrix strip or wedges are placed, the composite material is first placed around the post. Then, to avoid air bubbles, place additional increments of the composite while holding the matrix strip or wedges with finger pressure.

B-8.Curing

Refer to Light curing time and depth of cure in A-10

B-9. Preparation an abutment tooth

After curing the composite resin, prepare the abutment tooth.

[CAUTION]
Federal (U.S.A.) law restricts this device to sale by or on the order of a licensed dentist.

[WARRANTY]

Kuraray Noritake Dental Inc. will replace any product that is proved to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

[NOTE]

CLEARFIL, CLEARFIL PHOTO CORE and PANAVIA are trademarks of KURARAY CO.,LTD.

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