

Pet Microchip Reader

User Guide

Getting to know the reader.

This item is a wireless handheld microchip reader that can read ISO FDX-B coded radio frequency identification tags. It has very simple controls and a high-brightness OLED display which allows you to read the numbers even in bright sunlight.

With built-in memory, it is able to store up to 30 ID tag codes, which can be downloaded to a computer via a USB cable connection (Note: There must be a text field for the reader to input character in when downloading stored codes). Optimum performance delivers great animal traceability management.

Specifications

Power	Built-in Rechargeable Battery (Charge it for 2 hours before the first use)
Battery Life	4 hours (Work Continuously)
Operating Frequency	134.2KHz
Tag Compatibility	ISO FDX-B(ISO11784/85)
Interface	USB 2.0
Display	OLED
Memory	30 ID Numbers
Indicator	Audio Beeps and OLED
Read Distance	10-15cm (Glass Tag, vary on Chip Sizes & Types) 20-22cm (Ear Tag, vary on Chip Sizes and Types)
Language	English

Accessory	USB Cable *1
Operating Temperature	-10°C to 55°C
Storage Temperature	-30°C to 65°C
Relative Humidity	5% to 90% non-condensing
Package Dimensions	21.3cm*14.9cm*2.9cm
Weight	90g

Instructions:

1. Hold down the button for 3s to power up the reader and the following sign appears on the screen.
2. Hold the reader over the chip location. If the reader fails to detect any tags, it displays No Tag on the screen.
(Note: The reader will stop reading if no tags are detected within 20s)
3. If the reader detects a tag, it displays the number on the screen.
4. Press the button once to continue scanning tags.
5. To check scanned data, press the button twice while the screen shows the read sign.
6. If there are no further inputs, the scanned data will remain displayed for 5s and the reader will turn off in 60s.