### 5. ABOUT US

Driven by the passion for innovation, we at Dr Trust endeavour to provide our customers with the latest medical inventions with an objective to promote good health and wellness all around the world. All the medical devices and health monitors provided by Dr Trust are supported by accurate, latest and ground breaking technologies, innovated at our headquarters in NY, USA. All our products adhere to the most stringent CE and FDA guidelines and are strongly recommended by doctors and health practitioners. Our products are designed in the utmost exemplary ways to ensure that their accuracy and convenience are unrivalled. The ease of their use and operation makes them even more suitable for users of all age groups.

DrTrust strives to enhance the quality of lifestyle by providing with the most trusted and innovative health care and wellness products. Being a renowned global leader in health care products, Dr Trust ensures that our technically efficient team works dynamically and tirelessly to provide the best of the medical devices to our clients. The products that we have to offer are suitably designed for use at homes, laboratories and hospitals.

Our ground breaking solutions allow you to monitor your health in the easiest ways possible. In today's era when all of our lives are too hassled to handle, it becomes a bit difficult to pay attention to our health. But it has now become easier with the coming of the monitoring devices which can be conveniently used at homes and even on the go.

We bring to you a variety of best self medical devices, trusted and used by Doctors, medical professionals and home users all over the world.

# Dr Trust® **Health and Fitness Tracker** with Heart Rate-8001

Thank you for your purchase. Please read instructions manual carefully to use the Dr Trust Health and Fitness Tracker with Heart Rate-8001 correctly and retain it for future reference.

USER

INSTRUCTIONS

TABLE OF CONTENT	Dr Trust <sup>®</sup>	1. INTRODUCT	TION .	Dr Trust <sup>®</sup>
1. INTRODUCTION	3	Are you looking for a gear that you can wear constantly in order to monitor your fitness progress? With Dr Trust Health & Fitness Tracker with Heart Rate		
2. TECHNICAL SPECIFICATIONS	3	8001, you can track your activities and body composition easily. It works with both Android and Apple phones by connecting through an easy-to-use app that lets you know about the taken steps, burnt calories, walked distance,		
3. CLEANING AND MAINTENANCE	4	activity time, distance, % of your goal achieved, body composition and lots more. It has a touchscreen that displays, time/date with your body data.		
4. CUSTOMER SUPPORT	5	2. TECHNICAL SPECIFICATIONS		
		Main body Size	42*19*11.2mm	
5. ABOUT US	6	Wristband	204.6*16.8mm	
		Weight	31.8g	
		Display	0.96 'TFT	
		Data Memory	30 days	
		Walking distance	Max 999.99km	
		Calories burned	Max 9999.9kcal	
		Activity time	99h59min	
			on the fitness tracker case. fitness tracker into dish-wa ying machine.	shing machine, washing
2	www.drtrustusa.com		3	www.drtrustusa.com

Dr Trust<sup>®</sup>

# 4. CUSTOMER SUPPORT

Dr Trust<sup>®</sup>

## CONTACT ADDRESS

NURECA INC.USA

276 5th Avenue, Suite 704-397, New York (NY) - 10001, USA

TNDTA

USA

Dr Trust<sup>®</sup>

Corporate Office (Mumbai)

Nureca Private Limited #103/104, Orbit Plaza, New Prabhadevi Marg, Prabhadevi, Mumbai, Maharashtra - 400025

Contact us India: +91-7527013265 / +91-9356658436

Website: www.drtrustusa.com Email: customercare@nureca.com

Connect with us on social networks Facebook: @drtrust

Instagram: @drtrustisin Youtube: Nurecallsa

COPYRIGHT© 2020 DR TRUST, ALL RIGHTS RESERVED

Don't use detergent to clean the fitness tracker.

Don't charge the fitness tracker before it is totally dry.

or risk situation, otherwise, it may cause distractions.

temperature.

stop charging.

bearing.

liquid for a long time.

3, CLEANING AND MAINTENANCE Regular cleaning of fitness tracker is recommended to prevent the buildup of sweat and other particles. Take down the fitness tracker to clean with a disinfectant solution or rinse with water. Just allow it to dry completely before

Don't expose the fitness tracker under extreme high or low

Don't wear the fitness tracker in sauna bathroom or steam bathroom.

Don't soak the fitness tracker in hot spring, the sea or other chemical

During charging, if the fitness tracker's temperature is too high, please

Don't check the information on the fitness tracker screen during driving

Don't put the fitness tracker into fire as the battery may explode.

www.drtrustusa.com

www.drtrustusa.com

5