



Application Instructions

1. Locate MKO label at top of brace
2. Grip both sides at top of brace by spiral stays.
3. Pull the support over the foot and up to the knee.
4. Position knee buttress comfortably over the knee.
5. Adjust where needed making sure the brace has a comfortable fit. The brace should not cut off circulation nor slide down the limb.
6. To remove brace pull from spiral stays at bottom over ankle until brace is off limb. To avoid damaging the brace do not turn inside out.

To prevent moisture build up, store brace in an open dry area between uses.