

Lacing Wrist Brace



Figure 1



Figure 2



Figure 3



Figure 4

**NOT MADE WITH
natural rubber latex**

Please read the following instructions carefully and completely before use. Correct application is required to ensure proper function of the device.

Intended for single patient use.



WARNING: If your pain increases or persists or you experience any unusual reactions while using this product, consult your healthcare provider immediately.

APPLICATION INSTRUCTIONS

1. Unfasten straps and slide affected hand into brace (Figure 1).
2. Position palm stay firmly against the palm of the hand (Figure 2). Ensure the edge of the brace is below the knuckles (Figure 3).
3. Secure the thumb strap.
4. Pull and wrap the lacing straps around wrist. Secure when the desired fit and level of compression is reached (Figure 3).
5. Adjust straps as needed when swelling of the hand, wrist or thumb changes (Figure 4).

NOTE: Please refer to sew-in label on device for care instructions.

NOTICE: Every effort has been made to obtain the maximum compatibility of function, strength, durability and comfort with this product, however, there is no guarantee that injury will be prevented with its use.