Medix Biochemica

Product Manual Cat. No: #5000

DirectBlood Genotyping PCR Kit

Description

The kit contains all components necessary for rapid, sensitive, and reproducible real-time PCR detection of SNPs from EDTA blood samples without previous DNA isolation. Only target specific primers and probes have to be added.

The lyophilized 2x DirectBlood Genotyping PCR Mix includes an engineered, hot-start formulated DNA polymerase, optimized buffer components and ultrapure dNTPs. Additionally, the kit consists of a Rehydration Buffer for rehydration of the LyoCake.

Kit components

Component	S pack	M pack
DirectBlood Genotyping PCR Kit	1 x bag (containing 4 x LyoCakes)	5 x bags (each containing 4 x LyoCakes)
Rehydration Buffer	1 x 1.25 mL	5 x 1.25 mL

*Other pack sizes, bulk orders and customization are available upon request.

Storage and shipment

Transport with cool packs. The reagents should be stored at -20°C upon arrival. The reagents are stable until the expiration date if stored correctly. Protect all components from light.

Preparations before use

1. Rehydrate the Lyocake by adding exactly **218 \muL** of the included **Rehydration Buffer** onto the Lyocake, resulting in 250 μ L of ready-to-use 2x Master Mix.

2. Subsequently invert the closed tube a few times, briefly vortex and spin down the mixture before use. The tube should be placed on ice after rehydration.

3. The rehydrated 2x DirectBlood Genotyping PCR Mix is then ready to be used or stored at -20°C.

Recommended sample preparation

1. Dilute a given EDTA blood sample to 2% (v/v) by adding 10 μL of blood sample to 490 μL nuclease-free water.

2. Subsequently invert the closed tube a few times or briefly vortex the mixture before use. Samples can also be stored for a couple of weeks at 4°C, e.g., for re-testing purposes.

Optional step:

3. Heat the diluted EDTA blood sample to 80° C for ≥ 5 minutes to reduce potential infectivity of the blood sample.

4. Spin down the sample for 1–2 seconds but **do not centrifuge after this step anymore**.

Reaction Master Mix set-up

The recommended master mix set-up for a 20 μ L reaction volume is shown in the table below.

Reagent	Volume (µL)	Final concentration
2x DirectBlood Genotyping PCR Mix	10	1x
4x Primer/Probe Mix	5	1x
Diluted blood sample (2%) or control sample	5	0.5%

Please note, we recommend including one positive control and one negative control reaction in each PCR run.



myPOLS Biotec GmbH Byk-Gulden-Strasse 2 78467 Konstanz, Germany VAT no. DE294545185

www.medixbiochemica.com www.mypols.de mdx@medixbiochemica.com

Legal disclaimer

myPOLS Biotec part of Medix Biochemica group products meet their specifications if transported, stored, and used according to the instructions. Medix Biochemica's products may not be used or reproduced without Medix Biochemica's written permission.

Medix Biochemica

Product Manual Cat. No: #5000

Instrument and program set-up

Cycles	Steps	Temperature	Time
1	Initial denaturation	95°C	2 min
50	Denaturation	95°C	10 sec
	Annealing / Extension	60°C	40 sec

Technical information and support

This product is also available with already established SNP-assays. Please contact us for more information.

For technical enquiries or assay development support, please contact us via e-mail at: mdx@medixbiochemica.com.

Additional information and technical resources are available on our website at: info.medixbiochemica.com/resources.



myPOLS Biotec GmbH Byk-Gulden-Strasse 2 78467 Konstanz, Germany VAT no. DE294545185 www.medixbiochemica.com www.mypols.de mdx@medixbiochemica.com

Legal disclaimer

myPOLS Biotec part of Medix Biochemica group products meet their specifications if transported, stored, and used according to the instructions. Medix Biochemica's products may not be used or reproduced without Medix Biochemica's written permission.