

Stationary Massage Table

Assembly Instructions



Congratulations!

Congratulations on your purchase of a Master Stationary Massage Table. For years to come, you can enjoy relaxing massages on a table you can count on. With the utmost respect for our environment, our tables are made from the finest materials and workmanship, and will help enhance your healing touch.



Parts List

- a. Bolt screws 8mm x85mm- (8)
- b. Bolt screws 8mm x60mm- (4)
- c. Spring washer- (12)
- d. Soft-touch knobs- (8)
- e. Hexagon Wrench T30
- f. Bottom legs-(4)
- g. Top legs-(4)
- h. Upholstered table top
- i. Wooden accessory shelf
- j. Arm shelf
- k. Face cradle
- l. Face pillow
- m. Arm rests- (2)
- n. Neck bolster
- o. Half round bolster
- p. Spa music
- q. Disposable face pillow covers




1 Lay the upholstered table top (h) on the floor with the upholstery side down.

2 Leave the Top (g) and Bottom legs(f) screwed together. Before installing the four legs, make sure to match the shape of each leg with the pads under the upholstery table top (h).

3 To install the four legs, place the bolt (a) with a spring washer (c) from inside the table top towards the table leg. Tighten the bolt with a wrench (e). 2 bolts (a) are placed at each corner.

4 Place wooden accessory shelf (i) between the upholstered table top (h) with flat surface facing downward toward the underside frame. Line up the holds on the corner of the wooden shelf (i). Insert the bolt (b) with a spring washer (c) from inside the shelf toward the table leg. Tighten with wrench (e). One bolt is placed at each corner.



5 To adjust the height of the upholstered table top (h), unscrew the soft touch knobs (d), line up the desired holes with the 2 protruding screws from the top table leg (g), and carefully slide the bottom leg (f) onto the top leg.

6 Once in place, fasten securely with supplied Soft-Touch Knobs (d). Repeat this process for the other legs, and make sure that the protruding screws enter the same holes on all 4 legs so that all the table legs stand at equal heights.

7 Once this process has been completed carefully turn table over so all four legs are on the floor.

8 Attach supplied arm shelf (j) to face cradle support (k). Position the face cradle poles to slide into the 2 slots at the head of upholstered table top (h). Fasten the face pillow (l) onto face cradle (k) by lining up male to female Velcro fasteners.