Step 'N Rest Rolling Walker · A223

INDICATIONS FOR USE

This device is intended to provide ambulatory assistance for end users with moderate stability and balance issues. The operation of the device in sitting mode is intended to provide a resting surface during breaks in ambulation.

SAVE THESE INSTRUCTIONS. This sheet must accompany device with each use.

CAUTION! Read all instructions and warnings prior to assembly and use.

- If you do not understand these instructions, contact a healthcare professional before attempting to use.
- Inspect all parts for shipping damages before use. If there is shipping damage – **DO NOT** use. Contact the manufacturer for further instructions.
- If there are any problems with your rollator, **DO NOT** attempt to repair it yourself. Contact the manufacturer for any needed parts or repair instructions.
- Tires and brakes are non-durable components which will require maintenance through the life of the rollator.
- Brakes should be adjusted so that the wheels do not move while the brake is applied – if brakes are unable to be adjusted correctly, **DO NOT** continue to use the rollator, seek repair or replacement.
- Tires should be inspected regularly for damage and wear and replaced when the tires are excessively worn. Tires may wear faster if the device is used outdoors on rough surfaces.
- Compass Health Brands assumes no responsibilities for any damage or injury caused by improper installation, assembly, or use of this product.

INSTRUCTIONS FOR USE AND WARNINGS

- This rollator has a weight limit of 250 lbs. (113.4 kg). **ALWAYS** observe the weight limit on the labeling of your rollator.
- **ALWAYS INSPECT ROLLATOR PRIOR TO USE.** Properly inspecting and maintaining your rollator can reduce the risk of injury. Make sure that all parts are secure and are in good working order. **ALWAYS** check the fasteners, wheels, brakes and seat for tightness prior to use. Make sure that parking brakes are working, and the rollator wheels cannot turn with the brakes on. If rollator is not in good working order or fasteners cannot be tightened, **DISCONTINUE USE OF THE ROLLATOR.**
- All wheels and handles **MUST** be adjusted to the same height for safe use. **DO NOT** use the rollator if wheels or handles are out of alignment, or if hardware is not secure. As applicable to your rollator model, all hand knobs must be firmly tightened and height adjustment buttons must protrude fully through the height adjustment holes with an audible “click”.
- **DO NOT** use the rollator as a wheelchair.
- **ALWAYS** engage both parking brakes before sitting in the rollator. **DO NOT** move or scoot while seated in the rollator. Parking brakes **MUST** remain engaged at all times while seated.
- **ALWAYS** wear protective footwear while using the rollator. **ALWAYS** keep shoelaces, loose clothing and/or accessories clear of the wheels and handle bars during use.
- **ALWAYS** maintain good posture while using the rollator. When walking with the rollator, keep your feet in line with the rear wheels, and **DO NOT** let the rollator get too far in front of you.
- When sitting on the rollator, keep both feet in contact with the floor and **DO NOT** lean out of the seat.
- **DO NOT** use the rollator to pull yourself up from a chair or other seated position. Use the arms of the chair to support your weight as you stand.
- **DO NOT** attempt to reach objects that are out of your immediate reach while seated on the rollator. **DO NOT** lean forward, backward or to either side while seated on the rollator.
- **DO NOT** use the backrest bar to support your weight. Leaning back while seated on the rollator may cause the device to tip.
- **ALWAYS** have all wheels in contact with the floor/ground while using this device. Failure to maintain wheel contact could result in improper balance of the rollator and cause injury or damage.
- **DO NOT** sit on the rollator if it is parked on a slope. Only use the seat on level ground.
- **ALWAYS** make sure the rollator is fully opened and the folding bar is locked in the down position before use.
- **DO NOT** hang anything from the rollator handles or frame. Items may only be carried in the provided pouch/basket.
- **DO NOT** climb, kneel, or stand on the rollator.
- This is not a toy. **DO NOT** let children play on rollator.
- **BE AWARE** of your surroundings when operating the rollator. Look for hazards and avoid them. Avoid streets and surfaces with water, sand, gravel, dirt, leaves, and other debris. Wet, slick, uneven and/or rough surfaces may impair traction and contribute to possible accidents. Loose cords and unsecured rugs may move suddenly and cause a loss of balance. **DO NOT** use a rollator in mud, ice, or puddles.
- **NEVER** use near steps, sloped driveways/hills with a steep incline (never greater than 15 degrees). **DO NOT** walk backwards down a slope.
- **DO NOT** use stairs or escalators while using the rollator. **DO NOT** use near swimming pool areas.
- **ALWAYS** take additional care and move slowly when moving from carpeted to hard surface floors or cracks/separations in walkway surfaces; avoid sharp bumps, drainage grates, and sudden surface changes. Rollator may abruptly stop if a wheel becomes wedged.
- **DO NOT** use the hand brake as a primary stopping tool. It may be used if needed but is primarily designed to maintain a stopped stance.
- **DO NOT** use while under the influence of mind-altering drugs, alcohol or while suffering from dizziness.
- **DO NOT** store in direct sunlight or outdoors. Exposure to extreme temperatures (heat, cold, moisture, etc.) may make surfaces (seat, hand grips, etc.) hot to the touch, weak or brittle. Ensure hand grips **DO NOT** twist and skin contact surfaces are cool to the touch before use.
- **FAILURE TO USE COMMON SENSE AND HEED THE ABOVE WARNINGS FURTHER INCREASES RISK OF SERIOUS INJURY. USE AT YOUR OWN RISK AND WITH APPROPRIATE AND SERIOUS ATTENTION TO SAFE OPERATION. USE CAUTION.**

CLEANING & MAINTENANCE

- Wipe rollator frame with a clean, soft cloth at least once a week. All metal parts can be cleaned with auto wax or similar products. If exposed to moisture, dry immediately.
- Wipe hand grips with cold water and a clean towel.
- **NEVER** use ammonia-based cleaning solutions or soap on hand grips. This may cause them to loosen or slide off.
- **NEVER** use organic solvents, abrasive cleaners or scouring pads on any part of the rollator.
- Check all the parts and movable connections before each use to ensure that the rollator is safe to use.
- Check wheels/casters for rips or excessive wear before use.

COMPASSHEALTH 42-A22300_02 © 2019
Manufactured for / Fabricado por: Compass Health Brands Corp.
Middleburg Heights, OH 44130 • www.compasshealthbrands.com
Phone / Teléfono 800.376.7263 • Made in China / Hecho en China

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BACK

COMPASSHEALTH

CONTACT: 781.871.2710
FILE NAME: A223_Warning_Insert_02
FACTORY: Heibei HealthPlus
DATE: March 15, 2019

PRINTING INFORMATION

MATERIAL: N/A
FINISH: N/A
SIZE INCHES: 8.5"(L) x 11"(H)
SIZE MM: 215.9(L) x 279.4(H)

COLOR LEGEND: 4 Color

 Cyan
 Magenta
 Yellow
 Black

ARTWORK STATUS

READER ONLY
 FINAL ARTWORK